

LIVE THE GOOD LIGHT™



2017 PRODUCT AVAILABILITY GUIDE

LENS COMPARISON CHART	Transitions® Signature™	Transitions® Vantage™	Transitions® XTRActive™	Transitions® DRIVEWEAR™
Lens Description	Fully clear indoors, optimal tint outdoors	Variable polarization for even crisper, sharper vision outdoors	More activation for added protection	Polarized and adjust with changing sunlight to enhance the driving experience
Lens Type	For your primary prescription eyeglasses	For your primary prescription eyeglasses	For your primary prescription eyeglasses	For your prescription or non-prescription sunglasses
Usage	Indoors and outdoors, all day	Indoors and outdoors, all day	Indoors and outdoors, all day	Outdoors and in the car, daytime only
Lens color options	Clear to Gray Clear to Brown Clear to Graphite Green	Clear to Gray	Clear to Gray Clear to Brown Clear to Graphite Green	Olive Green to Copper to Dark Reddish Brown
Indoor clarity	Fully clear	Hint of tint	Hint of tint	Sunwear, not recommended for indoor use
Outdoor darkness	Dark	Dark	Our darkest lens, even in high temperatures	Dark
Speed to darken	Fast to activate	Fast to activate	Fast to activate	Fast to activate
Behind-the-windshield activation	No	No	Yes, moderate activation	Yes, fully activates
Polarization	No	Variable polarization that adjusts to match outdoor glare levels	No	Yes
Blocks 100% of UVA and UVB Light	Yes	Yes	Yes	Yes



Reinforce the value of your patients' purchase with the Transitions® lenses Certificate of Authenticity!

Transitions, the swirl and XTRActive are registered trademarks, Transitions adaptive lenses, the eclipse, and Vantage, are trademarks of Transitions Optical, Inc. ©2017 Transitions Optical, Inc. Drivewear is a registered trademark of Younger Mfg. Co.

Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.



	Transitions Signature			Transitions XTRActive			Transitions VANTAGE	Sphere	Cyl to
	Gray	Brown	Graphite Green	Gray	Brown	Graphite Green	Gray		
1.50 STANDARD INDEX FINISHED SINGLE VISION									
CARL ZEISS VISION									
	FSV	GR						-4.00 to +4.00	-2.00
ESSILOR LENS									
	Essilor Lens FSV Crizal Avancé UV*	GR						-4.00 to +4.00	-2.00
	Essilor Lens FSV Crizal Alizé UV*	GR						-4.00 to +4.00	-2.00
	Essilor Lens FSV	GR	BR	GG	XA-GR			-4.00 to +4.00	-2.00
HOYA VISION CARE									
	HOYA Vision Care FSV	GR						-4.00 to +4.00	-2.00
SEIKO OPTICAL PRODUCTS									
	Seiko Ultra 1.50 Spherical FSV ECP AR	GR	BR					-4.00 to +3.00	-2.00
	Seiko Ultra 1.50 Spherical FSV	GR	BR					-4.00 to +3.00	-2.00
SIGNET ARMORLITE									
	Signet Armorlite Kodak FSV 1.50 Clean&CleAR*	GR						-4.00 to +4.00	-2.00
YOUNGER OPTICS									
	Younger Optics FSV AR Super Hydrophobic	GR	BR					-4.00 to +4.00	-2.00
	Younger Optics FSV	GR	BR		XA-GR			-4.00 to +4.00	-2.00
1.53 TRIVEX® FINISHED SINGLE VISION									
HOYA VISION CARE									
	HOYA Phoenix® FSV	GR						-6.00 to +2.00	-2.00
	HOYA Phoenix® FSV with Super HiVision EX3 AR	GR						-6.00 to +2.00	-2.00
X-CEL OPTICAL									
	X-Cel Optical Aris® FSV	GR						-6.00 to Plano	-2.00
	X-Cel Optical Aris® FSV AR	GR						-6.00 to Plano	-2.00
1.59 POLYCARBONATE FINISHED SINGLE VISION									
CARL ZEISS VISION									
	FSV Aspheric Polycarbonate	GR						-4.00 to +2.00	-2.00
ESSILOR LENS									
	Essilor FSV Spherical Airwear® Crizal Avancé UV*	GR						-6.00 to +4.00	-2.00
	Essilor FSV Spherical Airwear® Crizal Alizé UV*	GR						-6.00 to +4.00	-2.00
	Essilor FSV Spherical Airwear®	GR	BR	GG				-6.00 to +4.00	-2.00
	Essilor FSV Aspheric Airwear®	GR						-6.00 to +4.00	-2.00
	Essilor FSV Aspheric Airwear® Crizal Avance UV*			GG				-6.00 to +4.00	-1.00
SEIKO OPTICAL PRODUCTS									
	Seiko DiamondClear Aspheric FSV	GR	BR					-6.00 to +3.00	-2.00
	Seiko DiamondClear ECP AR Aspheric FSV	GR	BR					-6.00 to +3.00	-2.00
SIGNET ARMORLITE									
	Signet Armorlite Signia® FSV Spherical Poly	GR						-6.00 to +4.00	-2.00
	Signet Armorlite Kodak Spherical FSV Poly Clean&CleAR*	GR						-6.00 to +4.00	-2.00
YOUNGER OPTICS									
	Younger FSV Poly AR Super Hydrophobic	GR	BR					-6.00 to +4.00	-2.00
	Younger Optics FSV Poly	GR	BR		XA-GR			-6.00 to +4.00	-2.00
1.67 SUPER HIGH INDEX FINISHED SINGLE VISION									
ESSILOR LENS									
	Essilor FSV Aspheric Thin&Lite® 1.67 Crizal Avancé UV*	GR						-8.00 to +0.00	-2.00
SEIKO OPTICAL PRODUCTS									
	Seiko 1.67 Super SV Aspheric	GR						-8.00 to -2.00	-2.00
	Seiko 1.67 Super SV ECP AR Aspheric	GR	BR					-8.00 to -2.00	-2.00
SIGNET ARMORLITE									
	Signet Armorlite Signia® FSV Aspherical 1.67 KODAK Clean&CleAR*	GR						-8.00 to -2.00	-2.00
YOUNGER OPTICS									
	Younger Optics FSV 1.67 AR Super Hydrophobic	GR	BR					-2.00 to -8.00	-2.00
1.74 ULTRA HIGH INDEX FINISHED SINGLE VISION									
SEIKO OPTICAL PRODUCTS									
	Seiko 1.74 Aspheric FSV AR	GR						-10.00 to -3.00	-2.00



	Transitions Signature			Transitions XTRActive			Transitions VANTAGE	Sphere	Cyl to
	Gray	Brown	Graphite Green	Gray	Brown	Graphite Green	Gray		
1.50 STANDARD INDEX									
CARL ZEISS VISION									
SOLA SFSV	GR	BR						-6.00 to +5.00	-4.00
synchrony Access HD (75, 125 & 175)	GR	BR					V	-9.00 to +10.00	-6.00
synchrony SV HDC	GR	BR					V	-9.00 to +10.00	-6.00
synchrony SV HDC Curves	GR	BR						-3.50 to +6.00	-4.00
ZEISS 3D Single Vision	GR	BR					V	-6.00 to +6.00	-4.00
ZEISS DriveSafe Individual SV	GR	BR						-6.00 to +6.00	-4.00
ZEISS DriveSafe SV	GR	BR						-6.00 to +6.00	-4.00
ZEISS EnergizeMe SV	GR	BR						-7.00 to +5.00	-4.00
ZEISS Individual® SV	GR	BR					V	-6.00 to +6.00	-4.00
ZEISS Individual® SV Sport	GR	BR						-4.00 to +4.00	-4.00
ZEISS Individual® SV Wrap	GR	BR						-4.00 to +4.00	-4.00
ZEISS Spherical Single Vision	GR	BR					V	-9.00 to +7.00	-4.00
ZEISS SV DriveSafe Individual Sport	GR	BR						-4.00 to +4.00	-4.00
ZEISS SV DriveSafe Sport	GR	BR						-4.00 to +4.00	-4.00
ESSILOR LENS									
Essilor Lens SFSV	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-10.50 to +10.00	-4.00
Essilor Stylistic Wrap® Single Vision	GR	BR	GG	XA-GR	XA-BR	XA-GG		-3.00 to +3.00	-1.00
Eyezen®+	GR	BR	GG	XA-GR	XA-BR		V	-10.00 to +6.00	-6.00
HOYA VISION CARE									
HOYA SFSV 1.50	GR	BR						-10.00 to +9.00	-4.00
HOYA SFSV CR-39	GR	BR						-10.00 to +7.00	-4.00
HOYA SFSV 1.50				XA-GR				-8.00 to +9.00	-4.00
HOYA SV iQ 1.50	GR	BR		XA-GR				-8.00 to +9.00	-4.00
Sync BKS	GR	BR		XA-GR			V	-10.00 to +3.75	-4.00
Sync Wrap BKS	GR	BR		XA-GR			V	-5.00 to +4.00	-4.00
NIKON OPTICAL									
Nikon Infinite Single Vision 1.50	GR	BR	GG					-10.00 to + 8.00	-6.00
SEIKO OPTICAL PRODUCTS									
Seiko PL 1.50 Transitions	GR	BR					V	-8.00 to +6.00	-5.00
Seiko 1.50 Super Atoric SV	GR	BR					V	-5.00 to +4.00	-5.00
Seiko 1.50 Superior SV	GR	BR					V	-5.00 to +7.00	-5.00
SHAMIR INSIGHT INC.									
Shamir Attitude III Single Vision	GR	BR		XA-GR			V	-6.50 to +6.00	-6.00
Shamir Autograph® II - Single Vision®	GR	BR		XA-GR			V	-14.00 to +8.00	-6.00
Shamir Autograph® II - Single Vision Attitude®	GR	BR		XA-GR			V	-6.00 to +6.00	-6.00
Shamir Spectrum® Single Vision	GR	BR		XA-GR			V	-14.00 to +8.00	-6.00
Shamir Relax® Single Vision	GR	BR		XA-GR			V	-12.25 to +6.00	-6.00
Shamir Autograph® III Single Vision	GR	BR		XA-GR			V	-12.00 to +7.25	-6.00
Shamir SFSV 1.5	GR	BR					V	-10.00 to +7.50	-6.00
SIGNET ARMORLITE									
Signet Armorlite SFSV	GR	BR					V	-6.00 to +10.00	-4.00
Signetek Processed SV	GR	BR	GG	XA-GR			V	-8.00 to +6.00	-6.00
KODAK Digital SV	GR	BR	GG	XA-GR			V	-8.00 to +6.00	-6.00
VISION EASE LENS									
Vision Ease SFDS	GR	BR						-11.00 to +6.00	-4.00
YOUNGER OPTICS									
Younger Optics SFSV (72mm)	GR	BR						-9.00 to +6.00	-4.00
Younger Optics SFSV DIGITAL	GR	BR						-10.00 to +8.50	-4.00
Younger Optics SFSV (76mm)				XA-GR	XA-BR		V	-11.00 to +6.00	-4.00



	Transitions Signature			Transitions XTRActive			Transitions	Sphere	Cyl to
	Gray	Brown	Graphite Green	Gray	Brown	Graphite Green	VANTAGE Gray		
1.53 TRIVEX® SINGLE VISION									
CARL ZEISS VISION									
synchrony Access HD (75, 125 & 175)	GR	BR						-11.00 to +8.00	-4.00
synchrony SV HDC	GR	BR						-11.00 to +8.00	-4.00
synchrony SV HDC Curves	GR	BR						-3.50 to +4.00	-4.00
ZEISS 3D Single Vision	GR	BR						-7.00 to +7.00	-4.00
ZEISS DriveSafe Individual SV	GR	BR						-7.00 to +7.00	-4.00
ZEISS DriveSafe SV	GR	BR						-7.00 to +7.00	-4.00
ZEISS EnergizeMe SV	GR	BR						-7.00 to +5.00	-4.00
ZEISS Individual® SV	GR	BR						-7.00 to +7.00	-4.00
ZEISS Individual® SV Sport	GR	BR						-4.00 to +4.00	-4.00
ZEISS Individual® SV Wrap	GR	BR						-4.00 to +4.00	-4.00
ZEISS Spherical Single Vision	GR	BR						-7.00 to +7.00	-4.00
ZEISS SV DriveSafe Individual Sport	GR	BR						-4.00 to +4.00	-4.00
ZEISS SV DriveSafe Sport	GR	BR						-4.00 to +4.00	-4.00
ESSILOR LENS									
Essilor Lens SFSV TREXA™	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	-4.00
Essilor eyecode™ Single Vision TREXA™	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	-6.00
Essilor Fit™ Single Vision TREXA™	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	-6.00
Essilor Azio™ Single Vision TREXA™	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	-6.00
Essilor 360™ Single Vision TREXA™	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	-6.00
Eyezen™+ TREXA™	GR	BR						-10.00 to +6.00	-6.00
HOYA VISION CARE									
HOYA SFSV Phoenix™	GR	BR						-13.00 to +9.00	-4.00
HOYA SFSV Phoenix™				XA-GR	XA-BR			-10.00 to +6.00	-4.00
HOYA SFSV Phoenix™							V	-8.00 to +9.00	-4.00
HOYA SV iQ Phoenix™	GR	BR						-13.00 to +9.00	-4.00
HOYA SV iQ Phoenix™				XA-GR	XA-BR			-10.00 to +6.00	-4.00
HOYA SV iQ Phoenix™							V	-8.00 to +9.00	-4.00
Sync BKS	GR	BR		XA-GR	XA-BR		V	-10.00 to +6.00	-4.00
Sync Wrap BKS	GR	BR		XA-GR	XA-BR		V	-5.00 to +4.00	-4.00
SEIKO OPTICAL PRODUCTS									
Seiko PL 1.53/Trivex® Transitions	GR	BR		XA-GR			V	-10.00 to +7.00	-5.00
Seiko 1.53/Trivex® Super Atoric SV	GR	BR		XA-GR			V	-7.00 to +5.00	-5.00
Seiko 1.53/Trivex® Superior SV	GR	BR		XA-GR			V	-7.00 to +8.00	-5.00
SHAMIR INSIGHT INC.									
Shamir Attitude III Single Vision	GR	BR		XA-GR			V	-6.50 to +6.00	-6.50
Shamir Autograph® II - Single Vision®	GR	BR		XA-GR			V	-14.00 to +7.50	-6.50
Shamir Autograph® II - Single Vision Attitude™	GR	BR		XA-GR			V	-6.00 to +6.00	-6.50
Shamir Spectrum™ Single Vision	GR	BR		XA-GR			V	-14.00 to +7.50	-6.50
Shamir Relax™ Single Vision	GR	BR		XA-GR			V	-13.25 to +6.00	-6.50
Shamir Golf™ Single Vision™	GR	BR		XA-GR			V	-5.00 to +4.00	-6.50
Shamir Autograph® III Single Vision	GR	BR		XA-GR			V	-14.00 to +8.00	-6.50
Shamir SFSV DLC™ Trivex™	GR	BR		XA-GR			V	-14.00 to +7.50	-6.50
SIGNET ARMORLITE									
Signetek Processed SV	GR	BR		XA-GR			V	-10.00 to +6.00	-6.00
KODAK Digital SV	GR	BR		XA-GR			V	-10.00 to +6.00	-6.00
X-CEL OPTICAL									
X-Cel SFSV	GR	BR	GG					-10.00 to +8.50	-4.00
X-Cel SFSV				XA-GR	XA-BR		V	-9.00 to +6.50	-4.00
YOUNGER OPTICS									
Younger Optics Trilogy® Aspheric SV	GR	BR		XA-GR	XA-BR			-9.00 to +6.00	-4.00
Younger Optics Trilogy® Spherical SV	GR	BR		XA-GR	XA-BR			-8.00 to +7.50	-4.00
Younger Optics Trilogy® Spherical SV							V	-8.25 to +7.50	-4.00



	Transitions Signature			Transitions XTRActive			Transitions	Sphere	Cyl to
	Gray	Brown	Graphite Green	Gray	Brown	Graphite Green	VANTAGE Gray		
1.55 MID INDEX SINGLE VISION									
X-CEL OPTICAL									
1.55 Mid Index SFSV	GR							-11.00 to +9.00	-4.00
1.56 TRIVEX® SINGLE VISION									
NIKON OPTICAL									
Nikon Infinite Single Vision TRIVEX®	GR	BR	GG					-11.00 to + 8.00	-6.00
1.59 POLYCARBONATE SINGLE VISION									
CARL ZEISS VISION									
SOLA ASL® Aspheric SFSV	GR	BR		XA-GR				-10.00 to +6.00	-4.00
synchrony Access HD (75, 125 & 175)	GR	BR		XA-GR			V	-10.00 to +8.00	-6.00
synchrony SV HDC	GR	BR		XA-GR			V	-10.00 to +8.00	-6.00
synchrony SV HDC Curves	GR	BR						-3.50 to +4.00	-4.00
ZEISS 3D Single Vision	GR	BR		XA-GR			V	-10.00 to +8.00	-4.00
ZEISS DriveSafe Individual SV	GR	BR						-10.00 to +8.00	-4.00
ZEISS DriveSafe SV	GR	BR						-10.00 to +8.00	-4.00
ZEISS EnergizeMe SV	GR	BR						-10.00 to +6.00	-4.00
ZEISS Individual® SV	GR	BR		XA-GR			V	-10.00 to +8.00	-4.00
ZEISS Individual® SV Sport	GR	BR						-4.00 to +4.00	-4.00
ZEISS Individual® SV Wrap	GR	BR		XA-GR				-4.00 to +4.00	-4.00
ZEISS Spherical Single Vision	GR	BR					V	-12.00 to +6.50	-4.00
ZEISS SV DriveSafe Individual Sport	GR	BR						-4.00 to +4.00	-4.00
ZEISS SV DriveSafe Sport	GR	BR						-4.00 to +4.00	-4.00
ESSILOR LENS									
Essilor SFSV Spherical Airwear®	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-10.00 to +8.00	-4.00
Essilor SFSV Aspheric Airwear®	GR	BR	GG	XA-GR		XA-GG		-10.00 to +8.00	-4.00
Essilor SV 360° Aspheric Airwear®	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-10.00 to +8.00	-6.00
Essilor Eyecode® Aspheric SV Airwear®	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-10.00 to +8.00	-6.00
Essilor FIT® Aspheric SV Airwear®	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-10.00 to +8.00	-6.00
Essilor AZIO® Aspheric SV Airwear®	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-10.00 to +8.00	-6.00
Essilor Stylistic Wrap® Single Vision Airwear®	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-4.00 to +4.00	-4.00
Eyezen®+ Airwear®	GR	BR	GG	XA-GR	XA-BR		V	-10.00 to +8.00	-6.00
HOYA VISION CARE									
HOYA SV iQ Polycarbonate	GR	BR		XA-GR			V	-6.00 to +6.00	-4.00
Sync BKS	GR	BR		XA-GR			V	-13.00 to +5.25	-4.00
Sync Wrap BKS	GR	BR		XA-GR			V	-6.00 to +4.00	-4.00
NIKON OPTICAL									
Nikon Infinite Single Vision Polycarbonate	GR	BR	GG					-10.00 to +8.00	-6.00
SEIKO OPTICAL PRODUCTS									
Seiko PL 1.59 Transitions	GR	BR		XA-GR			V	-10.00 to +7.00	-5.00
Seiko Spectaclear SFSV	GR	BR						-12.00 to +6.00	-5.00
Seiko 1.59 Super Atoric SV	GR	BR		XA-GR			V	-8.00 to +5.00	-5.00
Seiko 1.59 Superior SV	GR	BR		XA-GR			V	-7.00 to +8.00	-5.00
SHAMIR INSIGHT INC.									
Shamir Attitude III Single Vision	GR	BR		XA-GR			V	-6.50 to +6.00	-7.00
Shamir Autograph® II Single Vision® Poly	GR	BR		XA-GR			V	-14.00 to +8.75	-7.00
Shamir Autograph® II Single Vision Attitude® Poly	GR	BR		XA-GR			V	-5.50 to +6.00	-7.00
Shamir Spectrum® SFSV PolyPlus®	GR	BR		XA-GR			V	-14.00 to +8.75	-7.00
Shamir Relax® SFSV PolyPlus®	GR	BR		XA-GR			V	-14.00 to +8.00	-7.00
Shamir SFSV Poly	GR	BR		XA-GR			V	-11.25 to +9.00	-7.00
Shamir Golf® Single Vision®	GR	BR		XA-GR			V	-6.00 to +5.25	-7.00
Shamir Autograph® III Single Vision	GR	BR		XA-GR			V	-14.00 to +8.75	-7.00
SIGNET ARMORLITE									
Signet Armorlite Signia® Poly SFSV	GR	BR						-10.00 to +8.00	-4.00
Signetek Processed SV	GR	BR	GG	XA-GR			V	-10.00 to +6.50	-6.00
KODAK Digital SV	GR	BR	GG	XA-GR			V	-10.00 to +6.50	-6.00
VISION EASE LENS									
Vision Ease Polycarbonate SFSV 78mm	GR	BR						-11.00 to +8.00	-6.00
YOUNGER OPTICS									
Younger Optics SV	GR	BR		XA-GR	XA-BR		V	-9.00 to +6.00	-4.00



	Transitions Signature			Transitions XTRActive			Transitions	Sphere	Cyl to
	Gray	Brown	Graphite Green	Gray	Brown	Graphite Green	VANTAGE Gray		
1.60 HIGH INDEX SINGLE VISION									
ESSILOR LENS									
Essilor SFSV Spherical Thin&Lite® 1.60	GR	BR						-13.00 to +8.00	-4.00
Essilor SFSV Aspheric Thin&Lite® 1.60	GR	BR						-13.00 to +10.00	-4.00
HOYA VISION CARE									
HOYA SFSV 1.60	GR	BR						-10.00 to +8.00	-4.00
HOYA Aspheric SV 1.60	GR	BR						-10.00 to +8.00	-4.00
HOYA SV iQ 1.60	GR	BR						-10.00 to +8.00	-4.00
Sync BKS	GR	BR						-10.00 to +6.00	-4.00
NIKON OPTICAL									
Nikon Infinite Single Vision 1.60	GR	BR	GG					-12.00 to + 8.00	-6.00
Nikon SeeMax Single Vision 1.60	GR	BR	GG					-10.00 to + 8.00	-4.00
SEIKO OPTICAL PRODUCTS									
Seiko PL 1.60 Transitions	GR	BR						-12.00 to +10.00	-5.00
Seiko 1.60 Super Atoric SV	GR	BR						-12.00 to +10.00	-5.00
Seiko 1.60 Superior SV	GR	BR						-8.50 to +10.00	-5.00
SHAMIR INSIGHT INC.									
Shamir Attitude III Single Vision	GR	BR		XA-GR				-6.50 to +6.00	-7.00
Shamir Autograph® II Single Vision® SuperLite® 1.6	GR	BR		XA-GR				-14.00 to +8.00	-7.00
Shamir Golf® Single Vision®	GR	BR		XA-GR				-6.00 to +5.25	-7.00
Shamir Autograph® II Single Vision Attitude® SuperLite® 1.6	GR	BR		XA-GR				-6.50 to +5.25	-7.00
Shamir Spectrum® SFSV Superlite® 1.6	GR	BR		XA-GR				-14.00 to +8.00	-7.00
Shamir Relax® SFSV SuperLite® 1.6	GR	BR		XA-GR				-15.00 to +9.00	-7.00
Shamir SFSV SuperLite® 1.6	GR	BR		XA-GR				-14.00 to +8.00	-7.00
Shamir Autograph® III Single Vision	GR	BR		XA-GR				-14.00 to +8.75	-7.00
SIGNET ARMORLITE									
Signet Armoritite Signia® 1.60 SFSV	GR	BR						-13.00 to +8.00	-4.00
KODAK Digital SV	GR	BR						-10.00 to +7.50	-6.00
YOUNGER OPTICS									
Younger Optics Single Vision 1.60	GR	BR						-11.00 to +8.00	-4.00
1.67 SUPER HIGH INDEX SINGLE VISION									
CARL ZEISS VISION									
SOLA ASL 1.67 Aspheric SFSV	GR	BR		XA-GR				-12.00 to +8.00	-6.00
synchrony Access HD (75, 125 & 175)	GR	BR		XA-GR				-12.00 to +8.00	-6.00
synchrony SV HDC	GR	BR		XA-GR				-12.00 to +8.00	-6.00
ZEISS 3D Single Vision	GR	BR		XA-GR				-12.00 to +8.00	-4.00
ZEISS DriveSafe Individual SV	GR	BR						-12.00 to +8.00	-6.00
ZEISS DriveSafe SV	GR	BR						-12.00 to +8.00	-6.00
ZEISS EnergizeMe SV	GR	BR						-12.00 to +8.00	-4.00
ZEISS Individual® SV	GR	BR		XA-GR				-12.00 to +8.00	-4.00
ESSILOR LENS									
Essilor SFSV Spherical Thin&Lite® 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-13.00 to +9.00	-4.00
Essilor SFSV Aspheric Thin&Lite® 1.67	GR	BR	GG					-16.00 to +8.00	-4.00
Essilor SV 360° Spherical Thin&Lite® 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-12.00 to +9.00	-6.00
Essilor Stylistic Wrap® Single Vision Thin&Lite® 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-5.00 to +5.00	-3.00
Essilor Eyecode® SV Thin&Lite® 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-12.00 to +9.00	-6.00
Essilor Fit® SV Thin&Lite® 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-12.00 to +9.00	-6.00
Essilor Azio® SV Thin&Lite® 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-12.00 to +9.00	-6.00
Eyezen®+ Thin&Lite® 1.67	GR	BR	GG	XA-GR	XA-BR			-14.00 to +9.00	-6.00



	Transitions Signature			Transitions XTRActive			Transitions	Sphere	Cyl to
	Gray	Brown	Graphite Green	Gray	Brown	Graphite Green	VANTAGE Gray		
1.67 SUPER HIGH INDEX SINGLE VISION (CONTINUED)									
HOYA VISION CARE									
HOYA SFSV 1.67	GR	BR						-10.00 to +8.00	-4.00
HOYA Aspheric SV 1.67	GR	BR						-15.00 to +10.00	-4.00
HOYA SV iQ 1.67	GR	BR		XA-GR				-10.00 to +8.00	-4.00
Sync BKS	GR	BR		XA-GR				-13.00 to +7.50	-4.00
Sync Wrap BKS	GR	BR		XA-GR				-6.00 to +4.00	-4.00
NIKON OPTICAL									
Nikon Aspheric SV 1.67	GR	BR	GG					-12.00 to + 8.00	-4.00
Nikon Infinite Single Vision 1.67	GR	BR	GG					-12.00 to + 8.00	-6.00
SEIKO OPTICAL PRODUCTS									
Seiko PL 1.67 Transitions	GR	BR		XA-GR				-12.00 to +10.00	-5.00
Seiko 1.67 Super SV Aspheric	GR	BR						-15.00 to +10.00	-6.00
Seiko 1.67 Super Atoric SV	GR	BR		XA-GR				-15.00 to +10.00	-5.00
Seiko 1.67 Superior SV	GR	BR		XA-GR				-12.00 to +10.00	-5.00
SHAMIR INSIGHT INC.									
Shamir Attitude III Single Vision	GR	BR		XA-GR				-6.50 to +6.00	-8.00
Shamir Autograph® II Single Vision® SuperLite® 1.67	GR	BR		XA-GR				-16.00 to +16.00	-8.00
Shamir Autograph® II Single Vision Attitude® SuperLite® 1.67	GR	BR		XA-GR				-6.00 to +6.00	-8.00
Shamir Spectrum® SFSV Superlite® 1.67	GR	BR		XA-GR				-16.00 to +16.00	-8.00
Shamir SFSV Superlite® 1.67	GR	BR		XA-GR				-16.00 to +16.25	-8.00
Shamir Relax®	GR	BR		XA-GR				-16.00 to +8.25	-8.00
Shamir Autograph® III Single Vision	GR	BR		XA-GR				-16.00 to +16.00	-8.00
SIGNET ARMORLITE									
Signet Armorlite Signia® 1.67 SFSV	GR	BR						-13.00 to +9.00	-4.00
Signetek Processed SV	GR	BR	GG	XA-GR				-11.00 to +7.50	-6.00
KODAK Digital SV	GR	BR	GG	XA-GR				-11.00 to +7.50	-6.00
YOUNGER OPTICS									
Younger Optics Single Vision 1.67 MR-10™	GR	BR		XA-GR	XA-BR			-11.00 to +7.00	-4.00
1.70 ULTRA HIGH INDEX SINGLE VISION									
VISION EASE LENS									
Vision Ease SFSV Thindex 1.70	GR	BR						-12.00 to +8.00	-4.00
Vision Ease ASFSV Thindex 1.70	GR	BR						-16.00 to +10.00	-4.00
1.74 ULTRA HIGH INDEX SINGLE VISION									
ESSILOR LENS									
Essilor SFSV Spherical Thin&Lite® 1.74	GR							-20.00 to +13.00	-4.00
Essilor SFSV Aspheric Thin&Lite® 1.74	GR							-20.00 to +13.00	-4.00
Essilor Eyecode® SV Thin&Lite® 1.74	GR							-16.00 to +13.00	-6.00
Essilor Fit® SV Thin&Lite® 1.74	GR							-16.00 to +13.00	-6.00
Essilor 360® SV Thin&Lite® 1.74	GR							-16.00 to +13.00	-6.00
Eyezen®+ Thin&Lite® 1.74	GR							-20.00 to +13.00	-6.00
NIKON OPTICAL									
Nikon Aspheric SV 1.74	GR	BR	GG					-12.00 to + 4.00	-4.00
Nikon Infinite Single Vision 1.74	GR	BR	GG					-16.00 to + 10.00	-6.00
SEIKO OPTICAL PRODUCTS									
Seiko PL 1.74 Transitions	GR	BR						-17.00 to +10.00	-5.00
Seiko 1.74 Superior SV	GR	BR						-12.50 to +10.00	-5.00
Seiko 1.74 MX	GR	BR						-15.00 to +10.00	-4.00
SHAMIR INSIGHT INC.									
Shamir Spectrum® Single Vision	GR	BR						-16.00 to +12.50	-9.00
Shamir Autograph® II Single Vision®	GR	BR						-20.00 to +15.50	-9.00
Shamir Autograph® III Single Vision	GR	BR						-20.00 to +13.50	-9.00
Shamir SFSV Superlite® 1.74	GR	BR						-18.00 to +8.50	-9.00
SIGNET ARMORLITE									
Signetek® Processed SV	GR	BR						-13.00 to +9.00	-6.00
KODAK Digital SV	GR	BR						-13.00 to +9.00	-6.00
YOUNGER OPTICS									
Younger Optics Single Vision 1.74	GR	BR						-18.00 to +12.00	-4.00

MULTIFOCAL

	Transitions Signature Graphite Green			Transitions XTRActive Graphite Green			Transitions VANTAGE Gray	Sphere	Add	Cyl to
	Gray	Brown		Gray	Brown					
1.50 STANDARD INDEX SFST-24										
X-CEL OPTICAL										
	X-Cel RD-24	GR	BR					-4.00 to +5.50	+0.75 to +4.00	-4.00
1.50 STANDARD INDEX SFST-28										
CARL ZEISS VISION										
	SOLA SFFT 28	GR	BR					-9.00 to +5.00	+1.00 to +3.00	-4.00
ESSILOR LENS										
	Essilor Lenses FT 28	GR	BR	GG	XA-GR			-10.00 to +6.00	+1.00 to +3.00	-4.00
	Essilor FT 28 Bifocal AB	GR	BR	GG				-10.00 to +6.00	+1.00 to +3.00	-4.00
HOYA VISION CARE										
	HOYA Vision Care SFST 28	GR	BR		XA-GR			-11.00 to +6.00	+1.00 to +3.00	-4.00
	HOYA Vision Care SFST 28 iQ	GR	BR					-8.00 to +6.00	+1.00 to +3.00	-4.00
	HOYA Vision Care SFST 28 iQ				XA-GR			-10.00 to +6.00	+1.00 to +3.00	-4.00
SIGNET ARMORLITE										
	Signet Armorlite SFFT-28	GR						-10.00 to +6.00	+1.00 to +3.00	-4.00
	Signetek® Processed Signia® 1.50 FT-28	GR						-10.00 to +6.00	+1.00 to +3.00	-6.00
X-CEL OPTICAL										
	X-Cel D-28 76mm	GR	BR		XA-GR	XA-BR		-9.00 to +5.00	+1.00 to +3.00	-4.00
	X-Cel D-28 78mm	GR	BR					-4.00 to +4.25	+0.75 to +5.00	-4.00
	X-Cel D-28 78mm				XA-GR			-4.00 to +4.25	+0.75 to +4.00	-4.00
	X-Cel D-28 78mm					XA-BR		-4.00 to +4.25	+3.25 to +4.00	-4.00
	X-Cel D-28 72mm	GR						-4.50 to +7.00	+1.00 to +3.00	-4.00
YOUNGER OPTICS										
	Younger Optics SFFT 28	GR	BR		XA-GR	XA-BR		-11.00 to +6.00	+1.00 to +3.00	-4.00
	Younger Optics SFFT 28 HI-ADDS	GR	BR					-5.00 to +6.00	+3.25 to +4.00	-4.00
1.50 STANDARD INDEX SFST-35										
CARL ZEISS VISION										
	SOLA SFFT 35	GR						-4.00 to +6.00	+1.00 to +3.00	-4.00
X-CEL OPTICAL										
	X-Cel D-35	GR	BR					-7.00 to +5.50	+0.75 to +4.00	-4.00
YOUNGER OPTICS										
	Younger Optics SFFT 35	GR	BR					-5.00 to +6.00	+0.75 to +4.00	-4.00
1.50 STANDARD INDEX TRIFOCALS										
CARL ZEISS VISION										
	SOLA SF 7x28	GR						-4.00 to +6.00	+1.50 to +3.00	-4.00
ESSILOR LENS										
	Essilor Lenses SF 7x28	GR	BR	GG				-5.00 to +6.00	+1.50 to +3.00	-4.00
SIGNET ARMORLITE										
	Signet Armorlight Signia® 1.50 HC SF Trifocal 7x28	GR						-5.00 to +6.00	+1.50 to +3.00	-4.00
X-CEL OPTICAL										
	X-Cel 7x28	GR	BR					-7.00 to +7.50	+1.50 to +4.00	-4.00
	X-Cel 8x35	GR	BR					-7.00 to +5.00	+1.50 to +4.00	-4.00
YOUNGER OPTICS										
	Younger Optics 7x28	GR	BR		XA-GR			-6.00 to +5.00	+1.50 to +4.00	-4.00
	Younger Optics 8x35	GR	BR					-6.00 to +5.00	+1.50 to +4.00	-4.00
1.59 INDEX POLYCARBONATE COMPOSITE SFFT 28										
ESSILOR LENS										
	Bifocal FT28	GR						-10.00 to +6.00	+1.00 to +3.00	-4.00
YOUNGER OPTICS										
	Younger Optics Polycarbonate Composite SFFT 28	GR	BR					-9.00 to +7.00	+1.00 to +3.00	-4.00
	Younger Optics Polycarbonate Composite SFFT 35	GR	BR					-9.00 to +7.00	+1.00 to +3.50	-4.00
	Younger Optics Polycarbonate Composite SFFT 7X28	GR	BR					-6.00 to +5.00	+1.50 to +4.00	-4.00

MULTIFOCAL		Transitions®			Transitions® XTRActive®			Transitions® VANTAGE®	Sphere	Add	Cyl to
		Gray	Brown	Graphite Green	Gray	Brown	Graphite Green	Gray			
1.53 TRIVEX® SF RD-24											
X-CEL OPTICAL											
	1.53 Trivex® SF RD-24	GR							-8.00 to +5.50	+1.00 to +3.50	-4.00
1.53 TRIVEX® SFST-28											
HOYA VISION CARE											
	HOYA Phoenix™ SFST 28	GR	BR						-10.00 to +6.00	+1.00 to +3.50	-4.00
	HOYA Phoenix™ ST 28 iQ	GR	BR						-9.00 to +7.00	+1.00 to +3.50	-4.00
X-CEL OPTICAL											
	1.53 Trivex® SFST-28	GR	BR						-10.00 to +5.50	+1.00 to +4.00	-4.00
YOUNGER OPTICS											
	Younger Optics Trilogy® FLAT TOP 28	GR	BR						-9.00 to +7.00	+1.00 to +3.50	-4.00
1.53 TRIVEX® DD-28											
X-CEL OPTICAL											
	1.53 Trivex® DD-28	GR							-4.50 to +2.50	1.00 to +3.00	-4.00
1.53 TRIVEX® SFST-35											
X-CEL OPTICAL											
	1.53 Trivex® SFST-35	GR							-8.00 to +5.50	+1.00 to +4.00	-4.00
1.53 TRIVEX® TRIFOCALS											
X-CEL OPTICAL											
	1.53 Trivex® SFST-7X28	GR	BR						-8.00 to +5.50	+1.50 to +3.50	-4.00
	1.53 Trivex® SFST-8X35	GR							-8.00 to +5.50	+1.50 to +3.00	-4.00
1.55 MID INDEX SFST-28											
X-CEL OPTICAL											
	1.55 Mid Index SFST-28	GR							-11.00 to +9.00	+1.00 to +4.00	-4.00
1.55 MID INDEX SFST-35											
X-CEL OPTICAL											
	1.55 Mid Index SFST-35	GR							-7.00 to +5.50	+1.00 to +4.00	-4.00
1.55 MID INDEX TRIFOCALS											
X-CEL OPTICAL											
	1.55 Mid Index SFST-7X28	GR							-11.00 to +9.00	+1.50 to +4.00	-4.00
	1.55 Mid Index SFST-8X35	GR							-7.00 to +5.50	+1.50 to +3.00	-4.00

MULTIFOCAL		Transitions Signature®			Transitions® XTRActive®			Transitions® VANTAGE®	Sphere	Add	Cyl to	MFH*
		Gray	Brown	Graphite Green	Gray	Brown	Graphite Green	Gray				
1.50 STANDARD INDEX BLENDED BIFOCAL (DIGITAL)												
SHAMIR INSIGHT INC.												
	Shamir Duo	GR	BR					V	-9.00 to +6.00	+0.75 to +4.00	-6.00	15
1.53 TRIVEX® BLENDED BIFOCAL (DIGITAL)												
SHAMIR INSIGHT INC.												
	Shamir Duo	GR	BR					V	-13.25 to +6.00	+0.75 to +4.00	-6.50	15
1.59 POLYCARBONATE BLENDED BIFOCAL (DIGITAL)												
SHAMIR INSIGHT INC.												
	Shamir Duo	GR	BR					V	-14.75 to +6.00	+0.75 to +4.00	-7.00	15
1.60 HIGH INDEX BLENDED BIFOCAL (DIGITAL)												
SHAMIR INSIGHT INC.												
	Shamir Duo	GR	BR						-15.00 to +6.00	+0.75 to +4.00	-7.00	15
1.67 SUPER HIGH INDEX BLENDED BIFOCAL (DIGITAL)												
SHAMIR INSIGHT INC.												
	Shamir Duo	GR	BR						-16.75 to +8.00	+0.75 to +4.00	-8.00	15
1.74 ULTRA HIGH INDEX BLENDED BIFOCAL (DIGITAL)												
SHAMIR INSIGHT INC.												
	Shamir Duo	GR	BR						-18.75 to +12.00	+0.75 to +4.00	-9.00	15

DIGITAL	Transitions Signature Graphite Green			Transitions XTRActive Graphite Green			Transitions VANTAGE Gray	Sphere	Add	Cyl to	MFH*
	Gray	Brown	Graphite Green	Gray	Brown	Graphite Green	Gray				
1.50 STANDARD INDEX (DIGITAL)											
YOUNGER OPTICS											
Younger Optics Camber™	GR	BR						-11.25 to +6.25		-4.00	
1.53 TRIVEX® (DIGITAL)											
YOUNGER OPTICS											
Younger Optics Camber™	GR	BR						-11.00 to +6.25		-4.00	
1.59 POLYCARBONATE (DIGITAL)											
YOUNGER OPTICS											
Younger Optics Camber™	GR	BR						-10.00 to +6.25		-4.00	
1.60 HIGH INDEX (DIGITAL)											
YOUNGER OPTICS											
Younger Optics Camber™	GR	BR						-11.50 to +7.25		-4.00	
1.67 SUPER HIGH INDEX (DIGITAL)											
YOUNGER OPTICS											
Younger Optics Camber™	GR	BR						-14.00 to +8.00		-4.00	
1.74 ULTRA HIGH INDEX (DIGITAL)											
YOUNGER OPTICS											
Younger Optics Camber™	GR	BR						-19.00 to +9.00		-4.00	

PROGRESSIVES	Transitions Signature Graphite Green			Transitions XTRActive Graphite Green			Transitions VANTAGE Gray	Sphere	Add	Cyl to	MFH*
	Gray	Brown	Graphite Green	Gray	Brown	Graphite Green	Gray				
1.50 STANDARD INDEX											
CARL ZEISS VISION											
AO Easy	GR	BR						-9.00 to +5.00	+0.75 to +3.50	-4.00	18
AO Easy HD	GR	BR					V	-9.00 to +4.50	+0.75 to +3.50	-4.00	18
SOLA VIP	GR	BR						-7.00 to +6.00	+1.00 to +3.00	-4.00	18
SOLAMAX®	GR	BR						-8.00 to +4.00	+1.00 to +3.00	-4.00	17
synchrony Curves HDC	GR	BR						-3.50 to +4.00	+0.75 to +3.50	-4.00	13
synchrony Easy Adapt HD	GR	BR					V	-10.00 to +6.00	+0.75 to +3.50	-6.00	18
synchrony Easy Adapt (One)	GR	BR						-9.00 to +6.00	+1.00 to +3.50	-4.00	18
synchrony Easy M HD	GR	BR					V	-10.00 to +6.00	+0.75 to +3.50	-6.00	15
synchrony Easy M (Compact)	GR							-9.00 to +5.00	+1.00 to +3.50	-4.00	15
synchrony Easy S HD	GR	BR					V	-10.00 to +6.00	+0.75 to +3.50	-6.00	13
synchrony Easy S (Ultra)	GR	BR						-9.00 to +5.00	+0.75 to +3.50	-6.00	13
synchrony Easy View HD	GR	BR					V	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
synchrony Easy View (Instictive)	GR	BR						-9.00 to +5.00	+1.00 to +3.50	-4.00	17
synchrony Easy Wear HD	GR	BR					V	-10.00 to +6.00	+0.75 to +3.50	-6.00	18
synchrony Performance HD13	GR	BR					V	-10.00 to +6.00	+0.75 to +3.50	-6.00	13
synchrony Performance HD15	GR	BR					V	-10.00 to +6.00	+0.75 to +3.50	-6.00	15
synchrony Performance HD17	GR	BR					V	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
synchrony Performance HD19	GR	BR					V	-10.00 to +6.00	+0.75 to +3.50	-6.00	19
synchrony Performance HDV	GR	BR					V	-10.00 to +6.00	+0.75 to +3.50	-6.00	13
synchrony Starter HD	GR	BR					V	-10.00 to +6.00	+0.50 to +2.00	-6.00	14
synchrony Ultra HD 13	GR	BR					V	-10.00 to +6.00	+0.75 to +3.50	-6.00	13
synchrony Ultra HD 15	GR	BR					V	-10.00 to +6.00	+0.75 to +3.50	-6.00	15
synchrony Ultra HD 17	GR	BR					V	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
synchrony Ultra HD 19	GR	BR					V	-10.00 to +6.00	+0.75 to +3.50	-6.00	19
synchrony Ultra HD 21	GR	BR					V	-10.00 to +6.00	+0.75 to +3.50	-6.00	21
synchrony Ultra HDC	GR	BR					V	-10.00 to +6.00	+0.75 to +3.50	-6.00	13
synchrony Ultra HDV	GR	BR					V	-10.00 to +6.00	+0.75 to +3.50	-6.00	13
synchrony Work and Go HD	GR	BR					V	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
synchrony Work and Office HD	GR	BR					V	-10.00 to +6.00	+0.75 to +3.50	-6.00	16
synchrony Work and Read HD	GR	BR					V	-10.00 to +6.00	+0.75 to +3.00	-6.00	16
Synchrony™	GR							-9.00 to +6.00	+0.75 to +3.50	-4.00	17
ZEISS Choice 13	GR	BR					V	-7.00 to +5.00	+0.75 to +3.50	-4.00	13
ZEISS Choice 15	GR	BR					V	-7.00 to +5.00	+0.75 to +3.50	-4.00	15
ZEISS Choice 17	GR	BR					V	-7.00 to +5.00	+0.75 to +3.50	-4.00	17
ZEISS Choice 19	GR	BR					V	-7.00 to +5.00	+0.75 to +3.50	-4.00	19

PROGRESSIVES

1.50 STANDARD INDEX (CONTINUED)

CARL ZEISS VISION

	Transitions Signature			Transitions XTRActive			Transitions VANTAGE	Sphere	Add	Cyl to	MFH
	Gray	Brown	Graphite Green	Gray	Brown	Graphite Green					
ZEISS Choice Plus Sport 13	GR	BR						-4.00 to +4.00	+0.75 to +3.50	-4.00	13
ZEISS Choice Plus Sport 17	GR	BR						-4.00 to +4.00	+0.75 to +3.50	-4.00	17
ZEISS Digital 1000	GR	BR					V	-7.00 to +5.00	+1.00	-4.00	13
ZEISS Digital 1250	GR	BR					V	-7.00 to +5.00	+1.25	-4.00	13
ZEISS Digital 500	GR	BR					V	-7.00 to +5.00	+0.50	-4.00	13
ZEISS Digital 750	GR	BR					V	-7.00 to +5.00	+0.75	-4.00	13
ZEISS Digital Lens	GR	BR					V	-7.00 to +5.00	+0.50 to +1.25	-4.00	13
ZEISS Digital Sport	GR	BR						-4.00 to +4.00	+0.50 to +1.25	-4.00	13
ZEISS Digital Sport 1000	GR	BR						-4.00 to +4.00	+1.00	-4.00	13
ZEISS Digital Sport 1250	GR	BR						-4.00 to +4.00	+1.25	-4.00	13
ZEISS Digital Sport 500	GR	BR						-4.00 to +4.00	+0.50	-4.00	13
ZEISS Digital Sport 750	GR	BR						-4.00 to +4.00	+0.75	-4.00	13
ZEISS DriveSafe Individual	GR	BR						-7.00 to +5.00	+0.75 to +3.50	-4.00	13
ZEISS DriveSafe Individual Sport	GR	BR						-4.00 to +4.00	+0.75 to +3.50	-4.00	13
ZEISS DriveSafe	GR	BR						-7.00 to +5.00	+0.75 to +3.50	-4.00	13
ZEISS DriveSafe Sport	GR	BR						-4.00 to +4.00	+0.75 to +3.50	-4.00	13
ZEISS GT2 [®]	GR	BR						-10.00 to +6.00	+1.00 to +3.00	-6.00	17
ZEISS GT2 [®] Short	GR	BR						-10.00 to +6.00	+1.00 to +3.00	-6.00	13
ZEISS Individual [®] 2	GR	BR					V	-7.00 to +5.00	+0.75 to +3.50	-4.00	13
ZEISS Individual [®] 2I	GR	BR					V	-7.00 to +5.00	+0.75 to +3.50	-4.00	13
ZEISS Individual [®] 2N	GR	BR					V	-7.00 to +5.00	+0.75 to +3.50	-4.00	13
ZEISS Individual [®] Sport	GR	BR						-4.00 to +4.00	+0.75 to +3.50	-4.00	13
ZEISS Individual [®] Wrap	GR	BR						-4.00 to +4.00	+0.75 to +3.50	-4.00	13
ZEISS Plus 2 Short Sport	GR	BR						-4.00 to +4.00	+0.75 to +3.50	-4.00	14
ZEISS Plus 2 Sport	GR	BR						-4.00 to +4.00	+0.75 to +3.50	-4.00	14
ZEISS Precision Plus	GR	BR					V	-7.00 to +5.00	+0.75 to +3.50	-4.00	13
ZEISS Precision Pure 13	GR	BR					V	-7.00 to +5.00	+0.75 to +3.50	-4.00	13
ZEISS Precision Pure 15	GR	BR					V	-7.00 to +5.00	+0.75 to +3.50	-4.00	15
ZEISS Precision Pure 17	GR	BR					V	-7.00 to +5.00	+0.75 to +3.50	-4.00	17
ZEISS Precision Pure 19	GR	BR					V	-7.00 to +5.00	+0.75 to +3.50	-4.00	19
ZEISS Precision Pure 21	GR	BR					V	-7.00 to +5.00	+0.75 to +3.50	-4.00	21
ZEISS Precision Superb	GR	BR					V	-7.00 to +5.00	+0.75 to +3.50	-4.00	13

ESSILOR LENS

Varilux S 4D [™]	GR	BR		XA-GR	XA-BR			-10.00 to +6.00	+0.75 to +4.00	-6.00	14**
Varilux S Fit [™]	GR	BR		XA-GR	XA-BR			-10.00 to +6.00	+0.75 to +4.00	-6.00	14**
Varilux S Design [™]	GR	BR		XA-GR	XA-BR			-10.00 to +6.00	+0.75 to +4.00	-6.00	14**
Varilux [®] Physio [®] W3+ eyecode	GR	BR		XA-GR				-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux [®] Physio [®] W3+ FIT	GR	BR		XA-GR				-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux [®] Physio [®] W3+	GR	BR		XA-GR				-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort [®] W2+ eyecode	GR	BR		XA-GR				-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort [®] W2+ FIT	GR	BR		XA-GR				-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort [®] W2+	GR	BR		XA-GR				-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio Enhanced eyecode [™]	GR	BR						-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio Enhanced FIT [™]	GR	BR						-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio Enhanced [™]	GR	BR						-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio DRx [™]	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx [™]	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-10.00 to +6.00	+0.75 to +3.00	-6.00	14
Varilux [®] Physio [®]	GR	BR						-10.00 to +6.00	+0.75 to +3.50 & +4.00	-4.00	17
Varilux Physio Short [™]	GR							-10.00 to +6.00	+0.75 to +3.00	-4.00	14
Varilux Comfort Enhanced FIT [™]	GR	BR		XA-GR				-10.00 to +7.00	+0.75 to +3.50	-6.00	14
Varilux Comfort Enhanced [™]	GR	BR		XA-GR				-10.00 to +7.00	+0.75 to +3.50	-6.00	14
Varilux Comfort DRx [™]	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx [™]	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.00	-6.00	14
Varilux Comfort [®]	GR	BR		XA-GR				-10.00 to +7.00	+0.75 to +3.50	-4.00	17
Varilux Comfort Short [™]	GR	BR		XA-GR				-10.00 to +7.00	+0.75 to +3.00	-4.00	14
Varilux [®] Stylistic [®] Wrap	GR	BR	GG	XA-GR	XA-BR	XA-GG		-3.00 to +3.00	+0.75 to +3.50	-1.00	17
DEFINITY [®] 3 PLUS	GR	BR						-10.00 to +6.00	+1.00 to +3.50	-6.00	15
DEFINITY [®] 3	GR	BR						-10.00 to +6.00	+1.00 to +3.50	-6.00	15
DEFINITY [®]	GR							-10.00 to +6.00	+1.00 to +3.00	-4.00	18
DEFINITY SHORT [™]	GR							-10.00 to +6.00	+1.00 to +3.00	-4.00	15
Essilor Accolade Freedom [™]	GR	BR						-9.00 to +6.00	+0.75 to +3.50	-4.00	15

**For Varilux S, at fitting height less than 17mm, ADD power is up to +3.50.

Transitions lenses are compatible with all Crizal No-Glare technology

TRANSITIONS EVERYDAY LENSES

*MFH-Minimum Fitting Height

PROGRESSIVES

1.50 STANDARD INDEX (CONTINUED)

	Transitions Signature			Transitions XTRActive			Transitions VANTAGE	Sphere	Add	Cyl to	MFH
	Gray	Brown	Graphite Green	Gray	Brown	Graphite Green					
ESSILOR											
Essilor Accolade®	GR	BR						-9.00 to +6.00	+0.75 to +3.50	-4.00	17
Essilor Ideal Advanced®	GR	BR	GG					-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced Wrap®	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-4.00 to +3.00	+0.75 to +3.50	-4.00	17
Essilor Ideal®	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Fit	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-10.00 to +6.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Computer	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-10.00 to +6.00	+0.75 to +3.50	-6.00	13
Essilor Ideal® Driving	GR	BR	GG	XA-GR	XA-BR		V	-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Ovation® Digital	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Ovation®	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-4.00	17
Essilor Natural® Digital	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Essilor Natural®	GR							-10.00 to +6.00	+0.75 to +3.50	-4.00	18
Adaptar® Digital	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.00	-6.00	14
Essilor SmallFit® Digital	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.00	-6.00	14
Essilor SmallFit®	GR							-10.00 to +6.00	+0.75 to +3.00	-4.00	14
HOYALUX VISION CARE											
HOYALUX Summit CD® 1.50	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-4.00	14
HOYALUX Summit CD® 1.50				XA-GR				-10.00 to +6.00	+1.00 to +3.00	-4.00	14
HOYALUX Summit CD® iQ 1.50	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-6.00	14
HOYALUX Summit CD® iQ 1.50				XA-GR				-10.00 to +6.00	+1.00 to +3.00	-6.00	14
HOYALUX Summit ECP® 1.50	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-4.00	18
HOYALUX Summit ECP® 1.50				XA-GR				-10.00 to +6.00	+1.00 to +3.00	-4.00	18
HOYALUX Summit ECP® iQ 1.50	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-6.00	18
HOYALUX Summit ECP® iQ 1.50				XA-GR				-10.00 to +6.00	+1.00 to +3.00	-6.00	18
HOYALUX iD LifeStyle®2 Clarity 1.50	GR	BR						-8.00 to +6.00	+0.75 to +3.50	-4.00	18
HOYALUX iD LifeStyle®2 Clarity 1.50				XA-GR				-8.00 to +6.00	+1.00 to +3.50	-4.00	18
HOYALUX iD LifeStyle®2 Harmony 1.50	GR	BR						-8.00 to +6.00	+0.75 to +3.50	-4.00	18
HOYALUX iD LifeStyle®2 Harmony 1.50				XA-GR				-8.00 to +6.00	+1.00 to +3.50	-4.00	18
HOYALUX iD LifeStyle®2 Clarity CD 1.50	GR	BR						-8.00 to +6.00	+0.75 to +3.50	-4.00	14
HOYALUX iD LifeStyle®2 Clarity CD 1.50				XA-GR				-8.00 to +6.00	+1.00 to +3.50	-4.00	14
HOYALUX iD LifeStyle®2 Harmony CD 1.50	GR	BR						-8.00 to +6.00	+0.75 to +3.50	-4.00	14
HOYALUX iD LifeStyle®2 Harmony CD 1.50				XA-GR				-8.00 to +6.00	+1.00 to +3.50	-4.00	14
HOYALUX ARRAY® 11, 13, 15, 17, VL 1.50	GR	BR		XA-GR			V	-10.00 to +3.75	+0.75 to +3.50	-4.00	14, 17, 19, 21
HOYALUX ARRAY® Wrap 11, 13, 15, 17, VL 1.50	GR	BR		XA-GR			V	-5.00 to +4.00	+0.75 to +3.50	-4.00	14, 17, 19, 21
Amplitude 1.50	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-4.00	19
Amplitude iQ 1.50	GR							-10.00 to +6.00	+0.75 to +3.50	-4.00	19
Amplitude HiDef 1.50	GR	BR						-8.00 to +6.00	+0.75 to +3.50	-4.00	18
Amplitude Mini 1.50	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-4.00	14
Amplitude Mini iQ 1.50	GR							-10.00 to +6.00	+0.75 to +3.50	-4.00	14
Amplitude Mini HiDef 1.50	GR	BR						-8.00 to +6.00	+0.75 to +3.50	-4.00	14
Amplitude BKS 1.50	GR	BR		XA-GR			V	-10.00 to +3.75	+0.75 to +3.50	-4.00	19
Amplitude Mini BKS 1.50	GR	BR		XA-GR			V	-10.00 to +3.75	+0.75 to +3.50	-4.00	14
HOYALUX GP Wide	GR	BR						-10.00 to +3.75	+0.75 to +3.50	-4.00	18
NIKON OPTICAL											
Nikon SeeMax Power 1.50	GR	BR	GG					-10.00 to + 8.00	+0.75 to 4.00	-6.00	13, 15, 17
Nikon Presio Power 1.50	GR	BR	GG	XA-GR			V	-10.00 to + 8.00	+0.75 to 4.00	-6.00	13, 15, 17
Nikon DigiLife 1.50	GR	BR	GG	XA-GR			V	-10.00 to + 8.00	+0.75 to 4.00	-6.00	13, 15, 17
Nikon Presio i Digital 1.50	GR	BR	GG	XA-GR			V	-7.00 to + 4.00	+0.75 to 4.00	-6.00	15, 17
Nikon Move 1.50	GR	BR	GG	XA-GR			V	-7.00 to + 4.00	+0.75 to 3.50	-4.00	16, 18

PROGRESSIVES

	Transitions Signature			Transitions XTRActive			Transitions VANTAGE	Sphere	Add	Cyl to	MFH
	Gray	Brown	Graphite Green	Gray	Brown	Graphite Green					
1.50 STANDARD INDEX (CONTINUED)											
SEIKO OPTICAL PRODUCTS											
Seiko 1.50 Superior	GR	BR					V	-5.00 to +4.00	+0.50 to +4.00	-5.00	12
Seiko 1.50 Prominent	GR	BR					V	-5.00 to +4.00	+0.50 to +3.50	-5.00	14
Seiko 1.50 Premier 2	GR	BR					V	-5.00 to +4.00	+0.50 to +3.50	-4.00	14
Seiko 1.50 Prestige II	GR	BR					V	-5.00 to +4.00	+0.50 to +3.50	-5.00	14
Seiko 1.50 Prestige II Wide	GR	BR					V	-5.00 to +4.00	+0.50 to +3.50	-5.00	14
Seiko 1.50 Succeed	GR	BR					V	-5.00 to +4.00	+0.50 to +3.50	-4.00	17
Seiko 1.50 Succeed Ws	GR	BR					V	-5.00 to +4.00	+0.50 to +3.50	-4.00	15
Seiko 1.50 Supercede 2	GR	BR					V	-5.00 to +4.00	+0.50 to +3.50	-4.00	14
Seiko 1.50 Supernal	GR	BR					V	-5.00 to +4.00	+0.50 to +3.50	-5.00	14
Seiko 1.50 Surmount	GR	BR					V	-5.00 to +4.00	+0.50 to +3.50	-5.00	14
Seiko 1.50 Surmount Ws	GR	BR					V	-5.00 to +4.00	+0.50 to +3.50	-5.00	14
Seiko 1.50 PCWide	GR	BR						-3.00 to +5.00	-1.00, -1.50, -2.00	-5.00	15
SHAMIR INSIGHT INC.											
Shamir Autograph® III	GR	BR	GG	XA-GR			V	-12.25 to +6.00	+0.75 to +4.00	-6.00	11, 13, 15, 18
Shamir Autograph® II	GR	BR	GG	XA-GR			V	-12.25 to +6.00	+0.75 to +4.00	-6.00	11, 13, 15, 18
Shamir Autograph® II Attitude®	GR	BR	GG	XA-GR			V	-5.00 to +4.00	+0.75 to +3.50	-6.00	15, 18
Shamir Spectrum™	GR	BR	GG	XA-GR			V	-12.25 to +6.00	+0.75 to +3.50	-6.00	14, 16, 18
Shamir Spectrum™+	GR	BR	GG	XA-GR			V	-12.25 to +6.00	+0.75 to +3.50	-6.00	14, 16, 18
Shamir Element™	GR	BR	GG				V	-7.00 to +4.00	+0.75 to +3.50	-6.00	16, 19
Shamir FirstPal	GR	BR	GG	XA-GR			V	-12.25 to +6.00	+0.75 to +1.50	-6.00	15
Shamir Computer™	GR	BR	GG	XA-GR			V	-6.00 to +8.00	+0.75 to +2.00	-6.00	14 Lower, 10 Upper
Shamir Workspace™	GR	BR	GG	XA-GR			V	-6.00 to +8.00	+0.75 to +2.50	-6.00	14 Lower, 10 Upper
Shamir InTouch™	GR	BR	GG	XA-GR			V	-12.25 to +6.00	+0.75 to +3.50	-6.00	15, 18
Shamir Attitude III Sport	GR	BR	GG	XA-GR			V	-6.50 to +5.25	+0.75 to +3.50	-6.00	18
Shamir Attitude III Fashion	GR	BR	GG	XA-GR			V	-6.50 to +5.25	+0.75 to +3.50	-6.00	15, 18
SIGNET ARMORLITE											
DirecTek™ Signia™ 1.50	GR	BR	GG					-8.00 to +4.00	+0.75 to +3.50	-6.00	16.5
DirecTek™ Short Signia™ 1.50	GR	BR	GG					-8.00 to +4.00	+0.75 to +3.50	-6.00	13.5
Navigator® FBS Signia™ 1.50	GR	BR	GG					-8.00 to +4.00	+0.75 to +3.50	-6.00	18
Navigator® Short FBS Signia™ 1.50	GR	BR	GG					-8.00 to +4.00	+0.75 to +3.50	-6.00	14
KODAK Precise™ Signia™ 1.50 SF	GR							-9.00 to +5.00	+1.00 to +3.00	-4.00	17
KODAK Precise™ Short Signia™ 1.50 SF	GR							-8.50 to +5.25	+1.00 to +3.00	-4.00	13
KODAK Precise™ Digital Signia™ 1.50	GR							-8.00 to +5.00	+1.00 to +3.00	-6.00	17
KODAK Precise™ Short Digital Signia™ 1.50	GR							-8.00 to +5.00	+1.00 to +3.00	-6.00	13
KODAK Precise™ PB Signia™ 1.50	GR	BR	GG				V	-8.00 to +4.00	+0.75 to +3.50	-6.00	17
KODAK Precise™ Short PB Signia™ 1.50	GR	BR	GG				V	-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique™ Signia™ 1.50	GR	BR	GG	XA-GR			V	-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique™ HD Signia™ 1.50	GR	BR	GG	XA-GR			V	-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique DSII™ Signia™ 1.50	GR	BR						-8.00 to +4.00	+0.75 to +3.50	-6.00	13
X-CEL OPTICAL											
X-Cel Freedom ID™	GR							-4.50 to +5.50	+1.00 to +3.50	-4.00	17
YOUNGER OPTICS											
Younger Image®	GR	BR		XA-GR				-8.00 to +6.00	+1.00 to +3.00	-4.00	18

Transitions lenses are compatible with all Crizal No-Glare technology

*MFH-Minimum Fitting Height
† The lenses only have near and intermediate zones

PROGRESSIVES

1.53 TRIVEX®

CARL ZEISS VISION

	Transitions Signature		Transitions XTRActive			Transitions VANTAGE	Sphere	Add	Cyl to	MFH
	Gray	Brown	Graphite Green	Gray	Brown	Graphite Green				
AO Easy HD	GR	BR					-9.00 to +5.00	+0.75 to +3.50	-4.00	18
synchrony Curves HDC	GR	BR					-4.00 to +4.00	+0.75 to +3.50	-4.00	13
synchrony Easy Adapt HD	GR	BR					-9.00 to +5.00	+0.75 to +3.50	-6.00	18
synchrony Easy M HD	GR	BR					-9.00 to +5.00	+0.75 to +3.50	-6.00	15
synchrony Easy S HD	GR	BR					-9.00 to +5.00	+0.75 to +3.50	-6.00	13
synchrony Easy View HD	GR	BR					-9.00 to +5.00	+0.75 to +3.50	-6.00	17
synchrony Easy Wear HD	GR	BR					-9.00 to +5.00	+0.75 to +3.50	-6.00	18
synchrony Performance HD13	GR	BR					-9.00 to +5.00	+0.75 to +3.50	-6.00	13
synchrony Performance HD15	GR	BR					-9.00 to +5.00	+0.75 to +3.50	-6.00	15
synchrony Performance HD17	GR	BR					-9.00 to +5.00	+0.75 to +3.50	-6.00	17
synchrony Performance HD19	GR	BR					-9.00 to +5.00	+0.75 to +3.50	-6.00	19
synchrony Performance HDV	GR	BR					-9.00 to +5.00	+0.75 to +3.50	-6.00	13
synchrony Starter HD	GR	BR					-9.00 to +5.00	+0.50 to +2.00	-6.00	14
synchrony Ultra HD 13	GR	BR					-9.00 to +5.00	+0.75 to +3.50	-6.00	13
synchrony Ultra HD 15	GR	BR					-9.00 to +5.00	+0.75 to +3.50	-6.00	15
synchrony Ultra HD 17	GR	BR					-9.00 to +5.00	+0.75 to +3.50	-6.00	17
synchrony Ultra HD 19	GR	BR					-9.00 to +5.00	+0.75 to +3.50	-6.00	19
synchrony Ultra HD 21	GR	BR					-9.00 to +5.00	+0.75 to +3.50	-6.00	21
synchrony Ultra HDC	GR	BR					-9.00 to +5.00	+0.75 to +3.50	-6.00	13
synchrony Ultra HDV	GR	BR					-9.00 to +5.00	+0.75 to +3.50	-6.00	13
synchrony Work and Go HD	GR	BR					-9.00 to +5.00	+0.75 to +3.50	-6.00	17
synchrony Work and Office HD	GR	BR					-9.00 to +5.00	+0.75 to +3.50	-6.00	16
synchrony Work and Read HD	GR	BR					-9.00 to +5.00	+0.75 to +3.00	-6.00	16
ZEISS Choice 13	GR	BR					-7.00 to +5.00	+0.75 to +3.50	-4.00	13
ZEISS Choice 15	GR	BR					-7.00 to +5.00	+0.75 to +3.50	-4.00	15
ZEISS Choice 17	GR	BR					-7.00 to +5.00	+0.75 to +3.50	-4.00	17
ZEISS Choice 19	GR	BR					-7.00 to +5.00	+0.75 to +3.50	-4.00	19
ZEISS Choice Plus Sport 13	GR	BR					-4.00 to +4.00	+0.75 to +3.50	-4.00	13
ZEISS Choice Plus Sport 17	GR	BR					-4.00 to +4.00	+0.75 to +3.50	-4.00	17
ZEISS Digital 1000	GR	BR					-7.00 to +5.00	+1.00	-4.00	13
ZEISS Digital 1250	GR	BR					-7.00 to +5.00	+1.25	-4.00	13
ZEISS Digital 500	GR	BR					-7.00 to +5.00	+0.50	-4.00	13
ZEISS Digital 750	GR	BR					-7.00 to +5.00	+0.75	-4.00	13
ZEISS Digital Lens	GR	BR					-7.00 to +5.00	+0.50 to +1.25	-4.00	13
ZEISS Digital Sport	GR	BR					-4.00 to +4.00	+0.50 to +1.25	-4.00	13
ZEISS Digital Sport 1000	GR	BR					-4.00 to +4.00	+1.00	-4.00	13
ZEISS Digital Sport 1250	GR	BR					-4.00 to +4.00	+1.25	-4.00	13
ZEISS Digital Sport 500	GR	BR					-4.00 to +4.00	+0.50	-4.00	13
ZEISS Digital Sport 750	GR	BR					-4.00 to +4.00	+0.75	-4.00	13
ZEISS DriveSafe Individual	GR	BR					-7.00 to +5.00	+0.75 to +3.50	-4.00	13
ZEISS DriveSafe Individual Sport	GR	BR					-4.00 to +4.00	+0.75 to +3.50	-4.00	13
ZEISS DriveSafe	GR	BR					-7.00 to +5.00	+0.75 to +3.50	-4.00	13
ZEISS DriveSafe Sport	GR	BR					-4.00 to +4.00	+0.75 to +3.50	-4.00	13
ZEISS Individual® 2	GR	BR					-7.00 to +5.00	+0.75 to +3.50	-4.00	13
ZEISS Individual® 2I	GR	BR					-7.00 to +5.00	+0.75 to +3.50	-4.00	13
ZEISS Individual® 2N	GR	BR					-7.00 to +5.00	+0.75 to +3.50	-4.00	13
ZEISS Individual® Sport	GR	BR					-4.00 to +4.00	+0.75 to +3.50	-4.00	13
ZEISS Individual® Wrap	GR	BR					-4.00 to +4.00	+1.00 to +2.50	-4.00	13
ZEISS Plus 2 Short Sport	GR	BR					-4.00 to +4.00	+0.75 to +3.50	-4.00	14
ZEISS Plus 2 Sport	GR	BR					-4.00 to +4.00	+0.75 to +3.50	-4.00	18
ZEISS Precision Plus	GR	BR					-7.00 to +5.00	+0.75 to +3.50	-4.00	13
ZEISS Precision Pure 13	GR	BR					-7.00 to +5.00	+0.75 to +3.50	-4.00	13
ZEISS Precision Pure 15	GR	BR					-7.00 to +5.00	+0.75 to +3.50	-4.00	15
ZEISS Precision Pure 17	GR	BR					-7.00 to +5.00	+0.75 to +3.50	-4.00	17
ZEISS Precision Pure 19	GR	BR					-7.00 to +5.00	+0.75 to +3.50	-4.00	19
ZEISS Precision Pure 21	GR	BR					-7.00 to +5.00	+0.75 to +3.50	-4.00	21
ZEISS Precision Superb	GR	BR					-7.00 to +5.00	+0.75 to +3.50	-4.00	13

PROGRESSIVES

1.53 TRIVEX® (CONTINUED)

ESSILOR LENS

	Transitions Signature			Transitions XTRActive			Transitions VANTAGE	Sphere	Add	Cyl to	MFH
	Gray	Brown	Graphite Green	Gray	Brown	Graphite Green					
Varilux S 4D*	GR							-10.00 to +6.00	+0.75 to +4.00	-6.00	14**
Varilux S Fit*	GR							-10.00 to +6.00	+0.75 to +4.00	-6.00	14**
Varilux S Design*	GR							-10.00 to +6.00	+0.75 to +4.00	-6.00	14**
Varilux* Physio* W3+ eyecode TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux* Physio* W3+ FIT TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux* Physio* W3+ TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio DRx* TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Varilux Physio Short DRx* TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.00	-6.00	14
Varilux Physio Enhanced eyecode* TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio Enhanced FIT* TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio Enhanced* TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux* Physio* TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+1.00 to +3.50 & +4.00	-4.00	17
Varilux Comfort* W2+ eyecode	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort* W2+ FIT	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort* W2+	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort DRx* TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx* TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.00	-6.00	14
DEFINITY* 3 PLUS TREXA*	GR	BR						-10.00 to +6.00	+1.00 to +3.50	-6.00	15
DEFINITY* 3 TREXA*	GR	BR						-10.00 to +6.00	+1.00 to +3.50	-6.00	15
DEFINITY* TREXA*	GR							-10.00 to +6.00	+1.00 to +3.50 & +4.00	-4.00	18
DEFINITY Short* TREXA*	GR							-10.00 to +6.00	+1.00 to +3.50 & +4.00	-4.00	15
Essilor Accolade Freedom* TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.50	-4.00	15
Essilor Accolade* TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.50	-4.00	17
Essilor Ideal Advanced* TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Essilor Ideal* TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Essilor Ideal* Fit TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Essilor Ideal* Sport TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	18
Essilor Ideal* Computer TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	13
Essilor Ideal* Driving TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Ovation* Digital TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Essilor Natural* Digital TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Essilor SmallFit* Digital TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.00	-6.00	14
Adaptar* Digital TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Adaptar* Short Digital TREXA*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +3.00	-6.00	14

HOYALUX VISION CARE

HOYALUX Summit CD Phoenix*	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-4.00	14
HOYALUX Summit CD Phoenix*				XA-GR			V	-10.00 to +6.00	+1.00 to +3.00	-4.00	14
HOYALUX Summit CD iQ Phoenix*	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-6.00	14
HOYALUX Summit CD iQ Phoenix*				XA-GR			V	-10.00 to +6.00	+1.00 to +3.00	-6.00	14
HOYALUX Summit ECP Phoenix*	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-4.00	18
HOYALUX Summit ECP Phoenix*				XA-GR			V	-10.00 to +6.00	+1.00 to +3.00	-4.00	18
HOYALUX Summit ECP iQ Phoenix*	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-6.00	18
HOYALUX iD LifeStyle CD Phoenix*	GR	BR		XA-GR			V	-8.00 to +6.00	+0.75 to +3.50	-4.00	14
HOYALUX iD LifeStyle2 Clarity Phoenix*	GR	BR		XA-GR			V	-8.00 to +6.00	+0.75 to +3.50	-4.00	18
HOYALUX iD LifeStyle2 Harmony Phoenix*	GR	BR		XA-GR			V	-8.00 to +6.00	+0.75 to +3.50	-4.00	18
HOYALUX iD LifeStyle2 Clarity CD Phoenix*	GR	BR		XA-GR			V	-8.00 to +6.00	+0.75 to +3.50	-4.00	14
HOYALUX iD LifeStyle2 Harmony CD Phoenix*	GR	BR		XA-GR			V	-8.00 to +6.00	+0.75 to +3.50	-4.00	14
HOYALUX ARRAY 11, 13, 15, 17, VL Phoenix*	GR	BR		XA-GR	XA-BR		V	-10.00 to +6.00	+0.75 to +3.50	-4.00	14, 17, 19, 21
HOYALUX ARRAY Wrap 11, 13, 15, 17, VL Phoenix*	GR	BR		XA-GR	XA-BR		V	-5.00 to +4.00	+0.75 to +3.50	-4.00	14, 17, 19, 21
Amplitude Phoenix*	GR							-10.00 to +6.00	+0.75 to +3.50	-4.00	19
Amplitude iQ Phoenix*	GR							-10.00 to +6.00	+0.75 to +3.50	-4.00	19
Amplitude HiDef Phoenix*	GR	BR						-8.00 to +6.00	+0.75 to +3.50	-4.00	18
Amplitude Mini Phoenix*	GR							-10.00 to +6.00	+0.75 to +3.50	-4.00	14
Amplitude Mini iQ Phoenix*	GR							-10.00 to +6.00	+0.75 to +3.50	-4.00	14
Amplitude Mini HiDef Phoenix*	GR	BR						-8.00 to +6.00	+0.75 to +3.50	-4.00	14
Amplitude BKS Phoenix*	GR	BR		XA-GR	XA-BR		V	-10.00 to +6.00	+0.75 to +3.50	-4.00	19
Amplitude Mini BKS Phoenix*	GR	BR		XA-GR	XA-BR		V	-10.00 to +6.00	+0.75 to +3.50	-4.00	14
HOYALUX GP Wide	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-4.00	18

**For Varilux S, at fitting height less than 17mm, ADD power is up to +3.50.

PROGRESSIVES

	Transitions Signature			Transitions XTRActive			Transitions VANTAGE	Sphere	Add	Cyl to	MFH
	Gray	Brown	Graphite Green	Gray	Brown	Graphite Green					
1.53 TRIVEX® (CONTINUED)											
SEIKO OPTICAL PRODUCTS											
Seiko 1.53/Trivex® Superior	GR	BR		XA-GR			V	-7.00 to +5.00	+0.50 to +4.00	-5.00	12
Seiko 1.53/Trivex® Prominent	GR	BR		XA-GR			V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
Seiko 1.53/Trivex® Premier 2	GR	BR		XA-GR			V	-7.00 to +5.00	+0.50 to +3.50	-4.00	14
Seiko 1.53/Trivex® Prestige II	GR	BR		XA-GR			V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
Seiko 1.53/Trivex® Prestige II Wide	GR	BR		XA-GR			V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
Seiko 1.53/Trivex® Succeed	GR	BR		XA-GR			V	-7.00 to +5.00	+0.50 to +3.50	-4.00	17
Seiko 1.53/Trivex® Succeed Ws	GR	BR		XA-GR			V	-7.00 to +5.00	+0.50 to +3.50	-4.00	15
Seiko 1.53/Trivex® Supercede 2	GR	BR		XA-GR			V	-7.00 to +5.00	+0.50 to +3.50	-4.00	14
Seiko 1.53/Trivex® Supernal	GR	BR		XA-GR			V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
Seiko 1.53/Trivex® Surmount	GR	BR		XA-GR			V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
Seiko 1.53/Trivex® Surmount Ws	GR	BR		XA-GR			V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
Seiko 1.53/Trivex® PCWide	GR	BR		XA-GR				-5.00 to +6.00	-1.00, -1.50, -2.00	-5.00	15
SHAMIR INSIGHT INC.											
Shamir Autograph® III	GR	BR		XA-GR			V	-13.25 to +7.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
Shamir Autograph® II	GR	BR		XA-GR			V	-13.25 to +6.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
Shamir Autograph® II Attitude®	GR	BR		XA-GR			V	-5.00 to +4.00	+0.75 to +3.50	-6.50	15, 18
Shamir Golf™	GR	BR		XA-GR			V	-5.00 to +4.00	+1.00 to +3.50	-6.50	19
Shamir Spectrum™	GR	BR		XA-GR			V	-13.25 to +6.00	+0.75 to +3.50	-6.50	14, 16, 18
Shamir Spectrum+™	GR	BR		XA-GR			V	-13.25 to +6.00	+0.75 to +3.50	-6.50	14, 16, 18
Shamir Element™	GR	BR						-13.25 to +6.00	+0.75 to +3.50	-6.50	16, 19
Shamir Computer™	GR	BR		XA-GR				-6.00 to +8.00	+0.75 to +2.00	-6.50	14 Lower, 10 Upper
Shamir Workspace™	GR	BR		XA-GR				-6.00 to +8.00	+0.75 to +2.50	-6.50	14 Lower, 10 Upper
Shamir InTouch™	GR	BR		XA-GR				-13.25 to +6.00	+0.75 to +3.50	-6.50	15, 18
Shamir Attitude III Sport	GR	BR		XA-GR			V	-6.50 to +6.00	+0.75 to +3.50	-6.50	18
Shamir Attitude III Fashion	GR	BR		XA-GR			V	-6.50 to +6.00	+0.75 to +3.50	-6.50	15, 18
SIGNET ARMORLITE											
Navigator® FBS Signia® Trivex®	GR	BR						-10.00 to +4.50	+0.75 to +3.50	-6.00	18
Navigator® Short FBS Signia® Trivex®	GR	BR						-10.00 to +4.50	+0.75 to +3.50	-6.00	14
KODAK Precise® PB Trivex®	GR	BR					V	-10.00 to +4.50	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB Trivex®	GR	BR					V	-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique Trivex®	GR	BR		XA-GR			V	-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique HD Trivex®	GR	BR		XA-GR			V	-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique DSII™ Signia® Trivex®	GR	BR						-10.00 to +4.50	+0.75 to +3.50	-6.00	13
X-CEL OPTICAL											
X-Cel Optical Aris™ Freedom ID™	GR			XA-GR			V	-4.50 to +5.50	+1.00 to +3.00	-4.00	17
YOUNGER OPTICS											
Younger Trilogy® Image®	GR			XA-GR				-8.00 to +8.00	+1.00 to +3.00	-4.00	18
1.56 TRIVEX®											
NIKON OPTICAL											
Nikon SeeMax Power TRIVEX®	GR	BR	GG					-10.00 to + 8.00	+0.75 to +4.00	-6.00	13, 15, 17
Nikon Presio Power TRIVEX®	GR	BR	GG					-10.00 to + 8.00	+0.75 to +4.00	-6.00	13, 15, 17
Nikon DigiLife TRIVEX®	GR	BR	GG					-10.00 to + 8.00	+0.75 to +4.00	-6.00	13, 15, 17

Transitions lenses are compatible with all Crizal No-Glare technology

*MFH-Minimum Fitting Height
† The lenses only have near and intermediate zones

PROGRESSIVES

1.59 POLYCARBONATE

CARL ZEISS VISION

	Transitions Signature			Transitions XTRActive			Transitions VANTAGE	Sphere	Add	Cyl to	MFH
	Gray	Brown	Graphite Green	Gray	Brown	Graphite Green					
AO Easy	GR							-9.00 to +5.00	+1.00 to +3.00	-4.00	18
AO Easy HD	GR	BR		XA-GR			V	-9.00 to +6.00	+0.75 to +3.50	-4.00	18
SOLA VIP/VIPGold*	GR							-10.00 to +4.50	+1.00 to +3.00	-4.00	18
SOLAMAX*	GR							-10.00 to +6.00	+1.00 to +3.00	-4.00	18
synchrony Curves HDC	GR	BR						-4.00 to +4.00	+0.75 to +3.50	-4.00	13
synchrony Easy Adapt HD	GR	BR		XA-GR			V	-9.00 to +6.00	+0.75 to +3.50	-6.00	18
synchrony Easy Adapt (One)	GR							-9.00 to +6.00	+1.00 to +3.00	-4.00	18
synchrony Easy M HD	GR	BR		XA-GR			V	-9.00 to +6.00	+0.75 to +3.50	-6.00	15
synchrony Easy S HD	GR	BR		XA-GR			V	-9.00 to +6.00	+0.75 to +3.50	-6.00	13
synchrony Easy S (Ultra)	GR	BR						-9.00 to +5.00	+1.00 to +3.00	-6.00	13
synchrony Easy View HD	GR	BR		XA-GR			V	-9.00 to +6.00	+0.75 to +3.50	-6.00	17
synchrony Easy View (Instinctive)	GR							-9.00 to +5.00	+1.00 to +3.00	-4.00	17
synchrony Easy Wear HD	GR	BR		XA-GR			V	-9.00 to +6.00	+0.75 to +3.50	-6.00	18
synchrony Performance HD13	GR	BR		XA-GR			V	-9.00 to +6.00	+0.75 to +3.50	-6.00	13
synchrony Performance HD15	GR	BR		XA-GR			V	-9.00 to +6.00	+0.75 to +3.50	-6.00	15
synchrony Performance HD17	GR	BR		XA-GR			V	-9.00 to +6.00	+0.75 to +3.50	-6.00	17
synchrony Performance HD19	GR	BR		XA-GR			V	-9.00 to +6.00	+0.75 to +3.50	-6.00	19
synchrony Performance HDV	GR	BR		XA-GR			V	-9.00 to +6.00	+0.75 to +3.50	-6.00	13
synchrony Starter HD	GR	BR		XA-GR			V	-9.00 to +6.00	+0.50 to +2.00	-6.00	14
synchrony Ultra HD 13	GR	BR		XA-GR			V	-9.00 to +6.00	+0.75 to +3.50	-6.00	13
synchrony Ultra HD 15	GR	BR		XA-GR			V	-9.00 to +6.00	+0.75 to +3.50	-6.00	15
synchrony Ultra HD 17	GR	BR		XA-GR			V	-9.00 to +6.00	+0.75 to +3.50	-6.00	17
synchrony Ultra HD 19	GR	BR		XA-GR			V	-9.00 to +6.00	+0.75 to +3.50	-6.00	19
synchrony Ultra HD 21	GR	BR		XA-GR			V	-9.00 to +6.00	+0.75 to +3.50	-6.00	21
synchrony Ultra HDC	GR	BR		XA-GR			V	-9.00 to +6.00	+0.75 to +3.50	-6.00	13
synchrony Ultra HDV	GR	BR		XA-GR			V	-9.00 to +6.00	+0.75 to +3.50	-6.00	13
synchrony Work and Go HD	GR	BR		XA-GR			V	-9.00 to +6.00	+0.75 to +3.50	-6.00	17
synchrony Work and Office HD	GR	BR		XA-GR			V	-9.00 to +6.00	+0.75 to +3.50	-6.00	16
synchrony Work and Read HD	GR	BR		XA-GR			V	-9.00 to +6.00	+0.75 to +3.00	-6.00	16
Synchrony-	GR							-9.00 to +6.00	+1.00 to +3.00	-4.00	17
ZEISS Choice 13	GR	BR		XA-GR			V	-10.00 to +6.00	+0.75 to +3.50	-4.00	13
ZEISS Choice 15	GR	BR		XA-GR			V	-10.00 to +6.00	+0.75 to +3.50	-4.00	15
ZEISS Choice 17	GR	BR		XA-GR			V	-10.00 to +6.00	+0.75 to +3.50	-4.00	17
ZEISS Choice 19	GR	BR		XA-GR			V	-10.00 to +6.00	+0.75 to +3.50	-4.00	19
ZEISS Choice Plus Sport 13	GR	BR						-4.00 to +4.00	+0.75 to +3.50	-4.00	13
ZEISS Choice Plus Sport 17	GR	BR						-4.00 to +4.00	+0.75 to +3.50	-4.00	17
ZEISS Digital 1000	GR	BR		XA-GR			V	-10.00 to +6.00	+1.00	-4.00	13
ZEISS Digital 1250	GR	BR		XA-GR			V	-10.00 to +6.00	+1.25	-4.00	13
ZEISS Digital 500	GR	BR		XA-GR			V	-10.00 to +6.00	+0.50	-4.00	13
ZEISS Digital 750	GR	BR		XA-GR			V	-10.00 to +6.00	+0.75	-4.00	13
ZEISS Digital Lens	GR	BR		XA-GR			V	-10.00 to +6.00	+0.50 to +1.25	-4.00	13
ZEISS Digital Sport	GR	BR		XA-GR			V	-10.00 to +6.00	+0.50 to +1.25	-4.00	13
ZEISS Digital Sport 1000	GR	BR		XA-GR			V	-10.00 to +6.00	+1.00	-4.00	13
ZEISS Digital Sport 1250	GR	BR		XA-GR			V	-10.00 to +6.00	+1.25	-4.00	13
ZEISS Digital Sport 500	GR	BR		XA-GR			V	-10.00 to +6.00	+0.50	-4.00	13
ZEISS Digital Sport 750	GR	BR		XA-GR			V	-10.00 to +6.00	+0.75	-4.00	13
ZEISS DriveSafe Individual	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-4.00	13
ZEISS DriveSafe Individual Sport	GR	BR						-4.00 to +4.00	+0.75 to +3.50	-4.00	13
ZEISS DriveSafe	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-4.00	13
ZEISS DriveSafe Sport	GR	BR						-4.00 to +4.00	+0.75 to +3.50	-4.00	13
ZEISS GT2*	GR	BR						-10.00 to +6.00	+1.00 to +3.00	-4.00	17
ZEISS GT2* Short	GR	BR						-10.00 to +6.00	+1.00 to +3.00	-4.00	13
ZEISS Individual* 2	GR	BR		XA-GR			V	-10.00 to +6.00	+0.75 to +3.50	-4.00	13
ZEISS Individual* 2I	GR	BR		XA-GR			V	-10.00 to +6.00	+0.75 to +3.50	-4.00	13
ZEISS Individual* 2N	GR	BR		XA-GR			V	-10.00 to +6.00	+0.75 to +3.50	-4.00	13
ZEISS Individual* Sport	GR	BR						-4.00 to +4.00	+0.75 to +3.50	-4.00	13
ZEISS Individual* Wrap	GR	BR						-4.00 to +4.00	+1.00 to +2.50	-4.00	13
ZEISS Plus 2 Short Sport	GR	BR						-4.00 to +4.00	+0.75 to +3.50	-4.00	14
ZEISS Plus 2 Sport	GR	BR						-4.00 to +4.00	+0.75 to +3.50	-4.00	18
ZEISS Precision Plus	GR	BR		XA-GR			V	-10.00 to +6.00	+0.75 to +3.50	-4.00	13
ZEISS Precision Pure 13	GR	BR		XA-GR			V	-10.00 to +6.00	+0.75 to +3.50	-4.00	13
ZEISS Precision Pure 15	GR	BR		XA-GR			V	-10.00 to +6.00	+0.75 to +3.50	-4.00	15
ZEISS Precision Pure 17	GR	BR		XA-GR			V	-10.00 to +6.00	+0.75 to +3.50	-4.00	17
ZEISS Precision Pure 19	GR	BR		XA-GR			V	-10.00 to +6.00	+0.75 to +3.50	-4.00	19
ZEISS Precision Pure 21	GR	BR		XA-GR			V	-10.00 to +6.00	+0.75 to +3.50	-4.00	21
ZEISS Precision Superb	GR	BR		XA-GR			V	-10.00 to +6.00	+0.75 to +3.50	-4.00	13

PROGRESSIVES

1.59 POLYCARBONATE (CONTINUED)

ESSILOR LENS

	Transitions Signature			Transitions XTRActive			Transitions VANTAGE	Sphere	Add	Cyl to	MFH
	Gray	Brown	Graphite Green	Gray	Brown	Graphite Green					
Varilux S 4D™ Airwear*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14**
Varilux S Fit™ Airwear*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14**
Varilux S Design™ Airwear*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14**
Varilux® Physio® W3+ eyecode Airwear*	GR	BR	GG	XA-GR	XA-BR			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Airwear*	GR	BR	GG	XA-GR	XA-BR			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Airwear*	GR	BR	GG	XA-GR	XA-BR			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® W2+ eyecode Airwear*	GR	BR	GG	XA-GR	XA-BR			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® W2+ FIT Airwear*	GR	BR	GG	XA-GR	XA-BR			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® W2+ Airwear*	GR	BR	GG	XA-GR	XA-BR			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio Enhanced eyecode™ Airwear*	GR	BR	GG	XA-GR				-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio Enhanced FIT™ Airwear*	GR	BR	GG	XA-GR				-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio Enhanced™ Airwear*	GR	BR	GG	XA-GR				-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio DRx™ Airwear*	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-10.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Airwear*	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Varilux® Physio™ Airwear*	GR	BR	GG	XA-GR				-10.00 to +6.00	+0.75 to +3.50 & +4.00	-4.00	17
Varilux Physio Short™ Airwear*	GR							-10.00 to +6.00	+0.75 to +3.00	-4.00	14
Varilux Comfort Enhanced FIT™ Airwear*	GR	BR	GG	XA-GR				-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort Enhanced™ Airwear*	GR	BR	GG	XA-GR				-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort DRx™ Airwear*	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Airwear*	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-10.00 to +6.00	+0.75 to +3.00	-6.00	14
Varilux Comfort™ Airwear*	GR	BR	GG	XA-GR				-10.00 to +6.00	+0.75 to +3.50	-4.00	17
Varilux Comfort Short™ Airwear*	GR	BR	GG	XA-GR				-10.00 to +6.00	+0.75 to +3.00	-4.00	14
Varilux® Stylistic™ Wrap Airwear*	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-4.00 to +4.00	+0.75 to +3.50	-2.00	17
DEFINITY™ 3 PLUS Airwear*	GR	BR						-10.00 to +6.00	+1.00 to +3.50	-6.00	15
DEFINITY™ 3 Airwear*	GR	BR						-10.00 to +6.00	+1.00 to +3.50	-6.00	15
DEFINITY™ Airwear*	GR							-7.00 to +4.00	+1.00 to +3.00	-4.00	18
DEFINITY SHORT™ Airwear*	GR							-7.00 to +4.00	+1.00 to +3.00	-4.00	15
Essilor Ideal Advanced™ Airwear*	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced Wrap™ Airwear*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-4.00 to +3.00	+0.75 to +3.00	-4.00	17
Essilor Ideal™ Airwear*	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal™ Fit Airwear*	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal™ Sport Airwear*	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-10.00 to +8.00	+0.75 to +3.00	-6.00	18
Essilor Ideal™ Computer Airwear*	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal™ Driving Airwear*	GR	BR	GG	XA-GR	XA-BR	XA-GG	V	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Accolade Freedom™ Airwear*	GR	BR		XA-GR			V	-10.00 to +6.00	+0.75 to +3.50	-4.00	15
Essilor Accolade™ Airwear*	GR	BR		XA-GR			V	-10.00 to +6.00	+0.75 to +3.50	-4.00	17
Ovation™ Digital Airwear*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Ovation™ Airwear*	GR	BR						-10.00 to +8.00	+0.75 to +3.50	-4.00	17
Essilor Natural™ Digital Airwear*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Natural™ Airwear*	GR							-10.00 to +8.00	+0.75 to +3.50	-4.00	18
Essilor SmallFit™ Digital Airwear*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor SmallFit™ Airwear*	GR	BR						-10.00 to +8.00	+0.75 to +3.00	-4.00	14
Adaptar™ Digital Airwear*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar™ Short Digital Airwear*	GR	BR	GG	XA-GR	XA-BR	XA-GG		-10.00 to +8.00	+0.75 to +3.00	-6.00	14

HOYA VISION CARE

HOYALUX Summit CD Poly	GR	BR						-6.00 to +6.00	+1.00 to +3.00	-4.00	14
HOYALUX Summit CD iQ Poly	GR	BR						-6.00 to +6.00	+1.00 to +3.00	-4.00	14
HOYALUX Summit ECP Poly	GR	BR						-6.00 to +6.00	+1.00 to +3.00	-4.00	18
HOYALUX Summit ECP iQ Poly	GR	BR						-6.00 to +6.00	+1.00 to +3.00	-4.00	18
HOYALUX iD LifeStyle2 Clarity Poly	GR							-6.00 to +6.00	+1.00 to +3.00	-4.00	18
HOYALUX iD LifeStyle2 Harmony Poly	GR							-6.00 to +6.00	+1.00 to +3.00	-4.00	18
HOYALUX iD LifeStyle2 Clarity CD Poly	GR							-6.00 to +6.00	+1.00 to +3.00	-4.00	14
HOYALUX iD LifeStyle2 Harmony CD Poly	GR							-6.00 to +6.00	+1.00 to +3.00	-4.00	14
HOYALUX ARRAY 11, 13, 15, 17, VL Poly	GR	BR		XA-GR			V	-13.00 to +5.25	+0.75 to +3.50	-4.00	14, 17, 19, 21
HOYALUX ARRAY Wrap 11, 13, 15, 17, VL Poly	GR	BR		XA-GR			V	-6.00 to +4.00	+0.75 to +3.50	-4.00	14, 17, 19, 21
Amplitude BKS Poly	GR	BR		XA-GR			V	-13.00 to +5.25	+0.75 to +3.50	-4.00	19
Amplitude Mini BKS Poly	GR	BR		XA-GR			V	-13.00 to +5.25	+0.75 to +3.50	-4.00	14
HOYALUX GP Wide	GR	BR						-13.00 to +5.25	+0.75 to +3.50	-4.00	18

**For Varilux S, at fitting height less than 17mm, ADD power is up to +3.50.

Transitions lenses are compatible with all Crizal No-Glare technology

TRANSITIONS EVERYDAY LENSES

*MFH-Minimum Fitting Height

PROGRESSIVES

	Transitions Signature			Transitions XTRActive			Transitions VANTAGE	Sphere	Add	Cyl to	MFH
	Gray	Brown	Graphite Green	Gray	Brown	Graphite Green					
1.59 POLYCARBONATE (CONTINUED)											
NIKON OPTICAL											
Nikon SeeMax Power Polycarbonate	GR	BR	GG				V	-10.00 to + 8.00	+0.75 to +4.00	-6.00	13, 15, 17
Nikon Presio Power Polycarbonate	GR	BR	GG	XA-GR			V	-10.00 to + 8.00	+0.75 to +4.00	-6.00	13, 15, 17
Nikon DigiLife Polycarbonate	GR	BR	GG	XA-GR			V	-10.00 to + 8.00	+0.75 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital Polycarbonate	GR	BR	GG	XA-GR			V	-10.00 to + 5.00	+0.75 to +4.00	-4.00	15, 17
Nikon Move Polycarbonate	GR	BR	GG	XA-GR			V	-10.00 to + 5.00	+0.75 to +4.00	-4.00	16, 18
SEIKO OPTICAL PRODUCTS											
Seiko 1.59 Superior	GR	BR		XA-GR			V	-7.00 to +5.00	+0.50 to +4.00	-5.00	12
Seiko 1.59 Prominent	GR	BR		XA-GR			V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
Seiko 1.59 Premier 2	GR	BR		XA-GR			V	-7.00 to +5.00	+0.50 to +3.50	-4.00	14
Seiko 1.59 Prestige II	GR	BR		XA-GR			V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
Seiko 1.59 Prestige II Wide	GR	BR		XA-GR			V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
Seiko 1.59 Succeed	GR	BR		XA-GR			V	-7.00 to +5.00	+0.50 to +3.50	-4.00	17
Seiko 1.59 Succeed Ws	GR	BR		XA-GR			V	-7.00 to +5.00	+0.50 to +3.50	-4.00	15
Seiko 1.59 Supercede 2	GR	BR		XA-GR			V	-7.00 to +5.00	+0.50 to +3.50	-4.00	14
Seiko 1.59 Supernal	GR	BR		XA-GR			V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
Seiko 1.59 Surmount	GR	BR		XA-GR			V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
Seiko 1.59 Surmount Ws	GR	BR		XA-GR			V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
Seiko 1.59 PCWide	GR	BR		XA-GR				-5.00 to +6.00	-1.00, -1.50, -2.00	-5.00	15
SHAMIR INSIGHT INC.											
Shamir Autograph® III	GR	BR	GG	XA-GR			V	-14.75 to +8.00	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® II	GR	BR	GG	XA-GR			V	-14.75 to +6.50	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® II Attitude®	GR	BR	GG	XA-GR			V	-6.50 to +6.00	+0.75 to +3.50	-7.00	15, 18
Shamir Golf®	GR	BR	GG	XA-GR			V	-5.00 to +4.00	+1.00 to +3.50	-7.00	19
Shamir Spectrum®	GR	BR	GG	XA-GR			V	-14.75 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
Shamir Spectrum+®	GR	BR	GG	XA-GR			V	-14.75 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
Shamir Element®	GR	BR	GG					-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19
Shamir Computer®	GR	BR	GG	XA-GR				-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper
Shamir Workspace®	GR	BR	GG	XA-GR				-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper
Shamir InTouch®	GR	BR	GG	XA-GR			V	-14.75 to +6.00	+0.75 to +3.50	-7.00	15, 18
Shamir Attitude III Sport	GR	BR	GG	XA-GR			V	-6.50 to +6.00	+0.75 to +3.50	-7.00	18
Shamir Attitude III Fashion	GR	BR	GG	XA-GR			V	-6.50 to +6.00	+0.75 to +3.50	-7.00	15, 18
SIGNET ARMORLITE											
DirecTek® Signia® Poly	GR	BR	GG					-10.00 to +5.00	+0.75 to +3.50	-6.00	16.5
DirecTek® Short Signia® Poly	GR	BR	GG					-10.00 to +5.00	+0.75 to +3.50	-6.00	13.5
Navigator® FBS Signia® 1.67 Poly	GR	BR	GG					-10.00 to +5.00	+0.75 to +3.50	-6.00	18
Navigator® Short FBS Signia® 1.67 Poly	GR	BR	GG					-10.00 to +5.00	+0.75 to +3.50	-6.00	14
KODAK Precise® Signia® Poly SF	GR							-10.00 to +5.00	+1.00 to +3.00	-4.00	17
KODAK Precise® Short Signia® Poly SF	GR							-10.00 to +5.00	+1.00 to +3.00	-4.00	13
KODAK Precise® Digital Signia® Poly	GR							-10.00 to +6.00	+1.00 to +3.00	-6.00	17
KODAK Precise® Short Digital Signia® Poly	GR							-10.00 to +6.00	+1.00 to +3.00	-6.00	13
KODAK Precise® PB Signia® Poly	GR	BR	GG				V	-10.00 to +5.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB Signia® Poly	GR	BR	GG				V	-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique® Signia® Poly	GR	BR	GG	XA-GR			V	-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique® HD Signia® Poly	GR	BR	GG	XA-GR			V	-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique DSII® Signia® Poly	GR	BR						-10.00 to +5.00	+0.75 to +3.50	-6.00	13
YOUNGER OPTICS											
Younger Image®	GR	BR		XA-GR			V	-9.00 to +7.00	+1.00 to +3.00	-4.00	18
Younger Adage®	GR							-8.00 to +6.00	+1.00 to +3.00	-4.00	13

Transitions lenses are compatible with all Crizal No-Glare technology

*MFH-Minimum Fitting Height
† The lenses only have near and intermediate zones

PROGRESSIVES

	Transitions Signature			Transitions XTRActive			Transitions	Sphere	Add	Cyl to	MFH
	Gray	Brown	Graphite Green	Gray	Brown	Graphite Green	VANTAGE Gray				
1.60 HIGH INDEX											
CARL ZEISS VISION											
ZEISS DriveSafe Individual Sport	GR	BR						-4.00 to +4.00	+0.75 to +3.50	-4.00	13
ZEISS DriveSafe Sport	GR	BR						-4.00 to +4.00	+0.75 to +3.50	-4.00	13
ZEISS Individual [®] 2	GR	BR						-10.00 to +6.00	+0.75 to +4.00	-4.00	13
ZEISS Individual [®] 2I	GR	BR						-10.00 to +6.00	+0.75 to +4.00	-4.00	13
ZEISS Individual [®] 2N	GR	BR						-10.00 to +6.00	+0.75 to +4.00	-4.00	13
ZEISS Individual [®] Sport	GR	BR						-4.00 to +4.00	+0.75 to +3.50	-4.00	13
ZEISS Plus 2 Short Sport	GR	BR						-4.00 to +4.00	+0.75 to +3.50	-4.00	14
ZEISS Plus 2 Sport	GR	BR						-4.00 to +4.00	+0.75 to +3.50	-4.00	18
ZEISS Precision Plus	GR	BR						-10.00 to +6.00	+0.75 to +4.00	-4.00	13
ZEISS Precision Pure 13	GR	BR						-10.00 to +6.00	+0.75 to +4.00	-4.00	13
ZEISS Precision Pure 15	GR	BR						-10.00 to +6.00	+0.75 to +4.00	-4.00	15
ZEISS Precision Pure 17	GR	BR						-10.00 to +6.00	+0.75 to +4.00	-4.00	17
ZEISS Precision Pure 19	GR	BR						-10.00 to +6.00	+0.75 to +4.00	-4.00	19
ZEISS Precision Pure 21	GR	BR						-10.00 to +6.00	+0.75 to +4.00	-4.00	21
ZEISS Precision Superb	GR	BR						-10.00 to +6.00	+0.75 to +4.00	-4.00	13
ESSILOR LENS											
Varilux [®] Physio [®] W3+ eyecode Thin&Lite [®] 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux [®] Physio [®] W3+ FIT Thin&Lite [®] 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux [®] Physio [®] W3+ Thin&Lite [®] 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort [®] W2+ eyecode Thin&Lite [®] 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Varilux Comfort [®] W2+ FIT Thin&Lite [®] 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Varilux Comfort [®] W2+ Thin&Lite [®] 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Varilux Physio Enhanced eyecode [™] Thin&Lite [®] 1.60	GR							-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio Enhanced FIT [™] Thin&Lite [®] 1.60	GR							-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio Enhanced [™] Thin&Lite [®] 1.60	GR							-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio DRx [™] Thin&Lite [®] 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx [™] Thin&Lite [®] 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux [®] Physio [®] Thin&Lite [®] 1.60	GR							-12.00 to +8.00	+0.75 to +3.50 & +4.00	-4.00	17
Varilux Comfort Enhanced Fit [™] Thin&Lite [®] 1.60	GR							-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Varilux Comfort Enhanced [™] Thin&Lite [®] 1.60	GR							-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Varilux Comfort DRx [™] Thin&Lite [®] 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx [™] Thin&Lite [®] 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux Comfort [®] Thin&Lite [®] 1.60	GR							-12.00 to +8.00	+0.75 to +3.50	-4.00	17
Varilux Comfort Short [®] Thin&Lite [®] 1.60	GR							-12.00 to +8.00	+0.75 to +3.00	-4.00	14
Essilor Ideal Advanced [™] Thin&Lite [®] 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal [™] Thin&Lite [®] 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal [™] Fit Thin&Lite [®] 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	18
Essilor Ideal [™] Sport Thin&Lite [®] 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	18
Essilor Ideal [™] Computer Thin&Lite [®] 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	13
Essilor Ideal [™] Driving Thin&Lite [®] 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Ovation [®] Digital Thin&Lite [®] 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Natural [®] Digital Thin&Lite [®] 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	17
SmallFit [®] Digital Thin&Lite [®] 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.00	-6.00	14
HOYA VISION CARE											
HOYALUX Summit CD 1.60	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-4.00	14
HOYALUX Summit CD iQ 1.60	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-6.00	14
HOYALUX Summit ECP 1.60	GR	BR						-10.00 to +8.00	+0.75 to +3.50	-4.00	18
HOYALUX Summit ECP iQ 1.60	GR	BR						-10.00 to +8.00	+0.75 to +3.50	-6.00	18
HOYALUX iD LifeStyle2 Clarity 1.60	GR	BR						-11.00 to +8.00	+0.75 to +3.50	-4.00	18
HOYALUX iD LifeStyle2 Harmony 1.60	GR	BR						-11.00 to +8.00	+0.75 to +3.50	-4.00	18
HOYALUX iD LifeStyle2 Clarity CD 1.60	GR	BR						-11.00 to +8.00	+0.75 to +3.50	-4.00	14
HOYALUX iD LifeStyle2 Harmony CD 1.60	GR	BR						-11.00 to +8.00	+0.75 to +3.50	-4.00	14
HOYALUX ARRAY 11, 13, 15, 17, VL 1.60	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-4.00	14, 17, 19, 21
Amplitude 1.60	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-4.00	19
Amplitude HiDef 1.60	GR	BR						-11.00 to +8.00	+0.75 to +3.50	-4.00	18
Amplitude Mini 1.60	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-4.00	14
Amplitude Mini HiDef 1.60	GR	BR						-11.00 to +8.00	+0.75 to +3.50	-4.00	14
Amplitude BKS	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-4.00	19
Amplitude Mini BKS	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-4.00	14
NIKON OPTICAL											
Nikon SeeMax Power 1.60	GR	BR	GG					-10.00 to + 8.00	+0.75 to +3.50	-6.00	13, 15, 17
Nikon Presio Power 1.60	GR	BR	GG					-12.00 to + 8.00	+0.75 to +3.50	-6.00	13, 15, 17
Nikon DigiLife 1.60	GR	BR	GG					-12.00 to + 8.00	+0.75 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital 1.60	GR	BR	GG					-10.00 to + 5.00	+0.75 to +4.00	-6.00	15, 17
Nikon Move 1.60	GR	BR	GG					-10.00 to + 5.00	+0.75 to +3.50	-4.00	16, 18

PROGRESSIVES

	Transitions Signature		Transitions XTRActive			Transitions VANTAGE	Sphere	Add	Cyl to	MFH
	Gray	Brown	Graphite Green	Gray	Brown	Graphite Green				
1.60 HIGH INDEX (CONTINUED)										
SEIKO OPTICAL PRODUCTS										
Seiko 1.60 Superior	GR	BR					-8.50 to +6.50	+0.50 to +4.00	-5.00	12
Seiko 1.60 Prominent	GR	BR					-8.50 to +6.50	+0.50 to +3.50	-5.00	14
Seiko 1.60 Premier 2	GR	BR					-8.50 to +6.50	+0.50 to +3.50	-4.00	14
Seiko 1.60 Prestige II	GR	BR					-8.50 to +6.50	+0.50 to +3.50	-5.00	14
Seiko 1.60 Prestige II Wide	GR	BR					-8.50 to +6.50	+0.50 to +3.50	-5.00	14
Seiko 1.60 Succeed	GR	BR					-8.50 to +6.50	+0.50 to +3.50	-4.00	17
Seiko 1.60 Succeed Ws	GR	BR					-8.50 to +6.50	+0.50 to +3.50	-4.00	15
Seiko 1.60 Supercede 2	GR	BR					-8.50 to +6.50	+0.50 to +3.50	-4.00	14
Seiko 1.60 Supernal	GR	BR					-8.50 to +6.50	+0.50 to +3.50	-5.00	14
Seiko 1.60 Surmount	GR	BR					-8.50 to +6.50	+0.50 to +3.50	-5.00	14
Seiko 1.60 Surmount Ws	GR	BR					-8.50 to +6.50	+0.50 to +3.50	-5.00	14
Seiko 1.60 PCWide	GR	BR					-6.50 to +7.50	-1.00, -1.50, -2.00	-5.00	15
SHAMIR INSIGHT INC.										
Shamir Autograph® III SuperLite™ 1.6	GR	BR			XA-GR		-15.00 to +8.00	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® II SuperLite™ 1.6	GR	BR			XA-GR		-15.00 to +6.00	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® II Attitude™ SuperLite™ 1.6	GR	BR			XA-GR		-6.50 to +5.25	+0.75 to +3.50	-7.00	15, 18
Shamir Golf™ SuperLite™ 1.6	GR	BR			XA-GR		-6.00 to +5.25	+1.00 to +3.50	-7.00	19
Shamir Spectrum™ SuperLite™ 1.6	GR	BR			XA-GR		-14.75 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
Shamir Spectrum+™ SuperLite™ 1.6	GR	BR			XA-GR		-14.75 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
Shamir Element™ SuperLite™ 1.6	GR	BR			XA-GR		-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19
Shamir FirstPal 1.6	GR	BR			XA-GR		-15.00 to +6.00	+0.75 to +1.50	-7.00	15
Shamir Computer™	GR	BR			XA-GR		-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper
Shamir Workspace™	GR	BR			XA-GR		-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper
Shamir InTouch™	GR	BR			XA-GR		-15.00 to +6.00	+0.75 to +3.50	-7.00	15, 18
Shamir Attitude III Sport	GR	BR			XA-GR		-6.50 to +6.00	+0.75 to +3.50	-7.00	18
Shamir Attitude III Fashion	GR	BR			XA-GR		-6.50 to +6.00	+0.75 to +3.50	-7.00	15, 18
SIGNET ARMORLITE										
KODAK Unique	GR	BR					-10.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD	GR	BR					-10.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK DSII™	GR	BR					-10.00 to +6.00	+0.75 to +3.50	-6.00	13
1.67 SUPER HIGH INDEX										
CARL ZEISS VISION										
AO Easy	GR						-10.00 to +6.00	+0.75 to +3.50	-4.00	18
AO Easy HD	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-4.00	18
synchrony Easy Adapt HD	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-6.00	18
synchrony Easy Adapt (One)	GR	BR					-12.00 to +8.00	+1.00 to +3.50	-4.00	18
synchrony Easy M HD	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-6.00	15
synchrony Easy S HD	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-6.00	13
synchrony Easy S (Ultra)	GR						-10.00 to +6.00	+0.75 to +3.50	-6.00	13
synchrony Easy View HD	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-6.00	17
synchrony Easy View (Instinctive)	GR						-10.00 to +6.00	+0.75 to +3.50	-4.00	17
synchrony Easy Wear HD	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-6.00	18
synchrony Performance HD13	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-6.00	13
synchrony Performance HD15	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-6.00	15
synchrony Performance HD17	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-6.00	17
synchrony Performance HD19	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-6.00	19
synchrony Performance HDV	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-6.00	13
synchrony Starter HD	GR	BR			XA-GR		-12.00 to +8.00	+0.50 to +2.00	-6.00	14
synchrony Ultra HD 13	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-6.00	13
synchrony Ultra HD 15	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-6.00	15
synchrony Ultra HD 17	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-6.00	17
synchrony Ultra HD 19	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-6.00	19
synchrony Ultra HD 21	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-6.00	21
synchrony Ultra HDC	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-6.00	13
synchrony Ultra HDV	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-6.00	13
synchrony Work and Go HD	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-6.00	17
synchrony Work and Office HD	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-6.00	16
synchrony Work and Read HD	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.00	-6.00	16
Synchrony™	GR						-12.00 to +8.00	+1.00 to +3.50	-4.00	17
ZEISS Choice 13	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-4.00	13
ZEISS Choice 15	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-4.00	15
ZEISS Choice 17	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-4.00	17
ZEISS Choice 19	GR	BR			XA-GR		-12.00 to +8.00	+0.75 to +3.50	-4.00	19

Transitions lenses are compatible with all Crizal No-Glare technology

*MFH-Minimum Fitting Height
† The lenses only have near and intermediate zones

PROGRESSIVES

1.67 SUPER HIGH INDEX (CONTINUED)

CARL ZEISS VISION

	Transitions Signature			Transitions XTRActive			Transitions VANTAGE	Sphere	Add	Cyl to	MFH
	Gray	Brown	Graphite Green	Gray	Brown	Graphite Green					
ZEISS Digital 1000	GR	BR		XA-GR				-12.00 to +8.00	+1.00	-4.00	13
ZEISS Digital 1250	GR	BR		XA-GR				-12.00 to +8.00	+1.25	-4.00	13
ZEISS Digital 500	GR	BR		XA-GR				-12.00 to +8.00	+0.50	-4.00	13
ZEISS Digital 750	GR	BR		XA-GR				-12.00 to +8.00	+0.75	-4.00	13
ZEISS Digital Lens	GR	BR		XA-GR				-12.00 to +8.00	+0.50 to +1.25	-4.00	13
ZEISS DriveSafe Individual	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-4.00	13
ZEISS DriveSafe	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-4.00	13
ZEISS EnergizeMe Digital Lens	GR	BR						-12.00 to +8.00	+0.65	-4.00	14
ZEISS EnergizeMe 14	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-4.00	14
ZEISS EnergizeMe 16	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-4.00	16
ZEISS EnergizeMe 18	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-4.00	18
ZEISS GT2 [®]	GR	BR						-16.00 to +6.00	+1.00 to +3.00	-4.00	17
ZEISS GT2 [®] Short	GR	BR						-16.00 to +6.00	+1.00 to +3.00	-4.00	13
ZEISS Individual [®] 2	GR	BR		XA-GR				-12.00 to +8.00	+0.75 to +3.50	-4.00	13
ZEISS Individual [®] 2I	GR	BR		XA-GR				-12.00 to +8.00	+0.75 to +3.50	-4.00	13
ZEISS Individual [®] 2N	GR	BR		XA-GR				-12.00 to +8.00	+0.75 to +3.50	-4.00	13
ZEISS Precision Plus	GR	BR		XA-GR				-12.00 to +8.00	+0.75 to +3.50	-4.00	13
ZEISS Precision Pure 13	GR	BR		XA-GR				-12.00 to +8.00	+0.75 to +3.50	-4.00	13
ZEISS Precision Pure 15	GR	BR		XA-GR				-12.00 to +8.00	+0.75 to +3.50	-4.00	15
ZEISS Precision Pure 17	GR	BR		XA-GR				-12.00 to +8.00	+0.75 to +3.50	-4.00	17
ZEISS Precision Pure 19	GR	BR		XA-GR				-12.00 to +8.00	+0.75 to +3.50	-4.00	19
ZEISS Precision Pure 21	GR	BR		XA-GR				-12.00 to +8.00	+0.75 to +3.50	-4.00	21
ZEISS Precision Superb	GR	BR		XA-GR				-12.00 to +8.00	+0.75 to +3.50	-4.00	13

ESSILOR LENS

Varilux S 4D [™] Thin&Lite [®] 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-12.00 to +6.00	+0.75 to +4.00	-6.00	14**
Varilux S Fit [™] Thin&Lite [®] 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-12.00 to +6.00	+0.75 to +4.00	-6.00	14**
Varilux S Design [™] Thin&Lite [®] 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-12.00 to +6.00	+0.75 to +4.00	-6.00	14**
Varilux [®] Physio [®] W3+ eyecode Thin&Lite [®] 1.67	GR	BR		XA-GR				-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux [®] Physio [®] W3+ FIT Thin&Lite [®] 1.67	GR	BR		XA-GR				-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux [®] Physio [®] W3+ Thin&Lite [®] 1.67	GR	BR		XA-GR				-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort [®] W2+ eyecode Thin&Lite [®] 1.67	GR	BR		XA-GR				-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Comfort [®] W2+ FIT Thin&Lite [®] 1.67	GR	BR		XA-GR				-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Comfort [®] W2+ Thin&Lite [®] 1.67	GR	BR		XA-GR				-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Physio Enhanced eyecode [™] Thin&Lite [®] 1.67	GR	BR						-12.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio Enhanced FIT [™] Thin&Lite [®] 1.67	GR	BR						-12.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio Enhanced [™] Thin&Lite [®] 1.67	GR	BR						-12.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio DRx [™] Thin&Lite [®] 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-12.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx [™] Thin&Lite [®] 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-12.00 to +9.00	+0.75 to +3.00	-6.00	14
Varilux Physio [™] Thin&Lite [®] 1.67	GR	BR						-12.00 to +9.00	+0.75 to +3.50 & +4.00	-4.00	17
Varilux Physio Short [™] Thin&Lite [®] 1.67	GR	BR						-12.00 to +9.00	+0.75 to +3.00	-4.00	14
Varilux Comfort Enhanced [™] Thin&Lite [®] 1.67	GR	BR						-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Comfort Enhanced Fit [™] Thin&Lite [®] 1.67	GR	BR						-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Comfort DRx [™] Thin&Lite [®] 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-12.00 to +9.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx [™] Thin&Lite [®] 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-12.00 to +9.00	+0.75 to +3.00	-6.00	14
Varilux Comfort [™] Thin&Lite [®] 1.67	GR	BR						-14.00 to +9.00	+0.75 to +3.50	-4.00	17
Varilux Comfort Short [™] Thin&Lite [®] 1.67	GR	BR						-14.00 to +9.00	+0.75 to +3.00	-4.00	14
Varilux [®] Stylistic [®] Wrap Thin&Lite [®] 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-5.00 to +5.00	+0.75 to +3.50	-3.00	17
DEFINITY [®] 3 Plus Thin&Lite [®] 1.67	GR	BR						-12.00 to +9.00	+1.00 to +3.50	-6.00	15
DEFINITY [®] 3 Thin&Lite [®] 1.67	GR	BR						-12.00 to +9.00	+1.00 to +3.50	-6.00	15
DEFINITY [®] Thin&Lite [®] 1.67	GR							-12.00 to +10.00	+1.00 to +3.50 & +4.00	-4.00	18
DEFINITY SHORT [™] Thin&Lite [®] 1.67	GR							-12.00 to +10.00	+1.00 to +3.50 & +4.00	-4.00	15
Essilor Accolade Freedom [™] Thin&Lite [®] 1.67	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-4.00	15
Essilor Accolade [™] Thin&Lite [®] 1.67	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-4.00	17
Essilor Ideal Advanced [™] Thin&Lite [®] 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-12.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced Wrap [™] Thin&Lite [®] 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-5.00 to +3.00	+0.75 to +3.50	-6.00	17
Essilor Ideal [®] Thin&Lite [®] 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-12.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal [®] Fit Thin&Lite [®] 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-12.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal [®] Sport Thin&Lite [®] 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-12.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal [®] Computer Thin&Lite [®] 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-12.00 to +9.00	+0.75 to +3.00	-6.00	13
Essilor Ideal [®] Driving Thin&Lite [®] 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-12.00 to +9.00	+0.75 to +3.00	-6.00	14
Ovation [®] Digital Thin&Lite [®] 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-12.00 to +9.00	+0.75 to +3.50	-6.00	17
Ovation [®] Thin&Lite [®] 1.67	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-4.00	17
Essilor Natural [®] Digital Thin&Lite [®] 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-12.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor SmallFit [®] Digital Thin&Lite [®] 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-12.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor SmallFit [®] Thin&Lite [®] 1.67	GR	BR						-12.00 to +8.00	+0.75 to +3.00	-4.00	14
Adaptar [®] Digital Thin&Lite [®] 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-12.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar [®] Short Digital Thin&Lite [®] 1.67	GR	BR	GG	XA-GR	XA-BR	XA-GG		-12.00 to +9.00	+0.75 to +3.00	-6.00	14

**For Varilux S, at fitting height less than 17mm, ADD power is up to +3.50.

Transitions lenses are compatible with all Crizal No-Glare technology

TRANSITIONS EVERYDAY LENSES

*MFH-Minimum Fitting Height

PROGRESSIVES

	Transitions Signature			Transitions XTRActive			Transitions VANTAGE	Sphere	Add	Cyl to	MFH
	Gray	Brown	Graphite Green	Gray	Brown	Graphite Green					
1.67 SUPER HIGH INDEX (CONTINUED)											
HOYA VISION CARE											
HOYALUX Summit CD 1.67	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-4.00	14
HOYALUX Summit CD iQ 1.67	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-4.00	14
HOYALUX Summit ECP 1.67	GR	BR						-9.75 to +8.00	+0.75 to +3.50	-4.00	18
HOYALUX Summit ECP iQ 1.67	GR	BR						-9.75 to +8.00	+0.75 to +3.50	-4.00	18
HOYALUX iD LifeStyle2 Clarity 1.67	GR	BR						-13.00 to +8.00	+0.75 to +3.50	-4.00	18
HOYALUX iD LifeStyle2 Harmony 1.67	GR	BR						-13.00 to +8.00	+0.75 to +3.50	-4.00	18
HOYALUX iD LifeStyle2 Clarity CD 1.67	GR	BR						-13.00 to +8.00	+0.75 to +3.50	-4.00	14
HOYALUX iD LifeStyle2 Harmony CD 1.67	GR	BR						-13.00 to +8.00	+0.75 to +3.50	-4.00	14
HOYALUX ARRAY 11, 13, 15, 17, VL 1.67	GR	BR		XA-GR				-13.00 to +7.50	+0.75 to +3.50	-4.00	14, 17, 19, 21
HOYALUX ARRAY Wrap 11, 13, 15, 17, VL 1.67	GR	BR		XA-GR				-6.00 to +4.00	+0.75 to +3.50	-4.00	14, 17, 19, 21
Amplitude 1.67	GR							-10.00 to +6.00	+0.75 to +3.50	-4.00	19
Amplitude iQ 1.67	GR							-10.00 to +6.00	+0.75 to +3.50	-4.00	19
Amplitude HiDef 1.67	GR	BR						-13.00 to +8.00	+0.75 to +3.50	-4.00	18
Amplitude Mini HiDef 1.67	GR	BR						-13.00 to +8.00	+0.75 to +3.50	-4.00	14
Amplitude BKS	GR	BR		XA-GR				-13.00 to +7.50	+0.75 to +3.50	-4.00	19
Amplitude Mini BKS	GR	BR		XA-GR				-13.00 to +7.50	+0.75 to +3.50	-4.00	14
NIKON OPTICAL											
Nikon SeeMax Power 1.67	GR	BR	GG					-10.00 to + 8.00	+0.75 to +3.50	-6.00	13, 15, 17
Nikon Presio Power 1.67	GR	BR	GG					-14.00 to + 8.00	+0.75 to +3.50	-6.00	13, 15, 17
Nikon DigiLife 1.67	GR	BR	GG					-14.00 to + 8.00	+0.75 to +3.50	-6.00	13, 15, 17
Nikon Presio i Digital 1.67	GR	BR	GG					-12.00 to + 8.00	+0.75 to +3.50	-6.00	15, 17
Nikon Move 1.67	GR	BR	GG					-12.00 to + 5.00	+0.75 to +3.50	-4.00	16, 18
SEIKO OPTICAL PRODUCTS											
Seiko 1.67 Superior	GR	BR		XA-GR				-10.50 to +6.50	+0.50 to +4.00	-5.00	12
Seiko 1.67 Prominent	GR	BR		XA-GR				-10.50 to +6.50	+0.50 to +3.50	-5.00	14
Seiko 1.67 Premier 2	GR	BR		XA-GR				-10.50 to +6.50	+0.50 to +3.50	-4.00	14
Seiko 1.67 Prestige II	GR	BR		XA-GR				-12.50 to +6.50	+0.50 to +3.50	-5.00	14
Seiko 1.67 Prestige II Wide	GR	BR		XA-GR				-12.50 to +6.50	+0.50 to +3.50	-5.00	14
Seiko 1.67 Succeed	GR	BR		XA-GR				-10.50 to +6.50	+0.50 to +3.50	-4.00	17
Seiko 1.67 Succeed Ws	GR	BR		XA-GR				-10.50 to +6.50	+0.50 to +3.50	-4.00	15
Seiko 1.67 Supercede 2	GR	BR		XA-GR				-10.50 to +6.50	+0.50 to +3.50	-4.00	14
Seiko 1.67 Proceed III Super Short	GR	BR						-13.00 to +8.00	+0.75 to +3.25	-4.00	16
Seiko 1.67 Proceed II Short	GR	BR						-13.00 to +8.00	+0.50 to +3.50	-4.00	18
Seiko 1.67 Supernal	GR	BR		XA-GR				-12.50 to +6.50	+0.50 to +3.50	-5.00	14
Seiko 1.67 Surmount	GR	BR		XA-GR				-12.50 to +6.50	+0.50 to +3.50	-5.00	14
Seiko 1.67 Surmount Ws	GR	BR		XA-GR				-12.50 to +6.50	+0.50 to +3.50	-5.00	14
Seiko 1.67 PCWide	GR	BR		XA-GR				-10.50 to +7.50	-1.00, -1.50, -2.00	-5.00	15
SHAMIR INSIGHT INC.											
Shamir Autograph® III SuperLite® 1.67	GR	BR	GG	XA-GR				-16.75 to +9.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Autograph® II SuperLite® 1.67	GR	BR	GG	XA-GR				-16.75 to +8.25	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Autograph® II Attitude® SuperLite® 1.67	GR	BR	GG	XA-GR				-6.00 to +5.00	+0.75 to +3.50	-8.00	15, 18
Shamir Spectrum® SuperLite® 1.67	GR	BR	GG	XA-GR				-16.75 to +8.00	+0.75 to +3.50	-8.00	14, 16, 18
Shamir Spectrum+® SuperLite® 1.67	GR	BR	GG	XA-GR				-16.75 to +8.00	+0.75 to +3.50	-8.00	14, 16, 18
Shamir Element® SuperLite® 1.67	GR	BR	GG					-12.50 to +8.00	+0.75 to +3.50	-8.00	16, 19
Shamir InTouch®	GR	BR	GG	XA-GR				-16.75 to +8.00	+0.75 to +3.50	-8.00	15, 18
Shamir Attitude III Sport	GR	BR	GG	XA-GR				-6.50 to +6.00	+0.75 to +3.50	-8.00	18
Shamir Attitude III Fashion	GR	BR	GG	XA-GR				-6.50 to +6.00	+0.75 to +3.50	-8.00	15, 18
SIGNET ARMORLITE											
DirectTek® Signia® 1.67	GR	BR	GG					-11.00 to +6.00	+0.75 to +3.50	-6.00	16.5
DirectTek® Short Signia® 1.67	GR	BR	GG					-11.00 to +6.00	+0.75 to +3.50	-6.00	13.5
Navigator® FBS Signia® 1.67	GR	BR	GG					-11.00 to +6.00	+0.75 to +3.50	-6.00	18
Navigator® Short FBS Signia® 1.67	GR	BR	GG					-11.00 to +6.00	+0.75 to +3.50	-6.00	14
KODAK Precise® Signia® 1.67 SF	GR							-13.00 to +8.00	+1.00 to +3.00	-4.00	17
KODAK Precise® Short Signia® 1.67 SF	GR							-12.50 to +7.00	+1.00 to +3.00	-4.00	13
KODAK Precise® Digital Signia® 1.67	GR							-12.50 to +7.00	+1.00 to +3.00	-6.00	17
KODAK Precise® Short Digital Signia® 1.67	GR							-12.50 to +7.00	+1.00 to +3.00	-6.00	13
KODAK Precise® PB Signia® 1.67	GR	BR	GG					-11.00 to +6.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB Signia® 1.67	GR	BR	GG					-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique Signia® 1.67	GR	BR	GG	XA-GR				-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD Signia® 1.67	GR	BR	GG	XA-GR				-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK DSII® Signia® 1.67	GR	BR						-11.00 to +6.00	+0.75 to +3.50	-6.00	13
YOUNGER OPTICS											
Younger Image® 1.67 MR-10™	GR							-10.00 to +8.50	+1.00 to +3.50	-4.00	18

Transitions lenses are compatible with all Crizal No-Glare technology

*MFH-Minimum Fitting Height
† The lenses only have near and intermediate zones

PROGRESSIVES

	Transitions Signature			Transitions XTRActive			Transitions VANTAGE	Sphere	Add	Cyl to	MFH
	Gray	Brown	Graphite Green	Gray	Brown	Graphite Green	Gray				
1.74 ULTRA HIGH INDEX											
ESSILOR LENS											
Varilux S 4D [™] Thin&Lite [®] 1.74	GR	BR						-14.00 to +13.00	+0.75 to +4.00	-6.00	14**
Varilux S Fit [™] Thin&Lite [®] 1.74	GR	BR						-14.00 to +13.00	+0.75 to +4.00	-6.00	14**
Varilux S Design [™] Thin&Lite [®] 1.74	GR	BR						-14.00 to +13.00	+0.75 to +4.00	-6.00	14**
Varilux [®] Physio [®] W3+ eyecode Thin&Lite [®] 1.74	GR							-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux [®] Physio [®] W3+ FIT Thin&Lite [®] 1.74	GR							-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux [®] Physio [®] W3+ Thin&Lite [®] 1.74	GR							-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio Enhanced [™] Thin&Lite [®] 1.74	GR							-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio Enhanced eyecode [™] Thin&Lite [®] 1.74	GR							-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio Enhanced FIT [™] Thin&Lite [®] 1.74	GR							-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio DRx [™] Thin&Lite [®] 1.74	GR							-14.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx [™] Thin&Lite [®] 1.74	GR							-14.00 to +13.00	+0.75 to +3.50	-6.00	14
Varilux [®] Physio [®] Thin&Lite [®] 1.74	GR							-18.00 to +9.00	+0.75 to +3.50 & +4.00	-4.00	17
Essilor Ideal Advanced [™] Thin&Lite [®] 1.74	GR							-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal [®] Thin&Lite [®] 1.74	GR							-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal [®] Fit Thin&Lite [®] 1.74	GR							-20.00 to +13.00	+0.75 to +3.50	-6.00	18
Essilor Ideal [®] Sport Thin&Lite [®] 1.74	GR							-20.00 to +13.00	+0.75 to +3.50	-6.00	18
Essilor Ideal [®] Computer Thin&Lite [®] 1.74	GR							-20.00 to +13.00	+0.75 to +3.50	-6.00	13
Essilor Ideal [®] Driving Thin&Lite [®] 1.74	GR							-20.00 to +13.00	+0.75 to +3.50	-6.00	14
NIKON OPTICAL											
Nikon SeeMax Power 1.74	GR	BR	GG					-12.00 to + 8.00	+0.75 to 3.50	-6.00	13, 15, 17
Nikon Presio Power 1.74	GR	BR	GG					-16.00 to +10.00	+0.75 to 3.50	-6.00	13, 15, 17
Nikon DigiLife 1.74	GR	BR	GG					-16.00 to + 8.00	+0.75 to 4.00	-6.00	13, 15, 17
Nikon Presio i Digital 1.74	GR	BR	GG					-15.00 to + 8.00	+0.75 to 3.50	-6.00	15, 17
SEIKO OPTICAL PRODUCTS											
Seiko 1.74 Superior	GR	BR						-12.50 to +6.50	+0.50 to +4.00	-5.00	12
Seiko 1.74 Prominent	GR	BR						-12.50 to +6.50	+0.50 to +3.50	-5.00	14
Seiko 1.74 Premier 2	GR	BR						-12.50 to +6.50	+0.50 to +3.50	-4.00	14
Seiko 1.74 Prestige II	GR	BR						-12.50 to +6.50	+0.50 to +3.50	-5.00	14
Seiko 1.74 Prestige II Wide	GR	BR						-12.50 to +6.50	+0.50 to +3.50	-5.00	14
Seiko 1.74 Succeed	GR	BR						-12.50 to +6.50	+0.50 to +3.50	-4.00	17
Seiko 1.74 Succeed Ws	GR	BR						-12.50 to +6.50	+0.50 to +3.50	-4.00	15
Seiko 1.74 Supercede 2	GR	BR						-12.50 to +6.50	+0.50 to +3.50	-4.00	14
Seiko 1.74 Supernal	GR	BR						-12.50 to +6.50	+0.50 to +3.50	-5.00	14
Seiko 1.74 Surmount	GR	BR						-12.50 to +6.50	+0.50 to +3.50	-5.00	14
Seiko 1.74 Surmount Ws	GR	BR						-12.50 to +6.50	+0.50 to +3.50	-5.00	14
Seiko 1.74 PCWide	GR	BR						-10.50 to +7.50	-1.00, -1.50, -2.00	-5.00	15
SHAMIR INSIGHT INC.											
Shamir Autograph [®] III SuperLite [™] 1.74	GR	BR						-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir Autograph [®] II SuperLite [™] 1.74	GR	BR						-16.00 to +15.50	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir Spectrum [™] SuperLite [™] 1.74	GR	BR						-18.75 to +14.75	+0.75 to +3.50	-9.00	14, 16, 18
Shamir Spectrum+ [™] SuperLite [™] 1.74	GR	BR						-18.75 to +14.75	+0.75 to +3.50	-9.00	14, 16, 18
Shamir InTouch [™]	GR	BR						-18.75 to +15.50	+0.75 to +3.50	-9.00	15, 18
SIGNET ARMORLITE											
KODAK Unique	GR	BR						-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD	GR	BR						-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK DSII [™]	GR	BR						-13.00 to +9.00	+0.75 to +3.50	-6.00	13

**For Varilux S, at fitting height less than 17mm, ADD power is up to +3.50.

Transitions lenses are compatible with all Crizal No-Glare technology

TRANSITIONS EVERYDAY LENSES

*MFH-Minimum Fitting Height



Transitions®
ADAPTIVE SUNWEAR

FIXED IS OUT. CHANGE IS IN.

OUTDOOR LIGHT CHANGES, YOUR SUNWEAR SHOULD TOO

The Transitions® line of adaptive sunglasses, helmet shields and snow goggle lenses automatically lighten and darken, from sunrise to sunset, for an all-day, visually comfortable experience. Our adaptive sunwear technologies make sports and outdoor activities more convenient by eliminating the need to switch between different pairs of sunwear throughout the day. They help reduce glare, increase contrast, and improve distance and depth perception.

www.Transitions.com/Sunwear



TRANSITIONS® DRIVEWEAR® SUN LENSES

For those who drive, Transitions® Drivewear® sun lenses provide a safer, more comfortable driving experience by combining polarization, color optimization and photochromic technology to enhance a driver's vision by adapting to and managing the ever-changing daytime light & weather conditions faced by drivers.



Overcast
In low light or overcast conditions, the lenses are a green/yellow color that provide high contrast and minimize glare.



Behind the windshield: Sunny
Behind the windshield, the lenses activate to a copper color enhancing color recognition and depth perception.



Outdoors: Bright light
In bright outdoor light, the lenses activate to a dark red-brown filtering excess light to provide maximum comfort.

	COLOR	MATERIALS	SPHERE	ADD	CYL TO	MFH*
Transitions® Drivewear® Adaptive Sun Lenses						
Transitions Drivewear SFSV HC	Polarized, Olive Green to Copper to Brown	1.50	-9.00 to +7.50		-4.00	
Transitions Drivewear® Image®	Polarized, Olive Green to Copper to Brown	1.50	-8.00 to +6.00	+1.00 to +3.00	-4.00	
Transitions Drivewear SF FT28 HC	Polarized, Olive Green to Copper to Brown	1.50	-8.00 to +6.00	+1.00 to +3.00	-4.00	18
Transitions Drivewear SFSV HC	Polarized, Olive Green to Copper to Brown	Polycarbonate	-8.00 to +6.00		-4.00	
Transitions Drivewear® Image®	Polarized, Olive Green to Copper to Brown	Polycarbonate	-9.00 to +7.00	+1.00 to +3.00	-4.00	18
Transitions Drivewear Finished Plano HC/Backside AR	Polarized, Olive Green to Copper to Brown	Trilogy®	Plano 6 base			
Transitions Drivewear SFSV HC	Polarized, Olive Green to Copper to Brown	Trilogy®	-7.00 to +7.50		-4.00	
Transitions® Drivewear® Adaptive Sun Lenses (DIGITAL)						
Camber®	Polarized, Olive Green to Copper to Brown	Polycarbonate	-10.00 to +6.25		-4.00	

For more information visit DriveWearLens.com

ADAPTIVE SUNGLASSES

OAKLEY® TRANSITIONS® SUNGLASSES

CLEAR BLACK IRIDIUM

Oakley offers both prescription and plano Oakley® Transitions® sun lenses and sunglasses, designed to activate to darker, neutral tints to allow athletes to maintain peak performance.

CLEAR
Clear in low light and at night.



DARK GRAY
Neutral, truer color perception in a wide variety of light conditions.

	COLOR	MATERIALS	SPHERE	ADD	CYL TO	MFH*
OAKLEY®						
Oakley® Transitions Clear Black Iridium® sunglasses	Light Gray Iridium to Dark Gray	Oakley Plutonite®	Plano			
Oakley® Transitions SV lenses	Clear to Gray, Clear to Brown, Emerald Iridium®, Black Iridium® and Gold Iridium®	Oakley Plutonite®	-3.25 to +2.00		-3.00	
Oakley® Transitions PAL lenses	Clear to Gray, Clear to Brown, Emerald Iridium®, Black Iridium® and Gold Iridium®	Oakley Plutonite®	-3.25 to +2.00	+1.00 to +3.00	-3.00	

Rx sun lenses available through the Oakley lab. Visit Oakley.com or call 800-403-7449 for more info.

NIKE® MAX TRANSITIONS® SUNGLASSES

GOLF TINT

Engineered for athletes, Nike Max Optics delivers an advanced, high-precision solution to the world of sports. Together with Transitions technology, these lenses adapt to different athletic environments and changing light conditions.

VIOLET
Improves contour recognition on the greens and increases ball pop.



DARK PURPLE
Specifically tuned to the game of golf, enhances critical details on the fairway and green.

OUTDOOR TINT

GREEN
Brightens the shadows and increases contrast.



DARK GREEN
Enhances the visual spectrum found in natural environments.

SPEED TINT

RED
Designed to eliminate road glare.



DARK RED
Allows calming red light to enter the eye which helps during endurance sports.

	COLOR	MATERIALS	SPHERE	ADD	CYL TO	MFH*
NIKE® VISION						
Nike® MAX Transitions® Golf Tint sunglasses	Violet to Dark Purple	Polycarbonate	Plano			
Nike® MAX Transitions® Outdoor Tint sunglasses	Green to Dark Green	Polycarbonate	Plano			
Nike® MAX Transitions® Speed Tint sunglasses	Red to Dark Red	Polycarbonate	Plano			

Plano sun lenses available in Nike athletic frames and as accessory lenses for select Nike frames. Contact your Marchon sales rep or visit NikeVision.com/Transitions for more information.

ESS CROSSBOW® PHOTOCROMIC™ EYESHIELD

The ESS Crossbow® Photochromic™ eyeshield combines ballistic standards with Transitions Optical's light adapting technology to redefine performance and utility for eye protection. ClearZone™ FlowCoat technology minimizes scratches and Tri-Tech™ Fit enhances stability and comfort for any outdoor activity.

CLEAR
Clear in low light and at night.



DARK GRAY
Neutral dark gray in sunny conditions for comfort and decreased glare.

	COLOR	MATERIALS	SPHERE	ADD	CYL TO	MFH*
ESS EYEWEAR						
ESS Crossbow Photochromic	Clear to Dark Gray Ballistics Grade Eyeshield	Ess Crossbow One	Plano			

For more information visit www.esseyepro.com/Crossbow_211_category.html.



ADAPTIVE SHIELDS

TRANSITIONS® ADAPTIVE FACE SHIELDS

Motorcycle and bicycle helmet face shields that help you take your riding experience to the next level. Unlike clear or tinted shields, these shields automatically adjust their level of darkness to optimize vision in sun and shade, day and night.

CLEAR
Clear in low light and at night.



DARK GRAY
Reduces light intensity without altering the color of objects.

	SHIELD COLOR	STYLE	AVAILABILITY
BELL® HELMETS			
Bell® Transitions® face shield	Clear to Dark Gray	Motorcycle Helmet	Comes standard on the Qualifier DLX helmet. Available as an accessory shield on Bell's ClickRelease® helmets: Star, Vortex, RS-1, Qualifier and Revolver. Available for the MX-9 Adventure helmet.
Bell® Transitions® face shield	Clear to Dark Gray	Motorcycle Helmet	
Bell® Transitions® cycling shield	Clear to Dark Gray	Bike Helmet	Comes on the Bell Star Pro bicycle helmet through Bell Helmet retailers and BellHelmets.com.
SHOEI® HELMETS			
SHOEI® Transitions® CWR-1 face shield	Clear to Dark Gray	Motorcycle Helmet	Available in CWR-1 Pinlock® face shields for RF-1200 helmets through SHOEI helmet dealers and Revzilla.com.
LAZER® HELMETS			
Lazer Lumino® Transitions® face shield	Clear to Dark Gray	Motorcycle Helmet	Comes standard on Monaco EVO helmet. Sold as an accessory shield for Lazer's Kite, Osprey and Kestrel helmets.
KLIM HELMETS			
KLIM® Transitions® face shield	Clear to Dark Gray	Motorcycle Helmet	Comes standard on TK1200 ILLumino Hi-V12 Helmet



ADAPTIVE GOGGLES

NEW! TRANSITIONS® ADAPTIVE GOGGLES

Eliminate the need to stop and swap with new Transitions® adaptive goggle lens technology, designed to adapt to changing sunlight during snow sports and motocross activities. Transitions adaptive goggle lenses automatically darken in bright sunlight and lighten in cloudy or snowy conditions, giving you a clear visual advantage.

YELLOW, PINK, AMBER, CLEAR
High contrast tints for improved visibility in overcast or low light.



DARK GRAY
Dark gray in sunny conditions for increased comfort.



	GOGGLE COLOR	STYLE	AVAILABILITY
DRAGON® ALLIANCE			
Dragon® Transitions® Boost snow goggle lens	Yellow to Dark Gray Amber to Dark Gray Pink to Dark Gray Clear to Dark Gray	Snow Goggle	Dragon X1 and X1S and NFXS snow goggle frames through select Dragon goggle dealers and on www.dragonalliance.com .
Dragon® Transitions® Boost MX goggle lens	Yellow to Dark Gray Clear to Dark Gray	MX Snow Goggle	Dragon NFXs MOTO goggle frames through select Dragon goggle dealers and on www.dragonalliance.com .