

Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Design](#) / **Blended Bifocals**

Shamir Insights Inc. — Blended Bifocals

																	
Blended Bifocal (digital)														Sphere	Add	Cyl to	MFH
1.50 Standard Index Blended Bifocal (digital)																	
Shamir Duo	GR	BR							GR					-9.00 to +6.00	+0.75 to +4.00	-6.00	15
1.53 Trivex® Blended Bifocal (digital)																	
Shamir Duo	GR	BR							GR					-13.25 to +6.00	+0.75 to +4.00	-6.50	15
1.59 Polycarbonate Blended Bifocal (digital)																	
Shamir Duo	GR	BR							GR					-14.75 to +6.00	+0.75 to +4.00	-7.00	15
1.60 High Index Blended Bifocal (digital)																	
Shamir Duo	GR	BR							GR					-15.00 to +6.00	+0.75 to +4.00	-7.00	15
1.67 Super High Index Blended Bifocal (digital)																	
Shamir Duo	GR	BR							GR					-16.75 to +8.00	+0.75 to +4.00	-8.00	15
1.74 Ultra High Index Blended Bifocal (digital)																	
Shamir Duo	GR	BR												-18.75 to +12.00	+0.75 to +4.00	-9.00	15

TRANSITIONS® LIGHT INTELLIGENT LENSES™



GR
Gray



BR
Brown



GG
Graphite
Green



S
Sapphire



AME
Amethyst



AMB
Amber



E
Emerald



R
Ruby



Available
in style
mirrors

Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.