

Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Design](#) / [Bifocals](#) / [Hoya Vision Care](#)

Hoya Vision Care — Bifocals

HOYA	Transitions® Gen S							Transitions® Signature GEN 8			Transitions® SFST		Transitions® XTRACTIVE®			Transitions® XTRACTIVE® POLARIZED®					
Bifocal																			Sphere	Add	Cyl to
1.50 Standard Index																					
HOYA Vision Care SFST 28									GR	BR									-11.00 to +6.00	+1.00 to +3.00	-4.00
HOYA Vision Care SFST 28 iQ									GR	BR									-8.00 to +6.00	+1.00 to +3.00	-4.00
1.53 TRIVEX® SFST-28																					
HOYA Vision Care SFST 28									GR	BR									-11.00 to +6.00	+1.00 to +3.00	-4.00
HOYA Vision Care SFST 28 iQ									GR	BR									-8.00 to +6.00	+1.00 to +3.00	-4.00
1.59 Polycarbonate SFST-28																					
HOYA Vision Care SFST 28									GR	BR									-11.00 to +6.00	+1.00 to +3.00	-4.00
HOYA Vision Care SFST 28									GR	BR									-8.00 to +6.00	+1.00 to +3.00	-4.00

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.