





Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Design](#) / [Bifocals](#) / [Younger Optics](#)

Younger Optics – Bifocals

																
Bifocal Signature™														Sphere	Add	Cyl to
1.50 Standard Index SFST-28																
Younger Optics SFFT 28								GR						-11.00 to +6.00	+1.00 to +3.00	-4.00

TRANSITIONS® LIGHT INTELLIGENT LENSES™


GR
Gray


BR
Brown


GG
Graphite
Green


S
Sapphire


AME
Amethyst


AMB
Amber


E
Emerald


R
Ruby


Available
in style
mirrors

Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.