

Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Index](#) / **1.50 Index**

1.50 Standard Index

	Transitions Gen [®] S								Transitions [®] XTRACTIVE [™]			Transitions [®] XTRACTIVE [™] POLARIZED [™]			
Finished Single Vision														Sphere	Cyl to
1.50 Standard Index FSV															
Essilor FSV Spherical	GR	BR							GR					-4.00 to +4.00	-2.00
Essilor FSV Aspheric	GR													-4.00 to +4.00	-2.00
Essilor FSV 1.50 Crizal [®] Sapphire [™] HR	GR								GR					-4.00 to +4.00	-2.00
Essilor FSV 1.50 Crizal [®] Rock [™]	GR								GR					-4.00 to +4.00	-2.00
Essilor FSV 1.50 Crizal [®] Easy Pro	GR								GR					-4.00 to +4.00	-2.00

	Transitions Gen [®] S								Transitions [®] XTRACTIVE [™]			Transitions [®] XTRACTIVE [™] POLARIZED [™]			
Semi-Finished Single Vision														Sphere	Cyl to
1.50 Standard Index SFSV															
Essilor Spherical Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.50 to +10.00	-4.00
Essilor Stylistic [®] Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	-6.00
Eyezen [®] +0	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	-6.00
Eyezen [®] +1	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	-6.00
Eyezen [®] +2	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	-6.00
Eyezen [®] +3	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	-6.00
Eyezen [®] +4	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	-6.00
Eyezen [®] Start	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +6.00	-6.00

	Transitions Gen [®] S								Transitions [®] XTRACTIVE [™]			Transitions [®] XTRACTIVE [™] POLARIZED [™]					
Progressives														Sphere	Add	Cyl to	MFH
1.50 Standard Index																	
Varilux [®] XR Track [™]	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux [®] XR Design [™]	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux [®] X 4D [™]	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux [®] X FIT [™]	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux [®] X Design [™]	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux [®] Physio [®] Extensee [™]	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux [®] Physio [®] Extensee [™] Track	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux [®] Physio [®] Extensee [™] Classic Edition	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux [®] Physio [®] W3+ Eyecode	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux [®] Physio [®] W3+ FIT	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux [®] Physio [®] W3+	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort [®] Max Fit	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort [®] Max	GR	BR	GG						GR	BR	GG			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx [™]	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx [™]	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx [™]	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx [™]	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.00	-6.00	14

Varilux Comfort®	GR	BR							GR	BR	GG			-10.00 to +7.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short™	GR	BR							GR	BR	GG			-10.00 to +7.00	+0.75 to +3.00	-6.00	14
Varilux® Stylistic™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ FIT	GR	BR	GG	S	AME	AMB	E	R						-10.00 to +7.50	+0.75 to +3.50	-4.00	17
Essilor Ideal Advanced™ Wrap	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-4.00	17
Essilor Ideal® FIT	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	18
Essilor Ideal® Driving	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	13
Essilor Ideal®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal® Short	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Natural® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.50	-6.00	15
Essilor Accolade™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Brillenz®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Brillenz® Short	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Ovation® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Ovation®	GR	BR												-10.00 to +6.00	+0.75 to +3.50	-4.00	17
Essilor SmallFit® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +7.50	+0.75 to +3.00	-6.00	14

HOYA	Transitions Gen [®] S								Transitions [®] XTRACTIVE [™]			Transitions [®] XTRACTIVE [™] POLARIZED [®]			
Finished Single Vision														Sphere	Cyl to
1.50 Standard Index FSV															
HOYA FSV	GR													-4.00 to +4.00	-2.00

HOYA	Transitions Gen [®] S								Transitions [®] XTRACTIVE [™]			Transitions [®] XTRACTIVE [™] POLARIZED [®]			
Semi-Finished Single Vision														Sphere	Cyl to
1.50 Standard Index SFSV															
HOYA SFSV 1.50	GR	BR							GR	BR				-8.00 to +9.00	-4.00
HOYA MySV 150	GR	BR							GR	BR		GR		-8.00 to +9.00	-4.00
HOYA SV iQ 1.50	GR	BR							GR	BR		GR		-8.00 to +9.00	-4.00
Sync III	GR	BR							GR	BR		GR		-10.00 to +6.00	-4.00

HOYA	Transitions Gen [®] S								Transitions [®] XTRACTIVE [™]			Transitions [®] XTRACTIVE [™] POLARIZED [®]					
Progressives														Sphere	Add	Cyl to	MFH
1.50 Standard Index																	
HOYALUX iD LifeStyle®4 Indoor 1.50	GR!	BR!												-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Outdoor 1.50	GR!	BR!												-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Urban 1.50	GR!	BR!												-8.00 to +6.00	+0.75 to +4.00	-6.00	14
Amp HD3 1.50	GR!	BR!												-8.00 to +6.00	+0.75 to +4.00	-6.00	18
Amp Mini HD3 1.50	GR!	BR!												-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX ARRAY™ 11, 13, 15, 17, VL 1.50	GR	BR							GR	BR		GR		-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL 1.50	GR	BR							GR	BR		GR		-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ 2 11, 13, 15, 17, VL 1.50	GR	BR							GR	BR		GR		-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ 2 Wrap 11, 13, 15, 17, VL 1.50	GR	BR							GR	BR		GR		-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS 1.50	GR	BR							GR	BR		GR		-10.00 to +3.75	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS 1.50	GR	BR							GR	BR		GR		-10.00 to +3.75	+0.75 to +4.00	-6.00	14

															
Semi-Finished Single Vision														Sphere	Cyl to
1.50 Standard Index SFSV															
Signetek Processed SV	GR	BR	GG	S	AME	AMB	E	R	GR					-8.00 to +6.00	-6.00
KODAK Digital SV	GR	BR	GG	S	AME	AMB	E	R	GR					-8.00 to +6.00	-6.00
KODAK PowerUp Lens (.40/.66/1.00D)	GR	BR												-8.00 to +6.00	-6.00
Crossbows Custom SV	GR	BR												-8.00 to +6.00	-6.00

																
Bifocal Signature™													Sphere	Add	Cyl to	
1.50 Standard Index SFST-28																
Signetek™ Processed 1.50 FT-28									GR				-10.00 to +6.00	+1.00 to +3.00	-6.00	

	Transitions Gen [®] S								Transitions [®] XTRACTIVE [™]			Transitions [®] XTRACTIVE [™] PC-ARIZEP [™]					
Progressives													Sphere	Add	Cyl to	MFH	
1.50 Standard Index																	
KODAK Precise [®] 1.50 SF	GR												-9.00 to +5.00	+1.00 to +3.00	-4.00	17	
KODAK Precise [®] Digital 1.50	GR												-8.00 to +5.00	+1.00 to +3.00	-6.00	17	
KODAK Precise [®] Plus 1.50	GR	BR	GG	S	AME	AMB	E	R	GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	17	
KODAK Precise [®] Short Plus 1.50	GR	BR	GG	S	AME	AMB	E	R	GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
KODAK Precise [®] PB 1.50	GR	BR							GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	17	
KODAK Precise [®] Short PB 1.50	GR	BR							GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
KODAK Unique [™] 1.50	GR	BR	GG	S	AME	AMB	E	R	GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
KODAK Unique [™] HD 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
KODAK Unique DRO [™] 1.50	GR	BR	GG	S	AME	AMB	E	R	GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
KODAK Unique DRO [™] HD 1.50	GR	BR	GG	S	AME	AMB	E	R	GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
KODAK Unique Infinite [™] 1.50	GR	BR	GG	S	AME	AMB	E	R	GR	BR		GR	-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
KODAK Unique Infinite [™] HD 1.50	GR	BR	GG	S	AME	AMB	E	R	GR	BR		GR	-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
KODAK Easy 14 1.50	GR	BR	GG	S	AME	AMB	E	R					-8.00 to +4.00	+0.75 to +3.50	-6.00	14	
KODAK Easy 18 1.50	GR	BR	GG	S	AME	AMB	E	R					-8.00 to +4.00	+0.75 to +3.50	-6.00	18	
DirecTek 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	16.5	
DirecTek Short 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	13.5	
Navigator FBS 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	18	
Navigator Short FBS 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	14	
Crossbows Custom U 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
Crossbows Custom V 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
Crossbows Basic 13 Progressive 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
Crossbows Basic 17 Progressive 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	17	
Crossbows Junior	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	13	

															
Semi-Finished Single Vision														Sphere	Cyl to
1.50 Standard Index SFSV															
TRN 1.50 SV (72mm) UC/HC DIA.72									GR	BR				Base 0.50, 2.25, 4.25, 5.25, 6.25, 8.25	
TRN 1.50 DS (76mm) UC/HC DIA.72	GR	BR							GR	BR		GR		Base 0.50, 8.75, 10.50	
TRN 1.50 DS (76mm) UC/HC DIA.76	GR	BR							GR	BR		GR		Base 2.00, 3.50, 4.25, 5.25, 6.25, 7.25	

	Transitions [®] Gen [®] S								Transitions [®] XTRACTIVE [®]			Transitions [®] XTRACTIVE [®] POLARIZED [®]		
Semi-Finished Single Vision													Sphere	Cyl to
1.50 Standard Index SFSV														
Unity SVx	GR	BR	GG	S	AME	AMB	E	R	GR			GR	-12.25 to +6.00	-6.00
Unity SVxtra	GR	BR	GG	S	AME	AMB	E	R	GR			GR	-12.25 to +6.00	-6.00
Unity SVxtreme	GR	BR	GG	S	AME	AMB	E	R	GR			GR	-6.00 to +6.00	-6.00
Unity Relieve 50	GR	BR	GG	S	AME	AMB	E	R					-12.25 to +6.00	-6.00
Unity Relieve 70	GR	BR	GG	S	AME	AMB	E	R					-12.25 to +6.00	-6.00

																	
Progressives													Sphere	Add	Cyl to	MFH	
1.50 Standard Index																	
Unity Via II	GR	BR	GG	S	AME	AMB	E	R	GR				GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Plus II	GR	BR	GG	S	AME	AMB	E	R	GR				GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite VR	GR	BR	GG	S	AME	AMB	E	R	GR				GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II	GR	BR	GG	S	AME	AMB	E	R	GR				GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II 12	GR	BR	GG	S	AME	AMB	E	R	GR				GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II 14	GR	BR	GG	S	AME	AMB	E	R	GR				GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	14 mm
Unity Via Elite II 16	GR	BR	GG	S	AME	AMB	E	R	GR				GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	16 mm
Unity Via Elite II 18	GR	BR	GG	S	AME	AMB	E	R	GR				GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	18 mm
Unity Via Mobile II	GR	BR	GG	S	AME	AMB	E	R	GR				GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Wrap II	GR	BR	GG	S	AME	AMB	E	R	GR				GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Ethos	GR	BR	GG	S	AME	AMB	E	R	GR				GR	-9.00 to +4.50	+0.75 to +3.50	up to -4.00	14mm
Ethos Plus	GR	BR	GG	S	AME	AMB	E	R	GR				GR	-9.00 to +4.50	+0.75 to +3.50	up to -4.00	14mm

Unity Lenses – Near Variable Focus

																		
Near Variable Focus														Sphere	Add	Cyl to	MFH	
1.50 Standard Index																		
Unity OfficePro 5	GR	BR	GG	S	AME	AMB	E	R						-7.00 to +6.00	+1.00 to +3.50	Up to -6	14 mm	
Unity Office Pro 10	GR	BR	GG	S	AME	AMB	E	R						-7.00 to +6.00	+1.00 to +3.50	Up to -6	14 mm	

YOUNGER OPTICS  The Optical Lens Innovators	Transitions[®] Gen#S									Transitions[®] XTRACTIVE[™]			Transitions[®] XTRACTIVE[™] POLARIZED [®]				
Finished Single Vision																Sphere	Cyl to
1.50 Standard Index FSV																	
Younger Optics FSV AR Super Hydrophobic	GR	BR													-4.00 to +4.00	-2.00	

																		
Semi-Finished Single Vision																	Sphere	Cyl to
1.50 Standard Index SFSV																		
Younger Optics SFSV (76mm)	GR	BR	GG	SA	AME	AMB	E	R	GR	BR							-9.00 to +6.00	-4.00
Younger Optics SFSV DIGITAL													GR				-10.00 to +8.50	-4.00

YOUNGER OPTICS The Optical Lens Innovators	Transitions Gen#S						Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED				
Digital													Sphere	Add
1.50 Standard Index														

