










Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Index](#) / **1.74 Index**

1.74 Ultra High Index

																
	Semi-Finished Single Vision											Sphere	Cyl to			
1.74 Ultra High Index Single Vision																
Essilor Spherical Single Vision Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	-4.00	
Essilor Eyecode™ SV Thin&Lite® 1.74	GR	BR												-16.00 to +13.00	-6.00	
Essilor FIT™ SV Thin&Lite® 1.74	GR	BR												-16.00 to +13.00	-6.00	
Essilor 360™ SV Thin&Lite® 1.74	GR	BR												-16.00 to +13.00	-6.00	
Eyezen® +0 Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	-6.00	
Eyezen® +1 Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	-6.00	
Eyezen® +2 Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	-6.00	
Eyezen® +3 Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	-6.00	
Eyezen® +4 Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	-6.00	
Eyezen® Start Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	-6.00	

																		
	Progressives											Sphere	Add	Cyl to	MFH			
1.74 Ultra High Index																		
Varilux® XR Track™	GR	BR												-14.00 to +13.00	+0.75 to +4.00	-6.00	NVB Measurement	
Varilux® XR Design™	GR	BR												-14.00 to +13.00	+0.75 to +4.00	-6.00	14	
Varilux® X 4D™ Thin&Lite® 1.74	GR	BR												-14.00 to +13.00	+0.75 to +4.00	-6.00	14	
Varilux® X FIT™ Thin&Lite® 1.74	GR	BR												-14.00 to +13.00	+0.75 to +4.00	-6.00	14	
Varilux® X Design™ Thin&Lite® 1.74	GR	BR												-14.00 to +13.00	+0.75 to +4.00	-6.00	14	
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.74	GR	BR												-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14	
Varilux® Physio® W3+ FIT Thin&Lite® 1.74	GR	BR												-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14	
Varilux® Physio® W3+ Thin&Lite® 1.74	GR	BR												-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14	
Varilux Physio DRx™ Thin&Lite® 1.74	GR	BR												-14.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	17	
Varilux Physio Short DRx™ Thin&Lite® 1.74	GR	BR												-14.00 to +13.00	+0.75 to +3.50	-6.00	14	
Varilux Comfort DRx™ Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	+0.75 to +3.50	-6.00	17	
Varilux Comfort Short DRx™ Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	+0.75 to +3.50	-6.00	17	
Essilor Ideal Advanced™ FIT Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	+0.75 to +3.50	-6.00	17	
Essilor Ideal Advanced™ Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	+0.75 to +3.50	-6.00	17	
Essilor Ideal® FIT Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	+0.75 to +3.50	-6.00	17	
Essilor Ideal® Sport Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	+0.75 to +3.50	-6.00	18	
Essilor Ideal® Computer Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	+0.75 to +3.50	-6.00	18	
Essilor Ideal® Driving Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	+0.75 to +3.50	-6.00	18	
Essilor Ideal® Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	+0.75 to +3.50	-6.00	13	
Ovation® Digital Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	+0.75 to +3.50	-6.00	14	
Essilor Accolade Freedom™ Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	+0.75 to +3.50	-6.00	15	
Essilor Accolade® Thin&Lite® 1.74	GR	BR												-14.00 to +9.00	+0.75 to +3.50	-6.00	17	

	Transitions Gen#S		Transitions XTRACTIVE		Transitions XTRACTIVE POLARIZED								
Blended Bifocal (digital)										Sphere	Add	Cyl to	MFH
1.74 Ultra High Index Blended Bifocal (digital)													
Shamir Duo	GR	BR								-18.75 to +12.00	+0.75 to +4.00	-9.00	15

	Transitions Gen#S		Transitions XTRACTIVE		Transitions XTRACTIVE POLARIZED								
Progressives										Sphere	Add	Cyl to	MFH
1.74 Ultra High Index													
Shamir Autograph Intelligence™	GR	BR								-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR								-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR								-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir GENESIS HD	GR	BR								-6.00 to +6.00	+0.75 to +4.00	-4.00	16
Shamir Spectrum+™	GR	BR								-18.75 to +11.50	+0.75 to +3.50	-9.00	14, 16, 18
Shamir InTouch™	GR	BR								-18.75 to +15.50	+0.75 to +3.50	-9.00	15, 18

	Transitions Gen#S		Transitions XTRACTIVE		Transitions XTRACTIVE POLARIZED								
Semi-Finished Single Vision										Sphere	Cyl to		
1.74 Ultra High Index SFSV													
Signetek™ Processed SV	GR	BR								-13.00 to +9.00	-6.00		
KODAK Digital SV	GR	BR								-13.00 to +9.00	-6.00		

	Transitions Gen#S		Transitions XTRACTIVE		Transitions XTRACTIVE POLARIZED								
Semi-Finished Single Vision										Sphere	Cyl to		
1.74 ultra High Index Single Vision													
TOG 1.74 LeafECO SV (75mm) UC/HC DIA.65	GR	BR								Base 11.50, 12.50, 13.50, 14.50, 15.50, 16.50			
TOG 1.74 LeafECO SV (75mm) UC/HC DIA.70	GR	BR								Base 0.50, 1.50, 2.25, 6.75, 8.50, 10.25			
TOG 1.74 LeafECO SV (75mm) UC/HC DIA.75	GR	BR								Base 3.50, 5.00			

	Transitions Gen#S		Transitions XTRACTIVE		Transitions XTRACTIVE POLARIZED								
Semi-Finished Single Vision										Sphere	Cyl to		
1.74 Ultra High Index Single Vision													
Younger Optics Single Vision 1.74	GR	BR	GG	SA	AME	AMB	E	R	GR	BR		-18.00 to +12.00	-4.00

	Transitions Gen#S		Transitions XTRACTIVE		Transitions XTRACTIVE POLARIZED								
Digital										Sphere	Add		
1.74 Ultra High Index													
Camber™	GR	BR							GR	BR		-19.00 to +9.00	-4.00

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, *Transitions Signature*, and *XTRActive* are registered trademarks, and *XTRActive Polarized*, the *Transitions* logo and *Transitions Light Intelligent Lenses* are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

Unless indicated otherwise, all registered trademarks and trademarks are the property of Essilor International and/or its subsidiaries in the United States and in other countries.

All other trademarks are the property of their respective owners.