

Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Suppliers](#) / **Unity Lenses**

Unity Lenses – Entire Product Line

unity	Transitions Gen#S								Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED		Sphere	Cyl to
	Semi-Finished Single Vision														
1.50 Standard Index															
Unity SVx	GR	BR							GR		GR	-12.25 to +6.00	-6.00		
Unity SVxtra	GR	BR							GR		GR	-12.25 to +6.00	-6.00		
Unity SVxtreme	GR	BR							GR		GR	-6.00 to +6.00	-6.00		
1.53 TRIVEX® Single Vision															
Unity SVx	GR	BR							GR			-14.00 to +6.00	-6.00		
Unity SVxtra	GR	BR							GR			-14.00 to +6.00	-6.00		
Unity SVxtreme	GR	BR							GR			-6.50 to +6.00	-6.00		
1.59 Polycarbonate Single Vision															
Unity SVx	GR	BR							GR		GR	-14.00 to +6.00	-6.00		
Unity SVxtra	GR	BR							GR		GR	-14.00 to +6.00	-6.00		
Unity SVxtreme	GR	BR							GR		GR	-6.50 to +6.00	-6.00		
1.60 High Index Single Vision															
Unity SVx	GR	BR										-14.00 to +6.00	-6.00		
Unity SVxtra	GR	BR										-14.00 to +6.00	-6.00		
Unity SVxtreme	GR	BR										-6.50 to +6.00	-6.00		
1.67 Super High Index Single Vision															
Unity SVx	GR	BR							GR		GR	-16.00 to +8.00	-6.00		
Unity SVxtra	GR	BR							GR		GR	-16.00 to +8.00	-6.00		
Unity SVxtreme	GR	BR							GR		GR	-6.00 to +6.00	-6.00		

unity	Transitions Gen#S								Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED		Sphere	Add	Cyl to	MFH
	Progressives																
1.50 Standard Index																	
Unity Via II	GR	BR							GR		GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12		
Unity Via Plus II	GR	BR							GR		GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12		
Unity Via Elite II	GR	BR							GR		GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12		
Unity Via Mobile II	GR	BR							GR		GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12		
Unity Via Wrap II	GR	BR							GR		GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12		
Ethos	GR	BR							GR		GR	-9.00 to +4.50	0.75 to +3.50	-4.00	14		
Ethos Plus	GR	BR							GR		GR	-9.00 to +4.50	0.75 to +3.50	-4.00	14		
1.53 TRIVEX®																	
Unity Via II	GR	BR							GR			-10.00 to +6.00	+0.50 to +4.50	-6.00	12		
Unity Via Plus II	GR	BR							GR			-10.00 to +6.00	+0.50 to +4.50	-6.00	12		
Unity Via Elite II	GR	BR							GR			-10.00 to +6.00	+0.50 to +4.50	-6.00	12		
Unity Via Mobile II	GR	BR							GR			-10.00 to +6.00	+0.50 to +4.50	-6.00	12		
Unity Via Wrap II	GR	BR							GR			-5.00 to +4.00	+0.50 to +4.50	-6.00	12		
Ethos	GR	BR							GR			-9.00 to +6.00	0.75 to +3.50	-4.00	14		
Ethos Plus	GR	BR							GR			-9.00 to +6.00	0.75 to +3.50	-4.00	14		
1.59 Polycarbonate																	
Unity Via II	GR	BR							GR		GR	-10.00 to +6.50	+0.50 to +4.50	-6.00	12		
Unity Via Plus II	GR	BR							GR		GR	-10.00 to +6.50	+0.50 to +4.50	-6.00	12		
Unity Via Elite II	GR	BR							GR		GR	-10.00 to +6.50	+0.50 to +4.50	-6.00	12		

Unity Via Mobile II	GR	BR							GR			GR	-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR							GR			GR	-6.50 to +6.00	+0.50 to +4.50	-6.00	12
Ethos	GR	BR							GR			GR	-9.00 to +6.00	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR							GR			GR	-9.00 to +6.00	0.75 to +3.50	-4.00	14
1.60 High Index																
Unity Via II	GR	BR											-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR											-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR											-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR											-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR											-6.50 to +5.25	+0.50 to +4.50	-6.00	12
Ethos	GR	BR											-9.00 to +5.00	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR											-9.00 to +5.00	0.75 to +3.50	-4.00	14
1.67 Super High Index																
Unity Via II	GR	BR							GR			GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR							GR			GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR							GR			GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR							GR			GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR							GR			GR	-6.00 to +5.00	+0.50 to +4.50	-6.00	12
Ethos	GR	BR							GR			GR	-12.00 to +8.00	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR							GR			GR	-12.00 to +8.00	0.75 to +3.50	-4.00	14

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.