

Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Suppliers](#) / [Younger Optics](#)

Younger Optics — Entire Product Line

 The Optical Lens Innovators	Transitions Gen [®] S								Transitions XTRACTIVE [®]			Transitions XTRACTIVE [®] POLARIZED [®]				
Finished Single Vision															Sphere	Cyl to
1.50 Standard Index FSV																
Younger Optics FSV AR Super Hydrophobic	GR	BR													-4.00 to +4.00	-2.00
1.59 Polycarbonate Finished Single Vision																
Younger FSV Poly AR Super Hydrophobic	GR	BR													-6.00 to +4.00	-2.00
Younger Optics FSV Poly									GR	BR					-6.00 to +4.00	-2.00

 The Optical Lens Innovators	Transitions Gen [®] S								Transitions XTRACTIVE [®]			Transitions XTRACTIVE [®] POLARIZED [®]				
Semi-Finished Single Vision															Sphere	Cyl to
1.50 Standard Index SFSV																
Younger Optics SFSV (76mm)	GR	BR	GG	SA	AME	AMB	E	R	GR	BR					-9.00 to +6.00	-4.00
Younger Optics SFSV DIGITAL												GR			-10.00 to +8.50	-4.00
1.53 TRIVEX® Single Vision																
Younger Optics Trilogy® Spherical SV	GR	BR							GR	BR		GR			-8.00 to +7.50	-4.00
1.59 Polycarbonate Single Vision																
Younger Optics SV	GR	BR	GG	SA	AME	AMB	E	R	GR	BR		GR			-9.00 to +6.00	-4.00
1.60 High Index Single Vision																
Younger Optics Single Vision 1.60	GR	BR	GG	SA	AME	AMB	E	R							-11.00 to +8.00	-4.00
1.67 Super High Index Single Vision																
Younger Optics Single Vision 1.67 MR-10™	GR	BR	GG	SA	AME	AMB	E	R	GR	BR		GR			-11.00 to +7.00	-4.00
1.74 Ultra High Index Single Vision																
Younger Optics Single Vision 1.74	GR	BR	GG	SA	AME	AMB	E	R	GR	BR					-18.00 to +12.00	-4.00

 The Optical Lens Innovators	Transitions Gen [®] S								Transitions XTRACTIVE [®]			Transitions XTRACTIVE [®] POLARIZED [®]				
Digital															Sphere	Add
Camber™									GR	BR					-11.25 to +6.25	-4.00
1.53 TRIVEX®																
Camber™	GR	BR							GR	BR		GR			-11.00 to +6.25	-4.00
1.59 Polycarbonate																
Camber™	GR	BR							GR	BR		GR			-10.00 to +6.25	-4.00
1.60 High Index																
Camber™	GR	BR													-11.50 to +7.25	-4.00
1.67 Super High Index																
Camber™	GR	BR							GR	BR		GR			-14.00 to +8.00	-4.00
1.74 Ultra High Index																
Camber™	GR	BR							GR	BR					-19.00 to +9.00	-4.00

 The Optical Lens Innovators	Transitions Gen [®] S								Transitions XTRACTIVE [®]			Transitions XTRACTIVE [®] POLARIZED [®]					
Bifocal Signature™															Sphere	Add	Cyl to
1.50 Standard Index SFST-28																	
Younger Optics SFFT 28									GR						-11.00 to +6.00	+1.00 to +3.00	-4.00

														
Progressives											Sphere	Add	Cyl to	MFH
1.53 TRIVEX®														
Younger Trilogy® Image®	GR										-8.00 to +8.00	+1.00 to +3.00	-4.00	18
1.59 index polycarbonate														
Younger Image®	GR	BR						GR	BR		-9.00 to +7.00	+1.00 to +3.00	-4.00	18
1.67 Super High Index														
Younger Image® 1.67 MR-10™	GR	BR									-10.00 to +8.50	+1.00 to +3.50	-4.00	18

TRANSITIONS® LIGHT INTELLIGENT LENSES™



GR
Gray



BR
Brown



GG
Graphite
Green



S
Sapphire



AME
Amethyst



AMB
Amber



E
Emerald



R
Ruby



Available
in style
mirrors

Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.