

Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Design](#) / [Progressives](#) / [Hoya Vision Care](#)

Hoya Vision Care — Progressives

HOYA	Transitions® gen#S								Transitions® XTRACTIVE®			Transitions® XTRACTIVE® POLARIZED®						
Progressives															Sphere	Add	Cyl to	MFH
1.50 Standard Index																		
HOYALUX iD LifeStyle®4 Indoor 1.50	GR!	BR!													-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Outdoor 1.50	GR!	BR!													-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Urban 1.50	GR!	BR!													-8.00 to +6.00	+0.75 to +4.00	-6.00	14
Amp HD3 1.50	GR!	BR!													-8.00 to +6.00	+0.75 to +4.00	-6.00	18
Amp Mini HD3 1.50	GR!	BR!													-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX ARRAY™ 11, 13, 15, 17, VL 1.50	GR	BR							GR	BR		GR			-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL 1.50	GR	BR							GR	BR		GR			-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ 2 11, 13, 15, 17, VL 1.50	GR	BR							GR	BR		GR			-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ 2 Wrap 11, 13, 15, 17, VL 1.50	GR	BR							GR	BR		GR			-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS 1.50	GR	BR							GR	BR		GR			-10.00 to +3.75	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS 1.50	GR	BR							GR	BR		GR			-10.00 to +3.75	+0.75 to +4.00	-6.00	14
Summit cd BKS 1.50	GR	BR							GR	BR		GR			-10.00 to +3.75	+0.75 to +4.00	-6.00	14
Summit ecp BKS 1.50	GR	BR							GR	BR		GR			-10.00 to +3.75	+0.75 to +4.00	-6.00	18
Hoya GP Wide 1.50	GR	BR							GR	BR		GR			-10.00 to +3.75	+0.75 to +4.00	-6.00	18
Tact 40 BKS 1.50	GR	BR							GR	BR		GR			-10.00 to +3.75	+0.75 to +4.00	-6.00	18
Tact 60 BKS 1.50	GR	BR							GR	BR		GR			-10.00 to +3.75	+0.75 to +4.00	-6.00	18
1.53 TRIVEX®																		
HOYALUX iD LifeStyle®4 Indoor Phoenix™	GR!	BR!													-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Outdoor Phoenix™	GR!	BR!													-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Urban Phoenix™	GR!	BR!													-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Indoor Phoenix™	GR!	BR!													-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Outdoor Phoenix™	GR!	BR!													-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Urban Phoenix™	GR!	BR!													-8.00 to +6.00	+0.75 to +4.00	-6.00	14
Amp HD3 Phoenix™	GR!	BR!													-8.00 to +6.00	+0.75 to +4.00	-6.00	18
Amp Mini HD3 Phoenix™	GR!	BR!													-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX ARRAY™ 11, 13, 15, 17, VL Phoenix™	GR	BR							GR	BR		GR			-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL Phoenix™	GR	BR							GR	BR		GR			-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13, 15, 17, VL Phoenix™	GR	BR							GR	BR		GR			-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11, 13, 15, 17, VL Phoenix™	GR	BR							GR	BR		GR			-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS Phoenix™	GR	BR							GR	BR		GR			-10.00 to +6.00	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS Phoenix™	GR	BR							GR	BR		GR			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Summit cd BKS Phoenix™	GR	BR							GR	BR		GR			-10.00 to +6.00	+0.75 to +4.00	-6.00	18
Summit ecp BKS Phoenix™	GR	BR							GR	BR		GR			-10.00 to +6.00	+0.75 to +4.00	-6.00	18
Hoya GP Wide Phoenix™	GR	BR							GR	BR		GR			-10.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 40 BKS Phoenix™	GR	BR							GR	BR		GR			-10.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 60 BKS Phoenix™	GR	BR							GR	BR		GR			-10.00 to +6.00	+0.75 to +4.00	-6.00	18
1.59 Polycarbonate																		
HOYALUX iD LifeStyle®4 Indoor Poly	GR!	BR!													-6.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Outdoor Poly	GR!	BR!													-6.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Urban Poly	GR!	BR!													-6.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Indoor Poly	GR!	BR!													-6.00 to +6.00	+0.75 to +4.00	-6.00	14

Tact 40 BKS 1.67	GR	BR							GR	BR		GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	18
Tact 60 BKS 1.67	GR	BR							GR	BR		GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	18
1.74 ultra High Index																
HOYALUX iD LifeStyle®4 Indoor 1.74	GR!	BR!											-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle®4 Outdoor 1.74	GR!	BR!											-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle®4 Urban 1.74	GR!	BR!											-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle 3 Indoor 1.74	GR!	BR!											-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle 3 Outdoor 1.74	GR!	BR!											-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle 3 Urban 1.74	GR!	BR!											-13.00 to +8.00	+0.75 to +3.50	-6.00	14
Amp HD3 1.74	GR!	BR!											-13.00 to +8.00	+0.75 to +3.50	-6.00	14
Amp Mini HD3 1.74	GR!	BR!											-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX ARRAY 11, 13, 15, 17, VL 1.74	GR	BR											-13.00 to +8.00	+0.75 to +3.50	-6.00	14, 17, 19, 21
HOYALUX ARRAY Wrap 11, 13, 15, 17, VL 1.74	GR	BR											-6.00 to +4.00	+0.75 to +3.50	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13, 15, 17, VL 1.74	GR	BR											-13.00 to +8.00	+0.75 to +3.50	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 WRAP 11, 13, 15, 17, VL 1.74	GR	BR											-6.00 to +4.00	+0.75 to +3.50	-6.00	14, 17, 19, 21
Amplitude BKS	GR	BR											-13.00 to +8.00	+0.75 to +3.50	-6.00	19
Amplitude Mini BKS	GR	BR											-13.00 to +8.00	+0.75 to +3.50	-6.00	14
Summit cd BKS 1.74	GR	BR											-13.00 to +8.00	+0.75 to +3.50	-6.00	14
Summit ecp BKS 1.74	GR	BR											-13.00 to +8.00	+0.75 to +3.50	-6.00	19
Hoya GP Wide 1.74	GR	BR											-13.00 to +8.00	+0.75 to +3.50	-6.00	18
Tact 40 BKS 1.74	GR	BR											-13.00 to +8.00	+0.75 to +3.50	-6.00	18
Tact 60 BKS 1.74	GR	BR											-13.00 to +8.00	+0.75 to +3.50	-6.00	18

TRANSITIONS® LIGHT INTELLIGENT LENSES™



GR
Gray



BR
Brown



GG
Graphite
Green



S
Sapphire



AME
Amethyst



AMB
Amber



E
Emerald



R
Ruby



Available
in style
mirrors

Transitions, *Transitions Signature*, and *XTRActive* are registered trademarks, and *XTRActive Polarized*, the *Transitions* logo and *Transitions Light Intelligent Lenses* are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.