


Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Design](#) / **Progressives**

All Suppliers – Progressives

	Transitions® Gen#S								Transitions® XTRACTIVE®			Transitions® PO-AIRLESS®				
Progressives													Sphere	Add	Cyl to	MFH
1.50 Standard Index																
Varilux® XR Track™	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Varilux® Liberty™ 3.0	GR	BR											-10.00 to +7.50	+0.75 to +1.50	-6.00	14
Varilux Comfort®	GR	BR							GR	BR	GG		-10.00 to +7.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short™	GR	BR							GR	BR	GG		-10.00 to +7.00	+0.75 to +3.00	-6.00	14
Varilux® Stylistic™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ FIT	GR	BR	GG	S	AME	AMB	E	R					-10.00 to +7.50	+0.75 to +3.50	-4.00	17
Essilor Ideal Advanced™ Wrap	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-4.00	17
Essilor Ideal® FIT	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	18
Essilor Ideal® Driving	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	13
Essilor Ideal®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal® Short	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Natural® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +7.50	+0.75 to +3.50	-6.00	15
Essilor Accolade™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Brillenz®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Brillenz® Short	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	14

Ovation® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Ovation®	GR	BR											-10.00 to +6.00	+0.75 to +3.50	-4.00	17
Essilor SmallFit® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +7.50	+0.75 to +3.00	-6.00	14

	Transitions Gen#S								Transitions XTRACTIVE®			Transitions XTRACTIVE® POLARIZED®				
Progressives													Sphere	Add	Cyl to	MFH
1.53 TRIVEX®																
Varilux® XR Track™	GR	BR											-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR											-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ TREXA™	GR	BR											-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ TREXA™	GR	BR											-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ TREXA™	GR	BR											-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode TREXA™	GR	BR							GR	BR			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT TREXA™	GR	BR							GR	BR			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ TREXA™	GR	BR							GR	BR			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit TREXA™	GR	BR							GR	BR			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max TREXA™	GR	BR							GR	BR			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™ TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Essilor Ideal Advanced™ FIT TREXA™	GR	BR	GG										-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	18
Essilor Ideal® Computer TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	13
Essilor Ideal® Driving TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal® TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal® Short TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Natural® Digital TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	15
Essilor Accolade® TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Digital TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Brillenz® TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Brillenz® Short TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Ovation® Digital TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor SmallFit® Digital TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	+0.75 to +3.00	-6.00	14

	Transitions Gen#S								Transitions XTRACTIVE®			Transitions XTRACTIVE® POLARIZED®				
Progressives													Sphere	Add	Cyl to	MFH
1.59 Polycarbonate																
Varilux® XR Track™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement

Varilux® XR Design™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Varilux Comfort DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux® Liberty™ 3.0	GR	BR											-10.00 to +8.00	+0.75 to +1.50	-6.00	14
Varilux Comfort® Airwear®	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short™ Airwear®	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +3.00	-6.00	14
Varilux® Stylistic™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ FIT	GR	BR	GG	S	AME	AMB	E	R					-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Wrap	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Essilor Ideal® FIT Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Sport Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Driving Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Ideal® Short Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Natural® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +8.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Ovation® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Ovation® Airwear®	GR	BR											-10.00 to +8.00	+0.75 to +3.50	-4.00	17
Essilor SmallFit® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +8.00	+0.75 to +3.00	-6.00	14

																
Progressives													Sphere	Add	Cyl to	MFH
1.60 High Index																
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Thin&Lite® 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Thin&Lite® 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14

	Transitions Gen [#] S							Transitions [®] XTRACTIVE [™]			Transitions [®] XTRACTIVE [™] POLARIZED [®]					
Progressives													Sphere	Add	Cyl to	MFH
1.67 Super High Index																
Varilux [®] XR Track [™]	GR	BR	GG	S	AME	AMB	E	R					-12.00 to +9.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux [®] XR Design [™]	GR	BR	GG	S	AME	AMB	E	R					-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux [®] X 4D [™] Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R					-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux [®] X FIT [™] Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R					-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux [®] X Design [™] Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R					-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux [®] Physio [®] W3+ Eyecode Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R					-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux [®] Physio [®] W3+ FIT Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R					-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux [®] Physio [®] W3+ Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R					-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort [®] Max Fit Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R					-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Comfort [®] Max Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R					-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx [™] Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R				GR	-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx [™] Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R				GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx [™] Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R				GR	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx [™] Thin&Lite [®] 1.67	GR	BR	GG	S	AME	AMB	E	R				GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Varilux [®] Liberty [™] 3.0	GR	BR											-14.00 to +9.00	+0.75 to +1.50	-6.00	14
Varilux Comfort [®] Thin&Lite [®] 1.67	GR	BR											-14.00 to +9.00	+0.75 to +3.00	-6.00	14

Varilux Comfort® Short™ Thin&Lite® 1.67	GR	BR														-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Varilux® Stylistic™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		GR			-5.00 to +5.00	+0.75 to +3.50	-3.00	17
Essilor Ideal Advanced™ FIT Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R								-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal Advanced Wrap™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		GR			-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		GR			-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		GR			-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Sport Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		GR			-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		GR			-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Driving Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		GR			-14.00 to +9.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		GR			-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor Ideal® Short	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		GR			-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor Natural® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG					-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG					-14.00 to +9.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG					-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG					-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG					-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Brillenz® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		GR			-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		GR			-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Ovation® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG					-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Ovation® Thin&Lite® 1.67	GR	BR														-12.00 to +8.00	+0.75 to +3.50	-4.00	17
Essilor SmallFit® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG					-14.00 to +9.00	+0.75 to +3.00	-6.00	14

																	
Progressives														Sphere	Add	Cyl to	MFH
1.74 Ultra High Index																	
Varilux® XR Track™	GR	BR												-14.00 to +13.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR												-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ Thin&Lite® 1.74	GR	BR												-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ Thin&Lite® 1.74	GR	BR												-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ Thin&Lite® 1.74	GR	BR												-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.74	GR	BR												-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Thin&Lite® 1.74	GR	BR												-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Thin&Lite® 1.74	GR	BR												-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio DRx™ Thin&Lite® 1.74	GR	BR												-14.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Thin&Lite® 1.74	GR	BR												-14.00 to +13.00	+0.75 to +3.50	-6.00	14
Varilux Comfort DRx™ Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ FIT Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Computer Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Driving Thin&Lite® 1.74	GR	BR												-20.00 to +13.00	+0.75 to +3.50	-6.00	18

	Transitions® Gen ⁸ S								Transitions® XTRACTIVE™			Transitions® XTRACTIVE™ POLARISEZ™						
Progressives														Sphere	Add	Cyl to	MFH	
1.50 Standard Index																		
Nikon SeeMax Ultimate	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		GR	-10.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17	
Nikon Presio Power	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		GR	-10.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17	
Nikon Presio i Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		GR	-10.00 to +8.00	+0.50 to +4.00	-6.00	15, 17	
Nikon Move	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		GR	-10.00 to +8.00	+0.50 to	-6.00	15, 17	


Shamir Intouch™	GR	BR						GR			GR	-13.25 to +6.00	+0.75 to +3.50	-6.50	14, 16, 18
Shamir Golf™	GR	BR						GR			GR	-5.00 to +4.00	+1.00 to +3.50	-6.50	19
Shamir Spectrum+™	GR	BR						GR			GR	-13.25 to +6.00	+0.75 to +3.50	-6.50	14, 16, 18
Shamir Element™	GR	BR						GR			GR	-13.25 to +6.00	+0.75 to +3.50	-6.50	16, 19
Shamir GENESIS HD	GR	BR						GR			GR	-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir Computer™	GR	BR						GR			GR	-6.00 to +8.00	+0.75 to +2.00	-6.50	14 Lower, 10 Upper
Shamir Workspace™	GR	BR						GR			GR	-6.00 to +8.00	+0.75 to +2.50	-6.50	14 Lower, 10 Upper
Shamir Attitude III Sport	GR	BR						GR			GR	-6.50 to +6.00	+0.75 to +3.50	-6.50	18
Shamir Attitude III Fashion	GR	BR						GR			GR	-6.50 to +6.00	+0.75 to +3.50	-6.50	15, 18
1.59 Polycarbonate															
Shamir Autograph Intelligence™	GR	BR						GR			GR	-14.75 to +7.50	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR						GR			GR	-14.75 to +8.00	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR						GR			GR	-14.75to +6.50	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Intouch™	GR	BR						GR			GR	-14.75 to +6.00	+0.75 to +3.50	-7.00	15, 18
Shamir Golf™	GR	BR						GR			GR	-6.00 to +5.25	+1.00 to +3.50	-7.00	19
Shamir Spectrum+™	GR	BR						GR			GR	-14.75 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
Shamir Element™	GR	BR						GR			GR	-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19
Shamir GENESIS HD	GR	BR						GR			GR	-10.00 to +6.00	+0.75 to +3.50	-7.00	16
Shamir Computer™	GR	BR						GR			GR	-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper
Shamir Workspace™	GR	BR						GR			GR	-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper
Shamir Attitude III Sport	GR	BR						GR			GR	-6.50 to +6.00	+0.75 to +3.50	-6.50	18
Shamir Attitude III Fashion	GR	BR						GR			GR	-6.50 to +6.00	+0.75 to +3.50	-6.50	15, 18
1.60 High Index															
Shamir Autograph Intelligence™	GR	BR						GR				-14.75 to +8.75	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR						GR				-15.00 to +8.00	+0.75 to +3.50	-7.00	15,18
Shamir Autograph® II+	GR	BR						GR				-15.00 to +6.00	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Golf™ SuperLite™ 1.6	GR	BR						GR				-6.00 to +5.25	+1.00 to +3.50	-7.00	19
Shamir Spectrum+™	GR	BR						GR				-15.00 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
Shamir Element™	GR	BR										-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19
Shamir GENESIS HD	GR	BR						GR				-6.00 to +6.00	+0.75 to +3.50	-7.00	16
Shamir Computer™	GR	BR						GR				-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper
Shamir Workspace™	GR	BR						GR				-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper
Shamir InTouch™	GR	BR						GR				-15.00 to +6.00	+0.75 to +3.50	-7.00	15, 18
Shamir Attitude III Sport	GR	BR						GR				-6.50 to +6.00	+0.75 to +3.50	-7.00	18
Shamir Attitude III Fashion	GR	BR						GR				-6.50 to +6.00	+0.75 to +3.50	-7.00	15, 18
1.67 Super High Index															
Shamir Autograph Intelligence™	GR	BR						GR			GR	-16.75 to +10.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR						GR			GR	-16.75 to +9.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR						GR			GR	-16.75 to +8.25	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Spectrum+™	GR	BR						GR			GR	-16.75 to +8.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Element™ SuperLite™ 1.67	GR	BR						GR			GR	-12.50 to +8.00	+0.75 to +3.50	-8.00	16, 19
Shamir GENESIS HD	GR	BR						GR			GR	-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir InTouch™	GR	BR						GR			GR	-16.75 to +8.00	+0.75 to +3.50	-8.00	15, 18
Shamir Attitude III Sport	GR	BR						GR			GR	-6.50 to +6.00	+0.75 to +3.50	-8.00	18

Shamir Attitude III Fashion	GR	BR							GR			GR	-6.50 to +6.00	+0.75 to +3.50	-8.00	15, 18
1.74 Ultra High Index																
Shamir Autograph Intelligence™	GR	BR											-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR											-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR											-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir GENESIS HD	GR	BR											-6.00 to +6.00	+0.75 to +4.00	-4.00	16
Shamir Spectrum+™	GR	BR											-18.75 to +11.50	+0.75 to +3.50	-9.00	14, 16, 18
Shamir InTouch™	GR	BR											-18.75 to +15.50	+0.75 to +3.50	-9.00	15, 18

	Transitions Gen#S								Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED					
Progressives													Sphere	Add	Cyl to	MFH	
1.50 Standard Index																	
KODAK Precise® 1.50 SF	GR												-9.00 to +5.00	+1.00 to +3.00	-4.00	17	
KODAK Precise® Digital 1.50	GR												-8.00 to +5.00	+1.00 to +3.00	-6.00	17	
KODAK Precise® Plus 1.50	GR	BR	GG	S	AME	AMB	E	R	GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	17	
KODAK Precise® Short Plus 1.50	GR	BR	GG	S	AME	AMB	E	R	GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
KODAK Precise® PB 1.50	GR	BR							GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	17	
KODAK Precise® Short PB 1.50	GR	BR							GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
KODAK Unique™ 1.50	GR	BR	GG	S	AME	AMB	E	R	GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
KODAK Unique™ HD 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
KODAK Unique DRO™ 1.50	GR	BR	GG	S	AME	AMB	E	R	GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
KODAK Unique DRO™ HD 1.50	GR	BR	GG	S	AME	AMB	E	R	GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
KODAK Unique Infinite™ 1.50	GR	BR	GG	S	AME	AMB	E	R	GR	BR		GR	-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
KODAK Unique Infinite™ HD 1.50	GR	BR	GG	S	AME	AMB	E	R	GR	BR		GR	-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
KODAK Easy 14 1.50	GR	BR	GG	S	AME	AMB	E	R					-8.00 to +4.00	+0.75 to +3.50	-6.00	14	
KODAK Easy 18 1.50	GR	BR	GG	S	AME	AMB	E	R					-8.00 to +4.00	+0.75 to +3.50	-6.00	18	
DirecTek 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	16.5	
DirecTek Short 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	13.5	
Navigator FBS 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	18	
Navigator Short FBS 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	14	
Crossbows Custom U 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
Crossbows Custom V 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
Crossbows Basic 13 Progressive 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
Crossbows Basic 17 Progressive 1.50	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	17	
Crossbows Junior	GR	BR											-8.00 to +4.00	+0.75 to +3.50	-6.00	13	
1.53 TRIVEX®																	
KODAK Precise® PB Trivex®	GR	BR							GR				-10.00 to +4.50	+0.75 to +3.50	-6.00	17	
KODAK Precise® Short PB Trivex®	GR	BR							GR				-10.00 to +4.50	+0.75 to +3.50	-6.00	13	
KODAK Precise® Plus Trivex®	GR	BR							GR				-10.00 to +4.50	+0.75 to +3.50	-6.00	17	
KODAK Precise® Short Plus Trivex®	GR	BR							GR				-10.00 to +4.50	+0.75 to +3.50	-6.00	13	
KODAK Unique Trivex®	GR	BR							GR				-10.00 to +4.50	+0.75 to +3.50	-6.00	13	
KODAK Unique HD Trivex®	GR	BR											-10.00 to +4.50	+0.75 to +3.50	-6.00	13	
KODAK Unique DRO™ Trivex®	GR	BR							GR				-10.00 to +4.50	+0.75 to +3.50	-6.00	13	
KODAK Unique DRO™ HD Trivex®	GR	BR							GR				-10.00 to +4.50	+0.75 to	-6.00	13	

																		+3.50			
KODAK Unique Infinite™ Trivex®	GR	BR							GR	BR								-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ HD Trivex®	GR	BR							GR	BR								-10.00 to +4.50	+0.75 to +3.50	-6.00	13
Crossbows Custom U Trivex®	GR	BR																-10.00 to +4.50	+0.75 to +3.50	-6.00	13
Crossbows Custom V Trivex®	GR	BR																-10.00 to +4.50	+0.75 to +3.50	-6.00	13
1.59 Polycarbonate																					
KODAK Precise® Poly SF	GR																	-10.00 to +5.00	+1.00 to +3.00	-4.00	17
KODAK Precise® Digital Poly	GR																	-10.00 to +6.00	+1.00 to +3.00	-6.00	17
KODAK Precise® Plus Poly	GR	BR	GG	S	AME	AMB	E	R	GR									-10.00 to +5.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus Poly	GR	BR	GG	S	AME	AMB	E	R	GR									-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ Poly	GR	BR	GG	S	AME	AMB	E	R	GR				GR					-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Precise® PB Poly	GR	BR	GG						GR				GR					-10.00 to +5.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB Poly	GR	BR	GG						GR				GR					-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique Poly	GR	BR	GG	S	AME	AMB	E	R	GR				GR					-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD Poly	GR	BR																-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD Poly	GR	BR	GG	S	AME	AMB	E	R	GR				GR					-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ Poly	GR	BR	GG	S	AME	AMB	E	R	GR	BR			GR					-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ HD Poly	GR	BR	GG	S	AME	AMB	E	R	GR	BR			GR					-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Easy 14 Poly	GR	BR	GG	S	AME	AMB	E	R										-10.00 to +5.00	+0.75 to +3.50	-6.00	14
KODAK Easy 18 Poly	GR	BR	GG	S	AME	AMB	E	R										-10.00 to +5.00	+0.75 to +3.50	-6.00	18
DirecTek Poly	GR	BR																-10.00 to +5.00	+0.75 to +3.50	-6.00	16.5
DirecTek Short Poly	GR	BR																-10.00 to +5.00	+0.75 to +3.50	-6.00	13.5
Navigator FBS Poly	GR	BR																-10.00 to +5.00	+0.75 to +3.50	-6.00	18
Navigator Short FBS Poly	GR	BR																-10.00 to +5.00	+0.75 to +3.50	-6.00	14
Crossbows Custom U Poly	GR	BR																-10.00 to +5.00	+0.75 to +3.50	-6.00	13
Crossbows Custom V Poly	GR	BR																-10.00 to +5.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 13 Progressive Poly	GR	BR																-10.00 to +5.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 17 Progressive Poly	GR	BR																-10.00 to +5.00	+0.75 to +3.50	-6.00	17
Crossbows Junior Poly	GR	BR																-10.00 to +5.00	+0.75 to +3.50	-6.00	13
1.60 High Index																					
KODAK Unique	GR	BR																-10.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD	GR	BR																-10.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ 1.60	GR	BR																-10.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.60	GR	BR																-10.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ 1.60	GR	BR																-10.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ HD 1.60	GR	BR																-10.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Custom U 1.60	GR	BR																-10.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Custom V 1.60	GR	BR																-10.00 to +6.00	+0.75 to +3.50	-6.00	13
1.67 Super High Index																					
KODAK Precise® PB 1.67	GR	BR	GG						GR									-11.00 to +6.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB 1.67	GR	BR	GG						GR									-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Precise® Plus Poly	GR	BR	GG	S	AME	AMB	E	R	GR									-11.00 to +6.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus Poly	GR	BR	GG	S	AME	AMB	E	R	GR									-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique	GR	BR	GG	S	AME	AMB	E	R	GR									-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD	GR	BR																-11.00 to +6.00	+0.75 to +3.50	-6.00	13

KODAK Unique DRO™ 1.67	GR	BR	GG	S	AME	AMB	E	R	GR					-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.67	GR	BR	GG	S	AME	AMB	E	R	GR					-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR		GR		-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ HD 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR		GR		-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Easy 14 Poly	GR	BR	GG	S	AME	AMB	E	R						-11.00 to +6.00	+0.75 to +3.50	-6.00	14
KODAK Easy 18 Poly	GR	BR	GG	S	AME	AMB	E	R						-11.00 to +6.00	+0.75 to +3.50	-6.00	18
DirecTek 1.67	GR	BR												-11.00 to +6.00	+0.75 to +3.50	-6.00	16.5
DirecTek Short 1.67	GR	BR												-11.00 to +6.00	+0.75 to +3.50	-6.00	13.5
Navigator FBS 1.67	GR	BR												-11.00 to +6.00	+0.75 to +3.50	-6.00	18
Navigator Short FBS 1.67	GR	BR												-11.00 to +6.00	+0.75 to +3.50	-6.00	14
Crossbows Custom U 1.67	GR	BR												-11.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Custom V 1.67	GR	BR												-11.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 13 Progressive 1.67	GR	BR												-11.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 17 Progressive 1.67	GR	BR												-11.00 to +6.00	+0.75 to +3.50	-6.00	17
Crossbows Junior 1.67	GR	BR												-11.00 to +6.00	+0.75 to +3.50	-6.00	13
1.74 Ultra High Index																	
KODAK Unique	GR	BR												-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD	GR	BR												-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™	GR	BR												-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.74	GR	BR												-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™	GR	BR												-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ HD 1.74	GR	BR												-13.00 to +9.00	+0.75 to +3.50	-6.00	13

			Transitions Gen ⁸ S							Transitions [®] XTRACTIVE [®]			Transitions [®] XTRACTIVE [®] POLARIZED [®]				
Progressives														Sphere	Add	Cyl to	MFH
1.50 Standard Index																	
Unity Via II			GR	BR						GR			GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Plus II			GR	BR						GR			GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Elite II			GR	BR						GR			GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Mobile II			GR	BR						GR			GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Wrap II			GR	BR						GR			GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Ethos			GR	BR						GR			GR	-9.00 to +4.50	0.75 to +3.50	-4.00	14
Ethos Plus			GR	BR						GR			GR	-9.00 to +4.50	0.75 to +3.50	-4.00	14
1.53 TRIVEX [®]																	
Unity Via II			GR	BR						GR				-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Plus II			GR	BR						GR				-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Elite II			GR	BR						GR				-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Mobile II			GR	BR						GR				-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Wrap II			GR	BR						GR				-5.00 to +4.00	+0.50 to +4.50	-6.00	12
Ethos			GR	BR						GR				-9.00 to +6.00	0.75 to +3.50	-4.00	14
Ethos Plus			GR	BR						GR				-9.00 to +6.00	0.75 to +3.50	-4.00	14
1.59 Polycarbonate																	
Unity Via II			GR	BR						GR			GR	-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Plus II			GR	BR						GR			GR	-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Elite II			GR	BR						GR			GR	-10.00 to +6.50	+0.50 to +4.50	-6.00	12

Unity Via Mobile II	GR	BR							GR		GR	-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR							GR		GR	-6.50 to +6.00	+0.50 to +4.50	-6.00	12
Ethos	GR	BR							GR		GR	-9.00 to +6.00	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR							GR		GR	-9.00 to +6.00	0.75 to +3.50	-4.00	14
1.60 High Index															
Unity Via II	GR	BR										-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR										-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR										-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR										-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR										-6.50 to +5.25	+0.50 to +4.50	-6.00	12
Ethos	GR	BR										-9.00 to +5.00	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR										-9.00 to +5.00	0.75 to +3.50	-4.00	14
1.67 Super High Index															
Unity Via II	GR	BR							GR		GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR							GR		GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR							GR		GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR							GR		GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR							GR		GR	-6.00 to +5.00	+0.50 to +4.50	-6.00	12
Ethos	GR	BR							GR		GR	-12.00 to +8.00	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR							GR		GR	-12.00 to +8.00	0.75 to +3.50	-4.00	14

YOUNGER OPTICS The Optical Lens Innovators	Transitions Gen[#]S								Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED						
Progressives															Sphere	Add	Cyl to	MFH
1.53 TRIVEX®																		
Younger Trilogy® Image®	GR														-8.00 to +8.00	+1.00 to +3.00	-4.00	18
1.59 index polycarbonate																		
Younger Image®	GR	BR							GR	BR					-9.00 to +7.00	+1.00 to +3.00	-4.00	18
1.67 Super High Index																		
Younger Image® 1.67 MR-10™	GR	BR													-10.00 to +8.50	+1.00 to +3.50	-4.00	18

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

Unless indicated otherwise, all registered trademarks and trademarks are the property of Essilor International and/or its subsidiaries in the United States and in other countries.

All other trademarks are the property of their respective owners.