










Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Transitions® Signature® GEN 8™](#) / [Shamir Insights Inc.](#)

Shamir Insights Inc. – Transitions® Signature® GEN 8™

												
Progressives									Sphere	Add	Cyl to	MFH
1.50 Standard Index												
	Shamir Autograph® III	GR	BR						-12.25 to +6.00	+0.75 to +4.00	-6.00	11, 13, 15, 18
	Shamir Autograph® II+	GR	BR						-12.25 to +6.00	+0.75 to +4.00	-6.00	11, 13, 15, 18
	Shamir Intouch™	GR	BR						-12.25 to +6.00	+0.75 to +3.50	-6.00	14, 16, 18
	Shamir Spectrum+™	GR	BR						-12.25 to +6.00	+0.75 to +3.50	-6.00	14, 16, 18
	Shamir Element™	GR	BR						-7.00 to +4.00	+0.75 to +3.50	-6.00	16, 19
	Shamir GENESIS HD	GR	BR						-6.00 to +6.00	+0.75 to +3.50	-4.00	16
	Shamir Computer™	GR	BR						-7.00 to +4.00	+0.75 to +2.00	-6.00	14 Lower, 10 Upper
	Shamir Workspace™	GR	BR						-6.00 to +6.00	+0.75 to +2.50	-6.00	14 Lower, 10 Upper
	Shamir Attitude III Sport	GR	BR						-6.50 to +5.25	+0.75 to +3.50	-6.00	18
	Shamir Attitude III Fashion	GR	BR						-6.50 to +5.25	+0.75 to +3.50	-6.00	15, 18
1.53 TRIVEX®												
	Shamir Autograph Intelligence™	GR	BR						-13.25 to +7.00	+0.75 to +4.00	-6.5	11, 13, 15, 18
	Shamir Autograph® III	GR	BR						-13.25 to +7.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
	Shamir Autograph® II+	GR	BR						-13.25 to +6.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
	Shamir Intouch™	GR	BR						-13.25 to +6.00	+0.75 to +3.50	-6.50	14, 16, 18
	Shamir Golf™	GR	BR						-5.00 to +4.00	+1.00 to +3.50	-6.50	19
	Shamir Spectrum+™	GR	BR						-13.25 to +6.00	+0.75 to +3.50	-6.50	14, 16, 18
	Shamir Element™	GR	BR						-13.25 to +6.00	+0.75 to +3.50	-6.50	16, 19
	Shamir GENESIS HD	GR	BR						-6.00 to +6.00	+0.75 to +3.50	-4.00	16
	Shamir Computer™	GR	BR						-6.00 to +8.00	+0.75 to +2.00	-6.50	14 Lower, 10 Upper
	Shamir Workspace™	GR	BR						-6.00 to +8.00	+0.75 to +2.50	-6.50	14 Lower, 10 Upper
	Shamir Attitude III Sport	GR	BR						-6.50 to +6.00	+0.75 to +3.50	-6.50	18
	Shamir Attitude III Fashion	GR	BR						-6.50 to +6.00	+0.75 to +3.50	-6.50	15, 18
1.59 Polycarbonate												
	Shamir Autograph Intelligence™	GR	BR						-14.75 to +7.50	+0.75 to +4.00	-7	11, 13, 15, 18
	Shamir Autograph® III	GR	BR						-14.75 to +8.00	+0.75 to +4.00	-7.00	11, 13, 15, 18
	Shamir Autograph® II+	GR	BR						-14.75to +6.50	+0.75 to +4.00	-7.00	11, 13, 15, 18
	Shamir Intouch™	GR	BR						-14.75 to +6.00	+0.75 to +3.50	-7.00	15, 18
	Shamir Golf™	GR	BR						-6.00 to +5.25	+1.00 to +3.50	-7.00	19
	Shamir Spectrum+™	GR	BR						-14.75 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
	Shamir Element™	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19
	Shamir GENESIS HD	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-7.00	16
	Shamir Computer™	GR	BR						-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper
	Shamir Workspace™	GR	BR						-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper
	Shamir Attitude III Sport	GR	BR						-6.50 to +6.00	+0.75 to +3.50	-6.50	18
	Shamir Attitude III Fashion	GR	BR						-6.50 to +6.00	+0.75 to +3.50	-6.50	15, 18

TRANSITIONS® LIGHT INTELLIGENT LENSES™



GR
Gray



BR
Brown



GG
Graphite
Green



S
Sapphire



AME
Amethyst



AMB
Amber



E
Emerald



R
Ruby



Available
in style
mirrors

Transitions, *Transitions Signature*, and *XTRActive* are registered trademarks, and *XTRActive Polarized*, the *Transitions* logo and *Transitions Light Intelligent Lenses* are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.