



















## Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Transitions® Signature® GEN 8™](#) / **Unity Lenses**

### Unity Lenses – Transitions® Signature® GEN 8™

									
Semi-Finished Single Vision								Sphere	Cyl to
<b>1.50 Standard Index</b>									
Unity SVx	GR	BR						-12.25 to +6.00	-6.00
Unity SVxtra	GR	BR						-12.25 to +6.00	-6.00
Unity SVxtreme	GR	BR						-6.00 to +6.00	-6.00
<b>1.53 TRIVEX® Single Vision</b>									
Unity SVx	GR	BR						-14.00 to +6.00	-6.00
Unity SVxtra	GR	BR						-14.00 to +6.00	-6.00
Unity SVxtreme	GR	BR						-6.50 to +6.00	-6.00
<b>1.59 Polycarbonate Single Vision</b>									
Unity SVx	GR	BR						-14.00 to +6.00	-6.00
Unity SVxtra	GR	BR						-14.00 to +6.00	-6.00
Unity SVxtreme	GR	BR						-6.50 to +6.00	-6.00
<b>1.60 High Index Single Vision</b>									
Unity SVx	GR	BR						-14.00 to +6.00	-6.00
Unity SVxtra	GR	BR						-14.00 to +6.00	-6.00
Unity SVxtreme	GR	BR						-6.50 to +6.00	-6.00
<b>1.67 Super High Index Single Vision</b>									
Unity SVx	GR	BR						-16.00 to +8.00	-6.00
Unity SVxtra	GR	BR						-16.00 to +8.00	-6.00
Unity SVxtreme	GR	BR						-6.00 to +6.00	-6.00

											
Progressives								Sphere	Add	Cyl to	MFH
<b>1.50 Standard Index</b>											
Unity Via II	GR	BR						-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR						-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR						-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR						-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR						-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Ethos	GR	BR						-9.00 to +4.50	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR						-9.00 to +4.50	0.75 to +3.50	-4.00	14
<b>1.53 TRIVEX®</b>											
Unity Via II	GR	BR						-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR						-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR						-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR						-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR						-5.00 to +4.00	+0.50 to +4.50	-6.00	12
Ethos	GR	BR						-9.00 to +6.00	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR						-9.00 to +6.00	0.75 to +3.50	-4.00	14
<b>1.59 Polycarbonate</b>											
Unity Via II	GR	BR						-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR						-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR						-10.00 to +6.50	+0.50 to +4.50	-6.00	12

