

# Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Design](#) / **Trifocals**

	Transitions <b>Gen#S</b>								Transitions Signature GEN8			Transitions SFT		Transitions XTRACTIVE™			Transitions XTRACTIVE™ POLARIZED®				
Trifocal																		<b>Sphere</b>	<b>Add</b>	<b>Cyl to</b>	
<b>1.50 Standard Index SF Trifocal 7x28</b>																					
Signet Armorlite SF Trifocal 7x28									GR									-5.00 to +6.00	+1.50 to +3.00	-4.00	
Signetek™ Processed Trifocal 1.59 7x28									GR									-5.00 to +6.00	+1.50 to +3.00	-6.00	
<b>1.59 Polycarbonate SF Trifocal 7x28</b>																					
Signetek™ Processed Trifocal 1.59 7x28									GR									-6.00 to +6.00	+1.50 to +3.50	-6.00	

**TRANSITIONS® LIGHT INTELLIGENT LENSES™**



Transitions, Transitions Signature, and XTRACTIVE are registered trademarks, and XTRACTIVE Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2026 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.