

Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Transitions® XTRActive®](#) / [Hoya Vision Care](#)

Hoya Vision Care — Transitions® XTRActive®

HOYA		Transitions® XTRACTIVE			
Semi-Finished Single Vision				Sphere	Cyl to
1.50 Standard Index					
HOYA SFSV 1.50	GR BR			-8.00 to +9.00	-4.00
HOYA MySV 150	GR BR			-8.00 to +9.00	-4.00
HOYA SV iQ 1.50	GR BR			-8.00 to +9.00	-4.00
Sync III	GR BR			-10.00 to +6.00	-4.00
1.53 TRIVEX® Single Vision					
HOYA SFSV Phoenix™	GR BR			-13.00 to +9.00	-4.00
HOYA MySV Phoenix™	GR BR			-13.00 to +9.00	-4.00
HOYA SV iQ Phoenix™	GR BR			-13.00 to +9.00	-4.00
Sync III				-10.00 to +6.00	-4.00
1.59 Polycarbonate Single Vision					
HOYA SFSV Polycarbonate	GR BR			-6.00 to +6.00	-4.00
HOYA MySV Polycarbonate	GR BR			-10.00 to +6.00	-4.00
HOYA SV iQ Polycarbonate	GR BR			-10.00 to +6.00	-4.00
Sync III	GR BR			-13.00 to +5.25	-4.00
1.67 Super High Index Single Vision					
HOYA SFSV 1.67	GR BR			-10.00 to +8.00	-4.00
HOYA MySV 1.67	GR BR			-10.00 to +8.00	-4.00
HOYA SV iQ 1.67	GR BR			-10.00 to +8.00	-4.00
Sync III	GR BR			-13.00 to +7.50	-4.00
1.74 Super High Index Single Vision					
HOYA SFSV 1.74	GR BR			-13.00 to +10.00	-4.00

HOYA		Transitions® XTRACTIVE					
Progressives				Sphere	Add	Cyl to	MFH
1.50 Standard Index							
HOYALUX ARRAY™ 11, 13, 15, 17, VL 1.50	GR BR			-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL 1.50	GR BR			-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ 2 11, 13, 15, 17, VL 1.50	GR BR			-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ 2 Wrap 11, 13, 15, 17, VL 1.50	GR BR			-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS 1.50	GR BR			-10.00 to +3.75	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS 1.50	GR BR			-10.00 to +3.75	+0.75 to +4.00	-6.00	14
Summit cd BKS 1.50	GR BR			-10.00 to +3.75	+0.75 to +4.00	-6.00	14
Summit ecp BKS 1.50	GR BR			-10.00 to +3.75	+0.75 to +4.00	-6.00	18
Hoya GP Wide 1.50	GR BR			-10.00 to +3.75	+0.75 to +4.00	-6.00	18
Tact 40 BKS 1.50	GR BR			-10.00 to +3.75	+0.75 to +4.00	-6.00	18
Tact 60 BKS 1.50	GR BR			-10.00 to +3.75	+0.75 to +4.00	-6.00	18
1.53 TRIVEX®							
HOYALUX ARRAY™ 11, 13, 15, 17, VL Phoenix™	GR BR			-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL Phoenix™	GR BR			-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13, 15, 17, VL Phoenix™	GR BR			-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11, 13, 15, 17, VL Phoenix™	GR BR			-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS Phoenix™	GR BR			-10.00 to +6.00	+0.75 to +4.00	-6.00	19

Amplitude Mini BKS Phoenix™	GR	BR		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Summit cd BKS Phoenix™	GR	BR		-10.00 to +6.00	+0.75 to +4.00	-6.00	18
Summit ecp BKS Phoenix™	GR	BR		-10.00 to +6.00	+0.75 to +4.00	-6.00	18
Hoya GP Wide Phoenix™	GR	BR		-10.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 40 BKS Phoenix™	GR	BR		-10.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 60 BKS Phoenix™	GR	BR		-10.00 to +6.00	+0.75 to +4.00	-6.00	18
1.59 Polycarbonate							
HOYALUX ARRAY™ 11, 13, 15, 17, VL Poly	GR	BR		-13.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL Poly	GR	BR		-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13, 15, 17, VL Poly	GR	BR		-13.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11, 13, 15, 17, VL Poly	GR	BR		-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS Poly	GR	BR		-13.00 to +6.00	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS Poly	GR	BR		-13.00 to +6.00	+0.75 to +4.00	-6.00	14
Summit cd BKS Poly	GR	BR		-13.00 to +6.00	+0.75 to +4.00	-6.00	14
Summit ecp BKS Poly	GR	BR		-13.00 to +6.00	+0.75 to +4.00	-6.00	18
Hoya GP Wide Poly	GR	BR		-13.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 40 BKS Poly	GR	BR		-13.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 60 BKS Poly	GR	BR		-13.00 to +6.00	+0.75 to +4.00	-6.00	18
1.67 Super High Index							
HOYALUX ARRAY 11, 13,15,17,VL 1.67	GR	BR		-13.00 to +7.50	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY Wrap 11,13,15,17,VL 1.67	GR	BR		-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13,15,17,VL 1.67	GR	BR		-13.00 to +7.50	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11,13,15,17,VL 1.67	GR	BR		-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS 1.67	GR	BR		-13.00 to +7.50	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS 1.67	GR	BR		-13.00 to +7.50	+0.75 to +4.00	-6.00	14
Summit cd BKS 1.67	GR	BR		-13.00 to +7.50	+0.75 to +4.00	-6.00	14
Summit ecp BKS 1.67	GR	BR		-13.00 to +7.50	+0.75 to +4.00	-6.00	18
Hoya GP Wide 1.67	GR	BR		-13.00 to +7.50	+0.75 to +4.00	-6.00	18
Tact 40 BKS 1.67	GR	BR		-13.00 to +7.50	+0.75 to +4.00	-6.00	18
Tact 60 BKS 1.67	GR	BR		-13.00 to +7.50	+0.75 to +4.00	-6.00	18

TRANSITIONS® LIGHT INTELLIGENT LENSES™



GR
Gray



BR
Brown



GG
Graphite
Green



S
Sapphire



AME
Amethyst



AMB
Amber



E
Emerald



R
Ruby



Available
in style
mirrors

Style Mirrors are available where gray and brown *Transitions® XTRActive®* are available. Specify *Transitions* lenses in style mirrors **(no substitutions)** with your lab to ensure authenticity.

Transitions, *Transitions Signature*, and *XTRActive* are registered trademarks, and *XTRActive Polarized*, the *Transitions* logo and *Transitions Light Intelligent Lenses* are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.