







Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Transitions® XTRActive®](#) / [Unity Lenses](#)

Unity Lenses – Transitions® XTRActive®

unity		Transitions XTRACTIVE®				
Semi-Finished Single Vision					Sphere	Cyl to
1.50 Standard Index						
	Unity SVx	GR			-12.25 to +6.00	-6.00
	Unity SVxtra	GR			-12.25 to +6.00	-6.00
	Unity SVxtreme	GR			-6.00 to +6.00	-6.00
1.53 TRIVEX® Single Vision						
	Unity SVx	GR			-14.00 to +6.00	-6.00
	Unity SVxtra	GR			-14.00 to +6.00	-6.00
	Unity SVxtreme	GR			-6.50 to +6.00	-6.00
1.59 Polycarbonate Single Vision						
	Unity SVx	GR			-14.00 to +6.00	-6.00
	Unity SVxtra	GR			-14.00 to +6.00	-6.00
	Unity SVxtreme	GR			-6.50 to +6.00	-6.00
1.67 Super High Index Single Vision						
	Unity SVx	GR			-16.00 to +8.00	-6.00
	Unity SVxtra	GR			-16.00 to +8.00	-6.00
	Unity SVxtreme	GR			-6.00 to +6.00	-6.00

unity		Transitions XTRACTIVE®						
Progressives					Sphere	Add	Cyl to	MFH
1.50 Standard Index								
	Unity V3	GR			-7.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity V3 Plus	GR			-13.25 to +6.50	+0.50 to +4.50	-6.00	12
	Unity Via Elite VR	GR			-13.25 to +6.50	+0.50 to +4.50	-6.00	12
	Unity V3 Elite	GR			-13.25 to +6.50	+0.50 to +4.50	-6.00	12
	Unity V3 Elite 12	GR			-13.25 to +6.50	+0.50 to +4.50	-6.00	12
	Unity V3 Elite 14	GR			-13.25 to +6.50	+0.50 to +4.50	-6.00	14
	Unity V3 Elite 16	GR			-13.25 to +6.50	+0.50 to +4.50	-6.00	16
	Unity V3 Elite 18	GR			-13.25 to +6.50	+0.50 to +4.50	-6.00	18
	Unity V3 Mobile	GR			-13.25 to +6.50	+0.50 to +4.50	-6.00	12
	Unity V3 Wrap	GR			-7.00 to +6.00	+0.50 to +4.50	-4.00	12
	Ethos	GR			-9.00 to +4.50	0.75 to +3.50	up to -4.00	14
	Ethos Plus	GR			-9.00 to +4.50	0.75 to +3.50	up to -4.00	14
1.53 TRIVEX®								
	Unity V3	GR			-10.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity V3 Plus	GR			-13.50 to +7.00	+0.50 to +4.50	-7.00	12
	Unity Via Elite VR	GR			-13.50 to +7.00	+0.50 to +4.50	-7.00	12
	Unity V3 Elite	GR			-13.50 to +7.00	+0.50 to +4.50	-7.00	12
	Unity V3 Elite 12	GR			-13.50 to +7.00	+0.50 to +4.50	-7.00	12
	Unity V3 Elite 14	GR			-13.50 to +7.00	+0.50 to +4.50	-7.00	14
	Unity V3 Elite 16	GR			-13.50 to +7.00	+0.50 to +4.50	-7.00	16
	Unity V3 Elite 18	GR			-13.50 to +7.00	+0.50 to +4.50	-7.00	18
	Unity V3 Mobile	GR			-13.50 to +7.00	+0.50 to +4.50	-7.00	12
	Unity V3 Wrap	GR			-7.00 to +7.00	+0.50 to +4.50	-6.00	12
	Ethos	GR			-9.00 to +6.00	0.75 to +3.50	up to -4.00	14
	Ethos Plus	GR			-9.00 to +6.00	0.75 to +3.50	up to -4.00	14

1.59 Polycarbonate							
Unity V3	GR			-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity V3 Plus	GR			-15.00 to +9.00	+0.50 to +4.50	-8.00	12
Unity Via Elite VR	GR			-15.00 to +9.00	+0.50 to +4.50	-8.00	12
Unity V3 Elite	GR			-15.00 to +9.00	+0.50 to +4.50	-8.00	12
Unity V3 Elite 12	GR			-15.00 to +9.00	+0.50 to +4.50	-8.00	12
Unity V3 Elite 14	GR			-15.00 to +9.00	+0.50 to +4.50	-8.00	14
Unity V3 Elite 16	GR			-15.00 to +9.00	+0.50 to +4.50	-8.00	16
Unity V3 Elite 18	GR			-15.00 to +9.00	+0.50 to +4.50	-8.00	18
Unity V3 Mobile	GR			-15.00 to +9.00	+0.50 to +4.50	-8.00	12
Unity V3 Wrap	GR			-6.50 to +6.00	+0.50 to +4.50	-5.00	12
Ethos	GR			-9.00 to +6.00	0.75 to +3.50	up to -4.00	14
Ethos Plus	GR			-9.00 to +6.00	0.75 to +3.50	up to -4.00	14
1.67 Super High Index							
Unity V3	GR			-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity V3 Plus	GR			-18.00 to +10.00	+0.50 to +4.50	-9.00	12
Unity Via Elite VR	GR			-18.00 to +10.00	+0.50 to +4.50	-9.00	12
Unity V3 Elite	GR			-18.00 to +10.00	+0.50 to +4.50	-9.00	12
Unity V3 Elite 12	GR			-18.00 to +10.00	+0.50 to +4.50	-9.00	12
Unity V3 Elite 14	GR			-18.00 to +10.00	+0.50 to +4.50	-9.00	14
Unity V3 Elite 16	GR			-18.00 to +10.00	+0.50 to +4.50	-9.00	16
Unity V3 Elite 18	GR			-18.00 to +10.00	+0.50 to +4.50	-9.00	18
Unity V3 Mobile	GR			-18.00 to +10.00	+0.50 to +4.50	-9.00	12
Unity V3 Wrap	GR			-8.00 to +8.00	+0.50 to +4.50	-7.00	12
Ethos	GR			-12.00 to +8.00	0.75 to +3.50	up to -4.00	14
Ethos Plus	GR			-12.00 to +8.00	0.75 to +3.50	up to -4.00	14

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2026 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.