










Transitions® Lenses Product Availability Guide





< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Transitions® XTRActive®](#) / [Younger Optics](#)



Younger Optics — Transitions® XTRActive®

	Transitions® XTRACTIVE®				
Finished Single Vision				Sphere	Cyl to
1.50 Standard Index					
Younger Optics FSV	GR	BR		-4.00 to +4.00	-2.00
1.59 Polycarbonate Finished Single Vision					
Younger Optics FSV Poly	GR	BR		-6.00 to +4.00	-2.00

	Transitions® XTRACTIVE®				
Semi-Finished Single Vision				Sphere	Cyl to
1.50 Standard Index					
Younger Optics SFSV (76mm)	GR	BR		-9.00 to +6.00	-4.00
1.53 TRIVEX® Single Vision					
Younger Optics Trilogy® Spherical SV	GR	BR		-8.00 to +7.50	-4.00
1.59 Polycarbonate Single Vision					
Younger Optics SV	GR	BR		-9.00 to +6.00	-4.00
1.67 Super High Index Single Vision					
Younger Optics Single Vision 1.67 MR-10™	GR	BR		-11.00 to +7.00	-4.00
1.74 Ultra High Index Single Vision					
Younger Optics Single Vision 1.74	GR	BR		-18.00 to +12.00	-4.00

	Transitions® XTRACTIVE®					
Bifocal				Sphere	Add	Cyl to
1.50 Standard Index SFST-28						
Younger Optics SFFT 28	GR			-11.00 to +6.00	+1.00 to +3.00	-4.00

	Transitions® XTRACTIVE®				
Digital				Sphere	Cyl to
1.50 Index					
Camber™	GR	BR		-11.25 to +6.25	-4.00
1.53 TRIVEX®					
Camber™	GR	BR		-11.00 to +6.25	-4.00
1.59 Polycarbonate					
Camber™	GR	BR		-10.00 to +6.25	-4.00
1.67 Super High Index					
Camber™	GR	BR		-14.00 to +8.00	-4.00
1.74 Ultra High Index					
Camber™	GR	BR		-19.00 to +9.00	-4.00

	Transitions® XTRACTIVE®						
Progressives				Sphere	Add	Cyl to	MFH

1.59 index polycarbonate								
	Younger Image®	GR	BR		-9.00 to +7.00	+1.00 to +3.00	-4.00	18

TRANSITIONS® LIGHT INTELLIGENT LENSES™

- 

**GR**  
Gray



**BR**  
Brown



**GG**  
Graphite  
Green



**S**  
Sapphire



**AME**  
Amethyst



**AMB**  
Amber



**E**  
Emerald



**R**  
Ruby



Available  
in style  
mirrors

Style Mirrors are available where gray and brown *Transitions® XTRActive®* are available. Specify *Transitions* lenses in style mirrors **(no substitutions)** with your lab to ensure authenticity.

*Transitions*, *Transitions Signature*, and *XTRActive* are registered trademarks, and *XTRActive Polarized*, the *Transitions* logo and *Transitions Light Intelligent Lenses* are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.