





Transitions® Lenses Product Availability Guide


< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Transitions® XTRActive®](#)

All Suppliers – Transitions® XTRActive®

 Finished Single Vision	Transitions® XTRACTIVE®				
				Sphere	Cyl to
1.50 Standard Index					
Essilor FSV Spherical	GR			-4.00 to +4.00	-2.00
Essilor FSV 1.50 Crizal® Sapphire™ HR	GR			-4.00 to +4.00	-2.00
Essilor FSV 1.50 Crizal® Rock™	GR			-4.00 to +4.00	-2.00
Essilor FSV 1.50 Crizal® Easy Pro	GR			-4.00 to +4.00	-2.00
1.59 Polycarbonate Finished Single Vision					
Essilor FSV Aspheric Airwear®	GR			-6.00 to +4.00	-2.00
Essilor FSV Airwear® Crizal® Sapphire™ HR	GR			-6.00 to +4.00	-2.00
Essilor FSV Airwear® Crizal® Rock™	GR			-6.00 to +4.00	-2.00
Essilor FSV Airwear® Crizal® Easy Pro	GR			-6.00 to +4.00	-2.00
1.67 Super High Index Finished Single Vision					
Essilor FSV Spherical Thin&Lite® 1.67				-6.00 to +4.00	-2.00

 Semi-Finished Single Vision	Transitions® XTRACTIVE®				
				Sphere	Cyl to
1.50 Standard Index					
Essilor Spherical Single Vision	GR	BR	GG	-10.50 to +10.00	-4.00
Essilor Stylistic® Single Vision	GR	BR	GG	-10.00 to +7.50	-6.00
Eyezen® +0	GR	BR	GG	-10.00 to +7.50	-6.00
Eyezen® +1	GR	BR	GG	-10.00 to +7.50	-6.00
Eyezen® +2	GR	BR	GG	-10.00 to +7.50	-6.00
Eyezen® +3	GR	BR	GG	-10.00 to +7.50	-6.00
Eyezen® +4	GR	BR	GG	-10.00 to +7.50	-6.00
Eyezen® Start	GR	BR	GG	-10.00 to +6.00	-6.00
1.53 TRIVEX® Single Vision					
Essilor Spherical Single Vision TREXA™	GR	BR		-10.00 to +6.00	-4.00
Essilor Eyecode™ Single Vision TREXA™	GR	BR		-10.00 to +7.50	-6.00
Essilor FIT™ Single Vision TREXA™	GR	BR		-10.00 to +7.50	-6.00
Essilor Azio™ Single Vision TREXA™	GR	BR		-10.00 to +7.50	-6.00
Essilor 360°™ Single Vision TREXA™	GR	BR		-10.00 to +7.50	-6.00
Eyezen® +0 TREXA™	GR	BR		-10.00 to +7.50	-6.00
Eyezen® +1 TREXA™	GR	BR		-10.00 to +7.50	-6.00
Eyezen® +2 TREXA™	GR	BR		-10.00 to +7.50	-6.00
Eyezen® +3 TREXA™	GR	BR		-10.00 to +7.50	-6.00
Eyezen® +4 TREXA™	GR	BR		-10.00 to +7.50	-6.00
Eyezen® Start	GR	BR		-10.00 to +6.00	-6.00
1.59 Polycarbonate Single Vision					
Essilor Spherical Single Vision Airwear®	GR	BR	GG	-10.00 to +8.00	-4.00
Essilor Aspheric Single Vision Airwear®	GR	BR	GG	-10.00 to +8.00	-4.00
Essilor Eyecode™ Aspheric SV Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Essilor FIT™ Aspheric SV Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Essilor AZIO™ Aspheric SV Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00



Essilor SV 360°™ Aspheric Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Essilor Stylistic® Single Vision Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Eyezen® +0 Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Eyezen® +1 Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Eyezen® +2 Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Eyezen® +3 Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Eyezen® +4 Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Eyezen® Start Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
1.67 Super High Index Single Vision					
Essilor SFSV Spherical Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-4.00
Essilor Eyecode™ SV Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Essilor FIT™ SV Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Essilor Azio™ SV Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Essilor SV 360°™ Spherical Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Essilor Stylistic® Single Vision Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Eyezen® +0 Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Eyezen® +1 Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Eyezen® +2 Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Eyezen® +3 Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Eyezen® +4 Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Eyezen® Start Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00




 Progressives				Sphere	Add	Cyl to	MFH
							
1.50 Standard Index							
Varilux® XR Track™	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Varilux Comfort®	GR	BR	GG	-10.00 to +7.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short™	GR	BR	GG	-10.00 to +7.00	+0.75 to +3.00	-6.00	14
Varilux® Stylistic™	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Wrap	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-4.00	17
Essilor Ideal® FIT	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	18
Essilor Ideal® Driving	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	13
Essilor Ideal®	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal® Short	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	14

Essilor Natural® Digital	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	15
Essilor Accolade™	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Digital	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Brillenz®	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Brillenz® Short	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Ovation® Digital	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor SmallFit® Digital	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.00	-6.00	14
1.53 TRIVEX®							
Varilux® Physio® W3+ Eyecode TREXA™	GR	BR		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT TREXA™	GR	BR		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ TREXA™	GR	BR		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit TREXA™	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max TREXA™	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™ TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Essilor Ideal Advanced™ TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	18
Essilor Ideal® Computer TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	13
Essilor Ideal® Driving TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal® TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal® Short TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Natural® Digital TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	15
Essilor Accolade® TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Digital TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Brillenz® TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Brillenz® Short TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Ovation® Digital TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor SmallFit® Digital TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.00	-6.00	14
1.59 Polycarbonate							
Varilux® XR Track™	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode™ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Varilux Comfort DRx™ Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux Comfort® Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short™ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.00	-6.00	14





Varilux® Stylistic™ Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Wrap	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Essilor Ideal® FIT Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Sport Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Driving Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Ideal® Short Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Natural® Digital Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Ovation® Digital Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor SmallFit® Digital Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
1.67 Super High Index							
Varilux® XR Track™	GR	BR	GG	-12.00 to +9.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR	GG	-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ Thin&Lite® 1.67	GR	BR	GG	-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ Thin&Lite® 1.67	GR	BR	GG	-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ Thin&Lite® 1.67	GR	BR	GG	-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Varilux® Stylistic™ Thin&Lite® 1.67	GR	BR	GG	-5.00 to +5.00	+0.75 to +3.50	-3.00	17
Essilor Ideal Advanced Wrap™ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Sport Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Driving Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor Ideal® Short	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor Natural® Digital Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Brillenz® Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14

Brillenz® Short Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Ovation® Digital Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor SmallFit® Digital Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14

HOYA		Transitions XTRACTIVE				
Semi-Finished Single Vision					Sphere	Cyl to
1.50 Standard Index						
HOYA SFSV 1.50	GR	BR			-8.00 to +9.00	-4.00
HOYA MySV 150	GR	BR			-8.00 to +9.00	-4.00
HOYA SV iQ 1.50	GR	BR			-8.00 to +9.00	-4.00
Sync III	GR	BR			-10.00 to +6.00	-4.00
1.53 TRIVEX® Single Vision						
HOYA SFSV Phoenix™	GR	BR			-13.00 to +9.00	-4.00
HOYA MySV Phoenix™	GR	BR			-13.00 to +9.00	-4.00
HOYA SV iQ Phoenix™	GR	BR			-13.00 to +9.00	-4.00
Sync III					-10.00 to +6.00	-4.00
1.59 Polycarbonate Single Vision						
HOYA SFSV Polycarbonate	GR	BR			-6.00 to +6.00	-4.00
HOYA MySV Polycarbonate	GR	BR			-10.00 to +6.00	-4.00
HOYA SV iQ Polycarbonate	GR	BR			-10.00 to +6.00	-4.00
Sync III	GR	BR			-13.00 to +5.25	-4.00
1.67 Super High Index Single Vision						
HOYA SFSV 1.67	GR	BR			-10.00 to +8.00	-4.00
HOYA MySV 1.67	GR	BR			-10.00 to +8.00	-4.00
HOYA SV iQ 1.67	GR	BR			-10.00 to +8.00	-4.00
Sync III	GR	BR			-13.00 to +7.50	-4.00
1.74 Super High Index Single Vision						
HOYA SFSV 1.74	GR	BR			-13.00 to +10.00	-4.00

HOYA		Transitions XTRACTIVE						
Progressives					Sphere	Add	Cyl to	MFH
1.50 Standard Index								
HOYALUX ARRAY™ 11, 13, 15, 17, VL 1.50	GR	BR			-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL 1.50	GR	BR			-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ 2 11, 13, 15, 17, VL 1.50	GR	BR			-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ 2 Wrap 11, 13, 15, 17, VL 1.50	GR	BR			-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS 1.50	GR	BR			-10.00 to +3.75	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS 1.50	GR	BR			-10.00 to +3.75	+0.75 to +4.00	-6.00	14
Summit cd BKS 1.50	GR	BR			-10.00 to +3.75	+0.75 to +4.00	-6.00	14
Summit ecp BKS 1.50	GR	BR			-10.00 to +3.75	+0.75 to +4.00	-6.00	18
Hoya GP Wide 1.50	GR	BR			-10.00 to +3.75	+0.75 to +4.00	-6.00	18
Tact 40 BKS 1.50	GR	BR			-10.00 to +3.75	+0.75 to +4.00	-6.00	18
Tact 60 BKS 1.50	GR	BR			-10.00 to +3.75	+0.75 to +4.00	-6.00	18
1.53 TRIVEX®								
HOYALUX ARRAY™ 11, 13, 15, 17, VL Phoenix™	GR	BR			-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL Phoenix™	GR	BR			-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13, 15, 17, VL Phoenix™	GR	BR			-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11, 13, 15, 17, VL Phoenix™	GR	BR			-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS Phoenix™	GR	BR			-10.00 to +6.00	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS Phoenix™	GR	BR			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Summit cd BKS Phoenix™	GR	BR			-10.00 to +6.00	+0.75 to +4.00	-6.00	18
Summit ecp BKS Phoenix™	GR	BR			-10.00 to +6.00	+0.75 to +4.00	-6.00	18




Hoya GP Wide Phoenix™	GR	BR		-10.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 40 BKS Phoenix™	GR	BR		-10.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 60 BKS Phoenix™	GR	BR		-10.00 to +6.00	+0.75 to +4.00	-6.00	18
1.59 Polycarbonate							
HOYALUX ARRAY™ 11, 13, 15, 17, VL Poly	GR	BR		-13.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL Poly	GR	BR		-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13, 15, 17, VL Poly	GR	BR		-13.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11, 13, 15, 17, VL Poly	GR	BR		-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS Poly	GR	BR		-13.00 to +6.00	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS Poly	GR	BR		-13.00 to +6.00	+0.75 to +4.00	-6.00	14
Summit cd BKS Poly	GR	BR		-13.00 to +6.00	+0.75 to +4.00	-6.00	14
Summit ecp BKS Poly	GR	BR		-13.00 to +6.00	+0.75 to +4.00	-6.00	18
Hoya GP Wide Poly	GR	BR		-13.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 40 BKS Poly	GR	BR		-13.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 60 BKS Poly	GR	BR		-13.00 to +6.00	+0.75 to +4.00	-6.00	18
1.67 Super High Index							
HOYALUX ARRAY 11, 13,15,17,VL 1.67	GR	BR		-13.00 to +7.50	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY Wrap 11,13,15,17,VL 1.67	GR	BR		-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13,15,17,VL 1.67	GR	BR		-13.00 to +7.50	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11,13,15,17,VL 1.67	GR	BR		-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS 1.67	GR	BR		-13.00 to +7.50	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS 1.67	GR	BR		-13.00 to +7.50	+0.75 to +4.00	-6.00	14
Summit cd BKS 1.67	GR	BR		-13.00 to +7.50	+0.75 to +4.00	-6.00	14
Summit ecp BKS 1.67	GR	BR		-13.00 to +7.50	+0.75 to +4.00	-6.00	18
Hoya GP Wide 1.67	GR	BR		-13.00 to +7.50	+0.75 to +4.00	-6.00	18
Tact 40 BKS 1.67	GR	BR		-13.00 to +7.50	+0.75 to +4.00	-6.00	18
Tact 60 BKS 1.67	GR	BR		-13.00 to +7.50	+0.75 to +4.00	-6.00	18



	Transitions® XTRACTIVE®				
Semi-Finished Single Vision				Sphere	Cyl to
1.50 Standard Index					
Nikon SeeMax Infinite	GR	BR	GG	-10.00 to +8.00	-6.00
Nikon Infinite Single Vision	GR	BR	GG	-10.00 to +8.00	-6.00
Nikon Smart Boost	GR	BR	GG	-10.00 to +8.00	-6.00
1.59 Polycarbonate Single Vision					
Nikon SeeMax Infinite	GR	BR	GG	-10.00 to +8.00	-6.00
Nikon Infinite Single Vision	GR	BR	GG	-10.00 to +8.00	-6.00
Nikon Smart Boost	GR	BR	GG	-10.00 to +8.00	-6.00
1.67 Super High Index Single Vision					
Nikon SeeMax Infinite	GR	BR	GG	-14.00 to +9.00	-6.00
Nikon Infinite Single Vision	GR	BR	GG	-14.00 to +9.00	-6.00
Nikon Smart Boost	GR	BR	GG	-14.00 to +9.00	-6.00




	Transitions® XTRACTIVE®						
Progressives				Sphere	Add	Cyl to	MFH
1.50 Standard index							
Nikon SeeMax Ultimate	GR	BR	GG	-10.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR	GG	-10.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR	GG	-10.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
Nikon Move	GR	BR	GG	-10.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
1.59 Polycarbonate							
Nikon SeeMax Ultimate	GR	BR	GG	-11.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR	GG	-11.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17

	Nikon Presio i Digital	GR	BR	GG	-11.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
	Nikon Move	GR	BR	GG	-11.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
1.67 Super High Index								
	Nikon SeeMax Ultimate	GR	BR	GG	-14.00 to +9.00	+0.50 to +4.00	-6.00	13, 15, 17
	Nikon Presio Power	GR	BR	GG	-14.00 to +9.00	+0.50 to +4.00	-6.00	13, 15, 17
	Nikon Presio i Digital	GR	BR	GG	-14.00 to +9.00	+0.50 to +4.00	-6.00	15, 17
	Nikon Move	GR	BR	GG	-14.00 to +9.00	+0.50 to +4.00	-6.00	15, 17

SEIKO		Transitions XTRACTIVE						
	Semi-Finished Single Vision				Sphere	Cyl to		
1.53 TRIVEX® Single Vision								
	Seiko PL 1.53/Trivex Transitions®	GR			-10.00 to +7.00	-5.00		
1.67 Super High Index Single Vision								
	Seiko PL 1.67 Transitions®	GR			-12.00 to +10.00	-5.00		

shamir		Transitions XTRACTIVE						
	Semi-Finished Single Vision				Sphere	Cyl to		
1.50 Standard Index								
	Shamir Autograph® III Single Vision	GR	BR		-12.00 to +7.25	-6.00		
	Shamir Relax™	GR	BR		-14.00 to +8.00	-6.00		
1.53 SFSV								
	Shamir Autograph® III Single Vision	GR	BR		-14.00 to +8.00	-6.50		
	Shamir Relax™	GR	BR		-14.00 to +8.00	-6.50		
1.59 SFSV								
	Shamir Autograph® III Single Vision	GR	BR		-14.00 to +8.75	-7.00		
	Shamir Relax™	GR	BR		-14.75 to +6.00	-7.00		
1.60 SFSV								
	Shamir Autograph® III Single Vision	GR			-14.00 to +8.75	-7.00		
	Shamir Relax™	GR			-14.00 to +8.75	-7.00		
1.67 SFSV								
	Shamir Autograph® III Single Vision	GR	BR		-16.00 to +16.00	-8.00		
	Shamir Relax™	GR	BR		-16.75 to +16.00	-8.00		






shamir		Transitions XTRACTIVE						
	Blended Bifocal (digital)				Sphere	Add	Cyl to	MFH
1.50 Standard Index Blended Bifocal (digital)								
	Shamir Duo	GR			-9.00 to +6.00	+0.75 to +4.00	-6.00	15
1.53 Trivex® Blended Bifocal (digital)								
	Shamir Duo	GR			-13.25 to +6.00	+0.75 to +4.00	-6.50	15
1.59 Polycarbonate Blended Bifocal (digital)								
	Shamir Duo	GR			-14.75 to +6.00	+0.75 to +4.00	-7.00	15
1.60 High Index Blended Bifocal (digital)								
	Shamir Duo	GR			-15.00 to +6.00	+0.75 to +4.00	-7.00	15
1.67 Super High Index Blended Bifocal (digital)								
	Shamir Duo	GR			-16.75 to +8.00	+0.75 to +4.00	-8.00	15

shamir		Transitions XTRACTIVE						
	Progressives				Sphere	Add	Cyl to	MFH
1.50 Standard Index								
	Shamir Autograph Intelligence™	GR			-12.25 to +6.50	+0.75 to +4.00	-6.00	11, 13, 15, 18

Shamir Autograph® III	GR			-12.25 to +6.00	+0.75 to +4.00	-6.00	11, 13, 15, 18
Shamir Autograph® II+	GR			-12.25 to +6.00	+0.75 to +4.00	-6.00	11, 13, 15, 18
Shamir Intouch™	GR			-12.25 to +6.00	+0.75 to +3.50	-6.00	14, 16, 18
Shamir Spectrum+™	GR			-12.25 to +6.00	+0.75 to +3.50	-6.00	14, 16, 18
Shamir Element™	GR			-7.00 to +4.00	+0.75 to +3.50	-6.00	16, 19
Shamir GENESIS HD	GR			-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir Computer™	GR			-7.00 to +4.00	+0.75 to +2.00	-6.00	14 Lower, 10 Upper
Shamir Workspace™	GR			-6.00 to +6.00	+0.75 to +2.50	-6.00	14 Lower, 10 Upper
Shamir Attitude III Sport	GR			-6.50 to +5.25	+0.75 to +3.50	-6.00	18
Shamir Attitude III Fashion	GR			-6.50 to +5.25	+0.75 to +3.50	-6.00	15, 18
1.53 TRIVEX®							
Shamir Autograph Intelligence™	GR			-13.25 to +7.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
Shamir Autograph® III	GR			-13.25 to +7.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
Shamir Autograph® II+	GR			-13.25 to +6.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
Shamir Intouch™	GR			-13.25 to +6.00	+0.75 to +3.50	-6.50	14, 16, 18
Shamir Golf™	GR			-5.00 to +4.00	+1.00 to +3.50	-6.50	19
Shamir Spectrum+™	GR			-13.25 to +6.00	+0.75 to +3.50	-6.50	14, 16, 18
Shamir Element™	GR			-13.25 to +6.00	+0.75 to +3.50	-6.50	16, 19
Shamir GENESIS HD	GR			-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir Computer™	GR			-6.00 to +8.00	+0.75 to +2.00	-6.50	14 Lower, 10 Upper
Shamir Workspace™	GR			-6.00 to +8.00	+0.75 to +2.50	-6.50	14 Lower, 10 Upper
Shamir Attitude III Sport	GR			-6.50 to +6.00	+0.75 to +3.50	-6.50	18
Shamir Attitude III Fashion	GR			-6.50 to +6.00	+0.75 to +3.50	-6.50	15, 18
1.59 Polycarbonate							
Shamir Autograph Intelligence™	GR			-14.75 to +7.50	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® III	GR			-14.75 to +8.00	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® II+	GR			-14.75 to +6.50	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Intouch™	GR			-14.75 to +6.00	+0.75 to +3.50	-7.00	15, 18
Shamir Golf™	GR			-6.00 to +5.25	+1.00 to +3.50	-7.00	19
Shamir Spectrum+™	GR			-14.75 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
Shamir Element™	GR			-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19
Shamir GENESIS HD	GR			-10.00 to +6.00	+0.75 to +3.50	-7.00	16
Shamir Computer™	GR			-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper
Shamir Workspace™	GR			-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper
Shamir Attitude III Sport	GR			-6.50 to +6.00	+0.75 to +3.50	-6.50	18
Shamir Attitude III Fashion	GR			-6.50 to +6.00	+0.75 to +3.50	-6.50	15, 18
1.60 High Index							
Shamir Autograph Intelligence™	GR			-14.75 to +8.75	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® III	GR			-15.00 to +8.00	+0.75 to +3.50	-7.00	15, 18
Shamir Autograph® II+	GR			-15.00 to +6.00	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Golf™ SuperLite™ 1.6	GR			-6.00 to +5.25	+1.00 to +3.50	-7.00	19
Shamir Spectrum+™	GR			-15.00 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
Shamir GENESIS HD	GR			-6.00 to +6.00	+0.75 to +3.50	-7.00	16
Shamir Computer™	GR			-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper
Shamir Workspace™	GR			-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper
Shamir InTouch™	GR			-15.00 to +6.00	+0.75 to +3.50	-7.00	15, 18
Shamir Attitude III Sport	GR			-6.50 to +6.00	+0.75 to +3.50	-7.00	18
Shamir Attitude III Fashion	GR			-6.50 to +6.00	+0.75 to +3.50	-7.00	15, 18
1.67 Super High Index							
Shamir Autograph Intelligence™	GR			-16.75 to +10.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Autograph® III	GR			-16.75 to +9.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Autograph® II+	GR			-16.75 to +8.25	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Spectrum+™	GR			-16.75 to +8.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Element™ SuperLite™ 1.67	GR			-12.50 to +8.00	+0.75 to +3.50	-8.00	16, 19
Shamir GENESIS HD	GR			-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir InTouch™	GR			-16.75 to +8.00	+0.75 to +3.50	-8.00	15, 18

Shamir Attitude III Sport	GR			-6.50 to +6.00	+0.75 to +3.50	-8.00	18
Shamir Attitude III Fashion	GR			-6.50 to +6.00	+0.75 to +3.50	-8.00	15, 18






						
Semi-Finished Single Vision					Sphere	Cyl to
1.50 Standard Index						
Signetek Processed SV	GR				-8.00 to +6.00	-6.00
KODAK Digital SV	GR				-8.00 to +6.00	-6.00
1.53 TRIVEX® Single Vision						
Signetek Processed SV	GR				-10.00 to +6.00	-6.00
KODAK Digital SV	GR				-10.00 to +6.00	-6.00
1.59 Polycarbonate Single Vision						
Signetek Processed SV	GR				-10.00 to +6.50	-6.00
KODAK Digital SV	GR				-10.00 to +6.50	-6.00
1.67 Super High Index Single Vision						
Signetek Processed SV	GR				-11.00 to +7.50	-6.00
KODAK Digital SV	GR				-11.00 to +7.50	-6.00

							
Bifocal Signature™					Sphere	Add	Cyl to
1.50 Standard Index SFST-28							
Signetek™ Processed 1.50 FT-28	GR				-10.00 to +6.00	+1.00 to +3.00	-6.00






								
Progressives					Sphere	Add	Cyl to	MFH
1.50 Standard Index								
KODAK Precise® Plus 1.50	GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus 1.50	GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Precise® PB 1.50	GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB 1.50	GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique™ 1.50	GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ 1.50	GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.50	GR				-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ 1.50	GR	BR			-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ HD 1.50	GR	BR			-8.00 to +4.00	+0.75 to +3.50	-6.00	13
1.53 TRIVEX®								
KODAK Precise® PB Trivex®	GR				-10.00 to +4.50	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB Trivex®	GR				-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Precise® Plus Trivex®	GR				-10.00 to +4.50	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus Trivex®	GR				-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique Trivex®	GR				-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ Trivex®	GR				-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD Trivex®	GR				-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ Trivex®	GR	BR			-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ HD Trivex®	GR	BR			-10.00 to +4.50	+0.75 to +3.50	-6.00	13
1.59 Polycarbonate								
KODAK Precise® Plus Poly	GR				-10.00 to +5.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus Poly	GR				-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ Poly	GR				-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Precise® PB Poly	GR				-10.00 to +5.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB Poly	GR				-10.00 to +5.00	+0.75 to +3.50	-6.00	13






KODAK Unique Poly	GR				-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD Poly	GR				-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ Poly	GR	BR			-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ HD Poly	GR	BR			-10.00 to +5.00	+0.75 to +3.50	-6.00	13
1.67 Super High Index								
KODAK Precise® PB 1.67	GR				-11.00 to +6.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB 1.67	GR				-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Precise® Plus Poly	GR				-11.00 to +6.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus Poly	GR				-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique	GR				-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ 1.67	GR				-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.67	GR				-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ 1.67	GR	BR			-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ HD 1.67	GR	BR			-11.00 to +6.00	+0.75 to +3.50	-6.00	13


					
Semi-Finished Single Vision				Sphere	Cyl to
1.50 Standard Index					
TRN 1.50 SV (72mm) UC/HC DIA.72	GR	BR		Base 0.50, 2.25, 4.25, 5.25, 6.25, 8.25	
TRN 1.50 DS (76mm) UC/HC DIA.72	GR	BR		Base 0.50, 8.75, 10.50	
TRN 1.50 DS (76mm) UC/HC DIA.76	GR	BR		Base 2.00, 3.50, 4.25, 5.25, 6.25, 7.25	
1.53 TRIVEX® Single Vision					
TOG 1.53 Trivex (75mm) SV UC/HC DIA.70	GR	BR		Base 1.00	
TOG 1.53 Trivex (75mm) SV UC/HC DIA.75	GR	BR		Base 2.50, 4.00, 5.00, 6.00, 7.00, 8.00	
1.59 Polycarbonate Single Vision					
TRN 1.59 Polycarbonate (79mm) SV HC DIA.79	GR	BR		Base 1.50, 3.25, 4.25, 4.75, 5.50, 7.00	
TRN 1.59 Polycarbonate (74mm) SV HC DIA.74	GR	BR		Base 7.50, 8.50	
1.60 High Index Single Vision					
TOG 1.60 SV (75mm) UC/HC DIA.65	GR	BR		Base 10.00, 11.00	
TOG 1.60 SV (75mm) UC/HC DIA.70	GR	BR		Base 9.00	
TOG 1.60 SV (75mm) UC/HC DIA.75	GR	BR		Base 1.00, 2.50, 4.00, 5.00, 6.00, 7.00, 8.00	
TOG 1.60 SV (80mm) UC/HC DIA.80	GR	BR		Base 1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.50	
1.67 Super High Index Single Vision					
TOG 1.67 SV (65mm) UC/HC DIA.65	GR	BR		Base 10.00, 11.00	
TOG 1.67 SV (70mm) UC/HC DIA.70	GR	BR		Base 9.00	
TOG 1.67 SV (75mm) UC/HC DIA.75	GR	BR		Base 1.00, 2.50, 4.00, 5.00, 6.00, 7.00, 8.00	
TOG 1.67 SV (80mm) UC/HC DIA.80	GR	BR		Base 1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.50	

					
Semi-Finished Single Vision				Sphere	Cyl to
1.50 Standard Index					
Unity SVx	GR			-12.25 to +6.00	-6.00
Unity SVxtra	GR			-12.25 to +6.00	-6.00
Unity SVxtreme	GR			-6.00 to +6.00	-6.00
1.53 TRIVEX® Single Vision					
Unity SVx	GR			-14.00 to +6.00	-6.00
Unity SVxtra	GR			-14.00 to +6.00	-6.00
Unity SVxtreme	GR			-6.50 to +6.00	-6.00




1.59 Polycarbonate Single Vision					
Unity SVx	GR			-14.00 to +6.00	-6.00
Unity SVxtra	GR			-14.00 to +6.00	-6.00
Unity SVxtreme	GR			-6.50 to +6.00	-6.00
1.67 Super High Index Single Vision					
Unity SVx	GR			-16.00 to +8.00	-6.00
Unity SVxtra	GR			-16.00 to +8.00	-6.00
Unity SVxtreme	GR			-6.00 to +6.00	-6.00



								
Progressives					Sphere	Add	Cyl to	MFH
1.50 Standard Index								
Unity Via II	GR				-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR				-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR				-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR				-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR				-7.00 to +6.00	+0.50 to +4.50	-6.00	12
1.53 TRIVEX®								
Unity Via II	GR				-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR				-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR				-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR				-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR				-5.00 to +4.00	+0.50 to +4.50	-6.00	12
1.59 Polycarbonate								
Unity Via II	GR				-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR				-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR				-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR				-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR				-6.50 to +6.00	+0.50 to +4.50	-6.00	12
1.67 Super High Index								
Unity Via II	GR				-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR				-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR				-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR				-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR				-6.00 to +5.00	+0.50 to +4.50	-6.00	12




						
Finished Single Vision					Sphere	Cyl to
1.50 Standard Index						
Younger Optics FSV	GR	BR			-4.00 to +4.00	-2.00
1.59 Polycarbonate Finished Single Vision						
Younger Optics FSV Poly	GR	BR			-6.00 to +4.00	-2.00

						
Semi-Finished Single Vision					Sphere	Cyl to
1.50 Standard Index						
Younger Optics SFSV (76mm)	GR	BR			-9.00 to +6.00	-4.00
1.53 TRIVEX® Single Vision						
Younger Optics Trilogy® Spherical SV	GR	BR			-8.00 to +7.50	-4.00
1.59 Polycarbonate Single Vision						
Younger Optics SV	GR	BR			-9.00 to +6.00	-4.00
1.67 Super High Index Single Vision						

Younger Optics Single Vision 1.67 MR-10™	GR	BR		-11.00 to +7.00	-4.00
1.74 Ultra High Index Single Vision					
Younger Optics Single Vision 1.74	GR	BR		-18.00 to +12.00	-4.00

YOUNGER OPTICS The Optical Lens Innovators	Transitions XTRACTIVE					
Bifocal				Sphere	Add	Cyl to
1.50 Standard Index SFST-28						
Younger Optics SFFT 28	GR			-11.00 to +6.00	+1.00 to +3.00	-4.00

YOUNGER OPTICS The Optical Lens Innovators	Transitions XTRACTIVE					
Digital				Sphere	Cyl to	
1.50 Index						
Camber™	GR	BR		-11.25 to +6.25	-4.00	
1.53 TRIVEX®						
Camber™	GR	BR		-11.00 to +6.25	-4.00	
1.59 Polycarbonate						
Camber™	GR	BR		-10.00 to +6.25	-4.00	
1.67 Super High Index						
Camber™	GR	BR		-14.00 to +8.00	-4.00	
1.74 Ultra High Index						
Camber™	GR	BR		-19.00 to +9.00	-4.00	

YOUNGER OPTICS The Optical Lens Innovators	Transitions XTRACTIVE						
Progressives				Sphere	Add	Cyl to	MFH
1.59 index polycarbonate							
Younger Image®	GR	BR		-9.00 to +7.00	+1.00 to +3.00	-4.00	18

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Style Mirrors are available where gray and brown *Transitions® XTRACTIVE®* are available. Specify *Transitions* lenses in style mirrors (**no substitutions**) with your lab to ensure authenticity.

Transitions, *Transitions Signature*, and *XTRACTIVE* are registered trademarks, and *XTRACTIVE Polarized*, the *Transitions* logo and *Transitions Light Intelligent Lenses* are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

Unless indicated otherwise, all registered trademarks and trademarks are the property of Essilor International and/or its subsidiaries in the United States and in other countries.

All other trademarks are the property of their respective owners.