









Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Design](#) / **Bifocals**

All Suppliers – Bifocals

	Transitions Gen#S	Transitions XTRACTIVE	Transitions XTRACTIVE POLARIZED			
Bifocal Signature™				Sphere	Add	Cyl to
1.50 Standard Index SFST-28						
Signetek™ Processed 1.50 FT-28		GR		-10.00 to +6.00	+1.00 to +3.00	-6.00

	Transitions Gen#S	Transitions XTRACTIVE	Transitions XTRACTIVE POLARIZED			
Bifocal Signature™				Sphere	Add	Cyl to
1.50 Standard Index SFST-28						
Younger Optics SFFT 28		GR		-11.00 to +6.00	+1.00 to +3.00	-4.00

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.