






Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Index](#) / **1.60 Index**

1.60 High Index

													Sphere	Cyl to
	Transitions gen#S								Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED		
Semi-Finished Single Vision														
1.60 High Index Single Vision														
Essilor Spherical Single Vision Thin&Lite® 1.60	GR	BR	GG	S	AME	AMB	E	R					-13.00 to +8.00	-4.00
Eyezen® +0 Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	-6.00
Eyezen® +1 Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	-6.00
Eyezen® +2 Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	-6.00
Eyezen® +3 Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	-6.00
Eyezen® +4 Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	-6.00
Eyezen® Start Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	-6.00

													Sphere	Add	Cyl to	MFH
	Transitions gen#S								Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED				
1.60 High Index																
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Thin&Lite® 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Thin&Lite® 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit Thin&Lite® 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max Thin&Lite® 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™ Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux Comfort® Thin&Lite® 1.60	GR												-12.00 to +8.00	+0.75 to +3.50	-4.00	17
Varilux Comfort Short™ Thin&Lite® 1.60	GR												-12.00 to +8.00	+0.75 to +3.00	-4.00	14
Essilor Ideal Advanced™ FIT Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Sport Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Driving Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	13
Essilor Ideal® Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Short Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Essilor Natural® Digital Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Ovation® Digital Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	17
SmallFit® Digital Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.00	-6.00	14

																
	Transitions gen#S								Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED				

1.60 SFSV													
Shamir Autograph® III Single Vision	GR	BR						GR				-14.00 to +8.75	-7.00
Shamir Relax™	GR	BR						GR				-14.00 to +8.75	-7.00

	Transitions Gen#S		Transitions XTRACTIVE		Transitions XTRACTIVE POLARIZED							
Blended Bifocal (digital)									Sphere	Add	Cyl to	MFH
1.60 High Index Blended Bifocal (digital)												
Shamir Duo	GR	BR					GR		-15.00 to +6.00	+0.75 to +4.00	-7.00	15

	Transitions Gen#S		Transitions XTRACTIVE		Transitions XTRACTIVE POLARIZED							
Progressives									Sphere	Add	Cyl to	MFH
1.60 High Index												
Shamir Autograph Intelligence™	GR	BR					GR		-14.75 to +8.75	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR					GR		-15.00 to +8.00	+0.75 to +3.50	-7.00	15, 18
Shamir Autograph® II+	GR	BR					GR		-15.00 to +6.00	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Golf™ SuperLite™ 1.6	GR	BR					GR		-6.00 to +5.25	+1.00 to +3.50	-7.00	19
Shamir Spectrum+™	GR	BR					GR		-15.00 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
Shamir Element™	GR	BR							-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19
Shamir GENESIS HD	GR	BR					GR		-6.00 to +6.00	+0.75 to +3.50	-7.00	16
Shamir Computer™	GR	BR					GR		-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper
Shamir Workspace™	GR	BR					GR		-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper
Shamir InTouch™	GR	BR					GR		-15.00 to +6.00	+0.75 to +3.50	-7.00	15, 18
Shamir Attitude III Sport	GR	BR					GR		-6.50 to +6.00	+0.75 to +3.50	-7.00	18
Shamir Attitude III Fashion	GR	BR					GR		-6.50 to +6.00	+0.75 to +3.50	-7.00	15, 18

	Transitions Gen#S		Transitions XTRACTIVE		Transitions XTRACTIVE POLARIZED					
Semi-Finished Single Vision									Sphere	Cyl to
1.60 High Index SFSV										
Signetek Processed SV	GR	BR							-10.00 to +7.50	-6.00
KODAK Digital SV	GR	BR							-10.00 to +7.50	-6.00
Crossbows Custom SV	GR	BR							-10.00 to +7.50	-6.00

	Transitions Gen#S		Transitions XTRACTIVE		Transitions XTRACTIVE POLARIZED							
Progressives									Sphere	Add	Cyl to	MFH
1.60 High Index												
KODAK Unique	GR	BR							-10.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD	GR	BR							-10.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ 1.60	GR	BR							-10.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.60	GR	BR							-10.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ 1.60	GR	BR							-10.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ HD 1.60	GR	BR							-10.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Custom U 1.60	GR	BR							-10.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Custom V 1.60	GR	BR							-10.00 to +6.00	+0.75 to +3.50	-6.00	13

	Transitions Gen#S		Transitions XTRACTIVE		Transitions XTRACTIVE POLARIZED					
Semi-Finished Single Vision									Sphere	Cyl to

Transitions, *Transitions Signature*, and *XTRActive* are registered trademarks, and *XTRActive Polarized*, the *Transitions* logo and *Transitions Light Intelligent Lenses* are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

Unless indicated otherwise, all registered trademarks and trademarks are the property of Essilor International and/or its subsidiaries in the United States and in other countries.

All other trademarks are the property of their respective owners.