

Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / **Transitions® Lenses PAG**


<

All Suppliers – All Products

	Transitions Gen [®] S								Transitions XTRACTIVE [®]			Transitions XTRACTIVE [®] POLARIZED [®]				
Finished Single Vision															Sphere	Cyl to
1.50 Standard Index FSV																
Essilor FSV Spherical	GR	BR							GR						-4.00 to +4.00	-2.00
Essilor FSV Aspheric	GR														-4.00 to +4.00	-2.00
Essilor FSV 1.50 Crizal® Sapphire™ HR	GR								GR						-4.00 to +4.00	-2.00
Essilor FSV 1.50 Crizal® Rock™	GR								GR						-4.00 to +4.00	-2.00
Essilor FSV 1.50 Crizal® Easy Pro	GR								GR						-4.00 to +4.00	-2.00
1.59 Polycarbonate Finished Single Vision																
Essilor FSV Spherical Airwear®	GR	BR													-6.00 to +4.00	-2.00
Essilor FSV Aspheric Airwear®	GR	BR							GR						-6.00 to +4.00	-2.00
Essilor FSV Airwear® Crizal® Sapphire™ HR	GR								GR						-6.00 to +4.00	-2.00
Essilor FSV Airwear® Crizal® Rock™	GR								GR						-6.00 to +4.00	-2.00
Essilor FSV Airwear® Crizal® Easy Pro	GR								GR						-6.00 to +4.00	-2.00

	Transitions Gen [®] S								Transitions XTRACTIVE [®]			Transitions XTRACTIVE [®] POLARIZED [®]				
Semi-Finished Single Vision															Sphere	Cyl to
1.50 Standard Index SFSV																
Essilor Spherical Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR			-10.50 to +10.00	-4.00
Essilor Stylistic® Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR			-10.00 to +7.50	-6.00
Eyezen® +0	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR			-10.00 to +7.50	-6.00
Eyezen® +1	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR			-10.00 to +7.50	-6.00
Eyezen® +2	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR			-10.00 to +7.50	-6.00
Eyezen® +3	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR			-10.00 to +7.50	-6.00
Eyezen® +4	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR			-10.00 to +7.50	-6.00
Eyezen® Start	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR			-10.00 to +6.00	-6.00
1.53 TRIVEX® Single Vision																
Essilor Spherical Single Vision TREXA™	GR	BR	GG						GR	BR					-10.00 to +6.00	-4.00
Essilor Eyecode™ Single Vision TREXA™	GR	BR	GG						GR	BR					-10.00 to +7.50	-6.00
Essilor FIT™ Single Vision TREXA™	GR	BR	GG						GR	BR					-10.00 to +7.50	-6.00
Essilor Azio™ Single Vision TREXA™	GR	BR	GG						GR	BR					-10.00 to +7.50	-6.00
Essilor 360°™ Single Vision TREXA™	GR	BR	GG						GR	BR					-10.00 to +7.50	-6.00
Eyezen® +0 TREXA™	GR	BR	GG						GR	BR					-10.00 to +7.50	-6.00
Eyezen® +1 TREXA™	GR	BR	GG						GR	BR					-10.00 to +7.50	-6.00
Eyezen® +2 TREXA™	GR	BR	GG						GR	BR					-10.00 to +7.50	-6.00
Eyezen® +3 TREXA™	GR	BR	GG						GR	BR					-10.00 to +7.50	-6.00
Eyezen® +4 TREXA™	GR	BR	GG						GR	BR					-10.00 to +7.50	-6.00
Eyezen® Start	GR	BR	GG						GR	BR					-10.00 to +6.00	-6.00
1.59 Polycarbonate Single Vision																
Essilor Spherical Single Vision Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR			-10.00 to +8.00	-4.00
Essilor Aspheric Single Vision Airwear®	GR	BR	GG						GR	BR	GG				-10.00 to +8.00	-4.00
Essilor Eyecode™ Aspheric SV Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR			-10.00 to +8.00	-6.00
Essilor FIT™ Aspheric SV Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR			-10.00 to +8.00	-6.00

Essilor AZIO™ Aspheric SV Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00
Essilor SV 360°™ Aspheric Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00
Essilor Stylistic® Single Vision Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00
Eyezen® +0 Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00
Eyezen® +1 Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00
Eyezen® +2 Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00
Eyezen® +3 Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00
Eyezen® +4 Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00
Eyezen® Start Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00
Eyezen®Kids	GR	BR	GG	S	AME	AMB	E	R					-10.00 to +6.00	-6.00
1.60 High Index Single Vision														
Essilor Spherical Single Vision Thin&Lite® 1.60	GR	BR	GG	S	AME	AMB	E	R					-13.00 to +8.00	-4.00
Eyezen® +0 Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	-6.00
Eyezen® +1 Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	-6.00
Eyezen® +2 Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	-6.00
Eyezen® +3 Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	-6.00
Eyezen® +4 Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	-6.00
Eyezen® Start Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	-6.00
1.67 Super High Index Single Vision														
Essilor Spherical Single Vision Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	-4.00
Essilor Aspheric Single Vision Thin&Lite® 1.67	GR	BR	GG										-16.00 to +8.00	-4.00
Essilor Eyecode™ SV Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	-6.00
Essilor FIT™ SV Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	-6.00
Essilor Azio™ SV Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	-6.00
Essilor SV 360°™ Spherical Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	-6.00
Essilor Stylistic® Single Vision Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	-6.00
Eyezen® +0 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	-6.00
Eyezen® +1 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	-6.00
Eyezen® +2 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	-6.00
Eyezen® +3 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	-6.00
Eyezen® +4 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	-6.00
Eyezen® Start Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	-6.00
1.74 Ultra High Index Single Vision														
Essilor Spherical Single Vision Thin&Lite® 1.74	GR	BR											-20.00 to +13.00	-4.00
Essilor Eyecode™ SV Thin&Lite® 1.74	GR	BR											-16.00 to +13.00	-6.00
Essilor FIT™ SV Thin&Lite® 1.74	GR	BR											-16.00 to +13.00	-6.00
Essilor 360°™ SV Thin&Lite® 1.74	GR	BR											-16.00 to +13.00	-6.00
Eyezen® +0 Thin&Lite® 1.74	GR	BR											-20.00 to +13.00	-6.00
Eyezen® +1 Thin&Lite® 1.74	GR	BR											-20.00 to +13.00	-6.00
Eyezen® +2 Thin&Lite® 1.74	GR	BR											-20.00 to +13.00	-6.00
Eyezen® +3 Thin&Lite® 1.74	GR	BR											-20.00 to +13.00	-6.00
Eyezen® +4 Thin&Lite® 1.74	GR	BR											-20.00 to +13.00	-6.00
Eyezen® Start Thin&Lite® 1.74	GR	BR											-20.00 to +13.00	-6.00

															
Progressives												Sphere	Add	Cyl to	MFH
1.50 Standard Index															
Varilux® XR Track™	GR	BR	GG						GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR	GG						GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™	GR	BR	GG						GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14

Varilux® X FIT™	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Varilux® Liberty™ 3.0	GR	BR											-10.00 to +7.50	+0.75 to +1.50	-6.00	14
Varilux Comfort®	GR	BR							GR	BR	GG		-10.00 to +7.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short™	GR	BR							GR	BR	GG		-10.00 to +7.00	+0.75 to +3.00	-6.00	14
Varilux® Stylistic™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ FIT	GR	BR	GG	S	AME	AMB	E	R					-10.00 to +7.50	+0.75 to +3.50	-4.00	17
Essilor Ideal Advanced™ Wrap	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-4.00	17
Essilor Ideal® FIT	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	18
Essilor Ideal® Driving	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	13
Essilor Ideal®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal® Short	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Natural® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +7.50	+0.75 to +3.50	-6.00	15
Essilor Accolade™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Brillenz®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Brillenz® Short	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Ovation® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Ovation®	GR	BR											-10.00 to +6.00	+0.75 to +3.50	-4.00	17
Essilor SmallFit® Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +7.50	+0.75 to +3.00	-6.00	14

	Transitions gen ⁸ S								Transitions XTRACTIVE [®]			Transitions XTRACTIVE [®] PC-ARISE™						
Progressives															Sphere	Add	Cyl to	MFH
1.53 TRIVEX®																		
Varilux® XR Track™	GR	BR													-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR													-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ TREXA™	GR	BR													-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ TREXA™	GR	BR													-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ TREXA™	GR	BR													-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode TREXA™	GR	BR							GR	BR					-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT TREXA™	GR	BR							GR	BR					-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ TREXA™	GR	BR							GR	BR					-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit TREXA™	GR	BR							GR	BR					-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max TREXA™	GR	BR							GR	BR					-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ TREXA™	GR	BR	GG						GR	BR					-10.00 to +7.50	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ TREXA™	GR	BR	GG						GR	BR					-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™ TREXA™	GR	BR	GG						GR	BR					-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ TREXA™	GR	BR	GG						GR	BR					-10.00 to +7.50	+0.75 to +3.00	-6.00	14

	Transitions Gen [®] S								Transitions [®] XTRACTIVE [™]			Transitions [®] XTRACTIVE [™] POLARIZED [®]				
Progressives													Sphere	Add	Cyl to	MFH
1.59 Polycarbonate																
Varilux [®] XR Track [™]	GR	BR	GG	S	AME	AMB	E	R			GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux [®] XR Design [™]	GR	BR	GG	S	AME	AMB	E	R			GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux [®] X 4D [™] Airwear [®]	GR	BR	GG	S	AME	AMB	E	R			GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux [®] X FIT [™] Airwear [®]	GR	BR	GG	S	AME	AMB	E	R			GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux [®] X Design [™] Airwear [®]	GR	BR	GG	S	AME	AMB	E	R			GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux [®] Physio [®] W3+ Eyecode [™] Airwear [®]	GR	BR	GG	S	AME	AMB	E	R			GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux [®] Physio [®] W3+ FIT Airwear [®]	GR	BR	GG	S	AME	AMB	E	R			GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux [®] Physio [®] W3+ Airwear [®]	GR	BR	GG	S	AME	AMB	E	R			GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort [®] Max Fit Airwear [®]	GR	BR	GG	S	AME	AMB	E	R			GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort [®] Max Airwear [®]	GR	BR	GG	S	AME	AMB	E	R			GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx [™] Airwear [®]	GR	BR	GG	S	AME	AMB	E	R			GG	GR	-10.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx [™] Airwear [®]	GR	BR	GG	S	AME	AMB	E	R			GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Varilux Comfort DRx [™] Airwear [®]	GR	BR	GG	S	AME	AMB	E	R			GG	GR	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx [™] Airwear [®]	GR	BR	GG	S	AME	AMB	E	R			GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux [®] Liberty [™] 3.0	GR	BR											-10.00 to +8.00	+0.75 to +1.50	-6.00	14
Varilux Comfort [®] Airwear [®]	GR	BR	GG								GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short [™] Airwear [®]	GR	BR	GG								GG		-10.00 to +6.00	+0.75 to +3.00	-6.00	14
Varilux [®] Stylistic [™] Airwear [®]	GR	BR	GG	S	AME	AMB	E	R			GG	GR	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced [™] FIT	GR	BR	GG	S	AME	AMB	E	R					-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced [™] Wrap	GR	BR	GG	S	AME	AMB	E	R			GG	GR	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced [™]	GR	BR	GG	S	AME	AMB	E	R			GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Essilor Ideal [®] FIT Airwear [®]	GR	BR	GG	S	AME	AMB	E	R			GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal [®] Sport Airwear [®]	GR	BR	GG	S	AME	AMB	E	R			GG	GR	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal [®] Computer Airwear [®]	GR	BR	GG	S	AME	AMB	E	R			GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	18
Essilor Ideal [®] Driving Airwear [®]	GR	BR	GG	S	AME	AMB	E	R			GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal [®] Airwear [®]	GR	BR	GG	S	AME	AMB	E	R			GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Ideal [®] Short Airwear [®]	GR	BR	GG	S	AME	AMB	E	R			GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Natural [®] Digital Airwear [®]	GR	BR	GG	S	AME	AMB	E	R			GG		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom [™] Airwear [®]	GR	BR	GG	S	AME	AMB	E	R			GG		-10.00 to +8.00	+0.75 to +3.50	-6.00	15



	Transitions® Gen [®] S								Transitions® XTRACTIVE®			Transitions® XTRACTIVE® POLARIZED®				
Progressives													Sphere	Add	Cyl to	MFH
1.67 Super High Index																
Varilux® XR Track™	GR	BR	GG	S	AME	AMB	E	R			GG		-12.00 to +9.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR	GG	S	AME	AMB	E	R			GG		-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R			GG		-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R			GG		-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R			GG		-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode Thin&Lite®	GR	BR	GG	S	AME	AMB	E	R			GG		-14.00 to +9.00	+0.75 to +3.50 &	-6.00	14

																		
Progressives													Sphere	Add	Cyl to	MFH		
1.74 Ultra High Index																		
Varilux® XR Track™	GR	BR											-14.00 to +13.00	+0.75 to +4.00	-6.00	NVB Measurement		
Varilux® XR Design™	GR	BR											-14.00 to +13.00	+0.75 to +4.00	-6.00	14		
Varilux® X 4D™ Thin&Lite® 1.74	GR	BR											-14.00 to +13.00	+0.75 to +4.00	-6.00	14		
Varilux® X FIT™ Thin&Lite® 1.74	GR	BR											-14.00 to +13.00	+0.75 to +4.00	-6.00	14		
Varilux® X Design™ Thin&Lite® 1.74	GR	BR											-14.00 to +13.00	+0.75 to +4.00	-6.00	14		
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.74	GR	BR											-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14		
Varilux® Physio® W3+ FIT Thin&Lite® 1.74	GR	BR											-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14		
Varilux® Physio® W3+ Thin&Lite® 1.74	GR	BR											-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14		
Varilux Physio DRx™ Thin&Lite® 1.74	GR	BR											-14.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	17		
Varilux Physio Short DRx™ Thin&Lite® 1.74	GR	BR											-14.00 to +13.00	+0.75 to +3.50	-6.00	14		
Varilux Comfort DRx™ Thin&Lite® 1.74	GR	BR											-20.00 to +13.00	+0.75 to +3.50	-6.00	17		

Varilux Comfort Short DRx™ Thin&Lite® 1.74	GR	BR								-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ FIT Thin&Lite® 1.74	GR	BR								-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Thin&Lite® 1.74	GR	BR								-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.74	GR	BR								-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport Thin&Lite® 1.74	GR	BR								-20.00 to +13.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Computer Thin&Lite® 1.74	GR	BR								-20.00 to +13.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Driving Thin&Lite® 1.74	GR	BR								-20.00 to +13.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Thin&Lite® 1.74	GR	BR								-20.00 to +13.00	+0.75 to +3.50	-6.00	13
Ovation® Digital Thin&Lite® 1.74	GR	BR								-20.00 to +13.00	+0.75 to +3.50	-6.00	14
Essilor Accolade Freedom™ Thin&Lite® 1.74	GR	BR								-20.00 to +13.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Thin&Lite® 1.74	GR	BR								-14.00 to +9.00	+0.75 to +3.50	-6.00	17

[illegible][illegible]




HOYA SFSV 1.74	GR!	BR!							GR	BR			-13.00 to +10.00	-4.00
HOYA SV IQ 1.74	GR!	BR!											-13.00 to +10.00	-4.00
Sync III	GR!	BR!											-13.00 to +8.00	-4.00

HOYA	Transitions Gen [®] S							Transitions XTRACTIVE [®]			Transitions XTRACTIVE [®] POLARIZED [®]				
Progressives												Sphere	Add	Cyl to	MFH
1.50 Standard Index															
HOYALUX iD LifeStyle [®] 4 Indoor 1.50	GR!	BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle [®] 4 Outdoor 1.50	GR!	BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle [®] 4 Urban 1.50	GR!	BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	14
Amp HD3 1.50	GR!	BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	18
Amp Mini HD3 1.50	GR!	BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX ARRAY™ 11, 13, 15, 17, VL 1.50	GR	BR						GR	BR		GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL 1.50	GR	BR						GR	BR		GR	-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ 2 11, 13, 15, 17, VL 1.50	GR	BR						GR	BR		GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ 2 Wrap 11, 13, 15, 17, VL 1.50	GR	BR						GR	BR		GR	-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS 1.50	GR	BR						GR	BR		GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS 1.50	GR	BR						GR	BR		GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	14
Summit cd BKS 1.50	GR	BR						GR	BR		GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	14
Summit ecp BKS 1.50	GR	BR						GR	BR		GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	18
Hoya GP Wide 1.50	GR	BR						GR	BR		GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	18
Tact 40 BKS 1.50	GR	BR						GR	BR		GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	18
Tact 60 BKS 1.50	GR	BR						GR	BR		GR	-10.00 to +3.75	+0.75 to +4.00	-6.00	18
1.53 TRIVEX[®]															
HOYALUX iD LifeStyle [®] 4 Indoor Phoenix™	GR!	BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle [®] 4 Outdoor Phoenix™	GR!	BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle [®] 4 Urban Phoenix™	GR!	BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle [®] 3 Indoor Phoenix™	GR!	BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle [®] 3 Outdoor Phoenix™	GR!	BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle [®] 3 Urban Phoenix™	GR!	BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	14
Amp HD3 Phoenix™	GR!	BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	18
Amp Mini HD3 Phoenix™	GR!	BR!										-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX ARRAY™ 11, 13, 15, 17, VL Phoenix™	GR	BR						GR	BR		GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL Phoenix™	GR	BR						GR	BR		GR	-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13, 15, 17, VL Phoenix™	GR	BR						GR	BR		GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11, 13, 15, 17, VL Phoenix™	GR	BR						GR	BR		GR	-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS Phoenix™	GR	BR						GR	BR		GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS Phoenix™	GR	BR						GR	BR		GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Summit cd BKS Phoenix™	GR	BR						GR	BR		GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	18
Summit ecp BKS Phoenix™	GR	BR						GR	BR		GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	18
Hoya GP Wide Phoenix™	GR	BR						GR	BR		GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 40 BKS Phoenix™	GR	BR						GR	BR		GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 60 BKS Phoenix™	GR	BR						GR	BR		GR	-10.00 to +6.00	+0.75 to +4.00	-6.00	18
1.59 Polycarbonate															
HOYALUX iD LifeStyle [®] 4 Indoor Poly	GR!	BR!										-6.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle [®] 4 Outdoor Poly	GR!	BR!										-6.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle [®] 4 Urban Poly	GR!	BR!										-6.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle [®] 3 Indoor Poly	GR!	BR!										-6.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle [®] 3 Outdoor Poly	GR!	BR!										-6.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle [®] 3 Urban Poly	GR!	BR!										-6.00 to +6.00	+0.75 to +4.00	-6.00	14

Amp HD3 Poly	GR!	BR!									-6.00 to +6.00	+0.75 to +4.00	-6.00	18
Amp Mini HD3 Poly	GR!	BR!									-6.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX ARRAY™ 11, 13, 15, 17, VL Poly	GR	BR					GR	BR		GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL Poly	GR	BR					GR	BR		GR	-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13, 15, 17, VL Poly	GR	BR					GR	BR		GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11, 13, 15, 17, VL Poly	GR	BR					GR	BR		GR	-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS Poly	GR	BR					GR	BR		GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS Poly	GR	BR					GR	BR		GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	14
Summit cd BKS Poly	GR	BR					GR	BR		GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	14
Summit ecp BKS Poly	GR	BR					GR	BR		GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	18
Hoya GP Wide Poly	GR	BR					GR	BR		GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 40 BKS Poly	GR	BR					GR	BR		GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 60 BKS Poly	GR	BR					GR	BR		GR	-13.00 to +6.00	+0.75 to +4.00	-6.00	18
1.60 High Index														
HOYALUX iD LifeStyle®4 Indoor 1.60	GR!	BR!									-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Outdoor 1.60	GR!	BR!									-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Urban 1.60	GR!	BR!									-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Indoor 1.60	GR!	BR!									-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Outdoor 1.60	GR!	BR!									-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Urban 1.60	GR!	BR!									-11.00 to +8.00	+0.75 to +4.00	-6.00	14
Amp HD3 1.60	GR!	BR!									-11.00 to +8.00	+0.75 to +4.00	-6.00	14
Amp Mini HD3 1.60	GR!	BR!									-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX ARRAY 11, 13, 15, 17, VL 1.60	GR	BR									-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY Wrap 11,13,15,17,VL 1.60	GR	BR									-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13, 15, 17, VL 1.60	GR	BR									-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11,13,15,17,VL 1.60	GR	BR									-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS 1.60	GR	BR									-11.00 to +8.00	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS 1.60	GR	BR									-11.00 to +8.00	+0.75 to +4.00	-6.00	14
Summit cd BKS 1.60	GR	BR									-11.00 to +8.00	+0.75 to +4.00	-6.00	14
Summit ecp BKS 1.60	GR	BR									-11.00 to +8.00	+0.75 to +4.00	-6.00	18
Hoya GP Wide 1.60	GR	BR									-11.00 to +8.00	+0.75 to +4.00	-6.00	18
Tact 40 BKS 1.60	GR	BR									-11.00 to +8.00	+0.75 to +4.00	-6.00	18
Tact 60 BKS 1.60	GR	BR									-11.00 to +8.00	+0.75 to +4.00	-6.00	18
1.67 Super High Index														
HOYALUX iD LifeStyle®4 Indoor 1.67	GR!	BR!									-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Outdoor 1.67	GR!	BR!									-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Urban 1.67	GR!	BR!									-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle 3 Indoor 1.67	GR!	BR!									-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle 3 Outdoor 1.67	GR!	BR!									-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle 3 Urban 1.67	GR!	BR!									-13.00 to +8.00	+0.75 to +4.00	-6.00	14
Amp HD3 1.67	GR!	BR!									-13.00 to +8.00	+0.75 to +4.00	-6.00	18
Amp Mini HD3 1.67	GR!	BR!									-13.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX ARRAY 11, 13,15,17,VL 1.67	GR	BR					GR	BR		GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY Wrap 11,13,15,17,VL 1.67	GR	BR					GR	BR		GR	-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13,15,17,VL 1.67	GR	BR					GR	BR		GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11,13,15,17,VL 1.67	GR	BR					GR	BR		GR	-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS 1.67	GR	BR					GR	BR		GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS 1.67	GR	BR					GR	BR		GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	14
Summit cd BKS 1.67	GR	BR					GR	BR		GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	14
Summit ecp BKS 1.67	GR	BR					GR	BR		GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	18
Hoya GP Wide 1.67	GR	BR					GR	BR		GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	18
Tact 40 BKS 1.67	GR	BR					GR	BR		GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	18
Tact 60 BKS 1.67	GR	BR					GR	BR		GR	-13.00 to +7.50	+0.75 to +4.00	-6.00	18
1.74 ultra High Index														

HOYALUX iD LifeStyle®4 Indoor 1.74	GR! BR!									-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle®4 Outdoor 1.74	GR! BR!									-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle®4 Urban 1.74	GR! BR!									-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle 3 Indoor 1.74	GR! BR!									-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle 3 Outdoor 1.74	GR! BR!									-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX iD LifeStyle 3 Urban 1.74	GR! BR!									-13.00 to +8.00	+0.75 to +3.50	-6.00	14
Amp HD3 1.74	GR! BR!									-13.00 to +8.00	+0.75 to +3.50	-6.00	14
Amp Mini HD3 1.74	GR! BR!									-13.00 to +8.00	+0.75 to +3.50	-6.00	14
HOYALUX ARRAY 11, 13, 15, 17, VL 1.74	GR BR									-13.00 to +8.00	+0.75 to +3.50	-6.00	14, 17, 19, 21
HOYALUX ARRAY Wrap 11, 13, 15, 17, VL 1.74	GR BR									-6.00 to +4.00	+0.75 to +3.50	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13, 15, 17, VL 1.74	GR BR									-13.00 to +8.00	+0.75 to +3.50	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 WRAP 11, 13, 15, 17, VL 1.74	GR BR									-6.00 to +4.00	+0.75 to +3.50	-6.00	14, 17, 19, 21
Amplitude BKS	GR BR									-13.00 to +8.00	+0.75 to +3.50	-6.00	19
Amplitude Mini BKS	GR BR									-13.00 to +8.00	+0.75 to +3.50	-6.00	14
Summit cd BKS 1.74	GR BR									-13.00 to +8.00	+0.75 to +3.50	-6.00	14
Summit ecp BKS 1.74	GR BR									-13.00 to +8.00	+0.75 to +3.50	-6.00	19
Hoya GP Wide 1.74	GR BR									-13.00 to +8.00	+0.75 to +3.50	-6.00	18
Tact 40 BKS 1.74	GR BR									-13.00 to +8.00	+0.75 to +3.50	-6.00	18
Tact 60 BKS 1.74	GR BR									-13.00 to +8.00	+0.75 to +3.50	-6.00	18

	Transitions Gen#S								Transitions XTRACTIVE®			Transitions XTRACTIVE® POLARIZED®			
Semi-Finished Single Vision													Sphere	Cyl to	
1.50 Standard Index Single Vision															
Nikon SeeMax Infinite	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
Nikon Infinite Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
Nikon Smart Boost	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
1.53 TRIVEX® Single Vision															
Nikon SeeMax Infinite	GR	BR	GG										-10.00 to +7.00	-6.00	
Nikon Infinite Single Vision	GR	BR	GG										-10.00 to +7.00	-6.00	
Nikon Smart Boost	GR	BR	GG										-10.00 to +7.00	-6.00	
1.59 Polycarbonate Single Vision															
Nikon SeeMax Infinite	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
Nikon Infinite Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
Nikon Smart Boost	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	-6.00	
1.60 High Index Single Vision															
Nikon SeeMax Infinite	GR	BR											-12.00 to +8.00	-6.00	
Nikon Infinite Single Vision	GR	BR											-12.00 to +8.00	-6.00	
Nikon Smart Boost	GR	BR											-12.00 to +8.00	-6.00	
1.67 Super High Index Single Vision															
Nikon SeeMax Infinite	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	-6.00	
Nikon Infinite Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	-6.00	
Nikon Smart Boost	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	-6.00	
1.74 ultra High Index Single Vision															
Nikon SeeMax Infinite	GR	BR											-16.00 to +10.00	-6.00	
Nikon Infinite Single Vision	GR	BR											-16.00 to +10.00	-6.00	
Nikon Smart Boost	GR	BR											-16.00 to +10.00	-6.00	

																		
Progressives													Sphere	Add	Cyl to	MFH		
1.50 Standard Index																		
Nikon SeeMax Ultimate	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17		
Nikon Presio Power	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17		

Nikon Presio i Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
Nikon Move	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
1.53 TRIVEX®																
Nikon SeeMax Ultimate	GR	BR	GG										-10.00 to +7.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR	GG										-10.00 to +7.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR	GG										-10.00 to +7.00	+0.50 to +4.00	-6.00	15, 17
Nikon Move	GR	BR	GG										-10.00 to +7.00	+0.50 to +4.00	-6.00	15, 17
1.59 Polycarbonate																
Nikon SeeMax Ultimate	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-11.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-11.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-11.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
Nikon Move	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-11.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
1.60 High Index																
Nikon SeeMax Ultimate	GR	BR											-12.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR											-12.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR											-12.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
Nikon Move	GR	BR											-12.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
1.67 Super High Index																
Nikon SeeMax Ultimate	GR	BR	GG	S	AME	AMB	E		GR	BR	GG	GR	-14.00 to +9.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR	GG	S	AME	AMB	E		GR	BR	GG	GR	-14.00 to +9.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR	GG	S	AME	AMB	E		GR	BR	GG	GR	-14.00 to +9.00	+0.50 to +4.00	-6.00	15, 17
Nikon Move	GR	BR	GG	S	AME	AMB	E		GR	BR	GG	GR	-14.00 to +9.00	+0.50 to +4.00	-6.00	15, 17
1.74 Ultra High Index																
Nikon SeeMax Ultimate	GR	BR											-16.00 to +10.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR											-16.00 to +10.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR											-16.00 to +10.00	+0.50 to +4.00	-6.00	15, 17

SEIKO							
Semi-Finished Single Vision						Sphere	Cyl to
1.53 TRIVEX® SFSV							
Seiko PL 1.53/Trivex Transitions®						GR	-10.00 to +7.00 -5.00
1.67 Super High Index Single Vision							
Seiko PL 1.67 Transitions®						GR	-12.00 to +10.00 -5.00

							
Finished Single Vision						Sphere	Cyl to
1.59 Polycarbonate FSV							
Shamir FSV with Glacier PLUS™ AR	GR						-6.00 to +4.00 -2.00

							
Semi-Finished Single Vision						Sphere	Cyl to
1.50 SFSV							
Shamir Autograph® III Single Vision	GR	BR				GR BR	-12.00 to +7.25 -6.00
Shamir Relax™	GR	BR				GR BR	-14.00 to +8.00 -6.00
1.53 SFSV							
Shamir Autograph® III Single Vision	GR	BR				GR BR	-14.00 to +8.00 -6.50
Shamir Relax™	GR	BR				GR BR	-14.00 to +8.00 -6.50
1.59 SFSV							
Shamir Autograph® III Single Vision	GR	BR				GR BR	-14.00 to +8.75 -7.00
Shamir Relax™	GR	BR				GR BR	-14.75 to +6.00 -7.00
1.60 SFSV							

Shamir Autograph® III Single Vision	GR	BR							GR				-14.00 to +8.75	-7.00
Shamir Relax™	GR	BR							GR				-14.00 to +8.75	-7.00
1.67 SFSV														
Shamir Autograph® III Single Vision	GR	BR							GR	BR		GR	-16.00 to +16.00	-8.00
Shamir Relax™	GR	BR							GR	BR		GR	-16.75 to +16.00	-8.00
1.74 SFSV														
Shamir Autograph® III Single Vision	GR	GR											-20.00 to +13.50	-9.00

																		
Blended Bifocal (digital)															Sphere	Add	Cyl to	MFH
1.50 Standard Index Blended Bifocal (digital)																		
Shamir Duo	GR	BR							GR						-9.00 to +6.00	+0.75 to +4.00	-6.00	15
1.53 Trivex® Blended Bifocal (digital)																		
Shamir Duo	GR	BR							GR						-13.25 to +6.00	+0.75 to +4.00	-6.50	15
1.59 Polycarbonate Blended Bifocal (digital)																		
Shamir Duo	GR	BR							GR						-14.75 to +6.00	+0.75 to +4.00	-7.00	15
1.60 High Index Blended Bifocal (digital)																		
Shamir Duo	GR	BR							GR						-15.00 to +6.00	+0.75 to +4.00	-7.00	15
1.67 Super High Index Blended Bifocal (digital)																		
Shamir Duo	GR	BR							GR						-16.75 to +8.00	+0.75 to +4.00	-8.00	15
1.74 Ultra High Index Blended Bifocal (digital)																		
Shamir Duo	GR	BR													-18.75 to +12.00	+0.75 to +4.00	-9.00	15

																		
Progressives															Sphere	Add	Cyl to	MFH
1.50 Standard Index																		
Shamir Autograph Intelligence™	GR	BR							GR			GR			-12.25 to +6.50	+0.75 to +4.00	-6.00	11, 13,15,18
Shamir Autograph® III	GR	BR							GR			GR			-12.25 to +6.00	+0.75 to +4.00	-6.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR							GR			GR			-12.25 to +6.00	+0.75 to +4.00	-6.00	11, 13, 15, 18
Shamir Intouch™	GR	BR							GR			GR			-12.25 to +6.00	+0.75 to +3.50	-6.00	14, 16, 18
Shamir Spectrum+™	GR	BR							GR			GR			-12.25 to +6.00	+0.75 to +3.50	-6.00	14, 16, 18
Shamir Element™	GR	BR							GR			GR			-7.00 to +4.00	+0.75 to +3.50	-6.00	16, 19
Shamir GENESIS HD	GR	BR							GR			GR			-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir Computer™	GR	BR							GR			GR			-7.00 to +4.00	+0.75 to +2.00	-6.00	14 Lower, 10 Upper
Shamir Workspace™	GR	BR							GR			GR			-6.00 to +6.00	+0.75 to +2.50	-6.00	14 Lower, 10 Upper
Shamir Attitude III Sport	GR	BR							GR			GR			-6.50 to +5.25	+0.75 to +3.50	-6.00	18
Shamir Attitude III Fashion	GR	BR							GR			GR			-6.50 to +5.25	+0.75 to +3.50	-6.00	15, 18
1.53 TRIVEX®																		
Shamir Autograph Intelligence™	GR	BR							GR			GR			-13.25 to +7.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
Shamir Autograph® III	GR	BR							GR			GR			-13.25 to +7.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
Shamir Autograph® II+	GR	BR							GR			GR			-13.25 to +6.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
Shamir Intouch™	GR	BR							GR			GR			-13.25 to +6.00	+0.75 to +3.50	-6.50	14, 16, 18
Shamir Golf™	GR	BR							GR			GR			-5.00 to +4.00	+1.00 to +3.50	-6.50	19
Shamir Spectrum+™	GR	BR							GR			GR			-13.25 to +6.00	+0.75 to +3.50	-6.50	14, 16, 18
Shamir Element™	GR	BR							GR			GR			-13.25 to +6.00	+0.75 to +3.50	-6.50	16, 19
Shamir GENESIS HD	GR	BR							GR			GR			-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir Computer™	GR	BR							GR			GR			-6.00 to +8.00	+0.75 to +2.00	-6.50	14 Lower, 10 Upper
Shamir Workspace™	GR	BR							GR			GR			-6.00 to +8.00	+0.75 to +2.50	-6.50	14 Lower, 10 Upper
Shamir Attitude III Sport	GR	BR							GR			GR			-6.50 to +6.00	+0.75 to +3.50	-6.50	18
Shamir Attitude III Fashion	GR	BR							GR			GR			-6.50 to +6.00	+0.75 to +3.50	-6.50	15, 18

1.59 Polycarbonate														
Shamir Autograph Intelligence™	GR	BR						GR		GR	-14.75 to +7.50	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR						GR		GR	-14.75 to +8.00	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR						GR		GR	-14.75to +6.50	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Intouch™	GR	BR						GR		GR	-14.75 to +6.00	+0.75 to +3.50	-7.00	15, 18
Shamir Golf™	GR	BR						GR		GR	-6.00 to +5.25	+1.00 to +3.50	-7.00	19
Shamir Spectrum+™	GR	BR						GR		GR	-14.75 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
Shamir Element™	GR	BR						GR		GR	-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19
Shamir GENESIS HD	GR	BR						GR		GR	-10.00 to +6.00	+0.75 to +3.50	-7.00	16
Shamir Computer™	GR	BR						GR		GR	-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper
Shamir Workspace™	GR	BR						GR		GR	-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper
Shamir Attitude III Sport	GR	BR						GR		GR	-6.50 to +6.00	+0.75 to +3.50	-6.50	18
Shamir Attitude III Fashion	GR	BR						GR		GR	-6.50 to +6.00	+0.75 to +3.50	-6.50	15, 18
1.60 High Index														
Shamir Autograph Intelligence™	GR	BR						GR			-14.75 to +8.75	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR						GR			-15.00 to +8.00	+0.75 to +3.50	-7.00	15,18
Shamir Autograph® II+	GR	BR						GR			-15.00 to +6.00	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Golf™ SuperLite™ 1.6	GR	BR						GR			-6.00 to +5.25	+1.00 to +3.50	-7.00	19
Shamir Spectrum+™	GR	BR						GR			-15.00 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
Shamir Element™	GR	BR									-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19
Shamir GENESIS HD	GR	BR						GR			-6.00 to +6.00	+0.75 to +3.50	-7.00	16
Shamir Computer™	GR	BR						GR			-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper
Shamir Workspace™	GR	BR						GR			-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper
Shamir InTouch™	GR	BR						GR			-15.00 to +6.00	+0.75 to +3.50	-7.00	15, 18
Shamir Attitude III Sport	GR	BR						GR			-6.50 to +6.00	+0.75 to +3.50	-7.00	18
Shamir Attitude III Fashion	GR	BR						GR			-6.50 to +6.00	+0.75 to +3.50	-7.00	15, 18
1.67 Super High Index														
Shamir Autograph Intelligence™	GR	BR						GR		GR	-16.75 to +10.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR						GR		GR	-16.75 to +9.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR						GR		GR	-16.75 to +8.25	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Spectrum+™	GR	BR						GR		GR	-16.75 to +8.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Element™ SuperLite™ 1.67	GR	BR						GR		GR	-12.50 to +8.00	+0.75 to +3.50	-8.00	16, 19
Shamir GENESIS HD	GR	BR						GR		GR	-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir InTouch™	GR	BR						GR		GR	-16.75 to +8.00	+0.75 to +3.50	-8.00	15, 18
Shamir Attitude III Sport	GR	BR						GR		GR	-6.50 to +6.00	+0.75 to +3.50	-8.00	18
Shamir Attitude III Fashion	GR	BR						GR		GR	-6.50 to +6.00	+0.75 to +3.50	-8.00	15, 18
1.74 Ultra High Index														
Shamir Autograph Intelligence™	GR	BR									-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR									-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR									-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir GENESIS HD	GR	BR									-6.00 to +6.00	+0.75 to +4.00	-4.00	16
Shamir Spectrum+™	GR	BR									-18.75 to +11.50	+0.75 to +3.50	-9.00	14, 16, 18
Shamir InTouch™	GR	BR									-18.75 to +15.50	+0.75 to +3.50	-9.00	15, 18

 Transitions Gen#S Transitions XTRACTIVE Transitions XTRACTIVE POLARIZED														
Semi-Finished Single Vision													Sphere	Cyl to
1.50 Standard Index SFSV														
Signetek Processed SV	GR	BR	GG	S	AME	AMB	E	R	GR				-8.00 to +6.00	-6.00
KODAK Digital SV	GR	BR	GG	S	AME	AMB	E	R	GR				-8.00 to +6.00	-6.00
KODAK PowerUp Lens (.40/.66/1.00D)	GR	BR											-8.00 to +6.00	-6.00
Crossbows Custom SV	GR	BR											-8.00 to +6.00	-6.00
1.53 TRIVEX® Single Vision														
Signetek Processed SV	GR	BR	GG	S	AME	AMB	E	R	GR				-10.00 to +6.00	-6.00
KODAK Digital SV	GR	BR	GG	S	AME	AMB	E	R	GR				-10.00 to +6.00	-6.00

Crossbows Custom SV	GR	BR										-10.00 to +6.00	-6.00
1.59 Polycarbonate Single Vision													
Signetek Processed SV	GR	BR	GG	S	AME	AMB	E	R	GR		GR	-10.00 to +6.50	-6.00
KODAK Digital SV	GR	BR	GG	S	AME	AMB	E	R	GR		GR	-10.00 to +6.50	-6.00
KODAK PowerUp Lens (.40/.66/1.00D)	GR	BR										-10.00 to +6.50	-6.00
Crossbows Custom SV	GR	BR										-10.00 to +6.50	-6.00
1.60 High Index Single Vision													
Signetek Processed SV	GR	BR										-10.00 to +7.50	-6.00
KODAK Digital SV	GR	BR										-10.00 to +7.50	-6.00
Crossbows Custom SV	GR	BR										-10.00 to +7.50	-6.00
1.67 Super High Index Single Vision													
Signetek Processed SV	GR	BR	GG	S	AME	AMB	E	R	GR			-11.00 to +7.50	-6.00
KODAK Digital SV	GR	BR	GG	S	AME	AMB	E	R	GR			-11.00 to +7.50	-6.00
KODAK PowerUp Lens (.40/.66/1.00D)	GR	BR										-11.00 to +6.00	-6.00
Crossbows Custom SV	GR	BR										-11.00 to +7.50	-6.00
1.74 Ultra High Index Single Vision													
Signetek™ Processed SV	GR	BR										-13.00 to +9.00	-6.00
KODAK Digital SV	GR	BR										-13.00 to +9.00	-6.00

																
Bifocal Signature™													Sphere	Add	Cyl to	
1.50 Standard Index SFST-28																
Signetek™ Processed 1.50 FT-28								GR					-10.00 to +6.00	+1.00 to +3.00	-6.00	

	Transitions Gen [®] S							Transitions [®] XTRACTIVE [®]			Transitions [®] XTRACTIVE [®] POLARIZED [®]						
Progressives														Sphere	Add	Cyl to	MFH
1.50 Standard Index																	
KODAK Precise® 1.50 SF	GR													-9.00 to +5.00	+1.00 to +3.00	-4.00	17
KODAK Precise® Digital 1.50	GR													-8.00 to +5.00	+1.00 to +3.00	-6.00	17
KODAK Precise® Plus 1.50	GR	BR	GG	S	AME	AMB	E	R	GR					-8.00 to +4.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus 1.50	GR	BR	GG	S	AME	AMB	E	R	GR					-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Precise® PB 1.50	GR	BR							GR					-8.00 to +4.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB 1.50	GR	BR							GR					-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique™ 1.50	GR	BR	GG	S	AME	AMB	E	R	GR					-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique™ HD 1.50	GR	BR												-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ 1.50	GR	BR	GG	S	AME	AMB	E	R	GR					-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.50	GR	BR	GG	S	AME	AMB	E	R	GR					-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ 1.50	GR	BR	GG	S	AME	AMB	E	R	GR	BR		GR		-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ HD 1.50	GR	BR	GG	S	AME	AMB	E	R	GR	BR		GR		-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Easy 14 1.50	GR	BR	GG	S	AME	AMB	E	R						-8.00 to +4.00	+0.75 to +3.50	-6.00	14
KODAK Easy 18 1.50	GR	BR	GG	S	AME	AMB	E	R						-8.00 to +4.00	+0.75 to +3.50	-6.00	18
DirecTek 1.50	GR	BR												-8.00 to +4.00	+0.75 to +3.50	-6.00	16.5
DirecTek Short 1.50	GR	BR												-8.00 to +4.00	+0.75 to +3.50	-6.00	13.5
Navigator FBS 1.50	GR	BR												-8.00 to +4.00	+0.75 to +3.50	-6.00	18
Navigator Short FBS 1.50	GR	BR												-8.00 to +4.00	+0.75 to +3.50	-6.00	14
Crossbows Custom U 1.50	GR	BR												-8.00 to +4.00	+0.75 to +3.50	-6.00	13
Crossbows Custom V 1.50	GR	BR												-8.00 to +4.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 13 Progressive 1.50	GR	BR												-8.00 to +4.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 17 Progressive 1.50	GR	BR												-8.00 to +4.00	+0.75 to +3.50	-6.00	17
Crossbows Junior	GR	BR												-8.00 to +4.00	+0.75 to +3.50	-6.00	13
1.53 TRIVEX®																	
KODAK Precise® PB Trivex®	GR	BR							GR					-10.00 to +4.50	+0.75 to +3.50	-6.00	17

KODAK Easy 18 Poly	GR	BR	GG	S	AME	AMB	E	R								-11.00 to +6.00	+0.75 to +3.50	-6.00	18
DirecTek 1.67	GR	BR														-11.00 to +6.00	+0.75 to +3.50	-6.00	16.5
DirecTek Short 1.67	GR	BR														-11.00 to +6.00	+0.75 to +3.50	-6.00	13.5
Navigator FBS 1.67	GR	BR														-11.00 to +6.00	+0.75 to +3.50	-6.00	18
Navigator Short FBS 1.67	GR	BR														-11.00 to +6.00	+0.75 to +3.50	-6.00	14
Crossbows Custom U 1.67	GR	BR														-11.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Custom V 1.67	GR	BR														-11.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 13 Progressive 1.67	GR	BR														-11.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 17 Progressive 1.67	GR	BR														-11.00 to +6.00	+0.75 to +3.50	-6.00	17
Crossbows Junior 1.67	GR	BR														-11.00 to +6.00	+0.75 to +3.50	-6.00	13
1.74 Ultra High Index																			
KODAK Unique	GR	BR														-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD	GR	BR														-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™	GR	BR														-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.74	GR	BR														-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™	GR	BR														-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ HD 1.74	GR	BR														-13.00 to +9.00	+0.75 to +3.50	-6.00	13

														
Semi-Finished Single Vision													Sphere	Cyl to
1.50 Standard Index SFSV														
TRN 1.50 SV (72mm) UC/HC DIA.72											GR	BR		Base 0.50, 2.25, 4.25, 5.25, 6.25, 8.25
TRN 1.50 DS (76mm) UC/HC DIA.72	GR	BR									GR	BR	GR	Base 0.50, 8.75, 10.50
TRN 1.50 DS (76mm) UC/HC DIA.76	GR	BR									GR	BR	GR	Base 2.00, 3.50, 4.25, 5.25, 6.25, 7.25
1.53 TRIVEX® Single Vision														
TOG 1.53 Trivex (75mm) SV UC/HC DIA.70	GR	BR									GR	BR		Base 1.00
TOG 1.53 Trivex (75mm) SV UC/HC DIA.75	GR	BR									GR	BR		Base 2.50, 4.00, 5.00, 6.00, 7.00, 8.00
TRN 1.53 Trivex (78mm) SV UC/HC DIA.74													GR	Base 0.50, 1.00
TRN 1.53 Trivex (78mm) SV UC/HC DIA.78													GR	Base 2.00, 3.50, 4.25, 5.25, 6.25, 7.25, 8.75
1.59 Polycarbonate Single Vision														
TRN 1.59 Polycarbonate (79mm) SV HC DIA.79	GR	BR									GR	BR	GR	Base 1.50, 3.25, 4.25, 4.75, 5.50, 7.00
TRN 1.59 Polycarbonate (74mm) SV HC DIA.74	GR	BR									GR	BR		Base 7.50, 8.50
1.60 High Index Single Vision														
TOG 1.60 SV (75mm) UC/HC DIA.65	GR	BR									GR	BR		Base 10.00, 11.00
TOG 1.60 SV (75mm) UC/HC DIA.70	GR	BR									GR	BR		Base 9.00
TOG 1.60 SV (75mm) UC/HC DIA.75	GR	BR									GR	BR		Base 1.00, 2.50, 4.00, 5.00, 6.00, 7.00, 8.00
TOG 1.60 SV (80mm) UC/HC DIA.80											GR	BR		Base 1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.50
TRN 1.60 SV (80mm) UC/HC DIA.76													GR	Base 8.50
TRN 1.60 SV (80mm) UC/HC DIA.80													GR	Base 1.75, 2.75, 3.75, 4.75, 6.00, 6.75
1.67 Super High Index Single Vision														
TOG 1.67 SV (65mm) UC/HC DIA.65	GR	BR									GR	BR		Base 10.00, 11.00
TOG 1.67 SV (70mm) UC/HC DIA.70	GR	BR									GR	BR		Base 9.00
TOG 1.67 SV (75mm) UC/HC DIA.75	GR	BR									GR	BR		Base 1.00, 2.50, 4.00, 5.00, 6.00, 7.00, 8.00
TOG 1.67 SV (80mm) UC/HC DIA.80											GR	BR		Base 1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.50
TRN 1.67 SV (80mm) UC/HC DIA.80													GR	Base 1.75, 2.75, 3.75, 4.75, 6.00

TRN 1.67 SV (80mm) UC/HC DIA.81														GR	Base 6.75, 8.00	
1.74 ultra High Index Single Vision																
TOG 1.74 LeafECO SV (75mm) UC/HC DIA.65	GR	BR													Base 11.50, 12.50, 13.50, 14.50, 15.50, 16.50	
TOG 1.74 LeafECO SV (75mm) UC/HC DIA.70	GR	BR													Base 0.50, 1.50, 2.25, 6.75, 8.50, 10.25	
TOG 1.74 LeafECO SV (75mm) UC/HC DIA.75	GR	BR													Base 3.50, 5.00	

																
Semi-Finished Single Vision															Sphere	Cyl to
1.50 Standard Index																
Unity SVx	GR	BR	GG	S	AME	AMB	E	R	GR				GR		-12.25 to +6.00	-6.00
Unity SVxtra	GR	BR	GG	S	AME	AMB	E	R	GR				GR		-12.25 to +6.00	-6.00
Unity SVxtreme	GR	BR	GG	S	AME	AMB	E	R	GR				GR		-6.00 to +6.00	-6.00
Unity Relieve 50	GR	BR	GG	S	AME	AMB	E	R							-12.25 to +6.00	-6.00
Unity Relieve 70	GR	BR	GG	S	AME	AMB	E	R							-12.25 to +6.00	-6.00
1.53 TRIVEX® Single Vision																
Unity SVx	GR	BR	GG	S	AME	AMB	E	R	GR						-14.00 to +6.00	up to -6.00
Unity SVxtra	GR	BR	GG	S	AME	AMB	E	R	GR						-14.00 to +6.00	up to -6.00
Unity SVxtreme	GR	BR	GG	S	AME	AMB	E	R	GR						-6.50 to +6.00	up to -6.00
Unity Relieve 50	GR	BR	GG	S	AME	AMB	E	R							-14.00 to +6.00	up to -6.00
Unity Relieve 70	GR	BR	GG	S	AME	AMB	E	R							-14.00 to +6.00	up to -6.00
1.59 Polycarbonate Single Vision																
Unity SVx	GR	BR	GG	S	AME	AMB	E	R	GR				GR		-14.00 to +6.00	up to -6.00
Unity SVxtra	GR	BR	GG	S	AME	AMB	E	R	GR				GR		-14.00 to +6.00	up to -6.00
Unity SVxtreme	GR	BR	GG	S	AME	AMB	E	R	GR				GR		-6.50 to +6.00	up to -6.00
Unity Relieve 50	GR	BR	GG	S	AME	AMB	E	R							-14.00 to +6.00	up to -6.00
Unity Relieve 70	GR	BR	GG	S	AME	AMB	E	R							-14.00 to +6.00	up to -6.00
1.60 High Index Single Vision																
Unity SVx	GR	BR													-14.00 to +6.00	up to -6.00
Unity SVxtra	GR	BR													-14.00 to +6.00	up to -6.00
Unity SVxtreme	GR	BR													-6.50 to +6.00	up to -6.00
1.67 Super High Index Single Vision																
Unity SVx	GR	BR	GG	S	AME	AMB	E	R	GR				GR		-16.00 to +8.00	up to -6.00
Unity SVxtra	GR	BR	GG	S	AME	AMB	E	R	GR				GR		-16.00 to +8.00	up to -6.00
Unity SVxtreme	GR	BR	GG	S	AME	AMB	E	R	GR				GR		-6.00 to +6.00	up to -6.00
Unity Relieve 50	GR	BR	GG	S	AME	AMB	E	R							-16.00 to +8.00	up to -6.00
Unity Relieve 70	GR	BR	GG	S	AME	AMB	E	R							-16.00 to +8.00	up to -6.00

																	
Progressives													Sphere	Add	Cyl to	MFH	
1.50 Standard Index																	
Unity Via II	GR	BR	GG	S	AME	AMB	E	R	GR				GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Plus II	GR	BR	GG	S	AME	AMB	E	R	GR				GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite VR	GR	BR	GG	S	AME	AMB	E	R	GR				GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II	GR	BR	GG	S	AME	AMB	E	R	GR				GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II 12	GR	BR	GG	S	AME	AMB	E	R	GR				GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II 14	GR	BR	GG	S	AME	AMB	E	R	GR				GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	14 mm
Unity Via Elite II 16	GR	BR	GG	S	AME	AMB	E	R	GR				GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	16 mm
Unity Via Elite II 18	GR	BR	GG	S	AME	AMB	E	R	GR				GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	18 mm
Unity Via Mobile II	GR	BR	GG	S	AME	AMB	E	R	GR				GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Wrap II	GR	BR	GG	S	AME	AMB	E	R	GR				GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Ethos	GR	BR	GG	S	AME	AMB	E	R	GR				GR	-9.00 to +4.50	+0.75 to +3.50	up to -4.00	14mm

Unity Lenses — Near Variable Focus

Camber™	GR	BR								GR	BR		GR	-14.00 to +8.00	-4.00
1.74 Ultra High Index															
Camber™	GR	BR								GR	BR			-19.00 to +9.00	-4.00

YOUNGER OPTICS The Optical Lens Innovators 	Transitions Gen#S								Transitions XTRACTIVE	Transitions XTRACTIVE POLARIZED			
Bifocal Signature™											Sphere	Add	Cyl to
1.50 Standard Index SFST-28													
Younger Optics SFFT 28									GR		-11.00 to +6.00	+1.00 to +3.00	-4.00

YOUNGER OPTICS The Optical Lens Innovators	Transitions Gen#S								Transitions XTRACTIVE	Transitions XTRACTIVE POLARIZED				
Progressives											Sphere	Add	Cyl to	MFH
1.53 TRIVEX®														
Younger Trilogy® Image®	GR										-8.00 to +8.00	+1.00 to +3.00	-4.00	18
1.59 index polycarbonate														
Younger Image®	GR	BR							GR	BR	-9.00 to +7.00	+1.00 to +3.00	-4.00	18
1.67 Super High Index														
Younger Image® 1.67 MR-10™	GR	BR									-10.00 to +8.50	+1.00 to +3.50	-4.00	18

TRANSITIONS® LIGHT INTELLIGENT LENSES™



















GR
Gray

BR
Brown

GG
Graphite
Green

S
Sapphire

AME
Amethyst

AMB
Amber

E
Emerald

R
Ruby

Available
in style
mirrors

Style Mirrors are available where gray and brown *Transitions® XTRActive®* are available. Specify *Transitions* lenses in style mirrors **(no substitutions)** with your lab to ensure authenticity.

Transitions, *Transitions Signature*, and *XTRActive* are registered trademarks, and *XTRActive Polarized*, the *Transitions* logo and *Transitions Light Intelligent Lenses* are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

Unless indicated otherwise, all registered trademarks and trademarks are the property of Essilor International and/or its subsidiaries in the United States and in other countries.

All other trademarks are the property of their respective owners.