

Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Design](#) / [Bifocals](#) / [Younger Optics](#)

Younger Optics — Transitions® Bifocals & Trifocals

						
Bifocals with Transitions® Technology				Sphere	Add	Cyl to
1.50 Standard Index SFST-28						
Younger Optics SFFT 28	GR	BR		-11.00 to +6.00	+1.00 to +3.00	-4.00
1.59 index polycarbonate composite						
Younger Optics Polycarbonate Composite SFFT 28	GR			-9.00 to +7.00	+1.00 to +3.00	-4.00

						
Trifocals with Transitions® Technology				Sphere	Add	Cyl to
1.50 Standard Index Trifocals						
Younger Optics 7x28	GR			-6.00 to +5.00	+1.50 to +4.00	-4.00

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2026 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.