





Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Transitions® GEN S™](#) / [Shamir Insights Inc.](#)

Shamir Insights Inc. – Transitions® GEN S™

											
Finished Single Vision										Sphere	Cyl to
1.59 FSV											
Shamir FSV with Glacier PLUS™ AR	GR									-6.00 to +4.00	-2.00

											
Semi-Finished Single Vision										Sphere	Cyl to
1.50 SFSV											
Shamir Autograph® III Single Vision	GR BR									-12.00 to +7.25	-6.00
Shamir Relax™	GR BR									-14.00 to +8.00	-6.00
1.53 SFSV											
Shamir Autograph® III Single Vision	GR BR									-14.00 to +8.00	-6.50
Shamir Relax™	GR BR									-14.00 to +8.00	-6.50
1.59 SFSV											
Shamir Autograph® III Single Vision	GR BR									-14.00 to +8.75	-7.00
Shamir Relax™	GR BR									-14.75 to +6.00	-7.00
1.60 SFSV											
Shamir Autograph® III Single Vision	GR BR									-14.00 to +8.75	-7.00
Shamir Relax™	GR BR									-14.00 to +8.75	-7.00
1.67 SFSV											
Shamir Autograph® III Single Vision	GR BR									-16.00 to +16.00	-8.00
Shamir Relax™	GR BR									-16.75 to +16.00	-8.00
1.74 SFSV											
Shamir Autograph® III Single Vision	GR GR									-20.00 to +13.50	-9.00

													
Blended Bifocal (digital)										Sphere	Add	Cyl to	MFH
1.50 Standard Index Blended Bifocal (digital)													
Shamir Duo	GR BR									-9.00 to +6.00	+0.75 to +4.00	-6.00	15
1.53 Trivex® Blended Bifocal (digital)													
Shamir Duo	GR BR									-13.25 to +6.00	+0.75 to +4.00	-6.50	15
1.59 Polycarbonate Blended Bifocal (digital)													
Shamir Duo	GR BR									-14.75 to +6.00	+0.75 to +4.00	-7.00	15
1.60 High Index Blended Bifocal (digital)													
Shamir Duo	GR BR									-15.00 to +6.00	+0.75 to +4.00	-7.00	15
1.67 Super High Index Blended Bifocal (digital)													
Shamir Duo	GR BR									-16.75 to +8.00	+0.75 to +4.00	-8.00	15
1.74 Ultra High Index Blended Bifocal (digital)													
Shamir Duo	GR BR									-18.75 to +12.00	+0.75 to +4.00	-9.00	15

													
Progressives										Sphere	Add	Cyl to	MFH

1.50 Standard Index													
Shamir Autograph Intelligence™	GR	BR								-12.25 to +6.50	+0.75 to +4.00	-6.00	11, 13,15,18
Shamir Autograph® III	GR	BR								-12.25 to +6.00	+0.75 to +4.00	-6.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR								-12.25 to +6.00	+0.75 to +4.00	-6.00	11, 13, 15, 18
Shamir Intouch™	GR	BR								-12.25 to +6.00	+0.75 to +3.50	-6.00	14, 16, 18
Shamir Spectrum+™	GR	BR								-12.25 to +6.00	+0.75 to +3.50	-6.00	14, 16, 18
Shamir Element™	GR	BR								-7.00 to +4.00	+0.75 to +3.50	-6.00	16, 19
Shamir GENESIS HD	GR	BR								-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir Computer™	GR	BR								-7.00 to +4.00	+0.75 to +2.00	-6.00	14 Lower, 10 Upper
Shamir Workspace™	GR	BR								-6.00 to +6.00	+0.75 to +2.50	-6.00	14 Lower, 10 Upper
Shamir Attitude III Sport	GR	BR								-6.50 to +5.25	+0.75 to +3.50	-6.00	18
Shamir Attitude III Fashion	GR	BR								-6.50 to +5.25	+0.75 to +3.50	-6.00	15, 18
1.53 TRIVEX®													
Shamir Autograph Intelligence™	GR	BR								-13.25 to +7.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
Shamir Autograph® III	GR	BR								-13.25 to +7.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
Shamir Autograph® II+	GR	BR								-13.25 to +6.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
Shamir Intouch™	GR	BR								-13.25 to +6.00	+0.75 to +3.50	-6.50	14, 16, 18
Shamir Golf™	GR	BR								-5.00 to +4.00	+1.00 to +3.50	-6.50	19
Shamir Spectrum+™	GR	BR								-13.25 to +6.00	+0.75 to +3.50	-6.50	14, 16, 18
Shamir Element™	GR	BR								-13.25 to +6.00	+0.75 to +3.50	-6.50	16, 19
Shamir GENESIS HD	GR	BR								-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir Computer™	GR	BR								-6.00 to +8.00	+0.75 to +2.00	-6.50	14 Lower, 10 Upper
Shamir Workspace™	GR	BR								-6.00 to +8.00	+0.75 to +2.50	-6.50	14 Lower, 10 Upper
Shamir Attitude III Sport	GR	BR								-6.50 to +6.00	+0.75 to +3.50	-6.50	18
Shamir Attitude III Fashion	GR	BR								-6.50 to +6.00	+0.75 to +3.50	-6.50	15, 18
1.59 Polycarbonate													
Shamir Autograph Intelligence™	GR	BR								-14.75 to +7.50	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR								-14.75 to +8.00	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR								-14.75to +6.50	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Intouch™	GR	BR								-14.75 to +6.00	+0.75 to +3.50	-7.00	15, 18
Shamir Golf™	GR	BR								-6.00 to +5.25	+1.00 to +3.50	-7.00	19
Shamir Spectrum+™	GR	BR								-14.75 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
Shamir Element™	GR	BR								-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19
Shamir GENESIS HD	GR	BR								-10.00 to +6.00	+0.75 to +3.50	-7.00	16
Shamir Computer™	GR	BR								-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper
Shamir Workspace™	GR	BR								-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper
Shamir Attitude III Sport	GR	BR								-6.50 to +6.00	+0.75 to +3.50	-6.50	18
Shamir Attitude III Fashion	GR	BR								-6.50 to +6.00	+0.75 to +3.50	-6.50	15, 18
1.60 High Index													
Shamir Autograph Intelligence™	GR	BR								-14.75 to +8.75	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR								-15.00 to +8.00	+0.75 to +3.50	-7.00	15,18
Shamir Autograph® II+	GR	BR								-15.00 to +6.00	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Golf™ SuperLite™ 1.6	GR	BR								-6.00 to +5.25	+1.00 to +3.50	-7.00	19
Shamir Spectrum+™	GR	BR								-15.00 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
Shamir Element™	GR	BR								-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19
Shamir GENESIS HD	GR	BR								-6.00 to +6.00	+0.75 to +3.50	-7.00	16
Shamir Computer™	GR	BR								-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper
Shamir Workspace™	GR	BR								-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper
Shamir InTouch™	GR	BR								-15.00 to +6.00	+0.75 to +3.50	-7.00	15, 18
Shamir Attitude III Sport	GR	BR								-6.50 to +6.00	+0.75 to +3.50	-7.00	18
Shamir Attitude III Fashion	GR	BR								-6.50 to +6.00	+0.75 to +3.50	-7.00	15, 18
1.67 Super High Index													
Shamir Autograph Intelligence™	GR	BR								-16.75 to +10.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR								-16.75 to +9.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR								-16.75 to +8.25	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Spectrum+™	GR	BR								-16.75 to +8.00	+0.75 to +4.00	-8.00	11, 13, 15, 18

Shamir Element™ SuperLite™ 1.67	GR	BR						-12.50 to +8.00	+0.75 to +3.50	-8.00	16, 19
Shamir GENESIS HD	GR	BR						-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir InTouch™	GR	BR						-16.75 to +8.00	+0.75 to +3.50	-8.00	15, 18
Shamir Attitude III Sport	GR	BR						-6.50 to +6.00	+0.75 to +3.50	-8.00	18
Shamir Attitude III Fashion	GR	BR						-6.50 to +6.00	+0.75 to +3.50	-8.00	15, 18
1.74 Ultra High Index											
Shamir Autograph Intelligence™	GR	BR						-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR						-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR						-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir GENESIS HD	GR	BR						-6.00 to +6.00	+0.75 to +4.00	-4.00	16
Shamir Spectrum+™	GR	BR						-18.75 to +11.50	+0.75 to +3.50	-9.00	14, 16, 18
Shamir InTouch™	GR	BR						-18.75 to +15.50	+0.75 to +3.50	-9.00	15, 18

TRANSITIONS® LIGHT INTELLIGENT LENSES™



GR
Gray



BR
Brown



GG
Graphite
Green



S
Sapphire



AME
Amethyst



AMB
Amber



E
Emerald



R
Ruby



Available
in style
mirrors

Transitions, *Transitions Signature*, and *XTRActive* are registered trademarks, and *XTRActive Polarized*, the *Transitions* logo and *Transitions Light Intelligent Lenses* are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.