





Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Transitions® Signature® GEN S™](#)

All Suppliers – Transitions® GEN S™










										
Finished Single Vision									Sphere	Cyl to
1.50 Standard Index										
Essilor FSV Spherical	GR	BR							-4.00 to +4.00	-2.00
Essilor FSV Aspheric	GR								-4.00 to +4.00	-2.00
Essilor FSV 1.50 Crizal® Sapphire™ HR	GR								-4.00 to +4.00	-2.00
Essilor FSV 1.50 Crizal® Rock™	GR								-4.00 to +4.00	-2.00
Essilor FSV 1.50 Crizal® Easy Pro	GR								-4.00 to +4.00	-2.00
1.59 Polycarbonate Finished Single Vision										
Essilor FSV Spherical Airwear®	GR	BR							-6.00 to +4.00	-2.00
Essilor FSV Aspheric Airwear®	GR	BR							-6.00 to +4.00	-2.00
Essilor FSV Airwear® Crizal® Sapphire™ HR	GR								-6.00 to +4.00	-2.00
Essilor FSV Airwear® Crizal® Rock™	GR								-6.00 to +4.00	-2.00
Essilor FSV Airwear® Crizal® Easy Pro	GR								-6.00 to +4.00	-2.00
1.67 Super High Index Finished Single Vision										
Essilor FSV Spherical Thin&Lite® 1.67									-6.00 to +4.00	-2.00

										
Semi-Finished Single Vision									Sphere	Cyl to
1.50 Standard Index										
Essilor Spherical Single Vision	GR	BR	GG	S	AME	AMB	E	R	-10.50 to +10.00	-4.00
Essilor Stylistic® Single Vision	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +7.50	-6.00
Eyezen® +0	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +7.50	-6.00
Eyezen® +1	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +7.50	-6.00
Eyezen® +2	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +7.50	-6.00
Eyezen® +3	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +7.50	-6.00
Eyezen® +4	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +7.50	-6.00
Eyezen® Start	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +6.00	-6.00
1.53 TRIVEX® Single Vision										
Essilor Spherical Single Vision TREXA™	GR	BR	GG						-10.00 to +6.00	-4.00
Essilor Eyecode™ Single Vision TREXA™	GR	BR	GG						-10.00 to +7.50	-6.00
Essilor FIT™ Single Vision TREXA™	GR	BR	GG						-10.00 to +7.50	-6.00
Essilor Azio™ Single Vision TREXA™	GR	BR	GG						-10.00 to +7.50	-6.00
Essilor 360°™ Single Vision TREXA™	GR	BR	GG						-10.00 to +7.50	-6.00
Eyezen® +0 TREXA™	GR	BR	GG						-10.00 to +7.50	-6.00
Eyezen® +1 TREXA™	GR	BR	GG						-10.00 to +7.50	-6.00
Eyezen® +2 TREXA™	GR	BR	GG						-10.00 to +7.50	-6.00
Eyezen® +3 TREXA™	GR	BR	GG						-10.00 to +7.50	-6.00
Eyezen® +4 TREXA™	GR	BR	GG						-10.00 to +7.50	-6.00
Eyezen® Start	GR	BR	GG						-10.00 to +6.00	-6.00
1.59 Polycarbonate Single Vision										
Essilor Spherical Single Vision Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	-4.00
Essilor Aspheric Single Vision Airwear®	GR	BR	GG						-10.00 to +8.00	-4.00
Essilor Eyecode™ Aspheric SV Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	-6.00
Essilor FIT™ Aspheric SV Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	-6.00










Essilor AZIO™ Aspheric SV Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	-6.00
Essilor SV 360™ Aspheric Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	-6.00
Essilor Stylistic® Single Vision Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	-6.00
Eyezen® +0 Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	-6.00
Eyezen® +1 Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	-6.00
Eyezen® +2 Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	-6.00
Eyezen® +3 Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	-6.00
Eyezen® +4 Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	-6.00
Eyezen® Start Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	-6.00
Eyezen®Kids	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +6.00	-6.00
1.60 High Index Single Vision										
Essilor Spherical Single Vision Thin&Lite® 1.60	GR	BR	GG	S	AME	AMB	E	R	-13.00 to +8.00	-4.00
Eyezen® +0 Thin&Lite® 1.60	GR	BR							-12.00 to +8.00	-6.00
Eyezen® +1 Thin&Lite® 1.60	GR	BR							-12.00 to +8.00	-6.00
Eyezen® +2 Thin&Lite® 1.60	GR	BR							-12.00 to +8.00	-6.00
Eyezen® +3 Thin&Lite® 1.60	GR	BR							-12.00 to +8.00	-6.00
Eyezen® +4 Thin&Lite® 1.60	GR	BR							-12.00 to +8.00	-6.00
Eyezen® Start Thin&Lite® 1.60	GR	BR							-12.00 to +8.00	-6.00
1.67 Super High Index Single Vision										
Essilor Spherical Single Vision Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	-4.00
Essilor Aspheric Single Vision Thin&Lite® 1.67	GR	BR	GG						-16.00 to +8.00	-4.00
Essilor Eyecode™ SV Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	-6.00
Essilor FIT™ SV Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	-6.00
Essilor Azio™ SV Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	-6.00
Essilor SV 360™ Spherical Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	-6.00
Essilor Stylistic® Single Vision Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	-6.00
Eyezen® +0 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	-6.00
Eyezen® +1 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	-6.00
Eyezen® +2 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	-6.00
Eyezen® +3 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	-6.00
Eyezen® +4 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	-6.00
Eyezen® Start Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	-6.00
1.74 Ultra High Index Single Vision										
Essilor Spherical Single Vision Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	-4.00
Essilor Eyecode™ SV Thin&Lite® 1.74	GR	BR							-16.00 to +13.00	-6.00
Essilor FIT™ SV Thin&Lite® 1.74	GR	BR							-16.00 to +13.00	-6.00
Essilor 360™ SV Thin&Lite® 1.74	GR	BR							-16.00 to +13.00	-6.00
Eyezen® +0 Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	-6.00
Eyezen® +1 Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	-6.00
Eyezen® +2 Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	-6.00
Eyezen® +3 Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	-6.00
Eyezen® +4 Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	-6.00
Eyezen® Start Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	-6.00











 Progressives									Sphere	Add	Cyl to	MFH
1.50 Standard Index												
Varilux® XR Track™	GR	BR	GG						-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR	GG						-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™	GR	BR	GG						-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™	GR	BR	GG						-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™	GR	BR	GG						-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode	GR	BR	GG						-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT	GR	BR	GG						-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14











Varilux® Physio® W3+	GR	BR	GG							-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit	GR	BR	GG							-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max	GR	BR	GG							-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +7.50	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Varilux® Liberty™ 3.0	GR	BR								-10.00 to +7.50	+0.75 to +1.50	-6.00	14
Varilux Comfort®	GR	BR								-10.00 to +7.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short™	GR	BR								-10.00 to +7.00	+0.75 to +3.00	-6.00	14
Varilux® Stylistic™	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ FIT	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +7.50	+0.75 to +3.50	-4.00	17
Essilor Ideal Advanced™ Wrap	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +7.50	+0.75 to +3.50	-4.00	17
Essilor Ideal® FIT	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +7.50	+0.75 to +3.50	-6.00	18
Essilor Ideal® Driving	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +7.50	+0.75 to +3.50	-6.00	13
Essilor Ideal®	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal® Short	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Natural® Digital	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +7.50	+0.75 to +3.50	-6.00	15
Essilor Accolade™	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Digital	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Brillenz®	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Brillenz® Short	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Ovation® Digital	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Ovation®	GR	BR								-10.00 to +6.00	+0.75 to +3.50	-4.00	17
Essilor SmallFit® Digital	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +7.50	+0.75 to +3.00	-6.00	14

												
Progressives									Sphere	Add	Cyl to	MFH
1.53 TREXA®												
	Varilux® XR Track™	GR	BR						-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement
	Varilux® XR Design™	GR	BR						-10.00 to +6.00	+0.75 to +4.00	-6.00	14
	Varilux® X 4D™ TREXA™	GR	BR						-10.00 to +6.00	+0.75 to +4.00	-6.00	14
	Varilux® X FIT™ TREXA™	GR	BR						-10.00 to +6.00	+0.75 to +4.00	-6.00	14
	Varilux® X Design™ TREXA™	GR	BR						-10.00 to +6.00	+0.75 to +4.00	-6.00	14
	Varilux® Physio® W3+ Eyecode TREXA™	GR	BR						-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
	Varilux® Physio® W3+ FIT TREXA™	GR	BR						-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
	Varilux® Physio® W3+ TREXA™	GR	BR						-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
	Varilux Comfort® Max Fit TREXA™	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-6.00	14
	Varilux Comfort® Max TREXA™	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-6.00	14
	Varilux Physio DRx™ TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50 & +4.00	-6.00	17
	Varilux Physio Short DRx™ TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.00	-6.00	14
	Varilux Comfort DRx™ TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50	-6.00	17
	Varilux Comfort Short DRx™ TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.00	-6.00	14
	Essilor Ideal Advanced™ FIT TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50	-6.00	17
	Essilor Ideal Advanced™ TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50	-6.00	17
	Essilor Ideal® FIT TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50	-6.00	17
	Essilor Ideal® Sport TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50	-6.00	18
	Essilor Ideal® Computer TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50	-6.00	13










Essilor Ideal® Driving TREXA™	GR	BR	GG							-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal® TREXA™	GR	BR	GG							-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal® Short TREXA™	GR	BR	GG							-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Natural® Digital TREXA™	GR	BR	GG							-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ TREXA™	GR	BR	GG							-10.00 to +7.50	+0.75 to +3.50	-6.00	15
Essilor Accolade® TREXA™	GR	BR	GG							-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Digital TREXA™	GR	BR	GG							-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital TREXA™	GR	BR	GG							-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Brillenz® TREXA™	GR	BR	GG							-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Brillenz® Short TREXA™	GR	BR	GG							-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Ovation® Digital TREXA™	GR	BR	GG							-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor SmallFit® Digital TREXA™	GR	BR	GG							-10.00 to +7.50	+0.75 to +3.00	-6.00	14










												
Progressives									Sphere	Add	Cyl to	MFH
1.59 Polycarbonate												
Varilux® XR Track™	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Varilux Comfort DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux® Liberty™ 3.0	GR	BR							-10.00 to +8.00	+0.75 to +1.50	-6.00	14
Varilux Comfort® Airwear®	GR	BR	GG						-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short™ Airwear®	GR	BR	GG						-10.00 to +6.00	+0.75 to +3.00	-6.00	14
Varilux® Stylistic™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ FIT	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Wrap	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Essilor Ideal® FIT Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Sport Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Driving Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Ideal® Short Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Natural® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Ovation® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Ovation® Airwear®	GR	BR							-10.00 to +8.00	+0.75 to +3.50	-4.00	17
Essilor SmallFit® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.75 to +3.00	-6.00	14

 Progressives									Sphere	Add	Cyl to	MFH	
													
1.60 High Index													
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.60	GR	BR	GG							-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Thin&Lite® 1.60	GR	BR	GG							-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Thin&Lite® 1.60	GR	BR	GG							-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit Thin&Lite® 1.60	GR	BR	GG							-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max Thin&Lite® 1.60	GR	BR	GG							-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Thin&Lite® 1.60	GR	BR								-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Thin&Lite® 1.60	GR	BR								-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™ Thin&Lite® 1.60	GR	BR								-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Thin&Lite® 1.60	GR	BR								-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux Comfort® Thin&Lite® 1.60	GR									-12.00 to +8.00	+0.75 to +3.50	-4.00	17
Varilux Comfort Short™ Thin&Lite® 1.60	GR									-12.00 to +8.00	+0.75 to +3.00	-4.00	14
Essilor Ideal Advanced™ FIT Thin&Lite® 1.60	GR	BR								-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Thin&Lite® 1.60	GR	BR								-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.60	GR	BR								-12.00 to +8.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Sport Thin&Lite® 1.60	GR	BR								-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Thin&Lite® 1.60	GR	BR								-12.00 to +8.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Driving Thin&Lite® 1.60	GR	BR								-12.00 to +8.00	+0.75 to +3.50	-6.00	13
Essilor Ideal® Thin&Lite® 1.60	GR	BR								-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Short Thin&Lite® 1.60	GR	BR								-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Essilor Natural® Digital Thin&Lite® 1.60	GR	BR								-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Thin&Lite® 1.60	GR	BR								-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Thin&Lite® 1.60	GR	BR								-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Ovation® Digital Thin&Lite® 1.60	GR	BR								-12.00 to +8.00	+0.75 to +3.50	-6.00	17
SmallFit® Digital Thin&Lite® 1.60	GR	BR								-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Thin&Lite® 1.60	GR	BR								-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Thin&Lite® 1.60	GR	BR								-12.00 to +8.00	+0.75 to +3.00	-6.00	14

 Progressives									Sphere	Add	Cyl to	MFH	
													
1.67 Super High Index													
Varilux® XR Track™	GR	BR	GG	S	AME	AMB	E	R		-12.00 to +9.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR	GG	S	AME	AMB	E	R		-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Varilux® Liberty™ 3.0	GR	BR								-14.00 to +9.00	+0.75 to +1.50	-6.00	14
Varilux Comfort® Thin&Lite® 1.67	GR	BR								-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Varilux Comfort® Short™ Thin&Lite® 1.67	GR	BR								-14.00 to +9.00	+0.75 to +3.50	-6.00	17

Varilux® Stylistic™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-5.00 to +5.00	+0.75 to +3.50	-3.00	17
Essilor Ideal Advanced™ FIT Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal Advanced Wrap™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Sport Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Driving Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor Ideal® Short	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor Natural® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Brillenz® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Ovation® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Ovation® Thin&Lite® 1.67	GR	BR							-12.00 to +8.00	+0.75 to +3.50	-4.00	17
Essilor SmallFit® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	+0.75 to +3.00	-6.00	14



												
Progressives									Sphere	Add	Cyl to	MFH
1.74 Ultra High Index												
Varilux® XR Track™	GR	BR							-14.00 to +13.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR							-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ Thin&Lite® 1.74	GR	BR							-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ Thin&Lite® 1.74	GR	BR							-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ Thin&Lite® 1.74	GR	BR							-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.74	GR	BR							-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Thin&Lite® 1.74	GR	BR							-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Thin&Lite® 1.74	GR	BR							-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio DRx™ Thin&Lite® 1.74	GR	BR							-14.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Thin&Lite® 1.74	GR	BR							-14.00 to +13.00	+0.75 to +3.50	-6.00	14
Varilux Comfort DRx™ Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ FIT Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Computer Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Driving Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	+0.75 to +3.50	-6.00	13
Ovation® Digital Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	+0.75 to +3.50	-6.00	14
Essilor Accolade Freedom™ Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Thin&Lite® 1.74	GR	BR							-14.00 to +9.00	+0.75 to +3.50	-6.00	17


										
Finished Single Vision									Sphere	Cyl to
1.50 Standard Index										
HOYA FSV	GR								-4.00 to +4.00	-2.00

Amplitude Mini BKS 1.50	GR	BR								-10.00 to +3.75	+0.75 to +4.00	-6.00	14
Summit cd BKS 1.50	GR	BR								-10.00 to +3.75	+0.75 to +4.00	-6.00	14
Summit ecp BKS 1.50	GR	BR								-10.00 to +3.75	+0.75 to +4.00	-6.00	18
Hoya GP Wide 1.50	GR	BR								-10.00 to +3.75	+0.75 to +4.00	-6.00	18
Tact 40 BKS 1.50	GR	BR								-10.00 to +3.75	+0.75 to +4.00	-6.00	18
Tact 60 BKS 1.50	GR	BR								-10.00 to +3.75	+0.75 to +4.00	-6.00	18
1.53 TRIVEX®													
HOYALUX iD LifeStyle®4 Indoor Phoenix™	GR!	BR!								-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Outdoor Phoenix™	GR!	BR!								-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Urban Phoenix™	GR!	BR!								-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Indoor Phoenix™	GR!	BR!								-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Outdoor Phoenix™	GR!	BR!								-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Urban Phoenix™	GR!	BR!								-8.00 to +6.00	+0.75 to +4.00	-6.00	14
Amp HD3 Phoenix™	GR!	BR!								-8.00 to +6.00	+0.75 to +4.00	-6.00	18
Amp Mini HD3 Phoenix™	GR!	BR!								-8.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX ARRAY™ 11, 13, 15, 17, VL Phoenix™	GR	BR								-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL Phoenix™	GR	BR								-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13, 15, 17, VL Phoenix™	GR	BR								-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11, 13, 15, 17, VL Phoenix™	GR	BR								-5.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS Phoenix™	GR	BR								-10.00 to +6.00	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS Phoenix™	GR	BR								-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Summit cd BKS Phoenix™	GR	BR								-10.00 to +6.00	+0.75 to +4.00	-6.00	18
Summit ecp BKS Phoenix™	GR	BR								-10.00 to +6.00	+0.75 to +4.00	-6.00	18
Hoya GP Wide Phoenix™	GR	BR								-10.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 40 BKS Phoenix™	GR	BR								-10.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 60 BKS Phoenix™	GR	BR								-10.00 to +6.00	+0.75 to +4.00	-6.00	18
1.59 Polycarbonate													
HOYALUX iD LifeStyle®4 Indoor Poly	GR!	BR!								-6.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Outdoor Poly	GR!	BR!								-6.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Urban Poly	GR!	BR!								-6.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Indoor Poly	GR!	BR!								-6.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Outdoor Poly	GR!	BR!								-6.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Urban Poly	GR!	BR!								-6.00 to +6.00	+0.75 to +4.00	-6.00	14
Amp HD3 Poly	GR!	BR!								-6.00 to +6.00	+0.75 to +4.00	-6.00	18
Amp Mini HD3 Poly	GR!	BR!								-6.00 to +6.00	+0.75 to +4.00	-6.00	14
HOYALUX ARRAY™ 11, 13, 15, 17, VL Poly	GR	BR								-13.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL Poly	GR	BR								-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13, 15, 17, VL Poly	GR	BR								-13.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11, 13, 15, 17, VL Poly	GR	BR								-6.00 to +4.00	+0.75 to +4.00	-6.00	14, 17, 19, 21
Amplitude BKS Poly	GR	BR								-13.00 to +6.00	+0.75 to +4.00	-6.00	19
Amplitude Mini BKS Poly	GR	BR								-13.00 to +6.00	+0.75 to +4.00	-6.00	14
Summit cd BKS Poly	GR	BR								-13.00 to +6.00	+0.75 to +4.00	-6.00	14
Summit ecp BKS Poly	GR	BR								-13.00 to +6.00	+0.75 to +4.00	-6.00	18
Hoya GP Wide Poly	GR	BR								-13.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 40 BKS Poly	GR	BR								-13.00 to +6.00	+0.75 to +4.00	-6.00	18
Tact 60 BKS Poly	GR	BR								-13.00 to +6.00	+0.75 to +4.00	-6.00	18
1.60 High Index													
HOYALUX iD LifeStyle®4 Indoor 1.60	GR!	BR!								-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Outdoor 1.60	GR!	BR!								-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®4 Urban 1.60	GR!	BR!								-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Indoor 1.60	GR!	BR!								-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Outdoor 1.60	GR!	BR!								-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX iD LifeStyle®3 Urban 1.60	GR!	BR!								-11.00 to +8.00	+0.75 to +4.00	-6.00	14
Amp HD3 1.60	GR!	BR!								-11.00 to +8.00	+0.75 to +4.00	-6.00	14
Amp Mini HD3 1.60	GR!	BR!								-11.00 to +8.00	+0.75 to +4.00	-6.00	14
HOYALUX ARRAY 11, 13, 15, 17, VL 1.60	GR	BR								-10.00 to +6.00	+0.75 to +4.00	-6.00	14, 17, 19, 21

Nikon Infinite Single Vision	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	-6.00
Nikon Smart Boost	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	-6.00
1.53 TRIVEX® Single Vision										
Nikon SeeMax Infinite	GR	BR	GG						-10.00 to +7.00	-6.00
Nikon Infinite Single Vision	GR	BR	GG						-10.00 to +7.00	-6.00
Nikon Smart Boost	GR	BR	GG						-10.00 to +7.00	-6.00
1.59 Polycarbonate Single Vision										
Nikon SeeMax Infinite	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	-6.00
Nikon Infinite Single Vision	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	-6.00
Nikon Smart Boost	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	-6.00
1.60 High Index Single Vision										
Nikon SeeMax Infinite	GR	BR							-12.00 to +8.00	-6.00
Nikon Infinite Single Vision	GR	BR							-12.00 to +8.00	-6.00
Nikon Smart Boost	GR	BR							-12.00 to +8.00	-6.00
1.67 Super High Index Single Vision										
Nikon SeeMax Infinite	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	-6.00
Nikon Infinite Single Vision	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	-6.00
Nikon Smart Boost	GR	BR	GG	S	AME	AMB	E	R	-14.00 to +9.00	-6.00
1.74 ultra High Index Single Vision										
Nikon SeeMax Infinite	GR	BR							-16.00 to +10.00	-6.00
Nikon Infinite Single Vision	GR	BR							-16.00 to +10.00	-6.00
Nikon Smart Boost	GR	BR							-16.00 to +10.00	-6.00

												
Progressives									Sphere	Add	Cyl to	MFH
1.50 Standard index												
Nikon SeeMax Ultimate	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
Nikon Move	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
1.53 TRIVEX®												
Nikon SeeMax Ultimate	GR	BR	GG						-10.00 to +7.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR	GG						-10.00 to +7.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR	GG						-10.00 to +7.00	+0.50 to +4.00	-6.00	15, 17
Nikon Move	GR	BR	GG						-10.00 to +7.00	+0.50 to +4.00	-6.00	15, 17
1.59 Polycarbonate												
Nikon SeeMax Ultimate	GR	BR	GG	S	AME	AMB	E	R	-11.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR	GG	S	AME	AMB	E	R	-11.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR	GG	S	AME	AMB	E	R	-11.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
Nikon Move	GR	BR	GG	S	AME	AMB	E	R	-11.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
1.60 High Index												
Nikon SeeMax Ultimate	GR	BR							-12.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR							-12.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR							-12.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
Nikon Move	GR	BR							-12.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
1.67 Super High Index												
Nikon SeeMax Ultimate	GR	BR	GG	S	AME	AMB	E		-14.00 to +9.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR	GG	S	AME	AMB	E		-14.00 to +9.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR	GG	S	AME	AMB	E		-14.00 to +9.00	+0.50 to +4.00	-6.00	15, 17
Nikon Move	GR	BR	GG	S	AME	AMB	E		-14.00 to +9.00	+0.50 to +4.00	-6.00	15, 17
1.74 Ultra High Index												
Nikon SeeMax Ultimate	GR	BR							-16.00 to +10.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR							-16.00 to +10.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR							-16.00 to +10.00	+0.50 to +4.00	-6.00	15, 17

										
Finished Single Vision									Sphere	Cyl to
1.59 FSV										
Shamir FSV with Glacier PLUS™ AR	GR								-6.00 to +4.00	-2.00


										
Semi-Finished Single Vision									Sphere	Cyl to
1.50 SFSV										
Shamir Autograph® III Single Vision	GR	BR							-12.00 to +7.25	-6.00
Shamir Relax™	GR	BR							-14.00 to +8.00	-6.00
1.53 SFSV										
Shamir Autograph® III Single Vision	GR	BR							-14.00 to +8.00	-6.50
Shamir Relax™	GR	BR							-14.00 to +8.00	-6.50
1.59 SFSV										
Shamir Autograph® III Single Vision	GR	BR							-14.00 to +8.75	-7.00
Shamir Relax™	GR	BR							-14.75 to +6.00	-7.00
1.60 SFSV										
Shamir Autograph® III Single Vision	GR	BR							-14.00 to +8.75	-7.00
Shamir Relax™	GR	BR							-14.00 to +8.75	-7.00
1.67 SFSV										
Shamir Autograph® III Single Vision	GR	BR							-16.00 to +16.00	-8.00
Shamir Relax™	GR	BR							-16.75 to +16.00	-8.00
1.74 SFSV										
Shamir Autograph® III Single Vision	GR	GR							-20.00 to +13.50	-9.00










												
Blended Bifocal (digital)									Sphere	Add	Cyl to	MFH
1.50 Standard Index Blended Bifocal (digital)												
Shamir Duo	GR	BR							-9.00 to +6.00	+0.75 to +4.00	-6.00	15
1.53 Trivex® Blended Bifocal (digital)												
Shamir Duo	GR	BR							-13.25 to +6.00	+0.75 to +4.00	-6.50	15
1.59 Polycarbonate Blended Bifocal (digital)												
Shamir Duo	GR	BR							-14.75 to +6.00	+0.75 to +4.00	-7.00	15
1.60 High Index Blended Bifocal (digital)												
Shamir Duo	GR	BR							-15.00 to +6.00	+0.75 to +4.00	-7.00	15
1.67 Super High Index Blended Bifocal (digital)												
Shamir Duo	GR	BR							-16.75 to +8.00	+0.75 to +4.00	-8.00	15
1.74 Ultra High Index Blended Bifocal (digital)												
Shamir Duo	GR	BR							-18.75 to +12.00	+0.75 to +4.00	-9.00	15

												
Progressives									Sphere	Add	Cyl to	MFH
1.50 Standard Index												
Shamir Autograph Intelligence™	GR	BR							-12.25 to +6.50	+0.75 to +4.00	-6.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR							-12.25 to +6.00	+0.75 to +4.00	-6.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR							-12.25 to +6.00	+0.75 to +4.00	-6.00	11, 13, 15, 18
Shamir Intouch™	GR	BR							-12.25 to +6.00	+0.75 to +3.50	-6.00	14, 16, 18
Shamir Spectrum+™	GR	BR							-12.25 to +6.00	+0.75 to +3.50	-6.00	14, 16, 18
Shamir Element™	GR	BR							-7.00 to +4.00	+0.75 to +3.50	-6.00	16, 19

Shamir GENESIS HD	GR	BR								-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir Computer™	GR	BR								-7.00 to +4.00	+0.75 to +2.00	-6.00	14 Lower, 10 Upper
Shamir Workspace™	GR	BR								-6.00 to +6.00	+0.75 to +2.50	-6.00	14 Lower, 10 Upper
Shamir Attitude III Sport	GR	BR								-6.50 to +5.25	+0.75 to +3.50	-6.00	18
Shamir Attitude III Fashion	GR	BR								-6.50 to +5.25	+0.75 to +3.50	-6.00	15, 18
1.53 TRIVEX®													
Shamir Autograph Intelligence™	GR	BR								-13.25 to +7.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
Shamir Autograph® III	GR	BR								-13.25 to +7.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
Shamir Autograph® II+	GR	BR								-13.25 to +6.00	+0.75 to +4.00	-6.50	11, 13, 15, 18
Shamir Intouch™	GR	BR								-13.25 to +6.00	+0.75 to +3.50	-6.50	14, 16, 18
Shamir Golf™	GR	BR								-5.00 to +4.00	+1.00 to +3.50	-6.50	19
Shamir Spectrum+™	GR	BR								-13.25 to +6.00	+0.75 to +3.50	-6.50	14, 16, 18
Shamir Element™	GR	BR								-13.25 to +6.00	+0.75 to +3.50	-6.50	16, 19
Shamir GENESIS HD	GR	BR								-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir Computer™	GR	BR								-6.00 to +8.00	+0.75 to +2.00	-6.50	14 Lower, 10 Upper
Shamir Workspace™	GR	BR								-6.00 to +8.00	+0.75 to +2.50	-6.50	14 Lower, 10 Upper
Shamir Attitude III Sport	GR	BR								-6.50 to +6.00	+0.75 to +3.50	-6.50	18
Shamir Attitude III Fashion	GR	BR								-6.50 to +6.00	+0.75 to +3.50	-6.50	15, 18
1.59 Polycarbonate													
Shamir Autograph Intelligence™	GR	BR								-14.75 to +7.50	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR								-14.75 to +8.00	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR								-14.75 to +6.50	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Intouch™	GR	BR								-14.75 to +6.00	+0.75 to +3.50	-7.00	15, 18
Shamir Golf™	GR	BR								-6.00 to +5.25	+1.00 to +3.50	-7.00	19
Shamir Spectrum+™	GR	BR								-14.75 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
Shamir Element™	GR	BR								-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19
Shamir GENESIS HD	GR	BR								-10.00 to +6.00	+0.75 to +3.50	-7.00	16
Shamir Computer™	GR	BR								-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper
Shamir Workspace™	GR	BR								-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper
Shamir Attitude III Sport	GR	BR								-6.50 to +6.00	+0.75 to +3.50	-6.50	18
Shamir Attitude III Fashion	GR	BR								-6.50 to +6.00	+0.75 to +3.50	-6.50	15, 18
1.60 High Index													
Shamir Autograph Intelligence™	GR	BR								-14.75 to +8.75	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR								-15.00 to +8.00	+0.75 to +3.50	-7.00	15, 18
Shamir Autograph® II+	GR	BR								-15.00 to +6.00	+0.75 to +4.00	-7.00	11, 13, 15, 18
Shamir Golf™ SuperLite™ 1.6	GR	BR								-6.00 to +5.25	+1.00 to +3.50	-7.00	19
Shamir Spectrum+™	GR	BR								-15.00 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
Shamir Element™	GR	BR								-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19
Shamir GENESIS HD	GR	BR								-6.00 to +6.00	+0.75 to +3.50	-7.00	16
Shamir Computer™	GR	BR								-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper
Shamir Workspace™	GR	BR								-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper
Shamir InTouch™	GR	BR								-15.00 to +6.00	+0.75 to +3.50	-7.00	15, 18
Shamir Attitude III Sport	GR	BR								-6.50 to +6.00	+0.75 to +3.50	-7.00	18
Shamir Attitude III Fashion	GR	BR								-6.50 to +6.00	+0.75 to +3.50	-7.00	15, 18
1.67 Super High Index													
Shamir Autograph Intelligence™	GR	BR								-16.75 to +10.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Autograph® III	GR	BR								-16.75 to +9.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR								-16.75 to +8.25	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Spectrum+™	GR	BR								-16.75 to +8.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Element™ SuperLite™ 1.67	GR	BR								-12.50 to +8.00	+0.75 to +3.50	-8.00	16, 19
Shamir GENESIS HD	GR	BR								-6.00 to +6.00	+0.75 to +3.50	-4.00	16
Shamir InTouch™	GR	BR								-16.75 to +8.00	+0.75 to +3.50	-8.00	15, 18
Shamir Attitude III Sport	GR	BR								-6.50 to +6.00	+0.75 to +3.50	-8.00	18
Shamir Attitude III Fashion	GR	BR								-6.50 to +6.00	+0.75 to +3.50	-8.00	15, 18
1.74 Ultra High Index													
Shamir Autograph Intelligence™	GR	BR								-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18

Shamir Autograph® III	GR	BR								-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir Autograph® II+	GR	BR								-18.75 to +12.00	+0.75 to +4.00	-9.00	11, 13, 15, 18
Shamir GENESIS HD	GR	BR								-6.00 to +6.00	+0.75 to +4.00	-4.00	16
Shamir Spectrum+™	GR	BR								-18.75 to +11.50	+0.75 to +3.50	-9.00	14, 16, 18
Shamir InTouch™	GR	BR								-18.75 to +15.50	+0.75 to +3.50	-9.00	15, 18

										
Semi-Finished Single Vision									Sphere	Cyl to
1.50 Standard Index										
Signetek Processed SV	GR	BR	GG	S	AME	AMB	E	R	-8.00 to +6.00	-6.00
KODAK Digital SV	GR	BR	GG	S	AME	AMB	E	R	-8.00 to +6.00	-6.00
KODAK PowerUp Lens (.40/.66/1.00D)	GR	BR							-8.00 to +6.00	-6.00
Crossbows Custom SV	GR	BR							-8.00 to +6.00	-6.00
1.53 TRIVEX® Single Vision										
Signetek Processed SV	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +6.00	-6.00
KODAK Digital SV	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +6.00	-6.00
Crossbows Custom SV	GR	BR							-10.00 to +6.00	-6.00
1.59 Polycarbonate Single Vision										
Signetek Processed SV	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +6.50	-6.00
KODAK Digital SV	GR	BR	GG	S	AME	AMB	E	R	-10.00 to +6.50	-6.00
KODAK PowerUp Lens (.40/.66/1.00D)	GR	BR							-10.00 to +6.50	-6.00
Crossbows Custom SV	GR	BR							-10.00 to +6.50	-6.00
1.60 High Index Single Vision										
Signetek Processed SV	GR	BR							-10.00 to +7.50	-6.00
KODAK Digital SV	GR	BR							-10.00 to +7.50	-6.00
Crossbows Custom SV	GR	BR							-10.00 to +7.50	-6.00
1.67 Super High Index Single Vision										
Signetek Processed SV	GR	BR	GG	S	AME	AMB	E	R	-11.00 to +7.50	-6.00
KODAK Digital SV	GR	BR	GG	S	AME	AMB	E	R	-11.00 to +7.50	-6.00
KODAK PowerUp Lens (.40/.66/1.00D)	GR	BR							-11.00 to +6.00	-6.00
Crossbows Custom SV	GR	BR							-11.00 to +7.50	-6.00
1.74 Ultra High Index Single Vision										
Signetek™ Processed SV	GR	BR							-13.00 to +9.00	-6.00
KODAK Digital SV	GR	BR							-13.00 to +9.00	-6.00

												
Progressives									Sphere	Add	Cyl to	MFH
1.50 Standard Index												
KODAK Precise® 1.50 SF	GR								-9.00 to +5.00	+1.00 to +3.00	-4.00	17
KODAK Precise® Digital 1.50	GR								-8.00 to +5.00	+1.00 to +3.00	-6.00	17
KODAK Precise® Plus 1.50	GR	BR	GG	S	AME	AMB	E	R	-8.00 to +4.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus 1.50	GR	BR	GG	S	AME	AMB	E	R	-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Precise® PB 1.50	GR	BR							-8.00 to +4.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB 1.50	GR	BR							-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique™ 1.50	GR	BR	GG	S	AME	AMB	E	R	-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique™ HD 1.50	GR	BR							-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ 1.50	GR	BR	GG	S	AME	AMB	E	R	-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.50	GR	BR	GG	S	AME	AMB	E	R	-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ 1.50	GR	BR	GG	S	AME	AMB	E	R	-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ HD 1.50	GR	BR	GG	S	AME	AMB	E	R	-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Easy 14 1.50	GR	BR	GG	S	AME	AMB	E	R	-8.00 to +4.00	+0.75 to +3.50	-6.00	14
KODAK Easy 18 1.50	GR	BR	GG	S	AME	AMB	E	R	-8.00 to +4.00	+0.75 to +3.50	-6.00	18
DirecTek 1.50	GR	BR							-8.00 to +4.00	+0.75 to +3.50	-6.00	16.5

DirecTek Short 1.50	GR	BR								-8.00 to +4.00	+0.75 to +3.50	-6.00	13.5
Navigator FBS 1.50	GR	BR								-8.00 to +4.00	+0.75 to +3.50	-6.00	18
Navigator Short FBS 1.50	GR	BR								-8.00 to +4.00	+0.75 to +3.50	-6.00	14
Crossbows Custom U 1.50	GR	BR								-8.00 to +4.00	+0.75 to +3.50	-6.00	13
Crossbows Custom V 1.50	GR	BR								-8.00 to +4.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 13 Progressive 1.50	GR	BR								-8.00 to +4.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 17 Progressive 1.50	GR	BR								-8.00 to +4.00	+0.75 to +3.50	-6.00	17
Crossbows Junior	GR	BR								-8.00 to +4.00	+0.75 to +3.50	-6.00	13
1.53 TRIVEX®													
KODAK Precise® PB Trivex®	GR	BR								-10.00 to +4.50	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB Trivex®	GR	BR								-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Precise® Plus Trivex®	GR	BR								-10.00 to +4.50	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus Trivex®	GR	BR								-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique Trivex®	GR	BR								-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique HD Trivex®	GR	BR								-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ Trivex®	GR	BR								-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD Trivex®	GR	BR								-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ Trivex®	GR	BR								-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ HD Trivex®	GR	BR								-10.00 to +4.50	+0.75 to +3.50	-6.00	13
Crossbows Custom U Trivex®	GR	BR								-10.00 to +4.50	+0.75 to +3.50	-6.00	13
Crossbows Custom V Trivex®	GR	BR								-10.00 to +4.50	+0.75 to +3.50	-6.00	13
1.59 Polycarbonate													
KODAK Precise® Poly SF	GR									-10.00 to +5.00	+1.00 to +3.00	-4.00	17
KODAK Precise® Digital Poly	GR									-10.00 to +6.00	+1.00 to +3.00	-6.00	17
KODAK Precise® Plus Poly	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +5.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus Poly	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ Poly	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Precise® PB Poly	GR	BR	GG							-10.00 to +5.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB Poly	GR	BR	GG							-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique Poly	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD Poly	GR	BR								-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD Poly	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ Poly	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ HD Poly	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Easy 14 Poly	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +5.00	+0.75 to +3.50	-6.00	14
KODAK Easy 18 Poly	GR	BR	GG	S	AME	AMB	E	R		-10.00 to +5.00	+0.75 to +3.50	-6.00	18
DirecTek Poly	GR	BR								-10.00 to +5.00	+0.75 to +3.50	-6.00	16.5
DirecTek Short Poly	GR	BR								-10.00 to +5.00	+0.75 to +3.50	-6.00	13.5
Navigator FBS Poly	GR	BR								-10.00 to +5.00	+0.75 to +3.50	-6.00	18
Navigator Short FBS Poly	GR	BR								-10.00 to +5.00	+0.75 to +3.50	-6.00	14
Crossbows Custom U Poly	GR	BR								-10.00 to +5.00	+0.75 to +3.50	-6.00	13
Crossbows Custom V Poly	GR	BR								-10.00 to +5.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 13 Progressive Poly	GR	BR								-10.00 to +5.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 17 Progressive Poly	GR	BR								-10.00 to +5.00	+0.75 to +3.50	-6.00	17
Crossbows Junior Poly	GR	BR								-10.00 to +5.00	+0.75 to +3.50	-6.00	13
1.60 High Index													
KODAK Unique	GR	BR								-10.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD	GR	BR								-10.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ 1.60	GR	BR								-10.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.60	GR	BR								-10.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ 1.60	GR	BR								-10.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ HD 1.60	GR	BR								-10.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Custom U 1.60	GR	BR								-10.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Custom V 1.60	GR	BR								-10.00 to +6.00	+0.75 to +3.50	-6.00	13
1.67 Super High Index													
KODAK Precise® PB 1.67	GR	BR	GG							-11.00 to +6.00	+0.75 to +3.50	-6.00	17

KODAK Precise® Short PB 1.67	GR	BR	GG							-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Precise® Plus Poly	GR	BR	GG	S	AME	AMB	E	R		-11.00 to +6.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus Poly	GR	BR	GG	S	AME	AMB	E	R		-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique	GR	BR	GG	S	AME	AMB	E	R		-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD	GR	BR								-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ 1.67	GR	BR	GG	S	AME	AMB	E	R		-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.67	GR	BR	GG	S	AME	AMB	E	R		-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ 1.67	GR	BR	GG	S	AME	AMB	E	R		-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ HD 1.67	GR	BR	GG	S	AME	AMB	E	R		-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Easy 14 Poly	GR	BR	GG	S	AME	AMB	E	R		-11.00 to +6.00	+0.75 to +3.50	-6.00	14
KODAK Easy 18 Poly	GR	BR	GG	S	AME	AMB	E	R		-11.00 to +6.00	+0.75 to +3.50	-6.00	18
DirceTek 1.67	GR	BR								-11.00 to +6.00	+0.75 to +3.50	-6.00	16.5
DirceTek Short 1.67	GR	BR								-11.00 to +6.00	+0.75 to +3.50	-6.00	13.5
Navigator FBS 1.67	GR	BR								-11.00 to +6.00	+0.75 to +3.50	-6.00	18
Navigator Short FBS 1.67	GR	BR								-11.00 to +6.00	+0.75 to +3.50	-6.00	14
Crossbows Custom U 1.67	GR	BR								-11.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Custom V 1.67	GR	BR								-11.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 13 Progressive 1.67	GR	BR								-11.00 to +6.00	+0.75 to +3.50	-6.00	13
Crossbows Basic 17 Progressive 1.67	GR	BR								-11.00 to +6.00	+0.75 to +3.50	-6.00	17
Crossbows Junior 1.67	GR	BR								-11.00 to +6.00	+0.75 to +3.50	-6.00	13
1.74 Ultra High Index													
KODAK Unique	GR	BR								-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD	GR	BR								-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™	GR	BR								-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.74	GR	BR								-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™	GR	BR								-13.00 to +9.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ HD 1.74	GR	BR								-13.00 to +9.00	+0.75 to +3.50	-6.00	13

										
Semi-Finished Single Vision									Sphere	Cyl to
1.50 Standard Index										
TRN 1.50 DS (76mm) UC/HC DIA.72	GR	BR							Base 0.50, 8.75, 10.50	
TRN 1.50 DS (76mm) UC/HC DIA.76	GR	BR							Base 2.00, 3.50, 4.25, 5.25, 6.25, 7.25	
1.53 TRIVEX® Single Vision										
TOG 1.53 Trivex (75mm) SV UC/HC DIA.70	GR	BR							Base 1.00	
TOG 1.53 Trivex (75mm) SV UC/HC DIA.75	GR	BR							Base 2.50, 4.00, 5.00, 6.00, 7.00, 8.00	
1.59 Polycarbonate Single Vision										
TRN 1.59 Polycarbonate (79mm) SV HC DIA.79	GR	BR							Base 1.50, 3.25, 4.25, 4.75, 5.50, 7.00	
TRN 1.59 Polycarbonate (74mm) SV HC DIA.74	GR	BR							Base 7.50, 8.50	
1.60 High Index Single Vision										
TOG 1.60 SV (75mm) UC/HC DIA.65	GR	BR							Base 10.00, 11.00	
TOG 1.60 SV (75mm) UC/HC DIA.70	GR	BR							Base 9.00	
TOG 1.60 SV (75mm) UC/HC DIA.75	GR	BR							Base 1.00, 2.50, 4.00, 5.00, 6.00, 7.00, 8.00	
1.67 Super High Index Single Vision										
TOG 1.67 SV (65mm) UC/HC DIA.65	GR	BR							Base 10.00, 11.00	
TOG 1.67 SV (70mm) UC/HC DIA.70	GR	BR							Base 9.00	
TOG 1.67 SV (75mm) UC/HC DIA.75	GR	BR							Base 1.00, 2.50, 4.00, 5.00, 6.00, 7.00, 8.00	
1.74 ultra High Index Single Vision										
TOG 1.74 LeafECO SV (75mm) UC/HC DIA.65	GR	BR							Base 11.50, 12.50, 13.50, 14.50, 15.50, 16.50	

TOG 1.74 LeafECO SV (75mm) UC/HC DIA.70	GR	BR									Base 0.50, 1.50, 2.25, 6.75, 8.50, 10.25
TOG 1.74 LeafECO SV (75mm) UC/HC DIA.75	GR	BR									Base 3.50, 5.00

											
Semi-Finished Single Vision										Sphere	Cyl to
1.50 Standard Index											
	Unity SVx	GR	BR							-12.25 to +6.00	-6.00
	Unity SVxtra	GR	BR							-12.25 to +6.00	-6.00
	Unity SVxtreme	GR	BR							-6.00 to +6.00	-6.00
1.53 TRIVEX® Single Vision											
	Unity SVx	GR	BR							-14.00 to +6.00	-6.00
	Unity SVxtra	GR	BR							-14.00 to +6.00	-6.00
	Unity SVxtreme	GR	BR							-6.50 to +6.00	-6.00
1.59 Polycarbonate Single Vision											
	Unity SVx	GR	BR							-14.00 to +6.00	-6.00
	Unity SVxtra	GR	BR							-14.00 to +6.00	-6.00
	Unity SVxtreme	GR	BR							-6.50 to +6.00	-6.00
1.60 High Index Single Vision											
	Unity SVx	GR	BR							-14.00 to +6.00	-6.00
	Unity SVxtra	GR	BR							-14.00 to +6.00	-6.00
	Unity SVxtreme	GR	BR							-6.50 to +6.00	-6.00
1.67 Super High Index Single Vision											
	Unity SVx	GR	BR							-16.00 to +8.00	-6.00
	Unity SVxtra	GR	BR							-16.00 to +8.00	-6.00
	Unity SVxtreme	GR	BR							-6.00 to +6.00	-6.00





													
Progressives										Sphere	Add	Cyl to	MFH
1.50 Standard Index													
	Unity Via II	GR	BR							-7.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Plus II	GR	BR							-7.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Elite II	GR	BR							-7.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Mobile II	GR	BR							-7.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Wrap II	GR	BR							-7.00 to +6.00	+0.50 to +4.50	-6.00	12
	Ethos	GR	BR							-9.00 to +4.50	0.75 to +3.50	-4.00	14
	Ethos Plus	GR	BR							-9.00 to +4.50	0.75 to +3.50	-4.00	14
1.53 TRIVEX®													
	Unity Via II	GR	BR							-10.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Plus II	GR	BR							-10.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Elite II	GR	BR							-10.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Mobile II	GR	BR							-10.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Wrap II	GR	BR							-5.00 to +4.00	+0.50 to +4.50	-6.00	12
	Ethos	GR	BR							-9.00 to +6.00	0.75 to +3.50	-4.00	14
	Ethos Plus	GR	BR							-9.00 to +6.00	0.75 to +3.50	-4.00	14
1.59 Polycarbonate													
	Unity Via II	GR	BR							-10.00 to +6.50	+0.50 to +4.50	-6.00	12
	Unity Via Plus II	GR	BR							-10.00 to +6.50	+0.50 to +4.50	-6.00	12
	Unity Via Elite II	GR	BR							-10.00 to +6.50	+0.50 to +4.50	-6.00	12
	Unity Via Mobile II	GR	BR							-10.00 to +6.50	+0.50 to +4.50	-6.00	12
	Unity Via Wrap II	GR	BR							-6.50 to +6.00	+0.50 to +4.50	-6.00	12
	Ethos	GR	BR							-9.00 to +6.00	0.75 to +3.50	-4.00	14
	Ethos Plus	GR	BR							-9.00 to +6.00	0.75 to +3.50	-4.00	14

1.60 High Index													
Unity Via II	GR	BR								-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR								-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR								-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR								-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR								-6.50 to +5.25	+0.50 to +4.50	-6.00	12
Ethos	GR	BR								-9.00 to +5.00	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR								-9.00 to +5.00	0.75 to +3.50	-4.00	14
1.67 Super High Index													
Unity Via II	GR	BR								-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR								-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR								-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR								-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR								-6.00 to +5.00	+0.50 to +4.50	-6.00	12
Ethos	GR	BR								-12.00 to +8.00	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR								-12.00 to +8.00	0.75 to +3.50	-4.00	14

YOUNGER OPTICS The Optical Lens Innovators										Transitions gen#S						
Finished Single Vision															Sphere	Cyl to
1.50 Standard Index																
Younger Optics FSV AR Super Hydrophobic	GR	BR													-4.00 to +4.00	-2.00
1.59 Polycarbonate Finished Single Vision																
Younger FSV Poly AR Super Hydrophobic	GR	BR													-6.00 to +4.00	-2.00

YOUNGER OPTICS The Optical Lens Innovators										Transitions gen#S						
Semi-Finished Single Vision															Sphere	Cyl to
1.50 Standard Index																
Younger Optics SFSV (72mm)	GR	BR	GG	SA	AME	AMB	E	R							-9.00 to +6.00	-4.00
1.53 TRIVEX® Single Vision																
Younger Optics Trilogy® Spherical SV	GR	BR													-8.00 to +7.50	-4.00
1.59 Polycarbonate Single Vision																
Younger Optics SV	GR	BR	GG	SA	AME	AMB	E	R							-9.00 to +6.00	-4.00
1.60 High Index Single Vision																
Younger Optics Single Vision 1.60	GR	BR	GG	SA	AME	AMB	E	R							-11.00 to +8.00	-4.00
1.67 Super High Index Single Vision																
Younger Optics Single Vision 1.67 MR-10™	GR	BR	GG	SA	AME	AMB	E	R							-11.00 to +7.00	-4.00
1.74 Ultra High Index Single Vision																
Younger Optics Single Vision 1.74	GR	BR	GG	SA	AME	AMB	E	R							-18.00 to +12.00	-4.00

YOUNGER OPTICS The Optical Lens Innovators										Transitions gen#S						
Digital															Sphere	Cyl to
1.53 TRIVEX®																
Camber™	GR	BR													-11.00 to +6.25	-4.00
1.59 Polycarbonate																
Camber™	GR	BR													-10.00 to +6.25	-4.00
1.60 High Index																
Camber™	GR	BR													-11.50 to +7.25	-4.00
1.67 Super High Index																
Camber™	GR	BR													-14.00 to +8.00	-4.00
1.74 Ultra High Index																
Camber™	GR	BR													-19.00 to +9.00	-4.00

													
Progressives										Sphere	Add	Cyl to	MFH
1.53 TRIVEX®													
Younger Trilogy® Image®		GR								-8.00 to +8.00	+1.00 to +3.00	-4.00	18
1.59 index polycarbonate													
Younger Image®		GR	BR							-9.00 to +7.00	+1.00 to +3.00	-4.00	18
1.67 Super High Index													
Younger Image® 1.67 MR-10™		GR	BR							-10.00 to +8.50	+1.00 to +3.50	-4.00	18

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

Unless indicated otherwise, all registered trademarks and trademarks are the property of Essilor International and/or its subsidiaries in the United States and in other countries.

All other trademarks are the property of their respective owners.