

Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Suppliers](#) / **Essilor Lenses**

Essilor Lenses — Entire Product Line

														
Finished Single Vision													Sphere	Cyl to
1.50 Standard Index FSV														
Essilor FSV Spherical	GR	BR							GR				-4.00 to +4.00	-2.00
Essilor FSV Aspheric	GR												-4.00 to +4.00	-2.00
Essilor FSV 1.50 Crizal® Sapphire™ HR	GR								GR				-4.00 to +4.00	-2.00
Essilor FSV 1.50 Crizal® Rock™	GR								GR				-4.00 to +4.00	-2.00
Essilor FSV 1.50 Crizal® Easy Pro	GR								GR				-4.00 to +4.00	-2.00
1.59 Polycarbonate Finished Single Vision														
Essilor FSV Spherical Airwear®	GR	BR											-6.00 to +4.00	-2.00
Essilor FSV Aspheric Airwear®	GR	BR							GR				-6.00 to +4.00	-2.00
Essilor FSV Airwear® Crizal® Sapphire™ HR	GR								GR				-6.00 to +4.00	-2.00
Essilor FSV Airwear® Crizal® Rock™	GR								GR				-6.00 to +4.00	-2.00
Essilor FSV Airwear® Crizal® Easy Pro	GR								GR				-6.00 to +4.00	-2.00

														
Semi-Finished Single Vision													Sphere	Cyl to
1.50 Standard Index SFSV														
Essilor Spherical Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.50 to +10.00	-4.00
Essilor Stylistic® Single Vision	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	-6.00
Eyezen® +0	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	-6.00
Eyezen® +1	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	-6.00
Eyezen® +2	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	-6.00
Eyezen® +3	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	-6.00
Eyezen® +4	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +7.50	-6.00
Eyezen® Start	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +6.00	-6.00
1.53 TRIVEX® Single Vision														
Essilor Spherical Single Vision TREXA™	GR	BR	GG						GR	BR			-10.00 to +6.00	-4.00
Essilor Eyecode™ Single Vision TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	-6.00
Essilor FIT™ Single Vision TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	-6.00
Essilor Azio™ Single Vision TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	-6.00
Essilor 360°™ Single Vision TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	-6.00
Eyezen® +0 TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	-6.00
Eyezen® +1 TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	-6.00
Eyezen® +2 TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	-6.00
Eyezen® +3 TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	-6.00
Eyezen® +4 TREXA™	GR	BR	GG						GR	BR			-10.00 to +7.50	-6.00
Eyezen® Start	GR	BR	GG						GR	BR			-10.00 to +6.00	-6.00
1.59 Polycarbonate Single Vision														
Essilor Spherical Single Vision Airwear®	GR	BR	GG	S	AME	AMB	E	R			GG	GR	-10.00 to	-4.00

										GR	BR					+8.00	
Essilor Aspheric Single Vision Airwear®	GR	BR	GG							GR	BR	GG				-10.00 to +8.00	-4.00
Essilor Eyecode™ Aspheric SV Airwear®	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-10.00 to +8.00	-6.00
Essilor FIT™ Aspheric SV Airwear®	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-10.00 to +8.00	-6.00
Essilor AZIO™ Aspheric SV Airwear®	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-10.00 to +8.00	-6.00
Essilor SV 360°™ Aspheric Airwear®	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-10.00 to +8.00	-6.00
Essilor Stylistic® Single Vision Airwear®	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-10.00 to +8.00	-6.00
Eyezen® +0 Airwear®	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-10.00 to +8.00	-6.00
Eyezen® +1 Airwear®	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-10.00 to +8.00	-6.00
Eyezen® +2 Airwear®	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-10.00 to +8.00	-6.00
Eyezen® +3 Airwear®	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-10.00 to +8.00	-6.00
Eyezen® +4 Airwear®	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-10.00 to +8.00	-6.00
Eyezen® Start Airwear®	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-10.00 to +8.00	-6.00
Eyezen®Kids	GR	BR	GG	S	AME	AMB	E	R								-10.00 to +6.00	-6.00
1.60 High Index Single Vision																	
Essilor Spherical Single Vision Thin&Lite® 1.60	GR	BR	GG	S	AME	AMB	E	R								-13.00 to +8.00	-4.00
Eyezen® +0 Thin&Lite® 1.60	GR	BR														-12.00 to +8.00	-6.00
Eyezen® +1 Thin&Lite® 1.60	GR	BR														-12.00 to +8.00	-6.00
Eyezen® +2 Thin&Lite® 1.60	GR	BR														-12.00 to +8.00	-6.00
Eyezen® +3 Thin&Lite® 1.60	GR	BR														-12.00 to +8.00	-6.00
Eyezen® +4 Thin&Lite® 1.60	GR	BR														-12.00 to +8.00	-6.00
Eyezen® Start Thin&Lite® 1.60	GR	BR														-12.00 to +8.00	-6.00
1.67 Super High Index Single Vision																	
Essilor Spherical Single Vision Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-14.00 to +9.00	-4.00
Essilor Aspheric Single Vision Thin&Lite® 1.67	GR	BR	GG													-16.00 to +8.00	-4.00
Essilor Eyecode™ SV Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-14.00 to +9.00	-6.00
Essilor FIT™ SV Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-14.00 to +9.00	-6.00
Essilor Azio™ SV Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-14.00 to +9.00	-6.00
Essilor SV 360°™ Spherical Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-14.00 to +9.00	-6.00
Essilor Stylistic® Single Vision Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-14.00 to +9.00	-6.00
Eyezen® +0 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-14.00 to +9.00	-6.00
Eyezen® +1 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-14.00 to +9.00	-6.00
Eyezen® +2 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-14.00 to +9.00	-6.00
Eyezen® +3 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-14.00 to +9.00	-6.00
Eyezen® +4 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-14.00 to +9.00	-6.00
Eyezen® Start Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R		GR	BR	GG		GR		-14.00 to +9.00	-6.00
1.74 Ultra High Index Single Vision																	
Essilor Spherical Single Vision Thin&Lite® 1.74	GR	BR														-20.00 to +13.00	-4.00
Essilor Eyecode™ SV Thin&Lite® 1.74	GR	BR														-16.00 to +13.00	-6.00
Essilor FIT™ SV Thin&Lite® 1.74	GR	BR														-16.00 to +13.00	-6.00
Essilor 360°™ SV Thin&Lite® 1.74	GR	BR														-16.00 to +13.00	-6.00
Eyezen® +0 Thin&Lite® 1.74	GR	BR														-20.00 to +13.00	-6.00
Eyezen® +1 Thin&Lite® 1.74	GR	BR														-20.00 to +13.00	-6.00
Eyezen® +2 Thin&Lite® 1.74	GR	BR														-20.00 to +13.00	-6.00
Eyezen® +3 Thin&Lite® 1.74	GR	BR														-20.00 to	-6.00

	Transitions Gen [#] S								Transitions XTRACTIVE [®]			Transitions XTRACTIVE [®] POLARIZED [®]					
Progressives														Sphere	Add	Cyl to	MFH
1.53 TRIVEX [®]																	
Varilux [®] XR Track™	GR	BR												-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux [®] XR Design™	GR	BR												-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux [®] X 4D™ TREXA™	GR	BR												-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux [®] X FIT™ TREXA™	GR	BR												-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux [®] X Design™ TREXA™	GR	BR												-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux [®] Physio [®] W3+ Eyecode TREXA™	GR	BR							GR	BR				-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux [®] Physio [®] W3+ FIT TREXA™	GR	BR							GR	BR				-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux [®] Physio [®] W3+ TREXA™	GR	BR							GR	BR				-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort [®] Max Fit TREXA™	GR	BR							GR	BR				-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort [®] Max TREXA™	GR	BR							GR	BR				-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™ TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Essilor Ideal Advanced™ FIT TREXA™	GR	BR	GG											-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal [®] FIT TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal [®] Sport TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	18
Essilor Ideal [®] Computer TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	13
Essilor Ideal [®] Driving TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal [®] TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal [®] Short TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Natural [®] Digital TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	15
Essilor Accolade [®] TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar [®] Digital TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar [®] Short Digital TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Brillenz [®] TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Brillenz [®] Short TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Ovation [®] Digital TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor SmallFit [®] Digital TREXA™	GR	BR	GG						GR	BR				-10.00 to +7.50	+0.75 to +3.00	-6.00	14

	Transitions Gen [#] S								Transitions XTRACTIVE [®]			Transitions XTRACTIVE [®] POLARIZED [®]					
Progressives														Sphere	Add	Cyl to	MFH
1.59 Polycarbonate																	
Varilux [®] XR Track™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux [®] XR Design™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux [®] X 4D™ Airwear [®]	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux [®] X FIT™ Airwear [®]	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG			-10.00 to +6.00	+0.75 to +4.00	-6.00	14

Varilux® X Design™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Varilux Comfort DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux® Liberty™ 3.0	GR	BR											-10.00 to +8.00	+0.75 to +1.50	-6.00	14
Varilux Comfort® Airwear®	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short™ Airwear®	GR	BR	GG						GR	BR	GG		-10.00 to +6.00	+0.75 to +3.00	-6.00	14
Varilux® Stylistic™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ FIT	GR	BR	GG	S	AME	AMB	E	R					-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Wrap	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Essilor Ideal® FIT Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Sport Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Driving Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Ideal® Short Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Natural® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +8.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Ovation® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Ovation® Airwear®	GR	BR											-10.00 to +8.00	+0.75 to +3.50	-4.00	17
Essilor SmallFit® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-10.00 to +8.00	+0.75 to +3.00	-6.00	14

	Transitions Gen[®]S							Transitions[®] XTRA[®]CTIVE[®]			Transitions[®] XTRA[®]CTIVE[®] POLARIZED[®]					
Progressives													Sphere	Add	Cyl to	MFH
1.60 High Index																
Varilux [®] Physio [®] W3+ Eyecode Thin&Lite [®] 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux [®] Physio [®] W3+ FIT Thin&Lite [®] 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux [®] Physio [®] W3+ Thin&Lite [®] 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort [®] Max Fit Thin&Lite [®] 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Varilux Comfort [®] Max Thin&Lite [®] 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Thin&Lite [®] 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17

	Transitions gen#S							Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED						
Progressives													Sphere	Add	Cyl to	MFH	
1.67 Super High Index																	
Varilux® XR Track™	GR	BR	GG	S	AME	AMB	E	R					-12.00 to +9.00	+0.75 to +4.00	-6.00	NVB Measurement	
Varilux® XR Design™	GR	BR	GG	S	AME	AMB	E	R					-12.00 to +9.00	+0.75 to +4.00	-6.00	14	
Varilux® X 4D™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R					-12.00 to +9.00	+0.75 to +4.00	-6.00	14	
Varilux® X FIT™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R					-12.00 to +9.00	+0.75 to +4.00	-6.00	14	
Varilux® X Design™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R					-12.00 to +9.00	+0.75 to +4.00	-6.00	14	
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R					-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14	
Varilux® Physio® W3+ FIT Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R					-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14	
Varilux® Physio® W3+ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R					-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14	
Varilux Comfort® Max Fit Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R					-14.00 to +9.00	+0.75 to +3.50	-6.00	14	
Varilux Comfort® Max Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R					-14.00 to +9.00	+0.75 to +3.50	-6.00	14	
Varilux Physio DRx™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R				GR	-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	17	
Varilux Physio Short DRx™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R				GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	14	
Varilux Comfort DRx™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R				GR	-14.00 to +9.00	+0.75 to +3.50	-6.00	17	
Varilux Comfort Short DRx™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R				GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	14	
Varilux® Liberty™ 3.0	GR	BR											-14.00 to +9.00	+0.75 to +1.50	-6.00	14	
Varilux Comfort® Thin&Lite® 1.67	GR	BR											-14.00 to +9.00	+0.75 to +3.00	-6.00	14	
Varilux Comfort® Short™ Thin&Lite® 1.67	GR	BR											-14.00 to +9.00	+0.75 to +3.50	-6.00	17	
Varilux® Stylistic™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R				GR	-5.00 to +5.00	+0.75 to +3.50	-3.00	17	
Essilor Ideal Advanced™ FIT Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R					-14.00 to +9.00	+0.75 to +3.00	-6.00	18	

Essilor Ideal Advanced Wrap™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Sport Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Driving Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor Ideal® Short	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor Natural® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-14.00 to +9.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Brillenz® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Ovation® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Ovation® Thin&Lite® 1.67	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-4.00	17
Essilor SmallFit® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG		-14.00 to +9.00	+0.75 to +3.00	-6.00	14

	Transitions Gen#S										Transitions XTRACTIVE	Transitions XTRACTIVE POLARIZED				
Progressives													Sphere	Add	Cyl to	MFH
1.74 Ultra High Index																
Varilux® XR Track™	GR	BR											-14.00 to +13.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR											-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ Thin&Lite® 1.74	GR	BR											-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ Thin&Lite® 1.74	GR	BR											-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ Thin&Lite® 1.74	GR	BR											-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.74	GR	BR											-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Thin&Lite® 1.74	GR	BR											-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Thin&Lite® 1.74	GR	BR											-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio DRx™ Thin&Lite® 1.74	GR	BR											-14.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Thin&Lite® 1.74	GR	BR											-14.00 to +13.00	+0.75 to +3.50	-6.00	14
Varilux Comfort DRx™ Thin&Lite® 1.74	GR	BR											-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Thin&Lite® 1.74	GR	BR											-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ FIT Thin&Lite® 1.74	GR	BR											-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Thin&Lite® 1.74	GR	BR											-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.74	GR	BR											-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport Thin&Lite® 1.74	GR	BR											-20.00 to +13.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Computer Thin&Lite® 1.74	GR	BR											-20.00 to +13.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Driving Thin&Lite® 1.74	GR	BR											-20.00 to +13.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Thin&Lite® 1.74	GR	BR											-20.00 to +13.00	+0.75 to +3.50	-6.00	13
Ovation® Digital Thin&Lite® 1.74	GR	BR											-20.00 to +13.00	+0.75 to +3.50	-6.00	14
Essilor Accolade Freedom™ Thin&Lite® 1.74	GR	BR											-20.00 to +13.00	+0.75 to +3.50	-6.00	15

