









Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Suppliers](#) / [Unity Lenses](#)

Unity Lenses – Entire Product Line

 Semi-Finished Single Vision	 Transitions Gen#S								 Transitions XTRACTIVE			 Transitions XTRACTIVE POLARIZED	Sphere	Cyl to	
	GR	BR	GG	S	AME	AMB	E	R	GR						
1.50 Standard Index															
Unity SVx	GR	BR	GG	S	AME	AMB	E	R	GR			GR	-12.25 to +6.00	-6.00	
Unity SVxtra	GR	BR	GG	S	AME	AMB	E	R	GR			GR	-12.25 to +6.00	-6.00	
Unity SVxtreme	GR	BR	GG	S	AME	AMB	E	R	GR			GR	-6.00 to +6.00	-6.00	
Unity Relieve 50	GR	BR	GG	S	AME	AMB	E	R					-12.25 to +6.00	-6.00	
Unity Relieve 70	GR	BR	GG	S	AME	AMB	E	R					-12.25 to +6.00	-6.00	
1.53 TRIVEX® Single Vision															
Unity SVx	GR	BR	GG	S	AME	AMB	E	R	GR				-14.00 to +6.00	up to -6.00	
Unity SVxtra	GR	BR	GG	S	AME	AMB	E	R	GR				-14.00 to +6.00	up to -6.00	
Unity SVxtreme	GR	BR	GG	S	AME	AMB	E	R	GR				-6.50 to +6.00	up to -6.00	
Unity Relieve 50	GR	BR	GG	S	AME	AMB	E	R					-14.00 to +6.00	up to -6.00	
Unity Relieve 70	GR	BR	GG	S	AME	AMB	E	R					-14.00 to +6.00	up to -6.00	
1.59 Polycarbonate Single Vision															
Unity SVx	GR	BR	GG	S	AME	AMB	E	R	GR			GR	-14.00 to +6.00	up to -6.00	
Unity SVxtra	GR	BR	GG	S	AME	AMB	E	R	GR			GR	-14.00 to +6.00	up to -6.00	
Unity SVxtreme	GR	BR	GG	S	AME	AMB	E	R	GR			GR	-6.50 to +6.00	up to -6.00	
Unity Relieve 50	GR	BR	GG	S	AME	AMB	E	R					-14.00 to +6.00	up to -6.00	
Unity Relieve 70	GR	BR	GG	S	AME	AMB	E	R					-14.00 to +6.00	up to -6.00	
1.60 High Index Single Vision															
Unity SVx	GR	BR											-14.00 to +6.00	up to -6.00	
Unity SVxtra	GR	BR											-14.00 to +6.00	up to -6.00	
Unity SVxtreme	GR	BR											-6.50 to +6.00	up to -6.00	
1.67 Super High Index Single Vision															
Unity SVx	GR	BR	GG	S	AME	AMB	E	R	GR			GR	-16.00 to +8.00	up to -6.00	
Unity SVxtra	GR	BR	GG	S	AME	AMB	E	R	GR			GR	-16.00 to +8.00	up to -6.00	
Unity SVxtreme	GR	BR	GG	S	AME	AMB	E	R	GR			GR	-6.00 to +6.00	up to -6.00	
Unity Relieve 50	GR	BR	GG	S	AME	AMB	E	R					-16.00 to +8.00	up to -6.00	
Unity Relieve 70	GR	BR	GG	S	AME	AMB	E	R					-16.00 to +8.00	up to -6.00	

 Progressives	 Transitions Gen#S								 Transitions XTRACTIVE			 Transitions XTRACTIVE POLARIZED	Sphere	Add	Cyl to	MFH
	GR	BR	GG	S	AME	AMB	E	R	GR							
1.50 Standard Index																
Unity Via II	GR	BR	GG	S	AME	AMB	E	R	GR			GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Plus II	GR	BR	GG	S	AME	AMB	E	R	GR			GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite VR	GR	BR	GG	S	AME	AMB	E	R	GR			GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II	GR	BR	GG	S	AME	AMB	E	R	GR			GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II 12	GR	BR	GG	S	AME	AMB	E	R	GR			GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II 14	GR	BR	GG	S	AME	AMB	E	R	GR			GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	14 mm
Unity Via Elite II 16	GR	BR	GG	S	AME	AMB	E	R	GR			GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	16 mm
Unity Via Elite II 18	GR	BR	GG	S	AME	AMB	E	R	GR			GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	18 mm
Unity Via Mobile II	GR	BR	GG	S	AME	AMB	E	R	GR			GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Wrap II	GR	BR	GG	S	AME	AMB	E	R	GR			GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Ethos	GR	BR	GG	S	AME	AMB	E	R	GR			GR	-9.00 to +4.50	+0.75 to +3.50	up to -4.00	14mm

unity	Transitions Gen#S								Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED		Sphere	Add	Cyl to	MFH	
	Near Variable Focus	GR	BR	GG	S	AME	AMB	E	R									
1.50 Standard Index																		
Unity OfficePro 5	GR	BR	GG	S	AME	AMB	E	R							-7.00 to +6.00	+1.00 to +3.50	Up to -6	14 mm
Unity Office Pro 10	GR	BR	GG	S	AME	AMB	E	R							-7.00 to +6.00	+1.00 to +3.50	Up to -6	14 mm
1.53 TRIVEX®																		
Unity OfficePro 5	GR	BR	GG	S	AME	AMB	E	R							-10.00 to +6.00	+1.00 to +3.50	Up to -6	14 mm
Unity Office Pro 10	GR	BR	GG	S	AME	AMB	E	R							-10.00 to +6.00	+1.00 to +3.50	Up to -6	14 mm
1.59 Polycarbonate																		
Unity OfficePro 5	GR	BR	GG	S	AME	AMB	E	R							-10.00 to +6.50	+1.00 to +3.50	Up to -6	14 mm
Unity Office Pro 10	GR	BR	GG	S	AME	AMB	E	R							-10.00 to +6.50	+1.00 to +3.50	Up to -6	14 mm
1.67 Super High Index																		
Unity OfficePro 5	GR	BR	GG	S	AME	AMB	E	R							-12.00 to +8.25	+1.00 to +3.50	Up to -6	14 mm
Unity Office Pro 10	GR	BR	GG	S	AME	AMB	E	R							-12.00 to +8.25	+1.00 to +3.50	Up to -6	14 mm

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.