


Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Design](#) / [Progressives](#) / [Essilor Lenses](#) / **1.60 Index**

Essilor Lenses — Progressives 1.60 Index

																
Progressives													Sphere	Add	Cyl to	MFH
1.60 High Index																
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Thin&Lite® 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Thin&Lite® 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit Thin&Lite® 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max Thin&Lite® 1.60	GR	BR	GG										-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™ Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux Comfort® Thin&Lite® 1.60	GR												-12.00 to +8.00	+0.75 to +3.50	-4.00	17
Varilux Comfort Short™ Thin&Lite® 1.60	GR												-12.00 to +8.00	+0.75 to +3.00	-4.00	14
Essilor Ideal Advanced™ FIT Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Sport Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Driving Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	13
Essilor Ideal® Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Short Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Essilor Natural® Digital Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Ovation® Digital Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.50	-6.00	17
SmallFit® Digital Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Thin&Lite® 1.60	GR	BR											-12.00 to +8.00	+0.75 to +3.00	-6.00	14

TRANSITIONS® LIGHT INTELLIGENT LENSES™



GR
Gray



BR
Brown



GG
Graphite
Green



S
Sapphire



AME
Amethyst



AMB
Amber



E
Emerald



R
Ruby



Available
in style
mirrors

Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

Unless indicated otherwise, all registered trademarks and trademarks are the property of Essilor International and/or its subsidiaries in the United States and in other countries.

All other trademarks are the property of their respective owners.