

Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Lens Design](#) / [SFSV](#) / [Unity Lenses](#)

Unity Lenses – Semi-Finished Single Vision

unity		Transitions Gen#S								Transitions XTRACTIVE			Transitions XTRACTIVE POLARIZED			
Semi-Finished Single Vision		GR	BR	GG	S	AME	AMB	E	R	GR	BR	GG	GR	BR	Sphere	Cyl to
1.50 Standard Index																
	Unity SVx	GR	BR							GR			GR	-12.25 to +6.00	-6.00	
	Unity SVxtra	GR	BR							GR			GR	-12.25 to +6.00	-6.00	
	Unity SVxtreme	GR	BR							GR			GR	-6.00 to +6.00	-6.00	
1.53 TRIVEX® Single Vision																
	Unity SVx	GR	BR							GR			GR	-14.00 to +6.00	-6.00	
	Unity SVxtra	GR	BR							GR			GR	-14.00 to +6.00	-6.00	
	Unity SVxtreme	GR	BR							GR			GR	-6.50 to +6.00	-6.00	
1.59 Polycarbonate Single Vision																
	Unity SVx	GR	BR							GR			GR	-14.00 to +6.00	-6.00	
	Unity SVxtra	GR	BR							GR			GR	-14.00 to +6.00	-6.00	
	Unity SVxtreme	GR	BR							GR			GR	-6.50 to +6.00	-6.00	
1.60 High Index Single Vision																
	Unity SVx	GR	BR											-14.00 to +6.00	-6.00	
	Unity SVxtra	GR	BR											-14.00 to +6.00	-6.00	
	Unity SVxtreme	GR	BR											-6.50 to +6.00	-6.00	
1.67 Super High Index Single Vision																
	Unity SVx	GR	BR							GR			GR	-16.00 to +8.00	-6.00	
	Unity SVxtra	GR	BR							GR			GR	-16.00 to +8.00	-6.00	
	Unity SVxtreme	GR	BR							GR			GR	-6.00 to +6.00	-6.00	

TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.