

Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Transitions® Signature® GEN 8™](#) / **Essilor Lenses**

Essilor Lenses – Transitions® Signature® GEN 8™

									
Finished Single Vision								Sphere	Cyl to
1.50 Standard Index Finished Single Vision									
Essilor FSV Spherical	GR	BR						-4.00 to +4.00	-2.00
Essilor FSV Aspheric	GR							-4.00 to +4.00	-2.00
Essilor FSV 1.50 Crizal® Sapphire™ HR	GR							-4.00 to +4.00	-2.00
Essilor FSV 1.50 Crizal® Rock™	GR							-4.00 to +4.00	-2.00
Essilor FSV 1.50 Crizal® Easy Pro	GR							-4.00 to +4.00	-2.00
1.59 Polycarbonate Finished Single Vision									
Essilor FSV Spherical Airwear®	GR	BR						-6.00 to +4.00	-2.00
Essilor FSV Aspheric Airwear®	GR	BR						-6.00 to +4.00	-2.00
Essilor FSV Airwear® Crizal® Sapphire™ HR	GR							-6.00 to +4.00	-2.00
Essilor FSV Airwear® Crizal® Rock™	GR							-6.00 to +4.00	-2.00
Essilor FSV Airwear® Crizal® Easy Pro	GR							-6.00 to +4.00	-2.00
1.67 Super High Index Finished Single Vision									
Essilor FSV Spherical Thin&Lite® 1.67	GR	BR						-6.00 to +4.00	-2.00

									
Semi-Finished Single Vision								Sphere	Cyl to
1.50 Standard Index Single Vision									
Essilor Spherical Single Vision	GR	BR	GG	S	AME	AMB	E	-10.50 to +10.00	-4.00
Essilor Stylistic® Single Vision	GR	BR	GG	S	AME	AMB	E	-10.00 to +7.50	-6.00
Eyezen® +0	GR	BR	GG	S	AME	AMB	E	-10.00 to +7.50	-6.00
Eyezen® +1	GR	BR	GG	S	AME	AMB	E	-10.00 to +7.50	-6.00
Eyezen® +2	GR	BR	GG	S	AME	AMB	E	-10.00 to +7.50	-6.00
Eyezen® +3	GR	BR	GG	S	AME	AMB	E	-10.00 to +7.50	-6.00
Eyezen® +4	GR	BR	GG	S	AME	AMB	E	-10.00 to +7.50	-6.00
Eyezen® Start	GR	BR	GG	S	AME	AMB	E	-10.00 to +6.00	-6.00
1.53 TRIVEX® Single Vision									
Essilor Spherical Single Vision TREXA™	GR	BR	GG					-10.00 to +6.00	-4.00
Essilor Eyecode™ Single Vision TREXA™	GR	BR	GG					-10.00 to +7.50	-6.00
Essilor FIT™ Single Vision TREXA™	GR	BR	GG					-10.00 to +7.50	-6.00
Essilor Azio™ Single Vision TREXA™	GR	BR	GG					-10.00 to +7.50	-6.00
Essilor 360°™ Single Vision TREXA™	GR	BR	GG					-10.00 to +7.50	-6.00
Eyezen® +0 TREXA™	GR	BR						-10.00 to +7.50	-6.00
Eyezen® +1 TREXA™	GR	BR						-10.00 to +7.50	-6.00
Eyezen® +2 TREXA™	GR	BR						-10.00 to +7.50	-6.00
Eyezen® +3 TREXA™	GR	BR						-10.00 to +7.50	-6.00
Eyezen® +4 TREXA™	GR	BR						-10.00 to +7.50	-6.00
Eyezen® Start	GR	BR						-10.00 to +6.00	-6.00
1.59 Polycarbonate Single Vision									
Essilor Spherical Single Vision Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	-4.00
Essilor Aspheric Single Vision Airwear®	GR	BR	GG					-10.00 to +8.00	-4.00
Essilor Eyecode™ Aspheric SV Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	-6.00
Essilor FIT™ Aspheric SV Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	-6.00

Essilor AZIO™ Aspheric SV Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	-6.00
Essilor SV 360°™ Aspheric Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	-6.00
Essilor Stylistic® Single Vision Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	-6.00
Eyezen® +0 Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	-6.00
Eyezen® +1 Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	-6.00
Eyezen® +2 Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	-6.00
Eyezen® +3 Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	-6.00
Eyezen® +4 Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	-6.00
Eyezen® Start Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	-6.00
Eyezen® Kids	GR	BR	GG					-10.00 to +6.00	-6.00
1.60 High Index Single Vision									
Essilor Spherical Single Vision Thin&Lite® 1.60	GR	BR						-13.00 to +8.00	-4.00
Eyezen® +0 Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	-6.00
Eyezen® +1 Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	-6.00
Eyezen® +2 Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	-6.00
Eyezen® +3 Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	-6.00
Eyezen® +4 Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	-6.00
Eyezen® Start Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	-6.00
1.67 Super High Index Single Vision									
Essilor Spherical Single Vision Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	-4.00
Essilor Aspheric Single Vision Thin&Lite® 1.67	GR	BR	GG					-16.00 to +8.00	-4.00
Essilor Eyecode™ SV Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	-6.00
Essilor FIT™ SV Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	-6.00
Essilor Azio™ SV Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	-6.00
Essilor SV 360°™ Spherical Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	-6.00
Essilor Stylistic® Single Vision Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	-6.00
Eyezen® +0 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	-6.00
Eyezen® +1 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	-6.00
Eyezen® +2 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	-6.00
Eyezen® +3 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	-6.00
Eyezen® +4 Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	-6.00
Eyezen® Start Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	-6.00
1.74 Ultra High Index Single Vision									
Essilor Spherical Single Vision Thin&Lite® 1.74	GR							-20.00 to +13.00	-4.00
Essilor Eyecode™ SV Thin&Lite® 1.74	GR							-16.00 to +13.00	-6.00
Essilor FIT™ SV Thin&Lite® 1.74	GR							-16.00 to +13.00	-6.00
Essilor 360°™ SV Thin&Lite® 1.74	GR							-16.00 to +13.00	-6.00
Eyezen® +0 Thin&Lite® 1.74	GR							-20.00 to +13.00	-6.00
Eyezen® +1 Thin&Lite® 1.74	GR							-20.00 to +13.00	-6.00
Eyezen® +2 Thin&Lite® 1.74	GR							-20.00 to +13.00	-6.00
Eyezen® +3 Thin&Lite® 1.74	GR							-20.00 to +13.00	-6.00
Eyezen® +4 Thin&Lite® 1.74	GR							-20.00 to +13.00	-6.00
Eyezen® Start Thin&Lite® 1.74	GR							-20.00 to +13.00	-6.00

											
Bifocal Signature™									Sphere	Add	Cyl to
1.50 Standard Index SFST-28											
Bifocal FT28		GR	BR	GG					-10.00 to +6.00	+1.00 to +3.00	-4.00
Essilor Bifocal AB™ FT28		GR	BR	GG					-10.00 to +6.00	+1.00 to +3.00	-4.00
1.59 index polycarbonate composite SFFT 28											
Bifocal FT28		GR	BR						-10.00 to +6.00	+1.00 to +3.00	-4.00
Essilor Bifocal AB™ FT28		GR	BR						-10.00 to +6.00	+1.00 to +3.00	-4.00
1.59 index polycarbonate composite SFFT 35											
Bifocal FT35		GR	BR						-10.00 to +6.00	+1.00 to +3.00	-4.00

Varilux® XR Track™	GR	BR						-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR						-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ TREXA™	GR	BR						-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ TREXA™	GR	BR						-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ TREXA™	GR	BR						-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode TREXA™	GR	BR						-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT TREXA™	GR	BR						-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ TREXA™	GR	BR						-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit TREXA™	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max TREXA™	GR	BR						-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™ TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Essilor Ideal Advanced™ FIT TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50	-6.00	18
Essilor Ideal® Computer TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50	-6.00	13
Essilor Ideal® Driving TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal® TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal® Short TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Natural® Digital TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50	-6.00	15
Essilor Accolade® TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Digital TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Brillenz® TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Brillenz® Short TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Ovation® Digital TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor SmallFit® Digital TREXA™	GR	BR	GG					-10.00 to +7.50	+0.75 to +3.00	-6.00	14

												
Progressives									Sphere	Add	Cyl to	MFH
1.59 Polycarbonate												
Varilux® XR Track™		GR	BR	GG	S	AME	AMB	E	-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™		GR	BR	GG	S	AME	AMB	E	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ Airwear®		GR	BR	GG	S	AME	AMB	E	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ Airwear®		GR	BR	GG	S	AME	AMB	E	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ Airwear®		GR	BR	GG	S	AME	AMB	E	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode™ Airwear®		GR	BR	GG	S	AME	AMB	E	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Airwear®		GR	BR	GG	S	AME	AMB	E	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Airwear®		GR	BR	GG	S	AME	AMB	E	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit Airwear®		GR	BR	GG	S	AME	AMB	E	-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max Airwear®		GR	BR	GG	S	AME	AMB	E	-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Airwear®		GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Airwear®		GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Varilux Comfort DRx™ Airwear®		GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Airwear®		GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux® Liberty™ 3.0		GR	BR						-10.00 to +8.00	+0.75 to +1.50	-6.00	14
Varilux Comfort® Airwear®		GR	BR	GG					-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short™ Airwear®		GR	BR	GG					-10.00 to +6.00	+0.75 to +3.00	-6.00	14

Varilux® Stylistic™ Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ FIT	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Wrap	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Essilor Ideal® FIT Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Sport Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Driving Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Ideal® Short Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Natural® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Ovation® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Ovation® Airwear®	GR	BR						-10.00 to +8.00	+0.75 to +3.50	-4.00	17
Essilor SmallFit® Digital Airwear®	GR	BR	GG	S	AME	AMB	E	-10.00 to +8.00	+0.75 to +3.00	-6.00	14

											
Progressives								Sphere	Add	Cyl to	MFH
1.60 High Index											
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.60	GR	BR	GG					-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Thin&Lite® 1.60	GR	BR	GG					-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Thin&Lite® 1.60	GR	BR	GG					-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit Thin&Lite® 1.60	GR	BR	GG					-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max Thin&Lite® 1.60	GR	BR	GG					-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™ Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux Comfort® Thin&Lite® 1.60	GR							-12.00 to +8.00	+0.75 to +3.50	-4.00	17
Varilux Comfort Short™ Thin&Lite® 1.60	GR							-12.00 to +8.00	+0.75 to +3.00	-4.00	14
Essilor Ideal Advanced™ FIT Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Sport Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Driving Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	13
Essilor Ideal® Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Short Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Essilor Natural® Digital Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	14
Ovation® Digital Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-6.00	17
SmallFit® Digital Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Thin&Lite® 1.60	GR	BR						-12.00 to +8.00	+0.75 to +3.00	-6.00	14

											
---	---	--	--	--	--	--	--	--	--	--	--

Progressives								Sphere	Add	Cyl to	MFH
1.67 Super High Index											
Varilux® XR Track™	GR	BR	GG	S	AME	AMB	E	-12.00 to +9.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR	GG	S	AME	AMB	E	-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Varilux® Liberty™ 3.0	GR	BR						-14.00 to +9.00	+0.75 to +1.50	-6.00	14
Varilux Comfort® Thin&Lite® 1.67	GR	BR						-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Varilux Comfort® Short™ Thin&Lite® 1.67	GR	BR						-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Varilux® Stylistic™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-5.00 to +5.00	+0.75 to +3.50	-3.00	17
Essilor Ideal Advanced™ FIT Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal Advanced Wrap™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Sport Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Driving Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor Ideal® Short	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor Natural® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Brillenz® Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Ovation® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Ovation® Thin&Lite® 1.67	GR	BR						-12.00 to +8.00	+0.75 to +3.50	-4.00	17
Essilor SmallFit® Digital Thin&Lite® 1.67	GR	BR	GG	S	AME	AMB	E	-14.00 to +9.00	+0.75 to +3.00	-6.00	14

											
Progressives								Sphere	Add	Cyl to	MFH
1.74 Ultra High Index											
Varilux® XR Track™	GR	BR						-14.00 to +13.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR						-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ Thin&Lite® 1.74	GR	BR						-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ Thin&Lite® 1.74	GR	BR						-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ Thin&Lite® 1.74	GR	BR						-14.00 to +13.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.74	GR	BR						-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Thin&Lite® 1.74	GR	BR						-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Thin&Lite® 1.74	GR	BR						-18.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio DRx™ Thin&Lite® 1.74	GR	BR						-14.00 to +13.00	+0.75 to +3.50 & +4.00	-6.00	17

Varilux Physio Short DRx™ Thin&Lite® 1.74	GR	BR							-14.00 to +13.00	+0.75 to +3.50	-6.00	14
Varilux Comfort DRx™ Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ FIT Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Computer Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Driving Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	+0.75 to +3.50	-6.00	13
Ovation® Digital Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	+0.75 to +3.50	-6.00	14
Essilor Accolade Freedom™ Thin&Lite® 1.74	GR	BR							-20.00 to +13.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Thin&Lite® 1.74	GR	BR							-14.00 to +9.00	+0.75 to +3.50	-6.00	17

TRANSITIONS® LIGHT INTELLIGENT LENSES™



GR
Gray



BR
Brown



GG
Graphite
Green



S
Sapphire



AME
Amethyst



AMB
Amber



E
Emerald



R
Ruby



Available
in style
mirrors

Transitions, *Transitions Signature*, and *XTRActive* are registered trademarks, and *XTRActive Polarized*, the *Transitions* logo and *Transitions Light Intelligent Lenses* are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

Unless indicated otherwise, all registered trademarks and trademarks are the property of Essilor International and/or its subsidiaries in the United States and in other countries.

All other trademarks are the property of their respective owners.