

## Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Transitions® Signature® GEN 8™](#) / **Signet Armorlite**

### Signet Armorlite – Transitions® Signature® GEN 8™

									
Bifocal Signature™							Sphere	Add	Cyl to
<b>1.50 Standard Index SFST-28</b>									
Signet Armorlite SFFT-28	GR						-10.00 to +6.00	+1.00 to +3.00	-4.00
Signetek™ Processed 1.50 FT-28	GR						-10.00 to +6.00	+1.00 to +3.00	-6.00
<b>1.59 Polycarbonate SFST-28</b>									
Signetek™ Processed 1.59 FT-28	GR BR						-8.25 to +6.00	+1.00 to +3.00	-6.00
<b>1.59 Polycarbonate SFST-35</b>									
Signetek™ Processed 1.59 FT-35	GR						-8.25 to +6.00	+1.00 to +3.00	-6.00

									
Trifocal Signature™							Sphere	Add	Cyl to
<b>1.50 Standard Index SF Trifocal 7x28</b>									
Signet Armorlite SF Trifocal 7x28	GR						-5.00 to +6.00	+1.50 to +3.50	-4.00
Signetek™ Processed Trifocal 1.50 7x28	GR						-5.00 to +6.00	+1.50 to +3.50	-6.00
<b>1.59 Polycarbonate SF Trifocal 7x28</b>									
Signetek™ Processed Trifocal 1.59 7x28	GR						-6.00 to +6.00	+1.50 to +3.50	-6.00

										
Progressives							Sphere	Add	Cyl to	MFH
<b>1.50 Standard Index</b>										
KODAK Precise® 1.50 SF	GR						-9.00 to +5.00	+1.00 to +3.00	-4.00	17
KODAK Precise® Digital 1.50	GR						-8.00 to +5.00	+1.00 to +3.00	-6.00	17
<b>1.59 Polycarbonate</b>										
KODAK Precise® Poly SF	GR						-10.00 to +5.00	+1.00 to +3.00	-4.00	17
KODAK Precise® Digital Poly	GR						-10.00 to +6.00	+1.00 to +3.00	-6.00	17

#### TRANSITIONS® LIGHT INTELLIGENT LENSES™



Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.