









Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / **Transitions® XTRActive® style mirrors**

All Suppliers — Transitions® XTRActive® style mirrors

	Transitions® XTRACTIVE™				
Semi-Finished Single Vision				Sphere	Cyl to
1.59 Polycarbonate Single Vision					
Essilor Spherical Single Vision Airwear®	GR	BR	GG	-10.00 to +8.00	-4.00
Essilor Aspheric Single Vision Airwear®	GR	BR	GG	-10.00 to +8.00	-4.00
Essilor Eyecode™ Aspheric SV Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Essilor FIT™ Aspheric SV Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Essilor AZIO™ Aspheric SV Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Essilor SV 360°™ Aspheric Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Essilor Stylistic® Single Vision Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Eyezen® +0 Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Eyezen® +1 Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Eyezen® +2 Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Eyezen® +3 Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Eyezen® +4 Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Eyezen® Start Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
1.67 Super High Index Single Vision					
Essilor SFSV Spherical Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-4.00
Essilor Eyecode™ SV Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Essilor FIT™ SV Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Essilor Azio™ SV Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Essilor SV 360°™ Spherical Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Essilor Stylistic® Single Vision Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Eyezen® +0 Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Eyezen® +1 Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Eyezen® +2 Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Eyezen® +3 Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Eyezen® +4 Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Eyezen® Start Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00

	Transitions® XTRACTIVE™						
Progressives				Sphere	Add	Cyl to	MFH
1.59 Polycarbonate							
Varilux® XR Track™	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode™ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14

Varilux® Physio® W3+ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Varilux Comfort DRx™ Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux Comfort® Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short™ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.00	-6.00	14
Varilux® Stylistic™ Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Wrap	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Essilor Ideal® FIT Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Sport Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Driving Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Ideal® Short Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Natural® Digital Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Ovation® Digital Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor SmallFit® Digital Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
1.67 Super High Index							
Varilux® XR Track™	GR	BR	GG	-12.00 to +9.00	+0.75 to +4.00	-6.00	NVB Measurement
Varilux® XR Design™	GR	BR	GG	-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® X 4D™ Thin&Lite® 1.67	GR	BR	GG	-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ Thin&Lite® 1.67	GR	BR	GG	-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ Thin&Lite® 1.67	GR	BR	GG	-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Varilux® Stylistic™ Thin&Lite® 1.67	GR	BR	GG	-5.00 to +5.00	+0.75 to +3.50	-3.00	17
Essilor Ideal Advanced Wrap™ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Sport Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Computer Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	18

Essilor Ideal® Driving Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor Ideal® Short	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor Natural® Digital Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Brillenz® Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Ovation® Digital Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor SmallFit® Digital Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14

TRANSITIONS® LIGHT INTELLIGENT LENSES™



GR
Gray



BR
Brown



GG
Graphite
Green



S
Sapphire



AME
Amethyst



AMB
Amber



E
Emerald



R
Ruby



Available
in style
mirrors

Style Mirrors are available where gray and brown *Transitions® XTRActive®* are available. Specify *Transitions* lenses in style mirrors **(no substitutions)** with your lab to ensure authenticity.

Transitions, *Transitions Signature*, and *XTRActive* are registered trademarks, and *XTRActive Polarized*, the *Transitions* logo and *Transitions Light Intelligent Lenses* are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

Unless indicated otherwise, all registered trademarks and trademarks are the property of Essilor International and/or its subsidiaries in the United States and in other countries.

All other trademarks are the property of their respective owners.