









## Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Transitions® XTRActive®](#) / **Signet Armorlite**

### Signet Armorlite – Transitions® XTRActive®

		Transitions® XTRACTIVE				
Semi-Finished Single Vision					<b>Sphere</b>	<b>Cyl to</b>
<b>1.50 Standard Index</b>						
	Signetek Processed SV	GR			-8.00 to +6.00	-6.00
	KODAK Digital SV	GR			-8.00 to +6.00	-6.00
<b>1.53 TRIVEX® Single Vision</b>						
	Signetek Processed SV	GR			-10.00 to +6.00	-6.00
	KODAK Digital SV	GR			-10.00 to +6.00	-6.00
<b>1.59 Polycarbonate Single Vision</b>						
	Signetek Processed SV	GR			-10.00 to +6.50	-6.00
	KODAK Digital SV	GR			-10.00 to +6.50	-6.00
<b>1.67 Super High Index Single Vision</b>						
	Signetek Processed SV	GR			-11.00 to +7.50	-6.00
	KODAK Digital SV	GR			-11.00 to +7.50	-6.00

		Transitions® XTRACTIVE					
Bifocal Signature™					<b>Sphere</b>	<b>Add</b>	<b>Cyl to</b>
<b>1.50 Standard Index SFST-28</b>							
	Signetek™ Processed 1.50 FT-28	GR			-10.00 to +6.00	+1.00 to +3.00	-6.00

		Transitions® XTRACTIVE						
Progressives					<b>Sphere</b>	<b>Add</b>	<b>Cyl to</b>	<b>MFH</b>
<b>1.50 Standard Index</b>								
	KODAK Precise® Plus 1.50	GR			-8.00 to +4.00	+0.75 to +3.50	-6.00	17
	KODAK Precise® Short Plus 1.50	GR			-8.00 to +4.00	+0.75 to +3.50	-6.00	13
	KODAK Precise® PB 1.50	GR			-8.00 to +4.00	+0.75 to +3.50	-6.00	17
	KODAK Precise® Short PB 1.50	GR			-8.00 to +4.00	+0.75 to +3.50	-6.00	13
	KODAK Unique™ 1.50	GR			-8.00 to +4.00	+0.75 to +3.50	-6.00	13
	KODAK Unique DRO™ 1.50	GR			-8.00 to +4.00	+0.75 to +3.50	-6.00	13
	KODAK Unique DRO™ HD 1.50	GR			-8.00 to +4.00	+0.75 to +3.50	-6.00	13
	KODAK Unique Infinite™ 1.50	GR	BR		-8.00 to +4.00	+0.75 to +3.50	-6.00	13
	KODAK Unique Infinite™ HD 1.50	GR	BR		-8.00 to +4.00	+0.75 to +3.50	-6.00	13
<b>1.53 TRIVEX®</b>								
	KODAK Precise® PB Trivex®	GR			-10.00 to +4.50	+0.75 to +3.50	-6.00	17
	KODAK Precise® Short PB Trivex®	GR			-10.00 to +4.50	+0.75 to +3.50	-6.00	13
	KODAK Precise® Plus Trivex®	GR			-10.00 to +4.50	+0.75 to +3.50	-6.00	17
	KODAK Precise® Short Plus Trivex®	GR			-10.00 to +4.50	+0.75 to +3.50	-6.00	13
	KODAK Unique Trivex®	GR			-10.00 to +4.50	+0.75 to +3.50	-6.00	13
	KODAK Unique DRO™ Trivex®	GR			-10.00 to +4.50	+0.75 to +3.50	-6.00	13
	KODAK Unique DRO™ HD Trivex®	GR			-10.00 to +4.50	+0.75 to +3.50	-6.00	13
	KODAK Unique Infinite™ Trivex®	GR	BR		-10.00 to +4.50	+0.75 to +3.50	-6.00	13
	KODAK Unique Infinite™ HD Trivex®	GR	BR		-10.00 to +4.50	+0.75 to +3.50	-6.00	13
<b>1.59 Polycarbonate</b>								

KODAK Precise® Plus Poly	GR			-10.00 to +5.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus Poly	GR			-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ Poly	GR			-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Precise® PB Poly	GR			-10.00 to +5.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB Poly	GR			-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique Poly	GR			-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD Poly	GR			-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ Poly	GR	BR		-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ HD Poly	GR	BR		-10.00 to +5.00	+0.75 to +3.50	-6.00	13
<b>1.67 Super High Index</b>							
KODAK Precise® PB 1.67	GR			-11.00 to +6.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB 1.67	GR			-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Precise® Plus Poly	GR			-11.00 to +6.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus Poly	GR			-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique	GR			-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ 1.67	GR			-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.67	GR			-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ 1.67	GR	BR		-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique Infinite™ HD 1.67	GR	BR		-11.00 to +6.00	+0.75 to +3.50	-6.00	13

**TRANSITIONS® LIGHT INTELLIGENT LENSES™**



**GR**  
Gray



**BR**  
Brown



**GG**  
Graphite  
Green



**S**  
Sapphire



**AME**  
Amethyst



**AMB**  
Amber



**E**  
Emerald



**R**  
Ruby



Available  
in style  
mirrors

Style Mirrors are available where gray and brown *Transitions® XTRActive®* are available. Specify *Transitions* lenses in style mirrors (**no substitutions**) with your lab to ensure authenticity.

*Transitions*, *Transitions Signature*, and *XTRActive* are registered trademarks, and *XTRActive Polarized*, the *Transitions* logo and *Transitions Light Intelligent Lenses* are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.