






## Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Transitions® XTRActive®](#) / **Unity Lenses**

### Unity Lenses — Transitions® XTRActive®

			
Semi-Finished Single Vision	  	Sphere	Cyl to
<b>1.50 Standard Index</b>			
Unity SVx	GR	-12.25 to +6.00	-6.00
Unity SVxtra	GR	-12.25 to +6.00	-6.00
Unity SVxtreme	GR	-6.00 to +6.00	-6.00
<b>1.53 TRIVEX® Single Vision</b>			
Unity SVx	GR	-14.00 to +6.00	up to -6.00
Unity SVxtra	GR	-14.00 to +6.00	up to -6.00
Unity SVxtreme	GR	-6.50 to +6.00	up to -6.00
<b>1.59 Polycarbonate Single Vision</b>			
Unity SVx	GR	-14.00 to +6.00	up to -6.00
Unity SVxtra	GR	-14.00 to +6.00	up to -6.00
Unity SVxtreme	GR	-6.50 to +6.00	up to -6.00
<b>1.67 Super High Index Single Vision</b>			
Unity SVx	GR	-16.00 to +8.00	up to -6.00
Unity SVxtra	GR	-16.00 to +8.00	up to -6.00
Unity SVxtreme	GR	-6.00 to +6.00	up to -6.00

					
Progressives	  	Sphere	Add	Cyl to	MFH
<b>1.50 Standard Index</b>					
Unity Via II	GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Plus II	GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite VR	GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II	GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II 12	GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II 14	GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	14 mm
Unity Via Elite II 16	GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	16 mm
Unity Via Elite II 18	GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	18 mm
Unity Via Mobile II	GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Wrap II	GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Ethos	GR	-9.00 to +4.50	+0.75 to +3.50	up to -4.00	14 mm
Ethos Plus	GR	-9.00 to +4.50	+0.75 to +3.50	up to -4.00	14 mm
<b>1.53 TRIVEX®</b>					
Unity Via II	GR	-10.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Plus II	GR	-10.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite VR	GR	-10.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II	GR	-10.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II 12	GR	-10.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II 14	GR	-10.00 to +6.00	+0.50 to +4.50	up to -6.00	14 mm
Unity Via Elite II 16	GR	-10.00 to +6.00	+0.50 to +4.50	up to -6.00	16 mm
Unity Via Elite II 18	GR	-10.00 to +6.00	+0.50 to +4.50	up to -6.00	18 mm
Unity Via Mobile II	GR	-10.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Wrap II	GR	-5.00 to +4.00	+0.50 to +4.50	up to -6.00	12 mm

	Ethos	GR			-9.00 to +6.00	+0.75 to +3.50	up to -4.00	14 mm
	Ethos Plus	GR			-9.00 to +6.00	+0.75 to +3.50	up to -4.00	14 mm
<b>1.59 Polycarbonate</b>								
	Unity Via II	GR			-10.00 to +6.50	+0.50 to +4.50	up to -6.00	12 mm
	Unity Via Plus II	GR			-10.00 to +6.50	+0.50 to +4.50	up to -6.00	12 mm
	Unity Via Elite VR	GR			-10.00 to +6.50	+0.50 to +4.50	up to -6.00	12 mm
	Unity Via Elite II	GR			-10.00 to +6.50	+0.50 to +4.50	up to -6.00	12 mm
	Unity Via Elite II 12	GR			-10.00 to +6.50	+0.50 to +4.50	up to -6.00	12 mm
	Unity Via Elite II 14	GR			-10.00 to +6.50	+0.50 to +4.50	up to -6.00	14 mm
	Unity Via Elite II 16	GR			-10.00 to +6.50	+0.50 to +4.50	up to -6.00	16 mm
	Unity Via Elite II 18	GR			-10.00 to +6.50	+0.50 to +4.50	up to -6.00	18 mm
	Unity Via Mobile II	GR			-10.00 to +6.50	+0.50 to +4.50	up to -6.00	12 mm
	Unity Via Wrap II	GR			-6.50 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
	Ethos	GR			-9.00 to +6.00	+0.75 to +3.50	up to -4.00	14 mm
	Ethos Plus	GR			-9.00 to +6.00	+0.75 to +3.50	up to -4.00	14 mm
<b>1.67 Super High Index</b>								
	Unity Via II	GR			-12.00 to +8.25	+0.50 to +4.50	up to -6.00	12 mm
	Unity Via Plus II	GR			-12.00 to +8.25	+0.50 to +4.50	up to -6.00	12 mm
	Unity Via Elite VR	GR			-12.00 to +8.25	+0.50 to +4.50	up to -6.00	12 mm
	Unity Via Elite II	GR			-12.00 to +8.25	+0.50 to +4.50	up to -6.00	12 mm
	Unity Via Elite II 12	GR			-12.00 to +8.25	+0.50 to +4.50	up to -6.00	12 mm
	Unity Via Elite II 14	GR			-12.00 to +8.25	+0.50 to +4.50	up to -6.00	14 mm
	Unity Via Elite II 16	GR			-12.00 to +8.25	+0.50 to +4.50	up to -6.00	16 mm
	Unity Via Elite II 18	GR			-12.00 to +8.25	+0.50 to +4.50	up to -6.00	18 mm
	Unity Via Mobile II	GR			-12.00 to +8.25	+0.50 to +4.50	up to -6.00	12 mm
	Unity Via Wrap II	GR			-6.00 to +5.00	+0.50 to +4.50	up to -6.00	12 mm
	Ethos	GR			-12.00 to +8.00	+0.75 to +3.50	up to -4.00	14 mm
	Ethos Plus	GR			-12.00 to +8.00	+0.75 to +3.50	up to -4.00	14 mm

**TRANSITIONS® LIGHT INTELLIGENT LENSES™**



**GR**  
Gray



**BR**  
Brown



**GG**  
Graphite  
Green



**S**  
Sapphire



**AME**  
Amethyst



**AMB**  
Amber



**E**  
Emerald



**R**  
Ruby



Available  
in style  
mirrors

*Transitions*, *Transitions Signature*, and *XTRActive* are registered trademarks, and *XTRActive Polarized*, the *Transitions* logo and *Transitions Light Intelligent Lenses* are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.