






Transitions® Lenses Product Availability Guide






< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Transitions® XTRActive® new generation](#) / [Younger Optics](#)

All Suppliers – Transitions® XTRActive® new generation

					
Finished Single Vision				Sphere	Cyl to
1.50 Standard Index					
Essilor Lens FSV Spheric	GR			-4.00 to +4.00	-2.00
1.59 Polycarbonate Finished Single Vision					
Essilor FSV Aspheric Airwear®	GR			-6.00 to +4.00	-2.00

					
Semi-Finished Single Vision				Sphere	Cyl to
1.50 Standard Index					
Essilor Lens SFSV	GR	BR	GG	-10.50 to +10.00	-4.00
Essilor Stylistic® Single Vision	GR	BR	GG	-10.00 to +7.50	-6.00
Eyezen® +0	GR	BR	GG	-10.00 to +7.50	-6.00
Eyezen® +1	GR	BR	GG	-10.00 to +7.50	-6.00
Eyezen® +2	GR	BR	GG	-10.00 to +7.50	-6.00
Eyezen® +3	GR	BR	GG	-10.00 to +7.50	-6.00
Eyezen® +4	GR	BR	GG	-10.00 to +7.50	-6.00
Eyezen® Start	GR	BR	GG	-10.00 to +6.00	-6.00
1.53 TRIVEX® Single Vision					
Essilor Lens SFSV TREXA™	GR	BR		-10.00 to +6.00	-4.00
Essilor Eyecode™ Single Vision TREXA™	GR	BR		-10.00 to +7.50	-6.00
Essilor FIT™ Single Vision TREXA™	GR	BR		-10.00 to +7.50	-6.00
Essilor Azio™ Single Vision TREXA™	GR	BR		-10.00 to +7.50	-6.00
Essilor 360°™ Single Vision TREXA™	GR	BR		-10.00 to +7.50	-6.00
Essilor Stylistic® Single Vision	GR	BR		-10.00 to +7.50	-6.00
Eyezen® +0 TREXA™	GR	BR		-10.00 to +7.50	-6.00
Eyezen® +1 TREXA™	GR	BR		-10.00 to +7.50	-6.00
Eyezen® +2 TREXA™	GR	BR		-10.00 to +7.50	-6.00
Eyezen® +3 TREXA™	GR	BR		-10.00 to +7.50	-6.00

Eyezen® +4 TREXA™	GR	BR		-10.00 to +7.50	-6.00
Eyezen® Start	GR	BR		-10.00 to +6.00	-6.00
1.59 Polycarbonate Single Vision					
Essilor SFSV Spherical Airwear®	GR	BR	GG	-10.00 to +8.00	-4.00
Essilor SFSV Aspheric Airwear®	GR	BR	GG	-10.00 to +8.00	-4.00
Essilor Eyecode™ Aspheric SV Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Essilor FIT™ Aspheric SV Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Essilor AZIO™ Aspheric SV Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Essilor SV 360°™ Aspheric Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Essilor Stylistic® Single Vision Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Eyezen® +0 Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Eyezen® +1 Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Eyezen® +2 Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Eyezen® +3 Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Eyezen® +4 Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
Eyezen® Start Airwear®	GR	BR	GG	-10.00 to +8.00	-6.00
1.67 Super High Index Single Vision					
Essilor SFSV Spherical Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-4.00
Essilor Eyecode™ SV Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Essilor FIT™ SV Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Essilor Azio™ SV Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Essilor SV 360°™ Spherical Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Essilor Stylistic® Single Vision Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Eyezen® +0 Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Eyezen® +1 Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Eyezen® +2 Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Eyezen® +3 Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Eyezen® +4 Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00
Eyezen® Start Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	-6.00






							
Progressives				Sphere	Add	Cyl to	MFH
1.50 Standard Index							

Varilux® X 4D™	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® W2+ Eyecode	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® W2+ FIT	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® W2+	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Varilux Comfort®	GR	BR	GG	-10.00 to +7.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short™	GR	BR	GG	-10.00 to +7.00	+0.75 to +3.00	-6.00	14
Varilux® Stylistic™	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Wrap	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-4.00	17
Essilor Ideal Advanced™ FIT	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal®	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	18
Essilor Ideal® Computer	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	13
Essilor Ideal® Driving	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal® Short	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Natural® Digital	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	15
Essilor Accolade™	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Digital	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Brillenz®	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Brillenz® Short	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Ovation® Digital	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor SmallFit® Digital	GR	BR	GG	-10.00 to +7.50	+0.75 to +3.00	-6.00	14
1.53 TRIVEX®							
Varilux® Physio® W3+ Eyecode TREXA™	GR	BR		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT TREXA™	GR	BR		-10.00 to +6.00	+0.75 to +3.50 &	-6.00	14



					+4.00		
Varilux® Physio® W3+ TREXA™	GR	BR		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit TREXA™	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max TREXA™	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® W2+ Eyecode	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® W2+ FIT	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® W2+	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™ TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Essilor Ideal Advanced™ TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Wrap TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ FIT TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal® FIT TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Ideal® TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	18
Essilor Ideal® Sport TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	13
Essilor Ideal® Computer TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal® Driving TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Ideal® Short	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Essilor Natural® Digital TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	15
Essilor Accolade® TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Digital TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.00	-6.00	14
Brillenz® TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Brillenz® Short TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	14
Ovation® Digital TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.50	-6.00	17
Essilor SmallFit® Digital TREXA™	GR	BR		-10.00 to +7.50	+0.75 to +3.00	-6.00	14
1.59 Polycarbonate							
Varilux® X 4D™ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode™ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	14

Varilux Comfort® W2+ Eyecode™ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® W2+ FIT Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® W2+ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Varilux Comfort DRx™ Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Varilux Comfort® Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short™ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.00	-6.00	14
Varilux® Stylistic™ Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced Wrap™ Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Essilor Ideal Advanced™ FIT Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Essilor Ideal® FIT Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Computer Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Driving Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Ideal® Short	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Natural® Digital Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Ovation® Digital Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor SmallFit® Digital Airwear®	GR	BR	GG	-10.00 to +8.00	+0.75 to +3.00	-6.00	14
1.67 Super High Index							
Varilux® X 4D™ Thin&Lite® 1.67	GR	BR	GG	-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ Thin&Lite® 1.67	GR	BR	GG	-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ Thin&Lite® 1.67	GR	BR	GG	-12.00 to +9.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14






Varilux® Physio® W3+ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® Max Fit	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® Max	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® W2+ Eyecode Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® W2+ FIT Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® W2+ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx™ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx™ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx™ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx™ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Varilux® Stylistic™ Thin&Lite® 1.67	GR	BR	GG	-5.00 to +5.00	+0.75 to +3.50	-3.00	17
Essilor Ideal Advanced™ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced Wrap™ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced™ FIT Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Computer Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Driving Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor Ideal® Short	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor Natural® Digital Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom™ Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	15
Essilor Accolade® Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Brillenz® Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Brillenz® Short Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14
Ovation® Digital Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor SmallFit® Digital Thin&Lite® 1.67	GR	BR	GG	-14.00 to +9.00	+0.75 to +3.00	-6.00	14

					
Semi-Finished Single Vision				Sphere	Cyl to
1.50 Standard Index					






HOYA SV 1.50	GR	BR		-8.00 to +9.00	-4.00
HOYA MySV 1.50	GR!	BR!		-8.00 to +9.00	-4.00
HOYA SV iQ 1.50	GR	BR		-8.00 to +9.00	-4.00
Sync III	GR	BR		-10.00 to +6.00	-4.00
1.53 TRIVEX® Single Vision					
HOYA SFSV Phoenix™	GR	BR		-10.00 to +6.00	-4.00
HOYA MySV Phoenix™	GR!	BR!		-10.00 to +9.00	-4.00
HOYA SV iQ Phoenix™	GR	BR		-10.00 to +6.00	-4.00
1.59 Polycarbonate Single Vision					
HOYA SFSV Polycarbonate	GR	BR		-6.00 to +6.00	-4.00
HOYA MySV Polycarbonate	GR!	BR!		-10.00 to +6.00	-4.00
HOYA SV iQ Polycarbonate	GR	BR		-10.00 to +6.00	-4.00
Sync III	GR	BR		-13.00 to +5.25	-4.00
1.67 Super High Index Single Vision					
HOYA MySV 1.67	GR!	BR!		-10.00 to +8.00	-4.00
HOYA SV iQ 1.67	GR	BR		-10.00 to +8.00	-4.00
Sync III	GR	BR		-13.00 to +7.50	-4.00
1.74 Super High Index Single Vision					
HOYA SV 1.74	GR	BR		-13.00 to +10.00	-4.00






HOYA	Transitions XTRACTIVE® NEW GENERATION			Sphere	Add	Cyl to	MFH
	Progressives						
1.50 Standard Index							
HOYALUX ARRAY™ 11, 13, 15, 17, VL 1.50	GR	BR		-10.00 to +3.75	+0.75 to +3.50	-4.00	14, 17, 19, 21
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL 1.50	GR	BR		-5.00 to +4.00	+0.75 to +3.50	-4.00	14, 17, 19, 21
HOYALUX ARRAY™ 2 11, 13, 15, 17, VL 1.50	GR	BR		-10.00 to +3.75	+0.75 to +3.50	-4.00	14, 17, 19, 21
HOYALUX ARRAY™ 2 Wrap 11, 13, 15, 17, VL 1.50	GR	BR		-5.00 to +4.00	+0.75 to +3.50	-4.00	14, 17, 19, 21
Amplitude BKS 1.50	GR	BR		-10.00 to +3.75	+0.75 to +3.50	-4.00	19
Amplitude Mini BKS 1.50	GR	BR		-10.00 to +3.75	+0.75 to +3.50	-4.00	14
Summit cd BKS 1.50	GR	BR		-10.00 to +3.75	+0.75 to +3.50	-4.00	14
Summit ecp BKS 1.50	GR	BR		-10.00 to +3.75	+0.75 to +3.50	-4.00	19
Hoya GP Wide 1.50	GR	BR		-10.00 to +3.75	+0.75 to +3.50	-4.00	18
Tact 40 BKS 1.50	GR	BR		-10.00 to +3.75	+0.75 to +3.50	-4.00	18
Tact 60 BKS 1.50	GR	BR		-10.00 to +3.75	+0.75 to +3.50	-4.00	18
1.53 TRIVEX®							
HOYALUX ARRAY™ 11, 13, 15, 17, VL Phoenix™	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-4.00	14, 17, 19, 21
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL Phoenix™	GR	BR		-5.00 to +4.00	+0.75 to +3.50	-4.00	14, 17, 19, 21
HOYALUX ARRAY™ 2 11, 13, 15, 17, VL Phoenix™	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-4.00	14, 17, 19, 21

HOYALUX ARRAY 2 Wrap 11, 13, 15, 17, VL Phoenix™	GR	BR		-5.00 to +4.00	+0.75 to +3.50	-4.00	14, 17, 19, 21
Amplitude BKS Phoenix™	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-4.00	19
Amplitude Mini BKS Phoenix™	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-4.00	14
Amplitude BKS Phoenix™	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-4.00	14
Amplitude Mini BKS Phoenix™	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-4.00	19
Summit cd BKS Phoenix™	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-4.00	18
Summit ecp BKS Phoenix™	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-4.00	18
Hoya GP Wide Phoenix™	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-4.00	18
Tact 40 BKS Phoenix™	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-4.00	18
Tact 60 BKS Phoenix™	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-4.00	18
1.59 Polycarbonate							
HOYALUX ARRAY™ 11, 13, 15, 17, VL Poly	GR	BR		-13.00 to +5.25	+0.75 to +3.50	-4.00	14, 17, 19, 21
HOYALUX ARRAY™ Wrap 11, 13, 15, 17, VL Poly	GR	BR		-6.00 to +4.00	+0.75 to +3.50	-4.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13, 15, 17, VL Poly	GR	BR		-13.00 to +5.25	+0.75 to +3.50	-4.00	14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11, 13, 15, 17, VL Poly	GR	BR		-6.00 to +4.00	+0.75 to +3.50	-4.00	14, 17, 19, 21
Amplitude BKS Poly	GR	BR		-13.00 to +5.25	+0.75 to +3.50	-4.00	19
Amplitude Mini BKS Poly	GR	BR		-13.00 to +5.25	+0.75 to +3.50	-4.00	14
Summit cd BKS Poly	GR	BR		-13.00 to +5.25	+0.75 to +3.50	-4.00	14
Summit ecp BKS Poly	GR	BR		-13.00 to +5.25	+0.75 to +3.50	-4.00	19
Hoya GP Wide Poly	GR	BR		-13.00 to +5.25	+0.75 to +3.50	-4.00	18
Tact 40 BKS Poly	GR	BR		-13.00 to +5.25	+0.75 to +3.50	-4.00	18
Tact 60 BKS Poly	GR	BR		-13.00 to +5.25	+0.75 to +3.50	-4.00	18
1.67 Super High Index							
HOYALUX ARRAY 11, 13,15,17,VL 1.67	GR	BR		-13.00 to +7.50	+0.75 to +3.50	-4.00	14, 17, 19, 21
HOYALUX ARRAY Wrap 11,13,15,17,VL 1.67	GR	BR		-13.00 to +7.50	+0.75 to +3.50	-4.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13,15,17,VL 1.67	GR	BR		-13.00 to +7.50	+0.75 to +3.50	-4.00	14, 17, 19, 21
HOYALUX ARRAY 2 Wrap 11,13,15,17,VL 1.67	GR	BR		-6.00 to +4.00	+0.75 to +3.50	-4.00	14, 17, 19, 21
Amplitude BKS 1.67	GR	BR		-13.00 to +7.50	+0.75 to +3.50	-4.00	19
Amplitude Mini BKS 1.67	GR	BR		-13.00 to +7.50	+0.75 to +3.50	-4.00	14
Summit cd BKS 1.67	GR	BR		-13.00 to +7.50	+0.75 to +3.50	-4.00	14
Summit ecp BKS 1.67	GR	BR		-13.00 to +7.50	+0.75 to +3.50	-4.00	19
Hoya GP Wide 1.67	GR	BR		-13.00 to +7.50	+0.75 to +3.50	-4.00	18
Tact 40 BKS 1.67	GR	BR		-13.00 to +7.50	+0.75 to +3.50	-4.00	18
Tact 60 BKS 1.67	GR	BR		-13.00 to +7.50	+0.75 to +3.50	-4.00	18






					
Semi-Finished Single Vision				Sphere	Cyl to
1.50 Standard Index					

Nikon SeeMax Single Vision	GR	BR	GG	-10.00 to + 8.00	-6.00
Nikon Infinite Single Vision	GR	BR	GG	-10.00 to + 8.00	-6.00
Nikon Smart Boost	GR	BR	GG	-10.00 to + 8.00	-6.00
1.59 Polycarbonate Single Vision					
Nikon SeeMax Single Vision	GR	BR	GG	-10.00 to +8.00	-6.00
Nikon Infinite Single Vision	GR	BR	GG	-10.00 to +8.00	-6.00
Nikon Smart Boost	GR	BR	GG	-10.00 to +8.00	-6.00
1.67 Super High Index Single Vision					
Nikon SeeMax Single Vision 1.67	GR	BR	GG	-12.00 to + 8.00	-4.00
Nikon Infinite Single Vision 1.67	GR	BR	GG	-12.00 to + 8.00	-6.00
Nikon Smart Boost	GR	BR	GG	-12.00 to + 8.00	-6.00

							
Progressives				Sphere	Add	Cyl to	MFH
1.50 Standard Index							
Nikon SeeMax Ultimate	GR	BR	GG	-10.00 to + 8.00	+0.75 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR	GG	-10.00 to + 8.00	+0.75 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR	GG	-7.00 to + 4.00	+0.75 to +4.00	-6.00	15, 17
Nikon Move	GR	BR	GG	-7.00 to + 4.00	+0.75 to +3.50	-4.00	16, 18
1.59 Polycarbonate							
Nikon SeeMax Ultimate	GR	BR	GG	-10.00 to + 8.00	+0.75 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR	GG	-10.00 to + 8.00	+0.75 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR	GG	-10.00 to + 5.00	+0.75 to +4.00	-4.00	15, 17
Nikon Move	GR	BR	GG	-10.00 to + 5.00	+0.75 to +4.00	-4.00	16, 18
1.67 Super High Index							
Nikon SeeMax Ultimate	GR	BR	GG	-10.00 to + 8.00	+0.75 to +3.50	-6.00	13, 15, 17
Nikon Presio Power	GR	BR	GG	-14.00 to + 8.00	+0.75 to +3.50	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR	GG	-12.00 to + 8.00	+0.75 to +3.50	-6.00	15, 17
Nikon Move	GR	BR	GG	-12.00 to + 5.00	+0.75 to +3.50	-4.00	16, 18






					
Semi-Finished Single Vision				Sphere	Cyl to
1.53 TRIVEX® Single Vision					
Seiko PL 1.53/Trivex Transitions®	GR			-10.00 to +7.00	-5.00
1.67 Super High Index Single Vision					
Seiko PL 1.67 Transitions®	GR			-12.00 to +10.00	-5.00

					
Semi-Finished Single Vision				Sphere	Cyl to
1.50 Standard Index					
Shamir Attitude III Single Vision	GR			-6.50 to +6.00	-6.00
Shamir Autograph® III Single Vision	GR			-12.00 to +7.25	-6.00
Shamir Relax™	GR			-14.00 to +8.00	-6.00
Shamir Relax™ 50	GR			-14.00 to +8.00	-6.00
Shamir Relax™ 80	GR			-14.00 to +8.00	-6.00
1.53 SFSV					
Shamir Attitude III Single Vision	GR			-6.50 to +6.00	-6.50
Shamir Autograph® III Single Vision	GR			-14.00 to +8.00	-6.50
Shamir Relax™	GR			-14.00 to +8.00	-6.50
Shamir Relax™ 50	GR			-14.00 to +8.00	-6.50
Shamir Relax™ 80	GR			-14.00 to +8.00	-6.50
1.59 SFSV					
Shamir Attitude III Single Vision	GR			-6.50 to +6.00	-7.00
Shamir Autograph® III Single Vision	GR			-14.00 to +8.75	-7.00
Shamir Relax™	GR			-14.75 to +6.00	-7.00
Shamir Relax™ 50	GR			-14.75 to +6.00	-7.00
Shamir Relax™ 80	GR			-14.75 to +6.00	-7.00
1.60 SFSV					
Shamir Attitude III Single Vision	GR			-6.50 to +6.00	-7.00
Shamir Autograph® III Single Vision	GR			-14.00 to +8.75	-7.00
Shamir Relax™	GR			-14.00 to +8.75	-7.00
Shamir Relax™ 50	GR			-14.00 to +8.75	-7.00
Shamir Relax™ 80	GR			-14.00 to +8.75	-7.00
1.67 SFSV					
Shamir Attitude III Single Vision	GR			-6.50 to +6.00	-8.00
Shamir Autograph® III Single Vision	GR			-16.00 to +16.00	-8.00
Shamir Relax™	GR			-16.75 to +16.00	-8.00
Shamir Relax™50	GR			-16.75 to +16.00	-8.00
Shamir Relax™80	GR			-16.75 to +16.00	-8.00






					
Semi-Finished Single Vision				Sphere	Cyl to
1.50 Standard Index					

Signet Armorlite SFSV	GR			-6.00 to +10.00	-4.00
Signetek Processed SV	GR			-8.00 to +6.00	-6.00
KODAK Digital SV	GR			-8.00 to +6.00	-6.00
1.53 TRIVEX® Single Vision					
Signetek Processed SV	GR			-6.00 to +10.00	-4.00
KODAK Digital SV	GR			-8.00 to +6.00	-6.00
1.59 Polycarbonate Single Vision					
Signet Armorlite SFSV	GR			-8.00 to +10.00	-4.00
Signetek Processed SV	GR			-10.00 to +6.50	-6.00
KODAK Digital SV	GR			-10.00 to +6.50	-6.00
1.67 Super High Index Single Vision					
Signet Armorlite SFSV	GR			-13.00 to +9.00	-4.00
Signetek Processed SV	GR			-11.00 to +7.50	-6.00
KODAK Digital SV	GR			-11.00 to +7.50	-6.00

						
Multifocal				Sphere	Add	Cyl to
1.50 Standard Index SFST-28						
Signetek™ Processed 1.50 FT-28	GR			-10.00 to +6.00	+1.00 to +3.00	-6.00






							
Progressives				Sphere	Add	Cyl to	MFH
1.50 Standard Index							
KODAK Precise® Plus 1.50	GR			-8.00 to +4.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short Plus 1.50	GR			-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Precise® PB 1.50	GR			-8.00 to +4.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB 1.50	GR			-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique™ 1.50	GR			-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ 1.50	GR			-8.00 to +4.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.50	GR			-8.00 to +4.00	+0.75 to +3.50	-6.00	13
1.53 TRIVEX®							
KODAK Precise® PB Trivex®	GR			-10.00 to +4.50	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB Trivex®	GR			-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique Trivex®	GR			-10.00 to +4.50	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ Trivex®	GR			-10.00 to +4.50	+0.75 to +3.50	-6.00	13



KODAK Unique DRO™ HD Trivex®	GR			-10.00 to +4.50	+0.75 to +3.50	-6.00	13
1.59 Polycarbonate							
KODAK Unique DRO™ Poly	GR			-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Precise® PB Poly	GR!			-10.00 to +5.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB Poly	GR			-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique Poly	GR!			-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD Poly	GR!			-10.00 to +5.00	+0.75 to +3.50	-6.00	13
1.67 Super High Index							
KODAK Precise® PB 1.67	GR			-11.00 to +6.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB 1.67	GR			-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique	GR			-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ 1.67	GR!			-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.67	GR!			-11.00 to +6.00	+0.75 to +3.50	-6.00	13

					
Finished Single Vision				Sphere	Cyl to
1.50 Standard Index					
TRN 1.50 SV UC/HC DIA.65	GR	BR		+0.25 to +4.00	-2.00
TRN 1.50 SV UC/HC DIA.70	GR	BR		+2.00 to -4.00	-2.00

					
Semi-Finished Single Vision				Sphere	Cyl to
1.50 Standard Index					
TRN 1.50 SV (72mm) UC/HC DIA.72	GR	BR		Base 0.50, 2.25, 4.25, 5.25, 6.25, 8.25	
TRN 1.50 DS (76mm) UC/HC DIA.72	GR	BR		Base 0.50, 8.75, 10.50	
TRN 1.50 DS (76mm) UC/HC DIA.76	GR	BR		Base 2.00, 3.50, 4.25, 5.25, 6.25, 7.25	
1.53 TRIVEX® Single Vision					
TOG 1.53 Trivex (75mm) SV UC/HC DIA.70	GR	BR		Base 1.00	
TOG 1.53 Trivex (75mm) SV UC/HC DIA.75	GR	BR		Base 2.50, 4.00, 5.00, 6.00, 7.00, 8.00	
1.60 High Index Single Vision					
TOG 1.60 SV (75mm) UC/HC DIA.65	GR	BR		Base 10.00, 11.00	
TOG 1.60 SV (75mm) UC/HC DIA.70	GR	BR		Base 9.00	
TOG 1.60 SV (75mm) UC/HC DIA.75	GR	BR		Base 1.00, 2.50, 4.00, 5.00, 6.00, 7.00, 8.00	






TOG 1.60 SV (80mm) UC/HC DIA.80	GR	BR		Base 1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.50
1.67 Super High Index Single Vision				
TOG 1.67 SV (75mm) UC/HC DIA.65	GR	BR		Base 10.00, 11.00
TOG 1.67 SV (75mm) UC/HC DIA.70	GR	BR		Base 9.00
TOG 1.67 SV (75mm) UC/HC DIA.75	GR	BR		Base 1.00, 2.50, 4.00, 5.00, 6.00, 7.00, 8.00
TOG 1.67 SV (80mm) UC/HC DIA.80	GR	BR		Base 1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.50





					
Semi-Finished Single Vision				Sphere	Cyl to
1.50 Standard Index					
Unity SVx	GR	BR		-12.25 to +6.00	-6.00
Unity SVxtra	GR	BR		-12.25 to +6.00	-6.00
Unity SVxtreme	GR	BR		-6.00 to +6.00	-6.00
1.53 TRIVEX® Single Vision					
Unity SVx	GR	BR		-14.00 to +6.00	-6.00
Unity SVxtra	GR	BR		-14.00 to +6.00	-6.00
Unity SVxtreme	GR	BR		-6.50 to +6.00	-6.00
1.59 Polycarbonate Single Vision					
Unity SVx	GR	BR		-14.00 to +6.00	-6.00
Unity SVxtra	GR	BR		-14.00 to +6.00	-6.00
Unity SVxtreme	GR	BR		-6.50 to +6.00	-6.00
1.67 Super High Index Single Vision					
Unity SVx	GR	BR		-16.00 to +8.00	-6.00
Unity SVxtra	GR	BR		-16.00 to +8.00	-6.00
Unity SVxtreme	GR	BR		-6.00 to +6.00	-6.00






							
Progressives				Sphere	Add	Cyl to	MFH
1.50 Standard Index							
Unity Via	GR	BR		-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Plus	GR	BR		-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Elite	GR	BR		-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Mobile	GR	BR		-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Wrap	GR	BR		-7.00 to +6.00	+0.50 to +4.50	-6.00	12






Unity Via II	GR	BR		-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR		-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR		-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR		-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR		-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity PLx	GR	BR		-7.00 to +6.00	+0.50 to +4.50	-6.00	13-19
Unity PLxtra	GR	BR		-7.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity PLxtreme	GR	BR		-7.00 to +6.00	+0.50 to +4.50	-6.00	15-18
Ethos	GR	BR		-9.00 to +4.50	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR		-9.00 to +4.50	0.75 to +3.50	-4.00	14
1.53 TRIVEX®							
Unity Via	GR	BR		-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Plus	GR	BR		-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Elite	GR	BR		-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Mobile	GR	BR		-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Wrap	GR	BR		-5.00 to +4.00	+0.50 to +4.50	-6.00	12
Unity Via II	GR	BR		-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR		-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR		-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR		-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR		-5.00 to +4.00	+0.50 to +4.50	-6.00	12
Unity PLx	GR	BR		-10.00 to +6.00	+0.50 to +4.50	-6.00	13-19
Unity PLxtra	GR	BR		-10.00 to +6.00	+0.50 to +4.50	-6.00	12
Unity PLxtreme	GR	BR		-10.00 to +6.00	+0.50 to +4.50	-6.00	15-18
Ethos	GR	BR		-9.00 to +6.00	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR		-9.00 to +6.00	0.75 to +3.50	-4.00	14
1.59 Polycarbonate							
Unity Via	GR	BR		-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Plus	GR	BR		-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Elite	GR	BR		-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Mobile	GR	BR		-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Wrap	GR	BR		-6.50 to +6.00	+0.50 to +4.50	-6.00	12
Unity Via II	GR	BR		-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR		-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR		-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR		-10.00 to +6.50	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR		-6.50 to +6.00	+0.50 to +4.50	-6.00	12
Unity PLx	GR	BR		-10.00 to +6.50	+0.50 to +4.50	-6.00	13-19
Unity PLxtra	GR	BR		-10.00 to +6.50	+0.50 to +4.50	-6.00	12

Unity PLxtreme	GR	BR		-10.00 to +6.50	+0.50 to +4.50	-6.00	15-18
Ethos	GR	BR		-9.00 to +6.00	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR		-9.00 to +6.00	0.75 to +3.50	-4.00	14
1.67 Super High Index							
Unity Via	GR	BR		-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Plus	GR	BR		-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Elite	GR	BR		-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Mobile	GR	BR		-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Wrap	GR	BR		-6.00 to +5.00	+0.50 to +4.50	-6.00	12
Unity Via II	GR	BR		-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Plus II	GR	BR		-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Elite II	GR	BR		-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Mobile II	GR	BR		-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity Via Wrap II	GR	BR		-6.00 to +5.00	+0.50 to +4.50	-6.00	12
Unity PLx	GR	BR		-12.00 to +8.25	+0.50 to +4.50	-6.00	13-19
Unity PLxtra	GR	BR		-12.00 to +8.25	+0.50 to +4.50	-6.00	12
Unity PLxtreme	GR	BR		-12.00 to +8.25	+0.50 to +4.50	-6.00	15-18
Ethos	GR	BR		-12.00 to +8.00	0.75 to +3.50	-4.00	14
Ethos Plus	GR	BR		-12.00 to +8.00	0.75 to +3.50	-4.00	14

					
Finished Single Vision				Sphere	Cyl to
1.50 Standard Index					
Younger Optics FSV	GR	BR		-4.00 to +4.00	-2.00
1.59 Polycarbonate Finished Single Vision					
Younger Optics FSV Poly	GR	BR		-6.00 to +4.00	-2.00

					
Semi-Finished Single Vision				Sphere	Cyl to
1.50 Standard Index					
Younger Optics SFSV (72mm)	GR	BR		-9.00 to +6.00	-4.00
1.59 Polycarbonate Single Vision					
Younger Optics SV	GR	BR		-9.00 to +6.00	-4.00
1.67 Super High Index Single Vision					
Younger Optics Single Vision 1.67 MR-10™	GR	BR		-11.00 to +7.00	-4.00
1.74 Ultra High Index Single Vision					
Younger Optics Single Vision 1.74	GR	BR		-18.00 to +12.00	-4.00

 The Optical Lens Innovators		 NEW GENERATION				
Digital					Sphere	Cyl to
1.50 Standard Index						
	Camber™	GR	BR		-11.00 to +6.25	-4.00
1.53 TRIVEX®						
	Camber™	GR	BR		-11.00 to +6.25	-4.00
1.59 Polycarbonate						
	Camber™	GR	BR		-10.00 to +6.25	-4.00
1.67 Super High Index						
	Camber™	GR	BR		-14.00 to +8.00	-4.00
1.74 Ultra High Index						
	Camber™	GR	BR		-19.00 to +9.00	-4.00

 The Optical Lens Innovators		 NEW GENERATION						
Progressives					Sphere	Add	Cyl to	MFH
1.59 index polycarbonate								
	Younger Image®	GR	BR		-9.00 to +7.00	+1.00 to +3.00	-4.00	18

TRANSITIONS® LIGHT INTELLIGENT LENSES™

- 
GR
Gray
- 
BR
Brown
- 
GG
Graphite Green
- 
S
Sapphire
- 
AME
Amethyst
- 
AMB
Amber
- 
E
Emerald
- 
Available
in style
mirrors
- !
- #

Available Q4 2022

Available Q1 2023

Transitions, Transitions Signature, XTRActive, and Vantage are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN 8 and Light Under Control are trademarks of Transitions Optical Limited. ©2022 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.