




Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® Lenses PAG](#) / [Transitions® XTRActive® Polarized™](#) / [Nikon Optical USA](#)

Nikon Optical USA – Transitions® XTRActive® Polarized™

				
Semi-Finished Single Vision		Sphere	Cyl to	
1.50 Standard Index				
Nikon SeeMax Infinite	GR	-10.00 to +8.00	-6.00	
Nikon Infinite Single Vision	GR	-10.00 to +8.00	-6.00	
Nikon Smart Boost	GR	-10.00 to +8.00	-6.00	
1.59 Polycarbonate Single Vision				
Nikon SeeMax Infinite	GR	-10.00 to +8.00	-6.00	
Nikon Infinite Single Vision	GR	-10.00 to +8.00	-6.00	
Nikon Smart Boost	GR	-10.00 to +8.00	-6.00	
1.67 Super High Index Single Vision				
Nikon SeeMax Infinite	GR	-14.00 to +9.00	-6.00	
Nikon Infinite Single Vision	GR	-14.00 to +9.00	-6.00	
Nikon Smart Boost	GR	-14.00 to +9.00	-6.00	

					
Progressives		Sphere	Add	Cyl to	MFH
1.50 Standard index					
Nikon SeeMax Ultimate	GR	-10.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	-10.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	-10.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
Nikon Move	GR	-10.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
1.59 Polycarbonate					
Nikon SeeMax Ultimate	GR	-11.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	-11.00 to +8.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	-11.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
Nikon Move	GR	-11.00 to +8.00	+0.50 to +4.00	-6.00	15, 17
1.67 Super High Index					
Nikon SeeMax Ultimate	GR	-14.00 to +9.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	-14.00 to +9.00	+0.50 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	-14.00 to +9.00	+0.50 to +4.00	-6.00	15, 17
Nikon Move	GR	-14.00 to +9.00	+0.50 to +4.00	-6.00	15, 17

TRANSITIONS® LIGHT INTELLIGENT LENSES™



All other trademarks are the property of their respective owners.