

Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® XTRActive® Polarized™](#) / **Thai Optical Group USA**

Thai Optical Group USA – Transitions® XTRActive® Polarized™

			
Semi-Finished Single Vision		Sphere	Cyl to
1.50 Standard Index			
TRN 1.50 DS (76mm) UC/HC DIA.72	GR	Base 0.50, 8.75, 10.50	
TRN 1.50 DS (76mm) UC/HC DIA.76	GR	Base 2.00, 3.50, 4.25, 5.25, 6.25, 7.25	
1.53 TRIVEX® Single Vision			
TRN 1.53 Trivex (78mm) SV UC/HC DIA.74	GR	Base 0.50, 1.00	
TRN 1.53 Trivex (78mm) SV UC/HC DIA.78	GR	Base 2.00, 3.50, 4.25, 5.25, 6.25, 7.25, 8.75	
1.59 Polycarbonate Single Vision			
TRN 1.59 Polycarbonate (79mm) SV HC DIA.79	GR	Base 1.50, 3.25, 4.25, 4.75, 5.50, 7.00	
1.60 High Index Single Vision			
TRN 1.60 SV (80mm) UC/HC DIA.76	GR	Base 8.50	
TRN 1.60 SV (80mm) UC/HC DIA.80	GR	Base 1.75, 2.75, 3.75, 4.75, 6.00, 6.75	
1.67 Super High Index Single Vision			
TRN 1.67 SV (80mm) UC/HC DIA.80	GR	Base 1.75, 2.75, 3.75, 4.75, 6.00	
TRN 1.67 SV (80mm) UC/HC DIA.81	GR	Base 6.75, 8.00	

TRANSITIONS® LIGHT INTELLIGENT LENSES™



GR
Gray



BR
Brown



GG
Graphite
Green



S
Sapphire



AME
Amethyst



AMB
Amber



E
Emerald



R
Ruby



Available
in style
mirrors

Transitions, Transitions Signature, and XTRActive are registered trademarks, and XTRActive Polarized, the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.