






Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® XTRActive® Polarized™](#) / **Unity Lenses**

Unity Lenses – Transitions® XTRActive® Polarized™

				
Semi-Finished Single Vision		Sphere	Cyl to	
1.50 Standard Index				
Unity SVx	GR	-12.25 to +6.00	-6.00	
Unity SVxtra	GR	-12.25 to +6.00	-6.00	
Unity SVxtreme	GR	-6.00 to +6.00	-6.00	
1.59 Polycarbonate Single Vision				
Unity SVx	GR	-14.00 to +6.00	up to -6.00	
Unity SVxtra	GR	-14.00 to +6.00	up to -6.00	
Unity SVxtreme	GR	-6.50 to +6.00	up to -6.00	
1.67 Super High Index Single Vision				
Unity SVx	GR	-16.00 to +8.00	up to -6.00	
Unity SVxtra	GR	-16.00 to +8.00	up to -6.00	
Unity SVxtreme	GR	-6.00 to +6.00	up to -6.00	

					
Progressives		Sphere	Add	Cyl to	MFH
1.50 Standard Index					
Unity Via II	GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Plus II	GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite VR	GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II	GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II 12	GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II 14	GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	14 mm
Unity Via Elite II 16	GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	16 mm
Unity Via Elite II 18	GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	18 mm
Unity Via Mobile II	GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Wrap II	GR	-7.00 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Ethos	GR	-9.00 to +4.50	+0.75 to +3.50	up to -4.00	14 mm
Ethos Plus	GR	-9.00 to +4.50	+0.75 to +3.50	up to -4.00	14 mm
1.59 Polycarbonate					
Unity Via II	GR	-10.00 to +6.50	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Plus II	GR	-10.00 to +6.50	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite VR	GR	-10.00 to +6.50	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II	GR	-10.00 to +6.50	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II 12	GR	-10.00 to +6.50	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II 14	GR	-10.00 to +6.50	+0.50 to +4.50	up to -6.00	14 mm
Unity Via Elite II 16	GR	-10.00 to +6.50	+0.50 to +4.50	up to -6.00	16 mm
Unity Via Elite II 18	GR	-10.00 to +6.50	+0.50 to +4.50	up to -6.00	18 mm
Unity Via Mobile II	GR	-10.00 to +6.50	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Wrap II	GR	-6.50 to +6.00	+0.50 to +4.50	up to -6.00	12 mm
Ethos	GR	-9.00 to +6.00	+0.75 to +3.50	up to -4.00	14 mm
Ethos Plus	GR	-9.00 to +6.00	+0.75 to +3.50	up to -4.00	14 mm
1.67 Super High Index					
Unity Via II	GR	-12.00 to +8.25	+0.50 to +4.50	up to -6.00	12 mm

Unity Via Plus II	GR	-12.00 to +8.25	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite VR	GR	-12.00 to +8.25	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II	GR	-12.00 to +8.25	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II 12	GR	-12.00 to +8.25	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Elite II 14	GR	-12.00 to +8.25	+0.50 to +4.50	up to -6.00	14 mm
Unity Via Elite II 16	GR	-12.00 to +8.25	+0.50 to +4.50	up to -6.00	16 mm
Unity Via Elite II 18	GR	-12.00 to +8.25	+0.50 to +4.50	up to -6.00	18 mm
Unity Via Mobile II	GR	-12.00 to +8.25	+0.50 to +4.50	up to -6.00	12 mm
Unity Via Wrap II	GR	-6.00 to +5.00	+0.50 to +4.50	up to -6.00	12 mm
Ethos	GR	-12.00 to +8.00	+0.75 to +3.50	up to -4.00	14 mm
Ethos Plus	GR	-12.00 to +8.00	+0.75 to +3.50	up to -4.00	14 mm

TRANSITIONS® LIGHT INTELLIGENT LENSES™



GR
Gray



BR
Brown



GG
Graphite
Green



S
Sapphire



AME
Amethyst



AMB
Amber



E
Emerald



R
Ruby



Available
in style
mirrors

Transitions, *Transitions Signature*, and *XTRActive* are registered trademarks, and *XTRActive Polarized*, the *Transitions* logo and *Transitions Light Intelligent Lenses* are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.