



Transitions® Lenses Product Availability Guide

< [Home](#) / [Resources](#) / [Transitions® XTRActive® Polarized™](#) / [Unity Lenses](#)

Unity Lenses – Transitions® XTRActive® Polarized™

					
Semi-Finished Single Vision			Sphere	Cyl to	
1.50 Standard Index					
	Unity SVx	GR	-12.25 to +6.00	-6.00	
	Unity SVxtra	GR	-12.25 to +6.00	-6.00	
	Unity SVxtreme	GR	-6.00 to +6.00	-6.00	
1.59 Polycarbonate Single Vision					
	Unity SVx	GR	-14.00 to +6.00	-6.00	
	Unity SVxtra	GR	-14.00 to +6.00	-6.00	
	Unity SVxtreme	GR	-6.50 to +6.00	-6.00	
1.67 Super High Index Single Vision					
	Unity SVx	GR	-16.00 to +8.00	-6.00	
	Unity SVxtra	GR	-16.00 to +8.00	-6.00	
	Unity SVxtreme	GR	-6.00 to +6.00	-6.00	

						
Progressives			Sphere	Add	Cyl to	MFH
1.50 Standard Index						
	Unity Via II	GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Plus II	GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Elite II	GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Mobile II	GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12
	Unity Via Wrap II	GR	-7.00 to +6.00	+0.50 to +4.50	-6.00	12
	Ethos	GR	-9.00 to +4.50	0.75 to +3.50	-4.00	14
	Ethos Plus	GR	-9.00 to +4.50	0.75 to +3.50	-4.00	14
1.59 Polycarbonate						
	Unity Via II	GR	-10.00 to +6.50	+0.50 to +4.50	-6.00	12
	Unity Via Plus II	GR	-10.00 to +6.50	+0.50 to +4.50	-6.00	12
	Unity Via Elite II	GR	-10.00 to +6.50	+0.50 to +4.50	-6.00	12
	Unity Via Mobile II	GR	-10.00 to +6.50	+0.50 to +4.50	-6.00	12
	Unity Via Wrap II	GR	-6.50 to +6.00	+0.50 to +4.50	-6.00	12
	Ethos	GR	-9.00 to +6.00	0.75 to +3.50	-4.00	14
	Ethos Plus	GR	-9.00 to +6.00	0.75 to +3.50	-4.00	14
1.67 Super High Index						
	Unity Via II	GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
	Unity Via Plus II	GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
	Unity Via Elite II	GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
	Unity Via Mobile II	GR	-12.00 to +8.25	+0.50 to +4.50	-6.00	12
	Unity Via Wrap II	GR	-6.00 to +5.00	+0.50 to +4.50	-6.00	12
	Ethos	GR	-12.00 to +8.00	0.75 to +3.50	-4.00	14
	Ethos Plus	GR	-12.00 to +8.00	0.75 to +3.50	-4.00	14



Transitions, *Transitions Signature*, and *XTRActive* are registered trademarks, and *XTRActive Polarized*, the *Transitions* logo and *Transitions Light Intelligent Lenses* are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. GEN S and GEN 8 are trademarks of Transitions Optical Limited. ©2025 Transitions Optical Limited. Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.

All other trademarks are the property of their respective owners.