
YOUR GUIDE TO DISPENSING SUCCESS!









Transitions®

ADAPTIVE LENSES

SPRING 2013

LOG-ON TO WWW.TRANSITIONSPRO.COM FOR THE LATEST INFORMATION

There is something for every patient.

Use this as your guide to helping your patient choose the right Transitions® brand everyday lenses to fit their life.





For the patient who wants:

- Full clarity indoors, just the right amount of tint outdoors
- Fastest to fade back to clear indoors out of all *Transitions* lenses

RECOMMEND: Transitions®

COLORS: Gray Brown



For the patient who wants:

- More activation for extra protection
- Activation behind the windshield of a car

RECOMMEND: Transitions® XTRActive®

COLOR: Gray



For the patient who wants:

- Crisper, sharper vision outdoors
- The only everyday lenses with variable polarization

RECOMMEND: Transitions® Vantage™

COLOR: Gray


TRANSITIONS LENS COMPARISON CHART

	Transitions	Transitions XTRActive	Transitions Vantage	Regular clear lenses
Lens Description	Fully clear indoors, optimal tint outdoors	Always active for extra protection	Variable polarization for even crisper, sharper vision outdoors	Clear
Lens type	For your primary prescription eyeglasses	For your primary prescription eyeglasses	For your primary prescription eyeglasses	For your primary prescription eyeglasses
Optimal Usage	Indoors and outdoors, all day	Indoors and outdoors, all day	Indoors and outdoors, all day	Indoors and at night
Lens color options	Clear to Gray Clear to Brown	Clear to Gray	Clear to Gray	Clear
Indoor clarity	Fully clear	Hint of tint	Hint of tint	Clear
Outdoor darkness	Dark	Our darkest lens, even in high temperatures	Dark	Clear
Speed to darken	Fast to activate	Fast to activate	Fast to activate	Doesn't activate
Behind-the-windshield activation	No	Yes, moderate activation	No	No
Polarization	No	No	Variable polarization that adjusts to match outdoor glare levels	No

Transitions lenses help your eyes adjust to changing light better than they would on their own.



All *Transitions* lenses block 100% UVA/UVB rays.



Transitions everyday adaptive lenses are right for any age.





Reinforce the value of your patients purchase with the Transitions lenses Certificate of Authenticity!

Your patients have chosen *Transitions* lenses to enhance how they see their world. Reinforce the value of their purchase with the Transitions lenses Certificate of Authenticity so they feel confident knowing they are getting the world's #1 recommended photochromic lens.

TOI-ID number:

FSV Transitions Product Availability

	Transitions® GRAY	Transitions® BROWN	Transitions® XTRActive® GRAY	SPHERE	CYL TO
1.50 STANDARD INDEX FSV					
Essilor Lens FSV Crizal AVANCE™	GR			-4.00 to +4.00	-2.00
Essilor Lens FSV Crizal® Alizé®	GR			-4.00 to +4.00	-2.00
Essilor Lens FSV	GR	BR	XA	-4.00 to +4.00	-2.00
SOLA FSV	GR			-4.00 to +4.00	-2.00
HOYA Vision Care FSV	GR			-4.00 to +4.00	-2.00
Seiko Ultra 1.50 Spherical FSV	GR			-4.00 to +3.00	-2.00
Signet Armorlite FSV	GR	BR		-4.00 to +4.00	-2.00
Signet Armorlite Kodak FSV 1.50 Clean'N'CleAR	GR			-4.00 to +4.00	-2.00
X-Cel Optical 1.50 FSV	GR	BR	XA	-4.00 to +4.00	-2.00
Younger Optics FSV AR Super Hydrophobic	GR	BR		-4.00 to +4.00	-2.00
Younger Optics FSV	GR	BR	XA	-4.00 to +4.00	-2.00
1.53 Trivex® FSV					
HOYA Phoenix™ FSV	GR			-6.00 to +2.00	-2.00
HOYA Phoenix™ FSV with Super HiVision EX3 AR	GR			-6.00 to +2.00	-2.00
X-Cel Optical Aris™ FSV	GR			Plano to -6.00 total power	-2.00
1.59 POLYCARBONATE FSV					
Essilor FSV Aspheric Airwear® Crizal AVANCE™	GR			-6.00 to +4.00	-2.00
Essilor FSV Spherical Airwear® Crizal AVANCE™	GR			-6.00 to +4.00	-2.00
Essilor FSV Aspheric Airwear® Crizal® Alizé®	GR			-6.00 to +4.00	-2.00
Essilor FSV Spherical Airwear® Crizal® Alizé®	GR			-6.00 to +4.00	-2.00
Essilor FSV Junior Aspheric Airwear with Trio EZ Blue	GR			-6.00 to +4.00	-2.00
Essilor FSV Aspheric Airwear®	GR			-6.00 to +4.00	-2.00
Essilor FSV Spherical Airwear®	GR			-6.00 to +4.00	-2.00
Seiko Diamond Clear Aspheric FSV	GR			-6.00 to +3.00	-2.00
Seiko Diamond Clear ECP AR Aspheric FSV	GR			-6.00 to +3.00	-2.00
Signet Armorlite Signia™ FSV Spherical Poly	GR			-6.00 to +4.00	-2.00
Signet Armorite Kodak Spherical FSV Poly Clean'N'CleAR	GR			-6.00 to +4.00	-2.00
SOLA FSV Aspheric Polycarbonate	GR			-4.00 to +2.00	-2.00
Younger FSV Poly AR Super Hydrophobic	GR	BR		-6.00 to +4.00	-2.00
Younger Optics FSV Poly	GR	BR	XA	-6.00 to +4.00	-2.00
1.60 HIGH INDEX FSV					
Younger Optics FSV 1.60 AR Super Hydrophobic	GR	BR		-6.00 to +6.00	-2.00
1.67 SUPER HIGH INDEX FSV					
Essilor FSV Aspheric Thin&Lite® 1.67 Crizal Avancé™	GR			-8.00 to +2.00	-2.00
Essilor FSV Aspheric Thin&Lite® 1.67 Crizal® Alizé®	GR			-8.00 to +2.00	-2.00
Essilor FSV Aspheric Thin&Lite® 1.67	GR			-8.00 to +2.00	-2.00
Seiko UltraThin 1.67 Aspheric FSV	GR			-8.00 to +2.00	-2.00
Signet Armorlite Signia™ FSV Aspherical 1.67***	GR			-8.00 to +2.00	-2.00
1.74 ULTRA HIGH INDEX FSV					
Seiko 1.74 Aspheric FSV	GR			-10.00 to +3.00	-2.00

SINGLE VISION

	Transitions® GRAY	Transitions® BROWN	Transitions® XTRActive® GRAY	Transitions® VANTAGE® GRAY	SPHERE	CYL TO
1.50 STANDARD INDEX SINGLE VISION						
Essilor Lens SFSV	GR	BR	XA	V	-10.00 to +6.00	-4.00
Essilor Wrap SV	GR	BR			-4.00 to +3.00	-4.00
HOYA Vison Care SFSV	GR	BR	XA	V	-10.00 to +6.00	-4.00
HOYA Vison Care SV IQ	GR	BR	XA	V	-10.00 to +6.00	-4.00
Seiko ULTRA 1.50 SF				V	-5.00 to +4.00	-4.00
Shamir Autograph® II - Single Vision®	GR	BR		V	-14.00 to +8.00	-6.00
Shamir Autograph® II - Single Vision Attitude™	GR	BR		V	-6.00 to +6.00	-6.00
Shamir Spectrum™ Single Vision	GR	BR		V	-14.00 to +8.00	-6.00
Shamir Relax™ Single Vision	GR	BR		V	-12.25 to +6.00	-6.00
Signet Armorlite SFSV	GR	BR		V*	-10.00 to +6.00	-4.00
Signetek Processed SV	GR	BR	XA	V*	-8.00 to +6.00	-4.00
SOLA SFSV	GR	BR			-9.00 to +5.00	-4.00
SOLA ATL HD Single Vision	GR	BR		V	-9.00 to +7.00	-4.00
CZV Spherical Single Vision				V	-9.00 to +5.00	-4.00
ZEISS Individual® SV	GR	BR		V	-6.00 to +6.00	-4.00
ZEISS 3D Single Vision	GR	BR		V	-6.00 to +6.00	-4.00
ZEISS Individual® SV Wrap	GR	BR			-4.00 to +4.00	-4.00
X-Cel SFSV	GR	BR			-10.00 to +8.50	-4.00
X-Cel SFSV			XA	V	-9.00 to +6.50	-4.00
Younger Optics SFSV	GR	BR	XA	V	-11.00 to +6.00	-4.00
1.53 Trivex® SINGLE VISION						
Augen High Definition™ SFSV	GR	BR			-10.25 to +6.00	-4.00
Essilor Lens SFSV TREXA™	GR	BR			-10.00 to +6.00	-4.00
ZEISS Individual® SV	GR				-7.00 to +7.00	-4.00
ZEISS 3D Single Vision	GR				-7.00 to +7.00	-4.00
ZEISS Individual® SV Wrap	GR				-4.00 to +4.00	-4.00
SOLA ATL HD Single Vision	GR				-11.00 to +7.00	-4.00
HOYA Phoenix™ SV	GR	BR	XA	V	-10.00 to +6.00	-6.00
HOYA Phoenix™ SV IQ	GR	BR	XA	V	-10.00 to +6.00	-6.00
Seiko ULTRA 1.53/Trivex SF				V	-7.00 to +5.00	-4.00
Shamir Autograph® II - Single Vision®	GR	BR	XA	V	-14.00 to +7.50	-6.50
Shamir Autograph® II - Single Vision Attitude™	GR	BR	XA	V	-6.00 to +6.00	-6.50
Shamir Spectrum™ Single Vision	GR	BR	XA	V	-14.00 to +7.50	-6.50
Shamir Relax™ Single Vision				V	-13.25 to +6.00	-6.50

*Available later this year. Check with your supplier for more details.

***Summer 2013

SINGLE VISION CONTINUED

	Transitions® GRAY	Transitions® BROWN	Transitions® XTRActive® GRAY	Transitions® VANTAGE® GRAY	SPHERE	CYL TO
1.53 Trivex® SINGLE VISION						
Signetek Processed SV	GR	BR	XA	V*	-10.00 to +6.00	-4.00
X-Cel Optical Aris™ SV Aspheric	GR	BR			-10.00 to +8.50	-4.00
X-Cel Optical Aris™ SV Spherical	GR	BR	XA	V	-9.00 to +6.50	-4.00
Younger Optics Trilogy® Aspheric SV	GR	BR			-9.00 to +6.00	-4.00
Younger Optics Trilogy® Spherical SV	GR	BR	XA	V	-8.00 to +7.50	-4.00
1.55 INDEX SINGLE VISION						
X-Cel High-X™ SFSV	GR	BR			-11.00 to +9.00	-4.00
1.59 POLYCARBONATE SINGLE VISION						
Essilor SFSV Aspheric Airwear®	GR	BR	XA		-10.00 to +8.00	-4.00
Essilor SV 360°™ Aspheric Airwear®	GR				-10.00 to +8.00	-4.00
Essilor SFSV Spherical Airwear®	GR	BR	XA		-10.00 to +6.00	-4.00
Essilor SFSV Spherical Airwear®				V	-10.00 to +8.00	-6.00
Essilor Wrap SV Airwear®	GR	BR			-4.00 to +3.00	-4.00
Essilor Eyecode™ Aspheric SV Airwear®	GR	BR			-10.00 to +8.00	-4.00
Essilor Fit™ SV Airwear®	GR	BR			-10.00 to +8.00	-4.00
Essilor Azio™ Aspheric SV Airwear®	GR	BR			-10.00 to +8.00	-4.00
HOYA Vison Care SV IQ Polycarbonate	GR	BR	XA		-10.00 to +6.50	-4.00
HOYA Vison Care SV IQ Polycarbonate				V	-8.00 to +6.00	-4.00
Optima RESOLUTION® ASFSV	GR				-8.00 to +8.00	-4.00
Seiko ULTRA 1.59 SF				V	-7.00 to +5.00	-4.00
Shamir Autograph® II Single Vision® Poly	GR	BR	XA		-14.00 to +8.00	-7.00
Shamir Autograph® II Single Vision® Poly				V	-14.00 to +8.75	-7.00
Shamir Autograph® II Single Vision Attitude™ Poly	GR	BR	XA		-6.00 to +6.00	-8.00
Shamir Autograph® II Single Vision Attitude™ Poly				V	-6.50 to +6.00	-7.00
Shamir Spectrum™ SFSV PolyPlus™	GR	BR	XA	V	-14.00 to +8.75	-7.00
Shamir Relax™ SFSV PolyPlus™	GR	BR		V	-14.00 to +8.00	-7.00
Shamir SFSV PolyPlus™	GR	BR	XA	V	-11.25 to +9.00	-6.00
Signet Armorlite Signia™ Poly SFSV	GR	BR			-10.00 to +6.00	-4.00
Signetek Processed SV	GR	BR	XA	V	-10.00 to +6.50	-4.00
CZV Spherical Single Vision				V*	-12.00 to +6.50	-4.00
ZEISS Individual® SV	GR	BR	XA	V*	-10.00 to +8.00	-4.00
ZEISS 3D Single Vision	GR	BR	XA	V*	-10.00 to +8.00	-4.00
ZEISS Individual® SV Wrap	GR	BR	XA		-4.00 to +4.00	-4.00
SOLA ASL® Aspheric SFSV	GR	BR	XA		-10.00 to +6.00	-4.00
SOLA Spazio™ Wrap SFSV	GR				-6.00 to +1.50	-2.00
SOLA ATL HD Single Vision	GR	BR	XA	V*	-9.00 to +8.00	-4.00
Younger Optics SV	GR	BR	XA	V	-9.00 to +6.00	-4.00
1.60 HIGH INDEX SINGLE VISION						
Essilor SFSV Aspheric Thin&Lite® 1.6	GR	BR			-13.00 to +10.00	-4.00
Essilor SFSV Spherical Thin&Lite® 1.6	GR	BR			-13.00 to +8.00	-4.00
HOYA SV 1.60	GR	BR			-10.00 to +8.00	-6.00
HOYA Aspheric 1.60	GR	BR			-10.00 to +6.00	-6.00
Shamir Autograph® II Single Vision® SuperLite™ 1.6	GR	BR			-14.00 to +8.00	-7.00
Shamir Autograph® II Single Vision Attitude™ SuperLite™ 1.6	GR	BR			-6.00 to +6.00	-7.00
Shamir Autograph® II Single Vision Attitude™ SuperLite™ 1.6	GR	BR			-6.50 to +6.00	-7.00
Shamir Spectrum™ SFSV SuperLite™ 1.6	GR	BR			-14.00 to +8.75	-7.00
Shamir Spectrum™ SFSV SuperLite™ 1.6	GR	BR			-14.00 to +8.00	-7.00
Shamir Relax™ SFSV SuperLite™ 1.6	GR	BR			-14.00 to +8.00	-7.00
Shamir SFSV SuperLite™ 1.6	GR	BR	XA		-11.25 to +9.00	-6.00
Signet Armorlite Signia™ 1.60 SFSV	GR	BR			-13.00 to +8.00	-4.00
Younger Optics Single Vision 1.60	GR	BR			-11.00 to +8.00	-4.00
1.67 SUPER HIGH INDEX SINGLE VISION						
Essilor SFSV Aspheric Thin&Lite® 1.67	GR	BR			-16.00 to +8.00	-4.00
Essilor SV 360°™ Spherical Thin&Lite® 1.67	GR	BR			-13.00 to +8.00	-4.00
Essilor Spherical Thin&Lite® 1.67	GR	BR			-13.00 to +9.00	-6.00
Essilor Wrap SV Thin&Lite® 1.67	GR	BR			-5.00 to +3.00	-4.00
Essilor Eyecode™ SV Thin&Lite® 1.67	GR	BR			-13.00 to +8.00	-4.00
Essilor Fit™ SV Thin&Lite® 1.67	GR	BR			-13.00 to +8.00	-4.00
Essilor Azio™ SV Thin&Lite® 1.67	GR	BR			-13.00 to +8.00	-4.00
HOYA Aspheric SV 1.67	GR	BR			-15.00 to +10.00	-6.00
Seiko 1.67 Super SV Spherical	GR	BR			-15.00 to +10.00	-6.00
Shamir Autograph® II Single Vision® SuperLite™ 1.67	GR	BR	XA		-16.00 to +16.00	-8.00
Shamir Autograph® II Single Vision Attitude™ SuperLite™ 1.67	GR	BR	XA		-6.00 to +6.00	-8.00
Shamir Spectrum™ SFSV SuperLite™ 1.67	GR	BR	XA		-16.00 to +16.00	-6.00
Shamir SFSV 1.67	GR	BR	XA		-12.00 to +16.25	-6.00
Signet Armorlite Signia™ 1.67 Aspheric SFSV	GR	BR			-16.00 to +8.00	-4.00
Signet Armorlite Signia™ 1.67 SFSV	GR	BR			-11.75 to +7.00	-4.00
Signetek Processed SV	GR	BR			-12.00 to +7.50	-4.00
CZV Spherical Single Vision	GR				-13.50 to +7.50	-4.00
ZEISS Individual® SV	GR	BR	XA		-12.00 to +8.00	-4.00
ZEISS 3D Single Vision	GR	BR	XA		-12.00 to +8.00	-4.00
SOLA ASL 1.67 Aspheric SFSV	GR		XA		-12.00 to +8.00	-6.00
SOLA ATL Single Vision	GR	BR	XA		-12.00 to +8.00	-4.00
Younger Optics Single Vision 1.67 MR-10™	GR	BR	XA		-12.00 to +7.00	-4.00
1.74 ULTRA HIGH INDEX SINGLE VISION						
Essilor Eyecode™ SV Thin&Lite 1.74	GR				-14.00 to +10.00	-6.00
Essilor Fit™ SV Thin&Lite 1.74	GR				-14.00 to +10.00	-6.00
Essilor 360™ SV Thin&Lite 1.74	GR				-14.00 to +10.00	-6.00
Nikon Aspheric SV 1.74	GR	BR			-12.00 to +4.00	-4.00
Signetek™ Processed SV	GR	BR			-13.00 to +9.00	-4.00

*Available later this year. Check with your supplier for more details.

MULTIFOCAL

MULTIFOCAL		Transitions® GRAY	Transitions® BROWN	Transitions® XTRActive® GRAY	SPHERE	ADD	CYL TO
1.50 STANDARD INDEX SF RD-24	X-Cel RD-24	GR	BR		-4.00 to +5.50	+0.75 to +4.00	-4.00
1.50 STANDARD INDEX SFST-28	Essilor Lenses SFST 28	GR	BR	XA	-10.00 to +6.00	+1.00 to +3.00	-4.00
	Essilor Bifocal AB ST 28	GR	BR		-10.00 to +6.00	+1.00 to +3.00	-4.00
	HOYA Vision Care SFFT 28	GR	BR	XA	-10.00 to +6.00	+1.00 to +3.00	-4.00
	SOLA SFFT 28	GR	BR		-9.00 to +5.00	+1.00 to +3.00	-4.00
	Signet Armorlite SFFT-28	GR	BR		-10.00 to +6.00	+1.00 to +3.00	-4.00
	X-Cel D-28 76 mm	GR	BR	XA	-9.00 to +5.00	+1.00 to +3.00	-4.00
	X-Cel D-28 78mm	GR	BR	XA	-4.00 to +4.25	+0.75 to +5.00	-4.00
	X-Cel D-28 72 mm	GR			4.50 to +7.00	+1.00 to +3.00	-4.00
	Younger Optics SFFT 28	GR	BR	XA	-11.00 to +6.00	+1.00 to +3.00	-4.00
	Younger Optics SFFT 28 HI-ADDS	GR	BR		-11.00 to +6.00	+3.25 to +4.00	-4.00
1.50 STANDARD INDEX SFST-35	Signetek™ Processed Signia™ 1.50 FT-35	GR	BR		-6.00 to +6.00	+0.75 to +4.00	-4.00
	SOLA SFFT 35	GR			-4.00 to +6.00	+1.00 to +3.00	-4.00
	X-Cel D-35	GR	BR		-7.00 to +5.50	+0.75 to +4.00	-4.00
	Younger Optics SFFT 35	GR	BR		-5.00 to +6.00	+0.75 to +4.00	-4.00
1.50 STANDARD INDEX TRIFOCALS	Essilor Lens SF 7x28	GR			-5.00 to +6.00	+1.50 to +3.00	-4.00
	SOLA SF 7x28	GR			-4.00 to +6.00	+1.50 to +3.00	-4.00
	Signet Armorlight Signia™ 1.50 HC SF Trifocal 7x28	GR			-5.00 to +6.00	+1.50 to +3.00	-4.00
	X-Cel 7x28	GR	BR		-7.00 to +7.50	+1.50 to +4.00	-4.00
	X-Cel 8x35	GR	BR		-7.00 to +5.00	+1.50 to +4.00	-4.00
	Younger Optics SF 7x28	GR	BR	XA	-5.00 to +6.00	+1.50 to +4.00	-4.00
	Younger Optics 8x35	GR			-5.00 to +6.00	+1.50 to +4.00	-4.00
1.53 Trivex® SF RD-24	X-Cel Optical Aris™ RD-24	GR			-4.50 to +5.50	+1.00 to +3.00	-4.00
1.53 Trivex® SFST-28	Augen Flat-Top 28	GR			-8.50 to +4.75	+1.00 to +3.00	-4.00
	HOYA Phoenix™ SFFT 28	GR	BR		-10.00 to +6.00	+1.00 to +3.50	-4.00
	X-Cel Optical Aris™ D-28	GR	BR		-8.00 to +5.50	+1.00 to +4.00	-4.00
	Younger Optics Trilogy® FLAT TOP 28	GR	BR		-9.00 to +7.00	+1.00 to +3.50	-4.00
1.53 Trivex® DD-28	X-Cel Optical Aris™ DD-28	GR			-4.50 to +2.00	+1.00 to +3.00	-4.00
1.53 Trivex® SFST-35	X-Cel Optical Aris™ D-35	GR			-4.50 to +5.50	+1.00 to +3.00	-4.00
1.53 Trivex® TRIFOCALS	X-Cel Optical Aris™ 7x28	GR	BR		-8.00 to +5.50	+1.50 to +3.50	-4.00
	X-Cel Optical Aris™ 8x35	GR			-4.50 to +5.50	+1.50 to +3.00	-4.00
1.55 INDEX SFST-28	X-Cel High-X™ D-28	GR	BR		-11.00 to +9.00	+1.00 to +4.00	-4.00
1.55 INDEX SFST-35	X-Cel High-X™ D-35	GR	BR		-12.00 to +9.00	+1.00 to +4.00	-4.00
1.55 INDEX TRIFOCALS	X-Cel High-X™ 7x28	GR	BR		-8.00 to +9.00	+1.50 to +4.00	-4.00
	X-Cel High-X™ 8x35	GR			-7.00 to +5.50	+1.50 to +3.00	-4.00

PROGRESSIVES 1.50

1.50 STANDARD INDEX PAL		Transitions® GRAY	Transitions® BROWN	Transitions® XTRActive® GRAY	Transitions® VANTAGE® GRAY	SPHERE	ADD	CYL TO	MFH**
CARL ZEISS VISION	AO Easy	GR	BR			-9.00 to +5.00	+0.75 to +3.50	-4.00	18
	AO Compact®	GR				-9.00 to +5.00	+1.00 to +3.50	-4.00	15
	Instinctive	GR	BR			-9.00 to +5.00	+0.75 to +3.50	-4.00	17
	SOLA Compact Ultra™	GR	BR			-9.00 to +5.00	+0.75 to +3.50	-6.00	13
	SOLAOne®	GR	BR			-9.00 to +6.00	+1.00 to +3.50	-4.00	18
	SOLAMAX®	GR				-8.00 to +4.00	+1.00 to +3.00	-4.00	18
	SOLA VIP	GR	BR			-7.00 to +6.00	+1.00 to +3.00	-4.00	18
	ZEISS GT2®	GR	BR			-10.00 to +6.00	+1.00 to +3.00	-6.00	17
	ZEISS GT2® Short	GR	BR			-10.00 to +6.00	+0.75 to +3.50	-6.00	13
	ZEISS Gradal® Top	GR	BR			-6.00 to +6.00	+1.00 to +3.00	-4.00	18
	ZEISS Individual® 2	GR	BR		V	-7.00 to +5.00	+0.75 to +3.50	-4.00	13-35
	ZEISS Individual® 2I	GR	BR		V	-7.00 to +5.00	+0.75 to +3.50	-4.00	13-35
	ZEISS Individual® 2N	GR	BR		V	-7.00 to +5.00	+0.75 to +3.50	-4.00	13-35
	ZEISS Individual®	GR	BR			-6.00 to +5.00	+0.75 to +3.00	-4.00	13-35
	ZEISS Individual® Wrap	GR	BR			-4.00 to +4.00	+1.00 to +2.50	-4.00	13-35
	ZEISS Choice Plus V	GR	BR		V	-7.00 to +5.00	+0.75 to +3.50	-4.00	13-35
	ZEISS Choice Plus 13	GR	BR		V	-7.00 to +5.00	+0.75 to +3.50	-4.00	13
	ZEISS Choice Plus 15	GR	BR		V	-7.00 to +5.00	+0.75 to +3.50	-4.00	15
	ZEISS Choice Plus 17	GR	BR		V	-7.00 to +5.00	+0.75 to +3.50	-4.00	17
	ZEISS Choice Plus 19	GR	BR		V	-7.00 to +5.00	+0.75 to +3.50	-4.00	19
	ZEISS Choice Plus 21	GR	BR		V	-7.00 to +5.00	+0.75 to +3.50	-4.00	21
	ZEISS Choice 13	GR	BR		V	-7.00 to +5.00	+0.75 to +3.50	-4.00	13
	ZEISS Choice 15	GR	BR		V	-7.00 to +5.00	+0.75 to +3.50	-4.00	15
	ZEISS Choice 17	GR	BR		V	-7.00 to +5.00	+0.75 to +3.50	-4.00	17
	ZEISS Choice 19	GR	BR		V	-7.00 to +5.00	+0.75 to +3.50	-4.00	19
	ZEISS GT2® 3DV	GR	BR		V	-6.00 to +5.00	+0.75 to +3.00	-4.00	13-35
	ZEISS GT2® 3D	GR	BR		V	-6.00 to +5.00	+0.75 to +3.00	-4.00	17
	ZEISS GT2® 3DV Wrap	GR	BR			-4.00 to +4.00	+1.00 to +2.50	-4.00	13-35
	ZEISS GT2® 3D Wrap	GR	BR			-4.00 to +4.00	+1.00 to +2.50	-4.00	17
	ZEISS GT2® 3D Short	GR	BR		V	-6.00 to +5.00	+0.75 to +3.00	-4.00	13
	ZEISS GT2® 3D Short Wrap	GR	BR			-4.00 to +4.00	+1.00 to +2.50	-4.00	13
	AO Easy HD	GR	BR		V	-9.00 to +4.50	+0.75 to +3.50	-4.00	18
	SOLA HDV	GR	BR		V	-9.00 to +4.50	+0.75 to +3.50	-4.00	13-35
	SOLA Instinctive Performance 13	GR	BR		V	-9.00 to +4.50	+0.75 to +3.50	-4.00	13
	SOLA Instinctive Performance 15	GR	BR		V	-9.00 to +4.50	+0.75 to +3.50	-4.00	15
	SOLA Instinctive Performance 17	GR	BR		V	-9.00 to +4.50	+0.75 to +3.50	-4.00	17
	SOLA Compact Ultra™ HD	GR	BR		V	-9.00 to +4.50	+0.75 to +3.50	-4.00	13
	SOLAOne® HD	GR	BR		V	-9.00 to +4.50	+0.75 to +3.50	-4.00	18

PROGRESSIVES 1.50 CONTINUED

1.50 STANDARD INDEX PAL		Transitions® GRAY	Transitions® BROWN	Transitions® XTRActive® GRAY	Transitions® VANTAGE® GRAY	SPHERE	ADD	CYL TO	MFH**
ESSILOR LENS	Varilux S Design™	GR	BR			-10.00 to +6.00	+0.75 to +4.00	-6.00	17
	Varilux S Design Short™	GR	BR			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
	Varilux S Fit™	GR	BR			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
	Varilux S 4D™	GR	BR			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
	Varilux Comfort®	GR	BR	XA		-10.00 to +7.00	+0.75 to +3.50	-4.00	17
	Varilux Comfort Enhanced™	GR	BR	XA		-10.00 to +7.00	+0.75 to +3.50 +075 to +3.00 for MFH below 17	-6.00	14
	Varilux Comfort DRx™	GR	BR	XA	V	-8.00 to +6.00	+0.75 to +3.50	-6.00	17
	Varilux Comfort Short™	GR	BR	XA		-10.00 to +7.00	+0.75 to +3.00	-4.00	14
	Varilux Comfort Short DRx™	GR	BR	XA	V	-8.00 to +6.00	+0.75 to +3.00	-6.00	14
	Varilux® Ellipse®	GR				-10.00 to +6.00	+0.75 to +3.00	-4.00	14
	Varilux Ellipse 360°™	GR				-10.00 to +6.00	+0.75 to +3.00	-6.00	13
	Varilux Ipseo IV™	GR				-10.00 to +6.00	+0.75 to +3.50	-6.00	14
	Varilux Ipseo IV eyecode™	GR				-10.00 to +6.00	+0.75 to +3.50	-6.00	customized
	Varilux Physio™	GR	BR			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-4.00	17
	Varilux Physio 360°™	GR	BR			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	17
	Varilux Physio Enhanced™	GR	BR			-10.00 to +6.00	+0.75 to +3.50 & +4.00 +0.75 to +3.00 for MFH below 17	-6.00	14
	Varilux Physio Enhanced eyecode™	GR	BR			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
	Varilux Physio Enhanced FIT™	GR	BR			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
	Varilux Physio Enhanced AZIO™	GR	BR			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
	Varilux Physio Enhanced INDIA™	GR	BR			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
	Varilux Physio DRx™	GR	BR	XA	V	-8.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	17
	Varilux Physio Short™	GR	BR			-10.00 to +6.00	+0.75 to +3.00	-4.00	14
	Varilux Physio Short DRx™	GR	BR	XA	V	-8.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
	DEFINITY®	GR	BR			-10.00 to +6.00	+1.00 to +3.00	-4.00	18
	DEFINITY SHORT™	GR	BR			-10.00 to +6.00	+1.00 to +3.00	-4.00	15
	Essilor Accolade™	GR	BR			-10.00 to +6.00	+0.75 to +3.50	-4.00	17
	Essilor Accolade™ Freedom™	GR	BR			-9.00 to +6.00	+0.75 to +3.50	-4.00	15
	Essilor Ideal® Advanced	GR	BR	XA	V	-8.00 to +6.00	+0.75 to +3.50	-6.00	17
	Essilor Ideal® Advanced Wrap	GR	BR			-4.00 to +3.00	+0.75 to +3.50	-3.50	customized
	Essilor Ideal®	GR	BR	XA	V	-8.00 to +6.00	+0.75 to +3.50	-6.00	17
	Essilor Ideal® Short	GR	BR	XA	V	-8.00 to +6.00	+0.75 to +3.00	-6.00	14
	Adaptar™	GR				-10.00 to +7.00	+0.75 to +3.50	-4.00	18
	Ovation™	GR	BR			-10.00 to +6.00	+0.75 to +3.50	-4.00	17
	Essilor Natural™	GR				-10.00 to +6.00	+0.75 to +3.50	-4.00	18
	Essilor SmallFit™	GR				-10.00 to +6.00	+0.75 to +3.00	-4.00	14
HOYA VISION CARE	HOYALUX Summit CD®	GR	BR	XA		-10.00 to +6.00	+1.00 to +3.50	-4.00	14
	HOYALUX Summit CD® IQ	GR	BR	XA		-10.00 to +6.00	+1.00 to +3.50	-4.00	14
	HOYALUX Summit ECP®	GR	BR	XA		-10.00 to +6.00	+1.00 to +3.50	-4.00	18
	HOYALUX Summit ECP® IQ	GR	BR	XA		-10.00 to +6.00	+1.00 to +3.50	-4.00	18
	HOYALUX GP Wide ™	GR	BR			-10.00 to +6.00	+1.00 to +3.50	-4.00	18
	HOYALUX iD LifeStyle®	GR	BR	XA		-8.00 to +6.00	+0.75 to +3.50	-4.00	18
	HOYALUX iD LifeStyle cd	GR	BR	XA		-8.00 to +6.00	+0.75 to +3.50	-4.00	14
NIKON OPTICAL	Nikon Presio i Short 1.50	GR				-7.00 to +5.00	+0.75 to +3.50	-4.00	16
	Nikon W Short 1.50	GR	BR			-7.00 to +5.00	+0.75 to +3.50	-4.00	15
	Nikon Presio Power 1.50	GR	BR			-6.00 to +5.00	+0.75 to +3.50	-4.00	15, 16, 17
	Nikon Presio i Digital 1.50	GR	BR	XA	V	-7.00 to +4.00	+0.75 to +4.00	-6.00	13, 15, 17
	Nikon Move 1.50	GR	BR	XA		-7.00 to +4.00	+0.75 to +3.50	-4.00	16, 18
SEIKO OPTICAL PRODUCTS	1.50 Prominent	GR	BR		V	-5.00 to +4.00	+0.50 to +3.50	-5.00	14
	1.50 Premier 2	GR	BR		V	-5.00 to +4.00	+0.50 to +3.50	-4.00	14
	1.50 Prestige II	GR	BR		V	-5.00 to +4.00	+0.50 to +3.50	-5.00	14
	1.50 Prestige II Wide	GR	BR		V	-5.00 to +4.00	+0.50 to +3.50	-5.00	14
	1.50 Succeed	GR	BR		V	-5.00 to +4.00	+0.50 to +3.50	-4.00	17
	1.50 Succeed Ws	GR	BR		V	-5.00 to +4.00	+0.50 to +3.50	-4.00	15
	1.50 Supercede 2	GR	BR		V	-5.00 to +4.00	+0.50 to +3.50	-4.00	14
	1.50 Supernal	GR	BR		V	-5.00 to +4.00	+0.50 to +3.50	-5.00	14
	1.50 Surmount	GR	BR		V	-5.00 to +4.00	+0.50 to +3.50	-5.00	14
	1.50 Surmount Ws	GR	BR		V	-5.00 to +4.00	+0.50 to +3.50	-5.00	14
SHAMIR INSIGHT INC.	Shamir Autograph® II	GR	BR		V	-12.25 to +6.00	+0.75 to +3.50	-6.00	11, 13, 15, 18
	Shamir Autograph® II Attitude®	GR	BR		V	-5.00 to +4.00	+0.75 to +3.50	-6.00	15, 18
	Shamir Spectrum™	GR	BR		V	-12.25 to +6.00	+0.75 to +3.50	-6.00	14, 16, 18
	Shamir Element™	GR	BR		V	-7.00 to +4.00	+0.75 to +3.50	-6.00	16, 19
	Shamir FirstPal				V	-12.25 to +6.00	+0.75 to +1.50	-6.00	15
	Shamir Creation®	GR	BR			-8.50 to +6.75	+0.75 to +3.50	-6.00	19
	Shamir Genesis™	GR				-8.50 to +5.75	+0.75 to +3.50	-6.00	19
	Shamir Piccolo®	GR				-8.50 to +6.00	+0.75 to +3.50	-6.00	16
SIGNET ARMORLITE	KODAK Concise® Signia™ 1.50	GR	BR			-9.00 to +5.00	+1.00 to +3.00	-4.00	14
	KODAK Precise® Signia™ 1.50	GR	BR			-9.00 to +6.00	+1.00 to +3.00	-4.00	17
	KODAK Precise® Short Signia™ 1.50	GR	BR			-8.50 to +5.25	+1.00 to +3.00	-4.00	13
	KODAK Concise® Digital Signia™ 1.50	GR	BR			-8.00 to +5.00	+1.00 to +3.00	-4.00	14
	KODAK Precise® Digital Signia™ 1.50	GR	BR			-8.00 to +5.00	+1.00 to +3.00	-4.00	17
	KODAK Precise® Short Digital Signia™ 1.50	GR	BR			-8.00 to +5.00	+1.00 to +3.00	-4.00	13
	KODAK Precise® PB Signia™ 1.50	GR	BR			-8.00 to +4.00	+0.75 to +3.50	-4.00	17
	KODAK Precise® Short PB Signia™ 1.50	GR	BR			-8.00 to +4.00	+0.75 to +3.50	-4.00	13
	KODAK Unique™ Signia™ 1.50	GR	BR	XA	V	-8.00 to +4.00	+0.75 to +3.50	-4.00	13
	KODAK Unique DS™ Signia™ 1.50	GR	BR			-6.00 to +4.00	+0.75 to +3.50	-4.00	13
X-CEL OPTICAL	X-Cel Freedom ID™	GR				-4.50 to +5.50	+1.00 to +3.50	-4.00	17
YOUNGER OPTICS	Younger Image®	GR	BR	XA		-8.00 to +6.00	+1.00 to +3.00	-4.00	18

**MFH-Minimum Fitting Height

1.53 TRIVEX® PAL		Transitions® GRAY	Transitions® BROWN	Transitions® XTRActive® GRAY	Transitions® VANTAGE® GRAY	SPHERE	ADD	CYL TO	MFH**
AUGEN	Augen High Definition® Trinity Progressive 13/17	GR	BR			-10.00 to +6.00	+1.00 to +3.00	-4.00	17
	Augen High Definition® Trinity Short Progressive 8/12	GR	BR			-10.00 to +4.75	+1.00 to +3.00	-4.00	12
	Augen HDRx™ Trinity 17	GR	BR			-7.50 to +6.75	+0.75 to +3.50	-4.00	17
	Augen HDRx™ Trinity 12	GR	BR			-7.50 to +6.75	+0.75 to +3.50	-4.00	12
	Augen HDRx™ Trinity Spacia	GR	BR			-7.50 to +6.75	+0.75 to +3.50	-4.00	22
CARL ZEISS VISION	ZEISS Individual® 2	GR				-7.00 to +5.00	+0.75 to +3.50	-4.00	13-35
	ZEISS Individual® 2I	GR				-7.00 to +5.00	+0.75 to +3.50	-4.00	13-35
	ZEISS Individual® 2N	GR				-7.00 to +5.00	+0.75 to +3.50	-4.00	13-35
	ZEISS Individual®	GR				-7.00 to +5.00	+0.75 to +3.50	-4.00	13-35
	ZEISS Individual® Wrap	GR				-4.00 to +4.00	+0.75 to +3.50	-4.00	13-35
	ZEISS Choice Plus V	GR				-7.00 to +5.00	+0.75 to +3.50	-4.00	13-35
	ZEISS Choice Plus 13	GR				-7.00 to +5.00	+0.75 to +3.50	-4.00	13
	ZEISS Choice Plus 15	GR				-7.00 to +5.00	+0.75 to +3.50	-4.00	15
	ZEISS Choice Plus 17	GR				-7.00 to +5.00	+0.75 to +3.50	-4.00	17
	ZEISS Choice Plus 19	GR				-7.00 to +5.00	+0.75 to +3.50	-4.00	19
	ZEISS Choice Plus 21	GR				-7.00 to +5.00	+0.75 to +3.50	-4.00	21
	ZEISS Choice 13	GR				-7.00 to +5.00	+0.75 to +3.50	-4.00	13
	ZEISS Choice 15	GR				-7.00 to +5.00	+0.75 to +3.50	-4.00	15
	ZEISS Choice 17	GR				-7.00 to +5.00	+0.75 to +3.50	-4.00	17
	ZEISS Choice 19	GR				-7.00 to +5.00	+0.75 to +3.50	-4.00	19
	ZEISS GT2® 3DV	GR				-7.00 to +5.00	+0.75 to +3.50	-4.00	13-35
	ZEISS GT2® 3D	GR				-7.00 to +5.00	+0.75 to +3.50	-4.00	17
	ZEISS GT2® 3D Wrap	GR				-4.00 to +4.00	+1.00 to +2.50	-4.00	17
	ZEISS GT2® 3D Short	GR				-7.00 to +5.00	+0.75 to +3.50	-4.00	13
	ZEISS GT2® 3D Short Wrap	GR				-4.00 to +4.00	+1.00 to +2.50	-4.00	13
	AO Easy HD	GR				-9.00 to +5.00	+0.75 to +3.50	-4.00	18
	SOLA HDV	GR				-9.00 to +5.00	+0.75 to +3.50	-4.00	13-35
	SOLA Instinctive Performance 13	GR				-9.00 to +5.00	+0.75 to +3.50	-4.00	13
	SOLA Instinctive Performance 15	GR				-9.00 to +5.00	+0.75 to +3.50	-4.00	15
	SOLA Instinctive Performance 17	GR				-9.00 to +5.00	+0.75 to +3.50	-4.00	17
	SOLA Compact Ultra™ HD	GR				-9.00 to +5.00	+0.75 to +3.50	-4.00	13
	SOLAOne® HD	GR				-9.00 to +5.00	+0.75 to +3.50	-4.00	18
ESSILOR LENS	Varilux Comfort DRx™ TREXA™	GR	BR			-10.00 to +6.00	+0.75 to +3.50	-6.00	17
	Varilux Comfort Short DRx™ TREXA™	GR	BR			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
	Varilux®Physio® TREXA™	GR				-10.00 to +6.00	+1.00 to +3.50 & 4.00	-4.00	17
	Varilux Physio Enhanced™ TREXA™	GR				-10.00 to +6.00 +0.75 to +3.00 for MFH below 17	+1.00 to +3.50 & 4.00	-6.00	14
	Varilux Physio Enhanced eyecode™ TREXA™	GR				-10.00 to +6.00	+1.00 to +3.50 & 4.00	-6.00	14
	Varilux Physio Enhanced FIT™ TREXA™	GR				-10.00 to +6.00	+1.00 to +3.50 & 4.00	-6.00	14
	Varilux Physio Enhanced AZIO™ TREXA™	GR				-10.00 to +6.00	+1.00 to +3.50 & 4.00	-6.00	14
	Varilux Physio Enhanced INDIA™ TREXA™	GR				-10.00 to +6.00	+1.00 to +3.50 & 4.00	-6.00	14
	Varilux Physio Short™ TREXA™	GR				-10.00 to +6.00	+0.75 to +3.00	-4.00	14
	DEFINITY® TREXA™	GR	BR			-10.00 to +6.00	+1.00 to +3.50 & 4.00	-4.00	18
	DEFINITY Short™ TREXA™	GR	BR			-10.00 to +6.00	+1.00 to +3.50 & 4.00	-4.00	15
	Essilor Accolade™ TREXA™	GR	BR			-10.00 to +6.00	+0.75 to 3.50	-4.00	17
	Essilor Accolade Freedom™ TREXA™	GR	BR			-10.00 to +6.00	+0.75 to 3.50	-4.00	15
	Essilor Ideal® Advanced TREXA™	GR	BR			-10.00 to +6.00	+0.75 to +3.50	-6.00	17
	Essilor Ideal® TREXA™	GR	BR			-10.00 to +6.00	+0.75 to +3.50	-6.00	17
	Essilor Ideal® Short TREXA™	GR	BR			-10.00 to +6.00	+0.75 to +3.00	-6.00	14
HOYA VISION CARE	HOYALUX Summit CD Phoenix	GR	BR	XA	V	-10.00 to +6.00	+1.00 to +3.00	-4.00	14
	HOYALUX Summit CD IQ Phoenix	GR	BR	XA	V	-10.00 to +6.00	+1.00 to +3.00	-4.00	14
	HOYALUX Summit ECP Phoenix	GR	BR	XA	V	-10.00 to +6.00	+1.00 to +3.00	-4.00	18
	HOYALUX Summit ECP IQ Phoenix	GR	BR	XA	V	-10.00 to +6.00	+1.00 to +3.00	-4.00	18
	HOYALUX GP Wide Phoenix	GR	BR			-10.00 to +6.00	+1.00 to +3.00	-4.00	18
	HOYALUX iD LifeStyle Phoenix	GR	BR	XA	V	-8.00 to +6.00	+0.75 to 3.50	-4.00	18
	HOYALUX iD LifeStyle cd Phoenix	GR	BR	XA	V	-8.00 to +6.00	+0.75 to 3.50	-4.00	14
SEIKO OPTICAL PRODUCTS	1.53/Trivex Prominent	GR	BR	XA	V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
	1.53/Trivex Premier 2	GR	BR	XA	V	-7.00 to +5.00	+0.50 to +3.50	-4.00	14
	1.53/Trivex Prestige II	GR	BR	XA	V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
	1.53/Trivex Prestige II Wide	GR	BR	XA	V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
	1.53/Trivex Succeed	GR	BR	XA	V	-7.00 to +5.00	+0.50 to +3.50	-4.00	17
	1.53/Trivex Succeed Ws	GR	BR	XA	V	-7.00 to +5.00	+0.50 to +3.50	-4.00	15
	1.53/Trivex Supercede 2	GR	BR	XA	V	-7.00 to +5.00	+0.50 to +3.50	-4.00	14
	1.53/Trivex Supernal	GR	BR	XA	V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
	1.53/Trivex Surmount	GR	BR	XA	V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
	1.53/Trivex Surmount Ws	GR	BR	XA	V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
SHAMIR INSIGHT INC.	Shamir Autograph® II	GR	BR	XA	V	-13.25 to +6.00	+0.75 to +3.50	-6.50	11, 13, 15, 18
	Shamir Autograph® II Attitude®	GR	BR	XA	V	-5.00 to +4.00	+0.75 to +3.50	-6.50	15, 18
	Shamir Spectrum™	GR	BR	XA	V	-13.25 to +6.00	+0.75 to +3.50	-6.50	14, 16, 18
	Shamir Element™	GR	BR	XA	V	-13.25 to +6.00	+0.75 to +3.50	-6.50	16, 19
SIGNET ARMORLITE	KODAK Precise® PB Trivex®	GR	BR		V	-10.00 to +4.50	+0.75 to +3.50	-4.00	17
	KODAK Precise® Short PB Trivex®	GR	BR		V	-10.00 to +4.50	+0.75 to +3.50	-4.00	13
	KODAK Unique Trivex®	GR	BR	XA	V	-10.00 to +4.50	+0.75 to +3.50	-4.00	13
X-CEL OPTICAL	X-Cel Optical Aris™ Freedom ID™	GR		XA	V	-4.50 to +5.50	+1.00 to +3.00	-4.00	17
YOUNGER OPTICS	Younger Trilogy® Image®	GR		XA		-8.00 to +7.00	+1.00 to +3.00	-4.00	18

PROGRESSIVES 1.59 CONTINUED

1.59 POLYCARBONATE PAL		Transitions® GRAY	Transitions® BROWN	Transitions® XTRActive® GRAY	Transitions® VANTAGE® GRAY	SPHERE	ADD	CYL TO	MFH**
CARL ZEISS VISION	AO Easy	GR				-9.00 to +5.00	+1.00 to +3.00	-4.00	18
	Instinctive	GR				-9.00 to +5.00	+1.00 to +3.00	-4.00	17
	SOLA Compact Ultra™	GR	BR			-9.00 to +5.00	+1.00 to +3.00	-6.00	13
	SOLAOne®	GR				-9.00 to +6.00	+1.00 to +3.00	-4.00	18
	SOLAMAX®	GR				-10.00 to +6.00	+1.00 to +3.00	-4.00	18
	SOLA VIP/VIPGold®	GR				-10.00 to +4.50	+1.00 to +3.00	-4.00	18
	ZEISS GT2®	GR	BR			-10.00 to +6.00	+1.00 to +3.00	-4.00	17
	ZEISS GT2® Short	GR	BR			-10.00 to +6.00	+1.00 to +3.00	-4.00	13
	ZEISS Gradal®Top	GR				-6.00 to +4.50	+1.00 to +3.00	-4.00	18
	ZEISS Individual® 2	GR	BR	XA	V*	-10.00 to +6.00	+0.75 to +3.50	-4.00	13-35
	ZEISS Individual® 2I	GR	BR	XA	V*	-10.00 to +6.00	+0.75 to +3.50	-4.00	13-35
	ZEISS Individual® 2N	GR	BR	XA	V*	-10.00 to +6.00	+0.75 to +3.50	-4.00	13-35
	ZEISS Individual®	GR	BR	XA		-10.00 to +6.00	+0.75 to +3.50	-4.00	13-35
	ZEISS Individual® Wrap	GR	BR			-4.00 to +4.00	+1.00 to +2.50	-4.00	13-35
	ZEISS Choice Plus V	GR	BR	XA	V*	-10.00 to +6.00	+0.75 to +3.50	-4.00	13-35
	ZEISS Choice Plus 13	GR	BR	XA	V*	-10.00 to +6.00	+0.75 to +3.50	-4.00	13
	ZEISS Choice Plus 15	GR	BR	XA	V*	-10.00 to +6.00	+0.75 to +3.50	-4.00	15
	ZEISS Choice Plus 17	GR	BR	XA	V*	-10.00 to +6.00	+0.75 to +3.50	-4.00	17
	ZEISS Choice Plus 19	GR	BR	XA	V*	-10.00 to +6.00	+0.75 to +3.50	-4.00	19
	ZEISS Choice Plus 21	GR	BR	XA	V*	-10.00 to +6.00	+0.75 to +3.50	-4.00	21
	ZEISS Choice 13	GR	BR	XA	V*	-10.00 to +6.00	+0.75 to +3.50	-4.00	13
	ZEISS Choice 15	GR	BR	XA	V*	-10.00 to +6.00	+0.75 to +3.50	-4.00	15
	ZEISS Choice 17	GR	BR	XA	V*	-10.00 to +6.00	+0.75 to +3.50	-4.00	17
	ZEISS Choice 19	GR	BR	XA	V*	-10.00 to +6.00	+0.75 to +3.50	-4.00	19
	ZEISS GT2® 3DV	GR	BR	XA	V*	-10.00 to +6.00	+0.75 to +3.50	-4.00	13-35
	ZEISS GT2® 3D	GR	BR	XA	V*	-10.00 to +6.00	+0.75 to +3.50	-4.00	17
	ZEISS GT2® 3D Wrap	GR	BR			-4.00 to +4.00	+1.00 to +2.50	-4.00	17
	ZEISS GT2® 3D Short	GR	BR	XA	V*	-10.00 to +6.00	+0.75 to +3.50	-4.00	13
	ZEISS GT2® 3D Short Wrap	GR	BR			-4.00 to +4.00	+1.00 to +2.50	-4.00	17
	SOLA HDV	GR	BR	XA	V*	-9.00 to +6.00	+0.75 to +3.50	-4.00	13-35
	SOLA Instinctive Performance 13	GR	BR	XA	V*	-9.00 to +6.00	+0.75 to +3.50	-4.00	13
	SOLA Instinctive Performance 15	GR	BR	XA	V*	-9.00 to +6.00	+0.75 to +3.50	-4.00	15
	SOLA Instinctive Performance 17	GR	BR	XA	V*	-9.00 to +6.00	+0.75 to +3.50	-4.00	17
	SOLA Compact Ultra™ HD	GR	BR	XA	V*	-9.00 to +6.00	+0.75 to +3.50	-4.00	13
	SOLAOne® HD	GR	BR	XA	V*	-9.00 to +6.00	+0.75 to +3.50	-4.00	18
	AO Easy HD	GR	BR	XA	V*	-9.00 to +6.00	+0.75 to +3.50	-4.00	18
ESSILOR LENS	Varilux S Design™ Airwear®	GR	BR			-10.00 to +6.00	+0.75 to +4.00	-6.00	17
	Varilux S Design Short™ Airwear®	GR	BR			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
	Varilux S Fit™ Airwear®	GR	BR			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
	Varilux S 4D™ Airwear®	GR	BR			-10.00 to +6.00	+0.75 to +4.00	-6.00	14
	Varilux Comfort® Airwear®	GR	BR	XA		-10.00 to +6.00	+0.75 to +3.50	-4.00	17
	Varilux Comfort Enhanced™ Airwear®	GR	BR	XA		-10.00 to +6.00	+0.75 to +3.50 & +4.00 +0.75 to +3.00 for MFH below 17	-6.00	14
	Varilux Comfort DRx™ Airwear®	GR	BR	XA	V	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
	Varilux Comfort Short™ Airwear®	GR	BR	XA		-10.00 to +6.00	+0.75 to +3.00	-4.00	14
	Varilux Comfort Short DRx™ Airwear®	GR	BR	XA	V	-10.00 to +6.00	+0.75 to +3.00	-6.00	14
	Varilux® Ellipse® Airwear®	GR	BR			-10.00 to +6.00	+0.75 to +3.00	-4.00	14
	Varilux Ellipse 360°™ Airwear®	GR	BR			-10.00 to +6.00	+0.75 to +3.00	-6.00	13
	Varilux Ipseo IV™ Airwear®	GR				-10.00 to +6.00	+0.75 to +3.50	-6.00	14
	Varilux Ipseo IV eyecode™ Airwear®	GR				-10.00 to +6.00	+0.75 to +3.50	-6.00	customized
	Varilux Physio™ Airwear®	GR	BR	XA		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-4.00	17
	Varilux Physio 360°™ Airwear®	GR	BR			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	17
	Varilux Physio Enhanced™Airwear®	GR	BR	XA		-10.00 to +6.00	+0.75 to +3.50 & +4.00 +0.75 to +3.00 for MFH below 17	-6.00	14
	Varilux Physio Enhanced eyecode™ Airwear®	GR	BR	XA		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
	Varilux Physio Enhanced FIT™ Airwear®	GR	BR	XA		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
	Varilux Physio Enhanced AZIO™ Airwear®	GR	BR	XA		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
	Varilux Physio Enhanced INDIA™ Airwear®	GR	BR	XA		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
	Varilux Physio DRx™ Airwear®	GR	BR	XA	V	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	17
	Varilux Physio™ Short Airwear®	GR	BR			-10.00 to +6.00	+0.75 to +3.00	-4.00	14
	Varilux Physio Short 360°™ Airwear®	GR	BR			-10.00 to +6.00	+0.75 to +3.00	-6.00	14
	Varilux Physio Short DRx™ Airwear®	GR	BR	XA	V	-10.00 to +6.00	+0.75 to +3.00 & +4.00	-6.00	14
	DEFINITY® Airwear®	GR	BR			-7.00 to +4.00	+1.00 to +3.00	-4.00	18
	DEFINITY SHORT™ Airwear®	GR	BR			-7.00 to +4.00	+1.00 to +3.00	-4.00	15
	Essilor Ideal® Advanced Airwear®	GR	BR	XA	V	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
	Essilor Ideal® Advanced Wrap Airwear®	GR	BR	XA		-4.00 to +3.00	+0.75 to +3.00	-4.00	customized
	Essilor Ideal® Airwear®	GR	BR	XA	V	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
	Essilor Ideal® Short Airwear®	GR	BR	XA	V	-10.00 to +6.00	+0.75 to +3.00	-6.00	14
	Essilor Accolade™ Airwear®	GR	BR	XA	V	-10.00 to +6.00	+0.75 to +3.50	-4.00	17
	Essilor Accolade™ Freedom™ Airwear®	GR	BR	XA	V	-10.00 to +6.00	+0.75 to +3.50	-4.00	15
	Ovation® Airwear®	GR	BR			-10.00 to +6.00	+0.75 to +3.50	-4.00	17
	Essilor Natural® Airwear®	GR				-10.00 to +6.00	+0.75 to +3.00	-4.00	18
	Essilor SmallFit™ Airwear®	GR	BR			-10.00 to +6.00	+0.75 to +3.00	-4.00	14
HOYA VISION CARE	HOYALUX Summit CD Poly Avail Spring 2013	GR	BR			-6.00 to +6.00	+1.00 to +3.00	-4.00	14
	HOYALUX Summit CD IQ Poly Avail Spring 2013	GR	BR			-6.00 to +6.00	+1.00 to +3.00	-4.00	14
	HOYALUX Summit ECP Polycarbonate	GR	BR			-6.00 to +6.00	+1.00 to +3.00	-4.00	18
	HOYALUX Summit ECP IQ Polycarbonate	GR	BR			-6.00 to +6.00	+1.00 to +3.00	-4.00	18
	HOYALUX iD LifeStyle Polycarbonate	GR				-10.00 to +6.00	+0.75 to 3.50	-4.00	18
	HOYALUX iD LifeStyle cd Polycarbonate	GR				-10.00 to +6.00	+0.75 to 3.50	-4.00	14
NIKON OPTICAL	Nikon Presio i Digital Polycarbonate	GR	BR	XA	V	-10.00 to +5.00	+0.75 to 3.50	-4.00	13, 15, 17
	Nikon Move Polycarbonate	GR	BR	XA		-10.00 to + 5.00	+0.75 to 3.50	-4.00	16, 18
OPTIMA INC.	Optima RESOLUTION® RESPONSE®	GR				-10.00 to +6.00	+1.00 to +3.00	-4.00	18
SEIKO OPTICAL PRODUCTS	1.59 Prominent	GR	BR	XA	V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
	1.59 Premier 2	GR	BR	XA	V	-7.00 to +5.00	+0.50 to +3.50	-4.00	14
	1.59 Prestige II	GR	BR	XA	V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
	1.59 Prestige II Wide	GR	BR	XA	V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
	1.59 Succeed	GR	BR	XA	V	-7.00 to +5.00	+0.50 to +3.50	-4.00	17
	1.59 Succeed Ws	GR	BR	XA	V	-7.00 to +5.00	+0.50 to +3.50	-4.00	15

*Available later this year. Check with your supplier for more details.

**MFH-Minimum Fitting Height

PROGRESSIVES 1.59

1.59 POLYCARBONATE PAL		Transitions® GRAY	Transitions® BROWN	Transitions® XTRActive® GRAY	Transitions® VANTAGE® GRAY	SPHERE	ADD	CYL TO	MFH**
	1.59 Supercede 2	GR	BR	XA	V	-7.00 to +5.00	+0.50 to +3.50	-4.00	14
	1.59 Supernal	GR	BR	XA	V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
	1.59 Surmount	GR	BR	XA	V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
	1.59 Surmount Ws	GR	BR	XA	V	-7.00 to +5.00	+0.50 to +3.50	-5.00	14
SHAMIR INSIGHT INC.	Shamir Autograph® II Polycarbonate	GR	BR	XA	V	-14.75 to +6.50	+0.75 to +3.50	-7.00	11, 13, 15, 18
	Shamir Autograph® II Attitude® Polycarbonate	GR	BR	XA	V	-6.50 to +5.25	+0.75 to +3.50	-7.00	15
	Shamir Autograph® II Attitude® Polycarbonate	GR	BR	XA	V	-6.50 to +6.00	+0.75 to +3.50	-7.00	18
	Shamir Spectrum™ PolyPlus™	GR	BR	XA		-14.75 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
	Shamir Element™ Polycarbonate	GR	BR	XA		-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19
	Shamir FirstPal Polycarbonate				V	-14.75 to +6.00	+0.75 to +1.50	-7.00	15
	Shamir Creation® PolyPlus™	GR				-12.00 to +7.50	+1.00 to +3.50	-6.00	19
	Shamir Piccolo® PolyPlus™	GR				-8.50 to +6.00	+1.00 to +3.00	-6.00	16
SIGNET ARMORLITE	KODAK Precise® Signia™ Poly	GR	BR			-10.00 to +5.00	+1.00 to +3.00	-4.00	17
	KODAK Precise® Short Signia™ Poly	GR	BR			-10.00 to +6.00	+1.00 to +3.00	-4.00	13
	KODAK Precise® Digital Signia™ Poly	GR	BR			-10.00 to +6.00	+1.00 to +3.00	-4.00	17
	KODAK Precise® Short Digital Signia™ Poly	GR	BR			-10.00 to +6.00	+1.00 to +3.00	-4.00	13
	KODAK Precise® PB Signia™ Poly	GR	BR		V	-10.00 to +5.00	+0.75 to +3.50	-4.00	17
	KODAK Precise® Short PB Signia™ Poly	GR	BR		V	-10.00 to +5.00	+0.75 to +3.50	-4.00	13
	KODAK Unique Signia™ Poly	GR	BR	XA	V	-10.00 to +5.00	+0.75 to +3.50	-4.00	13
	KODAK Unique DS Signia™ Poly	GR	BR			-10.00 to +5.00	+0.75 to +3.50	-4.00	13
YOUNGER OPTICS	Younger Image®	GR	BR	XA	V	-9.00 to +7.00	+1.00 to +3.00	-4.00	18
	Younger Adage®	GR	BR			-8.00 to +6.00	+1.00 to +3.00	-4.00	13

PROGRESSIVES 1.60

1.60 HIGH INDEX PAL		Transitions® GRAY	Transitions® BROWN	Transitions® XTRActive® GRAY	Transitions® VANTAGE® GRAY	SPHERE	ADD	CYL TO	MFH**
ESSILOR LENS	Varilux Comfort® Thin&Lite® 1.6	GR	BR			-12.00 to +8.00	+0.75 to +3.50	-4.00	17
	Varilux Comfort Enhanced™ Thin&Lite® 1.6	GR	BR			-12.00 to +8.00	+0.75 to +3.50 +075 to +3.00 for MFH below 17	-6.00	14
	Varilux Comfort DRx™ Thin&Lite® 1.6	GR	BR			-12.00 to +8.00	+0.75 to +3.50	-6.00	17
	Varilux Comfort Short™ Thin&Lite® 1.6	GR	BR			-12.00 to +8.00	+0.75 to +3.00	-4.00	14
	Varilux Comfort Short DRx™ Thin&Lite® 1.6	GR	BR			-12.00 to +8.00	+0.75 to +3.00	-6.00	14
	Varilux® Ellipse® Thin&Lite® 1.6	GR				-12.00 to +8.00	+0.75 to +3.00	-4.00	14
	Varilux Ellipse 360°™ Thin&Lite® 1.6	GR				-12.00 to +8.00	+0.75 to +3.00	-4.00	13
	Varilux Physio™ Thin&Lite® 1.6	GR	BR			-12.00 to +8.00	+0.75 to +3.50 & +4.00	-4.00	17
	Varilux Physio 360°™ Thin&Lite® 1.6	GR	BR			-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17
	Varilux Physio Enhanced™ Thin&Lite® 1.6	GR	BR			-12.00 to +8.00	+0.75 to +3.50 & +4.00 +0.75 to +3.00 for MFH below 17	-6.00	14
	Varilux Physio Enhanced eyecode™ Thin&Lite® 1.6	GR	BR			-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
	Varilux Physio Enhanced FIT™ Thin&Lite® 1.6	GR	BR			-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
	Varilux Physio Enhanced AZIO™ Thin&Lite® 1.6	GR	BR			-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
	Varilux Physio Enhanced INDIA™ Thin&Lite® 1.6	GR	BR			-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	14
	Varilux Physio DRx™ Thin&Lite® 1.6	GR	BR			-12.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17
	Varilux Physio Short™ Thin&Lite® 1.6	GR	BR			-12.00 to +8.00	+0.75 to +3.00	-4.00	14
	Varilux Physio Short DRx™ Thin&Lite® 1.6	GR	BR			-12.00 to +8.00	+0.75 to +3.00	-6.00	14
	Essilor Accolade™ Thin&Lite® 1.6	GR	BR			-12.00 to +6.00	+0.75 to +3.50	-4.00	17
	Essilor Accolade™Freedom™ Thin&Lite® 1.6	GR	BR			-12.00 to +6.00	+0.75 to +3.50	-4.00	15
	Essilor Ideal® Advanced Thin&Lite® 1.6	GR	BR			-10.00 to +6.00	+0.75 to +3.50	-6.00	17
	Essilor Ideal® Thin&Lite® 1.6	GR	BR			-10.00 to +6.00	+0.75 to +3.50	-6.00	17
	Essilor Ideal® Short Thin&Lite® 1.6	GR	BR			-10.00 to +6.00	+0.75 to +3.00	-6.00	14
HOYA VISION CARE	HOYALUX Summit CD 1.60	GR	BR			-13.00 to +6.00	+0.75 to +3.50	-4.00	14
	HOYALUX Summit CD IQ 1.60	GR	BR			-13.00 to +6.00	+0.75 to +3.50	-4.00	14
	HOYALUX Summit ECP 1.60	GR	BR			-13.00 to +8.00	+1.00 to +3.00,+3.50,+4.00	-4.00	18
	HOYALUX Summit ECP IQ 1.60	GR	BR			-13.00 to +8.00	+1.00 to +3.00,+3.50,+4.00	-4.00	18
	HOYALUX GP WIDE 1.60	GR				-10.00 to +6.00	+0.75 to +3.50	-4.00	18
	HOYALUX iD LifeStyle 1.60	GR	BR			-10.00 to +8.00	+0.75 to +3.50	-4.00	18
	HOYALUX iD LifeStyle cd 1.60	GR	BR			-10.00 to +8.00	+0.75 to +3.50	-4.00	14
NIKON OPTICAL	Nikon Presio i Short 1.60	GR	BR			-10.00 to +5.00	+0.75 to +3.50	-4.00	16
	Nikon Presio i 1.60	GR	BR			-10.00 to +5.00	+0.75 to +3.50	-4.00	18
	Nikon SeeMax Power 1.60	GR	BR			-10.00 to +5.00	+0.75 to +3.50	-4.00	15, 16, 17
	Nikon Presio Power 1.60	GR	BR			-10.00 to +8.00	+0.75 to +3.50	-4.00	15, 16, 17
	Nikon Presio i Digital 1.60	GR	BR			-10.00 to +5.00	+0.75 to +4.00	-6.00	13, 15, 17
	Nikon Move 1.60	GR	BR			-10.00 to +5.00	+0.75 to +3.50	-4.00	16, 18
SEIKO OPTICAL PRODUCTS	1.60 Prominent	GR	BR			-8.50 to +6.50	+0.50 to +3.50	-5.00	14
	1.60 Premier 2	GR	BR			-8.50 to +6.50	+0.50 to +3.50	-4.00	14
	1.60 Prestige II	GR	BR			-8.50 to +6.50	+0.50 to +3.50	-5.00	14
	1.60 Prestige II Wide	GR	BR			-8.50 to +6.50	+0.50 to +3.50	-5.00	14
	1.60 Succeed	GR	BR			-8.50 to +6.50	+0.50 to +3.50	-4.00	17
	1.60 Succeed Ws	GR	BR			-8.50 to +6.50	+0.50 to +3.50	-4.00	15
	1.60 Supercede 2	GR	BR			-8.50 to +6.50	+0.50 to +3.50	-4.00	14
	1.60 Supernal	GR	BR			-8.50 to +6.50	+0.50 to +3.50	-5.00	14
	1.60 Surmount	GR	BR			-8.50 to +6.50	+0.50 to +3.50	-5.00	14
	1.60 Surmount Ws	GR	BR			-8.50 to +6.50	+0.50 to +3.50	-5.00	14
SHAMIR INSIGHT INC.	Shamir Autograph® II SuperLite™ 1.6	GR	BR			-15.00 to +6.00	+0.75 to +3.50	-7.00	11, 13, 15, 18
	Shamir Autograph® II Attitude® SuperLite™ 1.6	GR	BR			-6.50 to +5.25	+0.75 to +3.50	-7.00	15
	Shamir Autograph® II Attitude® SuperLite™ 1.6	GR	BR			-6.50 to +6.00	+0.75 to +3.50	-7.00	18
	Shamir Spectrum™ SuperLite™ 1.6	GR	BR			-15.00 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18
	Shamir Element™ SuperLite™ 1.6	GR	BR			-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19
	Shamir FirstPal 1.6	GR	BR			-15.00 to +6.00	+0.75 to +1.50	-7.00	15
	Shamir Creation® SuperLite™ 1.6	GR				-12.00 to +8.00	+1.00 to +3.50	-6.00	19
	Shamir Genesis™ SuperLite™ 1.6	GR				-11.25 to +7.00	+1.00 to +3.50	-6.00	19
	Shamir Piccolo® SuperLite™ 1.6	GR				-11.25 to +7.00	+1.00 to +3.50	-6.00	16
SIGNET ARMORLITE	KODAK Unique	GR	BR			-10.00 to +6.00	+0.75 to +3.50	-4.00	13

PROGRESSIVES 1.67 CONTINUED

1.67 SUPER HIGH INDEX PAL		Transitions® GRAY	Transitions® BROWN	Transitions® XTRActive® GRAY	Transitions® VANTAGE® GRAY	SPHERE	ADD	CYL TO	MFH**
CARL ZEISS VISION	AO Easy	GR				-10.00 to +6.00	+0.75 to +3.50	-4.00	18
	Instinctive	GR				-10.00 to +6.00	+0.75 to +3.50	-4.00	17
	SOLA Compact Ultra™	GR				-10.00 to +6.00	+0.75 to +3.50	-6.00	13
	SOLAOne®	GR	BR			-12.00 to +8.00	+1.00 to +3.50	-4.00	18
	ZEISS GT2®	GR	BR			-16.00 to +6.00	+1.00 to +3.00	-4.00	17
	ZEISS GT2® Short	GR	BR			-16.00 to +6.00	+1.00 to +3.00	-4.00	13
	ZEISS Gradal®Top	GR				-16.00 to +6.00	+1.00 to +3.00	-4.00	18
	ZEISS Individual® 2	GR	BR	XA		-10.00 to +6.00	+0.75 to +3.00	-4.00	13-35
	ZEISS Individual® 2I	GR	BR	XA		-10.00 to +6.00	+0.75 to +3.00	-4.00	13-35
	ZEISS Individual® 2N	GR	BR	XA		-10.00 to +6.00	+0.75 to +3.00	-4.00	13-35
	ZEISS Individual®	GR	BR	XA		-10.00 to +6.00	+0.75 to +3.00	-4.00	13-35
	ZEISS Choice Plus V	GR	BR	XA		-12.00 to +8.00	+0.75 to +3.50	-4.00	13-35
	ZEISS Choice Plus 13	GR	BR	XA		-12.00 to +8.00	+0.75 to +3.50	-4.00	13
	ZEISS Choice Plus 15	GR	BR	XA		-12.00 to +8.00	+0.75 to +3.50	-4.00	15
	ZEISS Choice Plus 17	GR	BR	XA		-12.00 to +8.00	+0.75 to +3.50	-4.00	17
	ZEISS Choice Plus 19	GR	BR	XA		-12.00 to +8.00	+0.75 to +3.50	-4.00	19
	ZEISS Choice Plus 21	GR	BR	XA		-12.00 to +8.00	+0.75 to +3.50	-4.00	21
	ZEISS Choice 13	GR	BR	XA		-12.00 to +8.00	+0.75 to +3.50	-4.00	13
	ZEISS Choice 15	GR	BR	XA		-12.00 to +8.00	+0.75 to +3.50	-4.00	15
	ZEISS Choice 17	GR	BR	XA		-12.00 to +8.00	+0.75 to +3.50	-4.00	17
	ZEISS Choice 19	GR	BR	XA		-12.00 to +8.00	+0.75 to +3.50	-4.00	19
	ZEISS GT2® 3DV	GR	BR	XA		-10.00 to +6.00	+0.75 to +3.00	-4.00	13-35
	ZEISS GT2® 3D	GR	BR	XA		-10.00 to +6.00	+0.75 to +3.00	-4.00	17
	ZEISS GT2® 3D Short	GR	BR	XA		-10.00 to +6.00	+0.75 to +3.00	-4.00	13
	SOLA HDV	GR	BR	XA		-12.00 to +8.00	+0.75 to +3.50	-4.00	13-35
	SOLA Instinctive Performance 13	GR	BR	XA		-12.00 to +8.00	+0.75 to +3.50	-4.00	13
	SOLA Instinctive Performance 15	GR	BR	XA		-12.00 to +8.00	+0.75 to +3.50	-4.00	15
	SOLA Instinctive Performance 17	GR	BR	XA		-12.00 to +8.00	+0.75 to +3.50	-4.00	17
	SOLA Compact Ultra™ HD	GR	BR	XA		-12.00 to +8.00	+0.75 to +3.50	-4.00	13
	SOLAOne® HD	GR	BR	XA		-12.00 to +8.00	+0.75 to +3.50	-4.00	13
	AO Easy HD	GR	BR	XA		-12.00 to +8.00	+0.75 to +3.50	-4.00	18
ESSILOR LENS	Varilux S Design™ Thin&Lite® 1.67	GR	BR			-12.00 to +9.00	+0.75 to +4.00	-6.00	17
	Varilux S Design Short™ Thin&Lite® 1.67	GR	BR			-12.00 to +9.00	+0.75 to +3.50	-6.00	14
	Varilux S Fit™ Thin&Lite® 1.67	GR	BR			-12.00 to +9.00	+0.75 to +4.00	-6.00	14
	Varilux S 4D™ Thin&Lite® 1.67	GR	BR			-12.00 to +9.00	+0.75 to +4.00	-6.00	14
	Varilux Comfort® Thin&Lite® 1.67	GR	BR			-14.00 to +9.00	+0.75 to +3.50	-4.00	17
	Varilux Comfort Enhanced™ Thin&Lite® 1.67	GR	BR			-14.00 to +9.00	+0.75 to +3.50 +075 to +3.00 for MFH below 17	-6.00	14
	Varilux Comfort DRx™ Thin&Lite® 1.67	GR	BR			-12.00 to +9.00	+0.75 to +3.50	-6.00	17
	Varilux Comfort Short™ Thin&Lite® 1.67	GR	BR			-14.00 to +9.00	+0.75 to +3.00	-4.00	14
	Varilux Comfort Short DRx™ Thin&Lite® 1.67	GR	BR			-12.00 to +9.00	+0.75 to +3.00	-6.00	14
	Varilux® Ellipse® Thin&Lite® 1.67	GR				-12.00 to +8.00	+0.75 to +3.00	-4.00	14
	Varilux Ellipse 360°™ Thin&Lite® 1.67	GR				-12.00 to +8.00	+0.75 to +3.00	-6.00	13
	Varilux Ipseo IV™ Thin&Lite® 1.67	GR	BR			-12.00 to +9.00	+0.75 to +3.50	-6.00	14
	Varilux Ipseo IV eyecode™ Thin&Lite® 1.67	GR	BR			-12.00 to +9.00	+0.75 to +3.50	-6.00	customized
	Varilux Physio™ Thin&Lite® 1.67	GR	BR			-12.00 to +9.00	+0.75 to +3.50 & +4.00	-4.00	17
	Varilux Physio 360°™ Thin&Lite® 1.67	GR	BR			-12.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	17
	Varilux Physio Enhanced™ Thin&Lite® 1.67	GR	BR			-12.00 to +9.00	+0.75 to +3.50 & +4.00 +0.75 to +3.00 for MFH below 17	-6.00	14
	Varilux Physio Enhanced eyecode™ Thin&Lite® 1.67	GR	BR			-12.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
	Varilux Physio Enhanced FIT™ Thin&Lite® 1.67	GR	BR			-12.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
	Varilux Physio Enhanced AZIO™ Thin&Lite® 1.67	GR	BR			-12.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
	Varilux Physio Enhanced INDIA™ Thin&Lite® 1.67	GR	BR			-12.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
	Varilux Physio DRx™ Thin&Lite® 1.67	GR	BR			-12.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	17
	Varilux Physio Short™ Thin&Lite® 1.67	GR	BR			-12.00 to +9.00	+0.75 to +3.00	-4.00	14
	Varilux Physio Short 360°™ Thin&Lite® 1.67	GR	BR			-12.00 to +9.00	+0.75 to +3.00	-6.00	14
	Varilux Physio Short DRx™ Thin&Lite® 1.67	GR	BR			-12.00 to +9.00	+0.75 to +3.00	-6.00	14
	DEFINITY® Thin&Lite® 1.67	GR	BR			-12.00 to +10.00	+1.00 to +3.50 & +4.00	-4.00	18
	DEFINITY SHORT™ Thin&Lite® 1.67	GR	BR			-12.00 to +10.00	+1.00 to +3.50 & +4.00	-4.00	15
	Essilor Accolade™ Thin&Lite® 1.67	GR	BR			-12.00 to +8.00	+0.75 to +3.50	-4.00	17
	Essilor Accolade™ Freedom™ Thin&Lite® 1.67	GR	BR			-12.00 to +8.00	+0.75 to +3.50	-4.00	15
	Essilor Ideal® Advanced Thin&Lite® 1.67	GR	BR			-12.00 to +8.00	+0.75 to +3.50	-6.00	17
	Essilor Ideal® Advanced Wrap Thin&Lite® 1.67	GR	BR			-5.00 to +3.00	+0.75 to +3.00	-4.00	customized
	Essilor Ideal® Thin&Lite® 1.67	GR	BR			-12.00 to +8.00	+0.75 to +3.50	-6.00	17
	Essilor Ideal® Short Thin&Lite® 1.67	GR	BR			-12.00 to +8.00	+0.75 to +3.00	-6.00	14
	Ovation® Thin&Lite® 1.67	GR	BR			-12.00 to +8.00	+0.75 to +3.50	-4.00	17
	Essilor SmallFit™ Thin&Lite® 1.67	GR	BR			-12.00 to +8.00	+0.75 to +3.50	-4.00	14

**MFH-Minimum Fitting Height

PROGRESSIVES 1.67 CONTINUED

1.67 SUPER HIGH INDEX PAL		Transitions® GRAY	Transitions® BROWN	Transitions® XTRActive® GRAY	Transitions® VANTAGE® GRAY	SPHERE	ADD	CYL TO	MFH**
HOYA VISION CARE	HOYALUX SummitCD 1.67	GR	BR			-13.00 to +8.00	+1.00 to +3.00	-4.00	14
	HOYALUX Summit ECP 1.67	GR	BR			-13.00 to +8.00	+0.75 to +4.00	-4.00	18
	HOYALUX Summit ECP IQ 1.67	GR	BR			-13.00 to +8.00	+0.75 to +4.00	-4.00	18
	HOYALUX iD LifeStyle 1.67	GR	BR			-10.00 to +8.00	+0.75 to +3.50	-4.00	18
	HOYALUX iD LifeStyle cd 1.67	GR	BR			-10.00 to +8.00	+0.75 to +3.50	-4.00	14
NIKON OPTICAL	Nikon Presio i Short 1.67	GR				-15.00 to +8.00	+0.75 to +3.50	-4.00	16
	Nikon W Short 1.67	GR	BR			-9.00 to +5.00	+0.75 to +3.50	-4.00	15
	Nikon SeeMax Power 1.67	GR	BR			-10.00 to +5.00	+0.75 to +3.50	-4.00	15, 16, 17
	Nikon Presio Power 1.67	GR	BR			-10.00 to +8.00	+0.75 to +3.50	-4.00	15, 16, 17
	Nikon Presio i Digital 1.67	GR	BR			-12.00 to +8.00	+0.75 to +3.50	-6.00	13, 15, 17
	Nikon Move 1.67	GR	BR			-12.00 to +5.00	+0.75 to +3.50	-4.00	16, 18'
SEIKO OPTICAL PRODUCTS	1.67 Prominent	GR	BR	XA		-10.50 to +6.50	+0.50 to +3.50	-5.00	14
	1.67 Premier 2	GR	BR	XA		-10.50 to +6.50	+0.50 to +3.50	-4.00	14
	1.67 Prestige II	GR	BR	XA		-12.50 to +6.50	+0.50 to +3.50	-5.00	14
	1.67 Prestige II Wide	GR	BR	XA		-12.50 to +6.50	+0.50 to +3.50	-5.00	14
	1.67 Succeed	GR	BR	XA		-10.50 to +6.50	+0.50 to +3.50	-4.00	17
	1.67 Succeed Ws	GR	BR	XA		-10.50 to +6.50	+0.50 to +3.50	-4.00	15
	1.67 Supercede 2	GR	BR	XA		-10.50 to +6.50	+0.50 to +3.50	-4.00	14
	1.67 Proceed® III Super Short	GR	BR			-13.00 to +8.00	+0.75 to +3.25	-4.00	16
	1.67 Proceed®II Short	GR	BR			-13.00 to +8.00	+0.50 to +3.50	-4.00	18
	1.67 Supernal	GR	BR	XA		-12.50 to +6.50	+0.50 to +3.50	-5.00	14
	1.67 Surmount	GR	BR	XA		-12.50 to +6.50	+0.50 to +3.50	-5.00	14
	1.67 Surmount Ws	GR	BR	XA		-12.50 to +6.50	+0.50 to +3.50	-5.00	14
SHAMIR INSIGHT INC.	Shamir Autograph® II SuperLite™ 1.67	GR	BR	XA		-16.00 to +8.25	+0.75 to +3.50	-8.00	11, 13, 15, 18
	Shamir Autograph® II Attitude® SuperLite™ 1.67	GR	BR	XA		-6.50 to +5.25	+0.75 to +3.50	-8.00	15
	Shamir Autograph® II Attitude® SuperLite™ 1.67	GR	BR	XA		-6.50 to +5.25	+1.00 to +3.50	-8.00	18
	Shamir Spectrum™ SuperLite™ 1.67	GR	BR	XA		-16.00 to +8.00	+0.75 to +3.50	-8.00	14, 16, 18
	Shamir Element™ SuperLite™ 1.67	GR	BR	XA		-12.50 to +8.00	+0.75 to +3.50	-6.00	16, 19
	Shamir Creation® SuperLite™ 1.67	GR				-12.50 to +8.00	+1.00 to +3.50	-6.00	19
	Shamir Piccolo® SuperLite™ 1.67	GR				-12.50 to +6.00	+1.00 to +3.00	-6.00	16
SIGNET ARMORLITE	KODAK Precise® Signia™ 1.67	GR	BR			-13.00 to +8.00	+1.00 to +3.00	-4.00	17
	KODAK Precise® Short Signia™ 1.67	GR	BR			-12.50 to +7.00	+1.00 to +3.00	-4.00	13
	KODAK Precise® Digital Signia™ 1.67	GR	BR			-12.00 to +6.00	+0.75 to +3.50	-4.00	17
	KODAK Precise® Short Digital Signia™ 1.67	GR	BR			-12.00 to +6.00	+1.00 to +3.00	-4.00	13
	KODAK Precise® PB Signia™ 1.67	GR	BR			-12.00 to +6.00	+1.00 to +3.00	-4.00	17
	KODAK Precise® Short PB Signia™ 1.67	GR	BR			-12.00 to +6.00	+1.00 to +3.00	-4.00	13
	KODAK Unique Signia™ 1.67	GR	BR	XA		-12.00 to +6.00	+0.75 to +3.50	-4.00	13-18
	KODAK Unique DS Signia™ 1.67	GR	BR			-12.00 to +6.00	+0.75 to +3.50	-4.00	13-18
X-CEL OPTICAL	X-Cel High-X™ 1.67 Freedom ID™	GR				-12.50 to +6.00	+1.00 to +3.00	-4.00	17
YOUNGER OPTICS	Younger Image® 1.67 MR-10™	GR				-10.00 to +7.00	+1.00 to +3.50	-4.00	18

PROGRESSIVES 1.74


1.74 ULTRA HIGH INDEX PAL		Transitions® GRAY	Transitions® BROWN	Transitions® XTRActive® GRAY	Transitions® VANTAGE® GRAY	SPHERE	ADD	CYL TO	MFH**
ESSILOR LENS	Varilux Physio DRx™Thin&Lite® 1.74	GR				-20.00 to +9.00	+0.75 to +3.50, +4.00	-6.00	17
	Varilux Physio Short DRx™Thin&Lite® 1.74	GR				-20.00 to +9.00	+0.75 to +3.00	-6.00	14
NIKON OPTICAL	Nikon Presio i Digital 1.74	GR	BR			-15.00 to + 8.00	+0.75 to 3.50	-6.00	13, 15, 17
SIGNET ARMORLITE	KODAK Unique	GR	BR			-12.00 to +6.00	+0.75 to +3.50	-4.00	13

ADAPTIVE SUN LENSES FOR YOUR DRIVING NEEDS.

TRANSITIONS® DRIVEWEAR® SUN LENSES

DRIVE

For those who drive, Transitions® Drivewear® sun lenses provide a safer, more comfortable driving experience by combining polarization, color optimization and photochromic technology to enhance a driver's vision by adapting to and managing the ever-changing daytime light & weather conditions faced by drivers.



In low light or overcast conditions, the lenses are a green/yellow color that provide high contrast and minimize glare.



Behind the windshield, the lenses activate to a copper color enhancing color recognition and depth perception.



In bright outdoor light, the lenses activate to a dark red-brown filtering excess light to provide maximum comfort.

		COLOR	MATERIAL	SPHERE	ADD	CYL TO	MFH**
YOUNGER OPTICS	Younger **DriveWear**® FSV	Polarized, green/yellow to copper to reddish brown	1.50 std index	Plano 6 base			
	Younger **DriveWear**® SV	Polarized, green/yellow to copper to reddish brown	1.50 std index	-9.00 to +7.50		-4.00	
	Younger **DriveWear**® SFFT28	Polarized, green/yellow to copper to reddish brown	1.50 std index	-8.00 to +6.00	+1.00 to +3.00	-4.00	
	Younger **DriveWear**® Image®	Polarized, green/yellow to copper to reddish brown	1.50 std index	-8.00 to +6.00	+1.00 to +3.00	-4.00	18
	Younger **DriveWear**® SV	Polarized, green/yellow to copper to reddish brown	polycarbonate	-8.00 to +6.00		-4.00	
	Younger **DriveWear**® Image®	Polarized, green/yellow to copper to reddish brown	polycarbonate	-9.00 to +7.00	+1.00 to +3.00	-4.00	18
		For more information visit DriveWearLens.com					

RANGE OF PRODUCTS FOR YOUR RANGE OF PATIENTS.

GOLF · PEDAL · RANGE · RIDE · Run · TRAIL · WATER

The Transitions® line of adaptive sunglasses change depending on outdoor light, to enhance vision during outdoor sports and activities. They are color optimized and increase contrast, improve distance and depth perception.

www.Transitions.com/Sunwear

AUTUMN GOLD® TRANSITIONS® SUN LENSES

Specifically designed to enhance contrast and improve clarity of images in cloudy, foggy or very bright outdoor environments – perfect for hunters, target and skeet shooters and archers.

RANGE

YELLOW → AMBER

28% TINT — 82% TINT

Enhances vision in hazy or foggy early morning light conditions.

Makes greens more vibrant and enhances contrast of targets against the sky.

BELL® TRANSITIONS® FACE SHIELD

Motorcycle face shield that helps you take your riding experience to the next level. Unlike static clear or tinted shields, this shield automatically adjusts its level of darkness to optimize vision in sun and shade, day and night.

RIDE

CLEAR → DARK GREY

CLEAR — 82% TINT

Clear in low light or night conditions.

Reduces light intensity without altering the color of objects.

CALLAWAY® NEOX® TRANSITIONS® SUNGLASSES

Using technology specifically designed for golfers, Callaway® NEOX® Transitions® sunglasses counteract the problems encountered with typical non-prescription sunglasses, like distortion in both color and object sight and decreased depth of field, delivering increased depth perception, better distance vision and superior clarity.

GOLF

GREEN/YELLOW → DARK GREEN

54% TINT — 90% TINT

Provides contrast enhancement in low light conditions.

Improves depth perception, makes the contours of the ground more distinct and the ball easier to see.

**MFH-Minimum Fitting Height

DEFINITY FAIRWAY™ TRANSITIONS® SUN LENSES

This progressive sunwear product improves contrast to sharpen colors and contours on the golf course. It comes systematically coated with Crizal SunShield™ Mirrors. In a golf study, players preferred the DEFINITY progressive lens design 7:1 over a competitive PAL for superiority.

GOLF

LIGHT AMBER — DARK BROWN

58% TINT — 85% TINT

Enhances contrast and adds brightness.

Improves contrast and depth perception to help golfers see the contour of the greens, and better see the ball.

NIKE® MAX TRANSITIONS® SUNGLASSES

Engineered for athletes, Nike Max Optics delivers an advanced, high-precision solution to the world of sport. Together with Transitions technology, these lenses adapt to different athletic environments and changing light conditions.

GOLF TINT

GOLF

VIOLET — DARK PURPLE

50% TINT — 74% TINT

Improves contour recognition on the greens and increases ball pop.

Specifically tuned to the game of golf, enhances critical details on the fairway and green.

OUTDOOR TINT

TRAIL

GREEN — DARK GREEN

55% TINT — 85% TINT

Brightens the shadows and increases contrast.

Enhances the visual spectrum found in natural environments.

SPEED TINT

Run

RED — DARK RED

50% TINT — 85% TINT

Designed to eliminate road glare.

Allows calming red light to enter the eye which helps during endurance sports.

OAKLEY® TRANSITIONS® SUNGLASSES

GOLF TRAIL PEDAL

Oakley offers both prescription and non-prescription Oakley® Transitions® sun lenses and sunglasses. Available in multiple colors they all activate to darker, neutral tints to allow athletes to maintain peak performance.

CLEAR BLACK IRIDIUM

GREY — DARK GREY

46% TINT — 89% TINT

Sport performance lens for flat to bright light.

Neutral transmission for truer color perception across a wide variety of light conditions.

VR50

LIGHT BROWN — BROWN

50% TINT — 88% TINT

Sport performance lens for flat to bright light

Subtle tint that boosts contrast.

G40

ROSE — DARK ROSE

60% TINT — 91% TINT

Sport performance lens for medium to extremely bright light.

Muted rose tint for improved contrast and image resolution (visual clarity).

PERSIMMON

ORANGE — BROWN

70% TINT — 91% TINT

Sport performance lens for flat to bright light

Subtle tint that boosts contrast.

XPERIO® TRANSITIONS® SUN LENSES

TRAIL WATER

Combining polarization with photochromic technology makes these sun lenses ideal for water sports such as boating and fishing, or for other activities when spending time outdoors.

ASH GRAY

GRAY — DARK GRAY

61% TINT — 88% TINT

Transmits colors evenly.

Allows for true color recognition in bright light.

CARAMEL BROWN

BROWN — DARK BROWN

70% TINT — 90% TINT

Improves contrast and depth perception.

Makes green colors, like grass, pop.

ADAPTIVE SUNGLASSES

	COLOR	MATERIAL	SPHERE	ADD	CYL TO	MFH**
BELL TRANSITIONS® FACE SHIELD						
	Clear to dark gray					
	Available through Bell retailers. For more information or to purchase online, go to www.BellSports.com.					
CALLAWAY® EYEWEAR						
NEOX® Transitions® Adaptive Sunglasses	G22 Green to dark green	polycarbonate	Plano			
NEOX® Transitions® Adaptive Sun lenses SV	G22 Green to dark green	polycarbonate	Limited Rx in Autograph II Attitude SV			
NEOX® Transitions® Adaptive Sun lenses Progressives	G22 Green to dark green	polycarbonate	Limited Rx in Autograph II Attitude PAL			
	For more information, call 1-877-863-2759					
ESSILOR LENS						
DEFINITY FAIRWAY™ Transitions®	Amber to dark brown	polycarbonate	-7.00 to +4.00	+1.00 to +3.00	-4.00	18
DEFINITY FAIRWAY™ Transitions® SHORT	Amber to dark brown	polycarbonate	-7.00 to +4.00	+1.00 to +3.00	-4.00	15
	For more information visit DEFINITY.com					
KBco						
Xperio™ Transitions® SFSV	Polarized Ash Gray to dark gray	1.50 std index	-7.00 to +6.00		-4.00	
Xperio™ Transitions® SFSV	Polarized Caramel to dark brown	1.50 std index	-7.00 to +6.00		-4.00	
PRO **Xperio™** Transitions®	Polarized Ash Gray to dark gray	1.50 std index	-10.00 to +6.00	+1.50 to +3.50	-4.00	18
PRO **Xperio™** Transitions®	Polarized Caramel to dark brown	1.50 std index	-10.00 to +6.00	+1.50 to +3.50	-4.00	18
Xperio™ Transitions® SFSV	Polarized Ash Gray to dark gray	polycarbonate	-10.00 to +6.00		-3.00	
Xperio™ Transitions® SFSV	Polarized Caramel to dark brown	polycarbonate	available soon			
	For more information from KBco call 800-722-8776					
NIKE VISION						
Nike® Max Transitions® Golf Tint Adaptive Sunglasses	Violet to dark purple	polycarbonate	Plano	Rx available soon		
Nike® Max Transitions® Outdoor Tint Adaptive Sunglasses	Green to dark green	polycarbonate	Plano	Rx available soon		
Nike® Max Transitions® Speed Tint Adaptive Sunglasses	Red to dark red	polycarbonate	Plano	Rx available soon		
	For more information visit NikeVision.com					
OAKLEY®						
Oakley® Transitions® SV lenses	Gray and Brown	Oakley Plutonite®	-5.25 to +2.00		-3.00	
Oakley® Iridium® Transitions™ SV lenses	Black, Emerald and Gold Flash Tint	Oakley Plutonite®	-3.25 to +2.00		-3.00	
Oakley® Transitions® Adaptive Sunglasses	Clear to Black Iridium, VR50, G40, Persimmon	Oakley Plutonite®	Plano			
Oakley® Transitions® PAL lenses	Gray and Brown	Oakley Plutonite®	-5.25 to +2.00	+1.00 to +3.00	-3.00	
Oakley® Iridium® Transitions™ PAL lenses	Black, Emerald and Gold Flash Tint	Oakley Plutonite®	-3.25 to +2.00	+1.00 to +3.00	-3.00	
	For more information on plano options visit Oakley.com					
X-CEL OPTICAL						
Autumn Gold™ Transitions® SV Spherical	Yellow to amber	1.53 Trivex®	-9.50 to +6.50		-4.00	

**MFH-Minimum Fitting Height

Transitions

Transitions.com/PRO

www.facebook.com/transitionslenses

twitter.com/transitions

Transitions, the swirl and *XTRActive* are registered trademarks, *Life well lit, Transitions adaptive lenses, the eclipse, Transitions Healthy Sight for Life Fund,* and *Vantage,* are trademarks of Transitions Optical, Inc. ©2013 Transitions Optical, Inc.
Drivewear is a registered trademark of Younger Mfg. Co.

Photochromic performance and polarization are influenced by temperature, UV exposure and lens material. U13P004

03/13